


Sales Account Information	Distributor/Dealer Information
Name:	Order No.
Email:	Order Date:
Phone:	Company Name:



Chassis type	Front Wheel Drive
Max. user weight	136 kg
Speed	6, 10, 12.5 km/h
Power module	R-Net 120A with CJSM joystick
Wheels	Driving wheels 3.00 - 8 (14") Casters 2.80-2.50 - 4 (9")
Motor	320 Watts
Overall dimension (LxWxH)	1200 x 610 x 1120 mm
Seat width:	420 - 520 mm, 50 mm increment
Seat depth:	400 - 600 mm, 25 mm increment
Front seat height:	415 - 465 mm, 10 mm increment
Backrest height	510 - 670 mm
Armrest height	70 - 300 mm
Footrest to seat distance	400 - 600 mm

User's name:	
User's weight (kg):	
User's height (mm) :	
① Top of head (mm) :	
② Top of shoulders (mm) :	
③ Seat pan to elbow (mm) :	
④ Elbow to hand (mm) :	
⑤ Leg length (mm):	
⑥ Depth (mm):	
⑦ Width (mm):	
⑧ Chest width (mm):	
⑨ Chest height (mm) :	
⑩ Shoulder width (mm):	

Base Price: £21,575.00


Please make your selection with a check in the box.

◆ Standard feature ○ Alternative options (Free of charge)
+£ Additional charge option (only apply to purchase the whole chair)

** Purchase parts only, please refer to spare part price list*

◆ Standard feature ○ Alternative options (Free of charge) £ Additional charge	
Chassis: EVO (powered support wheels are included)	
Colour (battery covers side panels)	
<input type="checkbox"/> White	◆
<input type="checkbox"/> Orange	○
<input type="checkbox"/> Silver Black	○
Speed	
<input type="checkbox"/> 6 km/h * Lights are optional	◆
<input type="checkbox"/> 10 km/h * Lights required, please select light kit	○
<input type="checkbox"/> 12.5 km/h * Lights required, please select light kit	○

Batteries (lead acid)		<input type="checkbox"/>	
	<input type="checkbox"/> 2x 74 Ah Batteries		◆
	<i>Notes: Battery box dimension 263.5 x 177 x 243 mm (LxWxH)</i>		
Charger		<input type="checkbox"/> Charger 8A * Suggested for batteries 70 - 80 A	◆
Lights		<input type="checkbox"/> No lights	◆
	<input type="checkbox"/> Lights kit		£325
Driving Wheel and Caster			
Rim	<input type="checkbox"/> Black rims		◆
Colour	<input type="checkbox"/> Black tyre		◆
Front tyre (Driving wheel)	<input type="checkbox"/> Air filled		◆
	<input type="checkbox"/> PU tyre		£50
Rear tyre (Caster)	<input type="checkbox"/> Air filled		◆
	<input type="checkbox"/> PU tyre		£20
	<input type="checkbox"/> None		◆
	<input type="checkbox"/> Bumper wheels on rear		
Seat Mounting Height (mechanical pre set, seat pan to ground)			
	<input type="checkbox"/> 415 mm		◆
	<input type="checkbox"/> 425 mm		○
	<input type="checkbox"/> 435 mm		○
	<input type="checkbox"/> 445 mm		○
	<input type="checkbox"/> 455 mm		○
	<input type="checkbox"/> 465 mm		○
Seat Mounting (Lift & Tilt)			
	<input checked="" type="checkbox"/> Powered seat lift (400 mm lift)		◆
	<input checked="" type="checkbox"/> Powered seat tilt 45°		◆
Standing Function (multiple selection: one or more options available)			
	<input checked="" type="checkbox"/> Sit to stand		○
	<input checked="" type="checkbox"/> Lay to stand		◆
	<input type="checkbox"/> 45° lay to stand		○
	<input type="checkbox"/> Tilt to stand		○
Seat			
Seat width	<input type="checkbox"/> 420 mm		○
	<input type="checkbox"/> 470 mm		○
	<input type="checkbox"/> 520 mm		○
Seat depth	<input type="checkbox"/> 400 mm	<input type="checkbox"/> 425 mm	<input type="checkbox"/> 450 mm
	<input type="checkbox"/> 475 mm	<input type="checkbox"/> 500 mm	<input type="checkbox"/> 525 mm
	<input type="checkbox"/> 550 mm	<input type="checkbox"/> 575 mm	<input type="checkbox"/> 600 mm
			○
Seat Cushion			
Upholstery			
	<input type="checkbox"/> No cushion		(£80)
	<input type="checkbox"/> 3D mesh		◆
<i>* Seat cushion size is based on seat width and seat depth, which you selected above for seat.</i>			
Hip Supports			
	<input type="checkbox"/> None		◆
	<input type="checkbox"/> 3D mesh hip supports, set		£360
	<input type="checkbox"/> Dartex hip supports, set (water repellent)		£385

Backrest			
Panel size	<input type="checkbox"/>	Effective backrest width	Effective backrest panel height @low mounting
	<input type="checkbox"/>	410 mm	<input type="checkbox"/> 530mm <input type="checkbox"/> 600mm
	<input type="checkbox"/>	460 mm	<input type="checkbox"/> 530mm <input type="checkbox"/> 600mm
	<input type="checkbox"/>	510 mm	<input type="checkbox"/> 530mm <input type="checkbox"/> 600mm <input type="checkbox"/> 670mm
	<input type="checkbox"/>		<input type="radio"/>
	<input type="checkbox"/>		<input type="radio"/>
	<input type="checkbox"/>		<input type="radio"/>
Mounting	<input type="checkbox"/>	Low mounting	◆
	<input type="checkbox"/>	High mounting (+ 30 mm) (in comb. lateral supports)	<input type="radio"/>
Adjustment	<input checked="" type="checkbox"/>	Power backrest recline	◆
Backrest Cushion			
Upholstery			
	<input type="checkbox"/>	No cushion	(£55)
	<input type="checkbox"/>	3D mesh	◆
	<input type="checkbox"/>	Dartex (water repellent)	£45
* Backrest cushion size is based on the backrest width and backrest height which you selected above.			
Backrest Accessories			
Side Support wedges			
	<input type="checkbox"/>	None	◆
	<input type="checkbox"/>	Side wedges S (530 mm backrest)	£75
	<input type="checkbox"/>	Side wedges M (600 mm backrest)	£75
	<input type="checkbox"/>	Side wedges L (670 mm backrest)	£75
* Side support wedges can be combined with lateral supports if needed, but cannot be combined with side panels.			
Side Panels (deep contour backrest)			
	<input type="checkbox"/>	None	◆
	<input type="checkbox"/>	3D mesh deep contour side panels, set	£275
	<input type="checkbox"/>	Dartex deep contour side panels, set (water repellent)	£275
Lateral Supports			
	<input type="checkbox"/>	None	◆
	<input type="checkbox"/>	3D mesh lateral supports set	£420
	<input type="checkbox"/>	Dartex lateral supports, set (water repellent)	£450
* Side panels and lateral supports cannot be mounted on the backrest simultaneously. Select either of them if needed.			
Rail for Lateral Supports			
	<input type="checkbox"/>	None	◆
	<input type="checkbox"/>	Height adjustment rail for lateral support, set	£35
Lumbar Support			
	<input type="checkbox"/>	None	◆
	<input type="checkbox"/>	Foam set lumbar support	£50
Armrest			
Distance Between Armrests (armrest setting can be set larger than backrest size)			
	<input type="checkbox"/>	set for 410 mm backrest	<input type="radio"/>
	<input type="checkbox"/>	set for 460 mm backrest	<input type="radio"/>
	<input type="checkbox"/>	set for 510 mm backrest	<input type="radio"/>
Adjustment			
	<input type="checkbox"/>	Angle adjustable only	◆
	<input type="checkbox"/>	Angle adjustable and side swivel	£85
Upholstery			
	<input type="checkbox"/>	3D mesh	◆
	<input type="checkbox"/>	Dartex (water repellent)	£20
Armrest Length			
Armrest left	<input type="checkbox"/>	320mm	◆
	<input type="checkbox"/>	400mm	<input type="radio"/>
Armrest right	<input type="checkbox"/>	320mm	◆
	<input type="checkbox"/>	400mm	<input type="radio"/>

Armrest Support		
	<input type="checkbox"/> None	◆
	<input type="checkbox"/> Telescopic armrest support, set	£275
Upper Body Fixation		
Chest Support		
	<input type="checkbox"/> Chest belt with soft padding and clip	◆
	<input type="checkbox"/> 3D mesh adjustable chest support	£875
	<input type="checkbox"/> Dartex adjustable chest support (<i>water repellent</i>)	£895
Legrest, Central Mounted with Adjustable Biomechanical Length Compensation		
Footplate Size		
	<input type="checkbox"/> 150 mm width	○
	<input type="checkbox"/> 180 mm width	○
	<input type="checkbox"/> 1-piece footplate 392 x 240 mm (WxL)	○
Adjustment		
	<input checked="" type="checkbox"/> Power legrest adjustment	◆
Knee Supports		
	<input type="checkbox"/> 3D mesh knee supports	◆
	<input type="checkbox"/> Dartex knee supports (<i>water repellent</i>)	£20
Calf Supports		
	<input type="checkbox"/> None	◆
	<input type="checkbox"/> 3D mesh calf supports	£360
	<input type="checkbox"/> Dartex calf supports (<i>water repellent</i>)	£375
Legrest Accessories		
	<input type="checkbox"/> None	◆
	<input type="checkbox"/> Footplate hinge for shorter legrest (<i>footplate to cushion 250 ~ 330 mm</i>) <i>* for more details, contact your sales rep.</i>	
Controls (select side control or tray table control)		
Side Control		
Joystick	<input type="checkbox"/> CJSM	○
Mounting position	<input type="checkbox"/> Left side mount	○
	<input type="checkbox"/> Right side mount	○
Mounting Type	<input type="checkbox"/> Fixed mounting	◆
	<input type="checkbox"/> Parallel swing away	£225
Tray Table Control (<i>tray table control and chest support cannot be combined</i>)		
Mounting position	<input type="checkbox"/> Left side mount	○
	<input type="checkbox"/> Right side mount	○
Table size (with flip-up control)	<input type="checkbox"/> for Backrest 410mm, flip up CJSM	£1,625
	<input type="checkbox"/> for Backrest 460mm, flip up CJSM	£1,625
	<input type="checkbox"/> for Backrest 510mm, flip up CJSM	£1,625
		
Attendant Control		
	<input type="checkbox"/> Attendant control, lever, CJSM backrest mount, left	£1,125
	<input type="checkbox"/> Attendant control, lever, CJSM backrest mount, right	£1,125
Table		
Table General for Special Controls e.g. head control, chip control...etc		
Mounting position	<input type="checkbox"/> Left mount, no cut out	○
	<input type="checkbox"/> Right mount, no cut out	○
Table size	<input type="checkbox"/> for backrest 410mm	£510
	<input type="checkbox"/> for backrest 460mm	£510
	<input type="checkbox"/> for backrest 510mm	£510
Headrest		
Headrest Model		
	<input type="checkbox"/> None	◆
	<input type="checkbox"/> Standard LTS headrest 3D mesh	£235
	<input type="checkbox"/> Standard LTS headrest Dartex (<i>water repellent</i>)	£250
Headrest Accessories		
	<input type="checkbox"/> None	◆
	<input type="checkbox"/> Headrest rail (<i>out of centre positioning</i>)	£35

Pelvic Positioning Belts		
	<input type="checkbox"/> Pelvic positioning belt, clip	◆
	<input type="checkbox"/> Pelvic positioning belt, auto retract	£25
Accessories		
 Backpack hooks	<input type="checkbox"/> None	◆
	<input type="checkbox"/> Push handles, set	£70
	<input type="checkbox"/> Backpack hooks, set	£35