

karma



MID Lectus

The Ultimate in Comfort and Manoeuvrability


The MID chassis offers the best manoeuvrability from Karma's Topline range, featuring 6-wheel independent suspension for unparalleled comfort and stability. The world-renowned Lectus seating system features biomechanical positioning and access to an extensive range of options and accessories.



MID Chassis

6-Wheel Independent Suspension

MID Lectus' suspension system evenly distributes and absorbs shocks from driving over uneven surfaces, allowing you to comfortably overcome terrain that would stop most other mid-wheel drive chassis in their tracks.



Lectus Seat



Relaxation Made Easy

The biomechanical power recline and biomechanical power elevating leg rest can each reach 177° angles, allowing the user to lay perfectly flat and relax or perform a supine transfer.

Power Positioning Function

45° Tilt-in-Space
171° Recline
177° Elevating Legrest
300 mm Lift

SPECIFICATION

Product Model	Nominal seat width (mm)		Nominal seat depth (mm)		Seat height w/Cushion (mm)		Backrest height w/Cushion (mm)		Power Tilt, Seat Plane Angle (X°)		Power Tilt, Seat Plane Angle (X°)		Power seat lift range (mm)		Power Tilt, Seat Plane Angle (X°)		Backrest recline, seat to back angle (X°)		Legrest angle, leg to seat (X°)	
MID Lectus (MID-LTS)	420-520		400-600		580-880		615-625		5-49		180-260		30		180-260		77-171		95-177	
	Drive range (km)	Max. speed (km)	Obstacle climbing (mm)	Min. turning radius w/o user (mm)	Ability to climb rated slope (X°)	Battery (Ah)	Electrical system	User weight (kg)	F/M/R Wheel (in)	Overall LxWxH (mm)										
	36	12	50	730	12	80*2	R-net	136	8/14/6	Min. 1165x615x690 Max. 2115x715x1750										

• KARMA reserves the right to modify information herein without further notice • Unit: mm – 10 mm – 0.5kg



KARMA Mobility Limited

Unit 6, Target Park, Shawbank Road,
Redditch, B98 8YN, UK
TEL:0845-630-3436 FAX: 0845-630-3736
Email: Info@karmamobility.co.uk

COLOUR OPTIONS

