



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

CATALYST 5

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

[Submit Form](#)

Effective: 6/17/2024

Catalyst 5 £2150
 300lb weight capacity

Frame

Seat Width

14" [50SWT014]	N/C
15" [50SWT015]	N/C
16" [50SWT016]	N/C
17" [50SWT017]	N/C
18" [50SWT018]	N/C
19" [50SWT019]	N/C
20" [50SWT020]	N/C
21" [50SWT021]	£483
22" [50SWT022]	£483

Seat Depth

14" [50SDE014]	N/C
15" [50SDE015]	N/C
16" [50SDE016]	N/C
17" [50SDE017]	N/C
18" [50SDE018]	N/C
19" [50SDE019]	£91
20" [50SDE020]	£604

Transit Option

Transit System (Weight Limit 300lbs) [50TRNOPT]	N/C
<i>Not available with 350 lb Weight Capacity</i>	
Omit [50TRNOP2]	N/C

Titanium Option

Titanium Frame [50FRMTI] £846
Includes hangers and base frame. Does not include Cross Braces.

Weight Capacity

300 lbs [50STD001]	N/C
350 lbs (Heavy Duty) [50HVD001]	£483
<i>Requires minimum 6" caster height</i>	

Cross Brace

Standard Cross Braces [50XBR001]	N/C
Seating System Cross Braces [50XBR003]	N/C
<i>Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.</i>	
Titanium Cross Braces [50XBR002]	£363
<i>Not available with 21" and 22" Seat Widths</i>	

Seat Pan

Removable Aluminum Seat Pan [50ACC004] £236
1/2" Shorter than Depth

Rear Frame Type

No selection needed with Titanium Frame

Standard (Curved) [50RFT001]	N/C
<i>Select from Standard or Angle Adjustable Back Post on chart on the next page</i>	
Depth Adjustable [50RFT002]	£200
<i>Select from depth adjustable back posts in chart on next page. Provides 2" forward and 2" rearward adjustment with minimum Seat Depth of 14" and maximum Seat Depth of 22". Available only with T-Arm or Flip Up Armrests.</i>	
Reclining Frame [50RFT003]	£1150
<i>Select from Reclining Back Post in chart on the next page. Not available with Dynamic Back, O2 Holder, Transit Option, 350 lb Weight Capacity, Swing Away Armrests, Flip Back Armrest, Height Adjustable Flip Back T-Arm, Back Upholstery or IV Pole.</i>	

Backrest

Backrest

Select Back Height in row corresponding to desired Back Post and Height Range. Depth Adjustable not available with Titanium Frame.

Back Post	Price	Height Range	Back Height				
Standard [50BOP002]	N/C	Short ^{3,4}	8.5"	9.5"	10.5"	11.5"	
		Medium ⁴	12.5" ⁶	13.5"	14.5"	15.5"	
		Tall	16.5"	17.5"	18.5"	19.5"	
Angle Adjustable w/ Push Handle [50BOP004]	N/C	Short ⁴	14"	15"	16"		
		Medium ⁴	16"	17"	18"		
		Tall	18"	19"	20"		
Stroller Handle Back Post (Depth Adjustable) [50BOP032]	N/C	Short	20"				
		Tall	24"				
Height Adjustable Straight with Push Handle (Depth Adjustable) [50BOP033] ^{2,9}	£73	Short ¹	13"	14"	15"	16"	17"
		Medium	16"	17"	18"	19"	20"
		Tall	20" ⁸	21" ⁵	22" ⁵	23" ⁵	24" ⁵
Reclining w/ Adjustable Handle [50BOP028] ^{1,7}	N/C	Short	22"				
		Tall	26"				
Half-Folding [50BOP029] ^{2,8}	£83	Short	12.5" ^{6,9}	13.5"	14.5"	15.5"	
		Tall	16.5"	17.5"	18.5"	19.5"	

- 1 Not available with Dynamic Back
- 2 Not available with Standard Back Upholstery
- 3 Not available with Tension Adjustable Back Upholstery
- 4 Not available with Standard Back Upholstery in combination with 21" and 22" Seat Widths
- 5 Not available with Tension Adjustable Back Upholstery in combination with 14" and 15" Seat Widths
- 6 Not available with Tension Adjustable Back Upholstery in combination with 22" Seat Width
- 7 Not available with O2 Holder, Transit Option, Heavy Duty Option, Swing Away Armrests, Flip Back Armrest, Height Adjustable Flip Back T-Arm or IV Pole
- 8 Not available with Transit Option, Heavy Duty Option, Flip Back Armrest, Angle Adjustable Locking Extendable Flip Up Armrest or Fold Down Push Handle
- 9 Not available with 8° Bend with Push Handle Back Post
- 10 Not available with Angle Adjustable Locking Extendable Flip Up Armrest
- 11 Not available with Tension Adjustable Back Upholstery in combination with Angle Adjustable Locking Extendable Flip Up Armrest

Reclining Back Options

- Select options below if Reclining Back is selected on chart above. Recline angle range is 0° - 65°.
- Standard Spring [50GSP001] _____ N/C
 - Heavy Duty Spring [50GSP002] _____ N/C

Standard Back Post

- Select option below only if a Standard Back Post option is selected on chart above
- Straight with Push Handle [50BCK004] _____ N/C
 - Straight without Push Handle [50BCK008] _____ N/C
 - 8 Degree Bend with Push Handle [50BOP003] _____ £37
Not available with Short Back Height Range
 - Straight with Fold Down Push Handle [50BCK006] _____ £158
Not available with 350 lb Weight Capacity
 - Caregiver Adjustable Height [50BOP034] _____ £196
Not available with Swing Away Armrests

Push Handle Options

- Not available with Stroller Handle Back Post, Straight Back Post with Fold Down Push Handle, Reclining Back Post or 350 lb Weight Capacity
- Low Profile Push Handle [50BCK007] _____ £67
2" shorter than standard push handles

Depth Adjustable Back Post Options

- Options below may be selected only if a Depth Adjustable back post option is selected on the chart above
- Stroller Handle Extensions [50BOP012] _____ £73
Not available with Height Adjustable Straight with Push Handle Back Post
 - Dynamic Back (E2398) [50BOP014] _____ £1932
Available only with T-Arms

Back Upholstery Options

- Not available with Reclining or Stroller Handle Back Post. If selecting both upholstery and Axiom Back, the upholstery will be installed.
- Standard [50BUP001] _____ N/C
 - Tension Adjustable with Standard Cover [50BUP002] _____ N/C
 - Tension Adjustable with Air Mesh Breathable Cover [50BUP007] _____ £182
 - Omit [50BUP003] _____ N/C

Frame Type / Footrest

Frame Type	PRO ELR Covers
Swing Away [50FRM001] _____ N/C	Available only on PRO ELR Pro ELR Gel Pad Covers (Pair) [50ELRPAD] _____ £304
Super Low - Swing Away [50FRM002] _____ N/C	
70° Fixed Front [50FRM003] _____ N/C	
80° Fixed Front [50FRM004] _____ N/C	

Footrest Hangers

Frame	Hanger	60°	70°	80°	90°	Retail
Swing Away	Extension Tube Mount	[50HAN001] ²	[50HAN002]	[50HAN003] ¹		N/C
	Front Mount ²		[50HAN016]	[50HAN017] ¹	[50HAN004] ¹	N/C
	Pro Elevating Leg Rests ³					[50HAN044] £483
	Omit Hangers					[50HAN046] N/C
Fixed	Extension Tube Mount		[50FRM003]	[50FRM004] ¹		N/C

¹ Not available with 7" or 8" casters

² Not available with Titanium Frame

³ Available only with Classic Hanger Release. Not available in combination with 24" wheels and Push to Lock Flush Mount Wheel Locks.

Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hangers and Fixed Front Frame measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length												
Extension Tube Hangers (60°, 70°, 80°)	Short (11.5" - 15")	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"				
	Medium (13" - 17")	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"				
	Long (15.5" - 19.5")	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"				
Front Mount Hangers ¹	Short Front Mount Hanger (70°, 80°, 90°)	Standard Selections									Requires 3" Extension			
		4" ²	4.5" ²	5"	5.5"	6"	6.5"					10"	10.5"	11"
	Medium Front Mount Hanger (70°, 80°, 90°)	7"	7.5"	8"	8.5"	9"	9.5"					11.5"	12"	12.5"
		4" ²	4.5" ²	5"	5.5"	6"	6.5"	7"	7.5"			12"	12.5"	13"
	Long Front Mount Hanger (70°, 80°)	8"	8.5"	9"	9.5"	10"	10.5"	11"	11.5"			13.5"	14"	14.5"
		5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"	9"		14"	14.5"	15"
		9.5"	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15.5"	16"	16.5"	
PRO ELR Low Pivot Setting	Pro ELR Low Pivot Setting Short (12" - 17.5")	12 - 17.5												
	Pro ELR Low Pivot Setting Long (15" - 21")	15 - 21												
PRO ELR High Pivot Setting	Pro ELR High Pivot Setting Short (10.5" - 15.5")	10.5 - 15.5												
	Pro ELR High Pivot Setting Long (13" - 19")	13 - 19												
Fixed Front Frame (70°, 80°)	Short (11.5" - 15")	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"					
	Medium (13" - 17")	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"				
	Long (15.5" - 19.5")	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"				

¹ Not available with Titanium Frame. See Minimum Front Seat Height chart below

² Available with 90° only

3" Extension

Available only with Front Mount Hangers

3" Extension [50EXT008] _____	£37
3" Extension (Provided as Parts) [50EXT018] _____	£37

Footplates

Composite [50FPL001] _____	N/C
<i>Not available with Front Mount Hanger</i>	
Composite Angle Adjustable [50FPL014] _____	N/C
<i>Not available with Front Mount Hanger or 14" Seat Width</i>	
Aluminum Angle Adjustable [50FPL002] _____	£115
Aluminum Locking Multi-Angle Adjustable [50FPL020] _____	£290
<i>Available only with Front Mount Hanger</i>	
One Piece Flip Up Angle Adjustable [50FPL009] _____	£272
<i>Not available with Front Mount Hanger or PRO ELR</i>	

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13.5"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

Hanger Release

If no selection is made Classic Hanger Release will be used

Classic [50REL001] _____	N/C
4-Way - In-Line [50REL002] _____	N/C
<i>Not available with Pro ELR</i>	
4-Way - Out [50REL003] _____	N/C

Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

Heel Loops [50ACC008]	N/C
Adjustable Heel Loops [50ACC029]	£25
Omit Heel Loops [50ACC011]	N/C

Calf Straps

Velcro Adjustable [50ACC002]	N/C
Padded Velcro Adjustable [50ACC003]	£67
Padded Airmesh Calf Strap [50ACC010]	£103

Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
Left	Left Limb Support	£170	Omit Left Hanger [50RLS004]	Smooth [50RLS006]	Air Mesh [50RLS007]	5 x 8 [50RLS008]	7 x 10 [50RLS009]	9 x 10 [50RLS010]	9 x 14 [50RLS011]	11 x 10 [50RLS012]	11 x 14 [50RLS013]
Right	Right Limb Support	£170	Omit Right Hanger [50RLS005]	Smooth [50RLS020]	Air Mesh [50RLS021]	5 x 8 [50RLS014]	7 x 10 [50RLS015]	9 x 10 [50RLS016]	9 x 14 [50RLS017]	11 x 10 [50RLS018]	11 x 14 [50RLS019]

Front Seat Height and Casters

Front Seat Height

Select Front Seat Height in row corresponding to desired caster size.

Note: Selections in first three Front Seat Height columns use Super Low Front Frame. Super Low Front Frame not available with Titanium Frame or Fixed Front Frame.

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
4"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
6"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"		
7"	17"	17.5"		18"	18.5"	19"	19.5"	20"	20.5"	21"				
8"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"						

Castors

Caster Type		Price	Caster Size				
			4"	5"	6"	7"	8"
	.75" Rollerblade	N/C	[50CAS004]				
	.75" Lighted Rollerblade	£81	[50CAS005]				
	1" Poly	N/C	[50CAS006]	[50CAS010]	[50CAS014]	[50CAS019]	[50CAS020]
	1" Poly Aluminum	£96	[50CAS007]	[50CAS011]			
	1" Pneumatic	£103			[50CAS017]		[50CAS021]
	1.5" Poly	N/C	[50CAS032]	[50CAS025]	[50CAS026]		
	1.5" Soft Roll Aluminum	£96	[50CAS009]	[50CAS013]	[50CAS016]		
	2" Poly	N/C			[50CAS024]		
	2" Pneumatic w/Foam Insert	£242					[50CAS023]

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump Standard Caster Housing	1.5"	2"	2"	2"	2.5"	2.5"	2.5"
Maximum Dump Performance Caster Housing	3"	3"	3"	3.5"	3.5"	4"	4"

Depending on configuration, chairs with Fixed Front Frame Type may achieve up to an additional .75"

Fork Stems

Best Fit [50CST006]	N/C
<i>Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height</i>	
Frog Legs II Suspension System [50CST014]	£424
<i>Lightweight shock absorbing fork. Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.</i>	
Frog Legs Suspension Fork [50CST013]	£344
<i>Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.</i>	
Standard Fork Stem [50CST001]	N/C
3/4" Longer Fork Stem [50CST002]	N/C
1 1/2" Longer Fork Stem [50CST003]	N/C

Caster Housing

Standard Housing [50CSH001]	N/C
Performance Housing [50CSH002]	£37
<i>Not available with 350 lb Weight Capacity</i>	

Wheel Locks

Wheel locks

Push to Lock [50WLK002]	N/C
Push to Lock Flush Mount [50WLK001]	£67
Pull to Lock [50WLK003]	N/C
Low Profile Scissor Lock [50WLK012]	£67
Drum Brake [50WLK015]	£532
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>	
Dual - Drum Brake with Push to Lock [50WLK016]	£653
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>	
Dual - Drum Brake with Pull to Lock [50WLK017]	£653
<i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Reclining Backrest, Straight without Push Handle or Straight with Fold Down Push Handle Back Posts, or Low Profile Push Handle. Axle is Quick Release.</i>	

Fork / Stem Options

Frog Legs Options

Only for Frog Legs

User Weight Under 80 lbs [50CST004]	N/C
User Weight 80-120 lbs [50CST008]	N/C
User Weight 120-160 lbs [50CST009]	N/C
User Weight 160-220 lbs [50CST010]	N/C
User Weight 220-260 lbs [50CST011]	N/C
User Weight 260-300 lbs [50CST012]	N/C

Caster Pin Locks

Caster Pin Locks [50CST005]	£194
<i>Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork</i>	

Wheel Lock Options

Extension Handles [50WLK004]	£49
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	
Grade Aids [50WLK009]	£145
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>	

Rear Seat Height and Wheels

Rear Seat Height

Select Rear Seat Height in row corresponding to Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" wheel sizes.

Wheel Size	Tire Profile	Rear Seat Height									
20"	Low	11.5" ³	12" ¹	12.5" ³	13" ¹	13.5" ³	14"	14.5" ³	15"	15.5" ³	16" ³
	Full	12" ³	12.5" ¹	13" ³	13.5" ¹	14" ³	14.5"	15" ³	15.5"	16" ³	16.5" ³
22"	Low	12.5" ³	13" ¹	13.5" ³	14" ¹	14.5" ³	15"	16"	16.5" ³	17" ³	
	Full	13" ³	13.5" ¹	14" ³	14.5" ¹	15" ³	15.5"	16" ³	16.5"	17" ³	17.5" ³
24"	Low	13.5" ³	14" ¹	14.5" ³	15" ¹	15.5" ³	16"	16.5" ³	17"	17.5" ³	18" ³
	Full	14" ³	14.5" ¹	15" ³	15.5" ¹	16" ³	16.5"	17" ³	17.5"	18" ³	18.5" ³
25"	Low	14" ³	14.5" ¹	15" ³	15.5" ¹	16" ³	16.5"	17" ³	17.5"	18" ³	18.5" ³
26"	Low		15" ^{1,2}	15.5" ³	16" ¹	16.5" ³	17"	17.5" ³	18"	18.5" ³	19" ³
	Full	15" ^{2,3}	15.5" ^{1,2}	16" ³	16.5" ¹	17" ³	17.5"	18" ³	18.5"	19" ³	19.5" ³

¹ Not available with Amputee Axle Plate

² Not available with Height Adjustable Flip Back T-Arm

³ Requires Universal Axle Plate

Axle Plates

No selection needed with One Arm Drive

One Arm Drive

Universal [50AXP009] _____ N/C
 Amputee [50AXP003] _____ £170
*Not available with Swing Away Armrests or Wheelchair Drive Ready Option.
 Required with Reclining Backrest.*

One Arm Drive - Quick Release [50DRV007] _____ £846
*Available only with 20", 22" and 24" spoke wheels and Quick Release Axles.
 Spoke wheel style is exclusive to OAD. Not available with 4 degrees of
 camber. Standard handrim will be Aluminum Anodized.*

Drive Side

Left Side Drive [50DRV002] _____ N/C
 Right Side Drive [50DRV003] _____ N/C

Drive Handrim

Aluminum Drive Handrim [50DRV004] _____ N/C
 Poly Coated Drive Handrim [50DRV005] _____ £194

Wheels

No selection needed with One Arm Drive

		Rear Wheel Size					
		Price	20"	22"	24"	25"	26"
Rear Wheel Type	Maxx Spoke	N/C	[50WH1120]	[50WH1122]		[50WH1125]	[50WH1126]
	Maxx Spoke	N/C			[50WH1124]		
	Mag	N/C	[50WHE120]	[50WHE122]	[50WHE124]		
	Superlight Spoke ¹	N/C			[50WHE224]	[50WHE225] ²	
	Spinergy Spox	£768		[50WHE422]	[50WHE424]	[50WHE425]	[50WHE426]
	Spinergy LX	£816		[50WHE522]	[50WHE524]	[50WHE525]	[50WHE526]

¹ Not available with 350 lb Weight Capacity

² Available only with Superlight and Plastic Coated Handrims

Spinergy Options

Select spoke color for Spinergy wheels

Black Spokes [50SSC001] _____ N/C
 Red Spokes [50SSC002] _____ Blue Spokes [50SSC005] _____
 White Spokes [50SSC003] _____ £91

Rear Tire Type

Full Profile

Street Pneumatic 1 3/8" [50TIR001] _____ N/C
Not available on 25" wheel size
 Full Poly 1 3/8" [50TIR002] _____ N/C
Not available on Superlight wheel or 25" or 26" wheel sizes
 Pneumatic with Airless Insert 1 3/8" [50TIR004] _____ £110
Not available on Superlight wheel or 25" wheel size
 Schwalbe HS-127 1 3/8" [50TIR017] _____ £110
Available only on 24" wheel size. Not available on Mag wheel.
 Schwalbe Downtown 1 3/8" [50TIR020] _____ £110

Low Profile

Low Poly 1" [50TIR003] _____ £110
Not available on Superlight wheel, 25" or 26" wheel sizes or 350 lb Weight Cap.
 Sentinel Solid Tire 1" [50TIR005] _____ £110
Not available on Superlight wheels
 High Pressure 1" [50TIR006] _____ N/C
Not available on Mag wheel
 Schwalbe One 1" [50TIR018] _____ £110
Not available on Mag wheel, 20", 22" or 26" wheel sizes
 Schwalbe Marathon Plus Evolution 1" [50TIR010] _____ £110
Not available on Mag wheel, 20" or 22" wheel sizes

Handrims

Aluminum Anodized [50HRM001] _____ N/C
 Superlight Aluminum Anodized [50HRM010] _____ N/C
*Available only on 24" and 25" wheel sizes. Available only on MAXX Spoke,
 Superlight or Spinergy wheels. Not available with One Arm Drive.*
 Plastic Coated [50HRM002] _____ £110
 Projection (8 Vertical Projections) [50HRM003] _____ £285
*Available only with full profile tires in combination with 22", 24", and 26"
 MAXX Spoke wheels, or 22" and 24" Mag wheels. Not available with One Arm
 Drive.*
 Omit Handrim [50HRM012] _____ N/C

Ergonomic Handrims

Not available with 20" or 26" wheels

Ergo Lite Regular Grip [50HRM005] _____ £272
 Ergo Standard Regular Grip [50HRM008] _____ £272
 Ergo Lite Extra Grip [50HRM006] _____ £272
 Ergo Standard Extra Grip [50HRM009] _____ £272
 Ergo Lite Omit Grip [50HRM004] _____ £272
 Ergo Standard Omit Grip [50HRM007] _____ £272

Handrim Options

Handrim Mount

If no selection is made, Standard Mount will be used.

Close Mount [50TAB001] _____ N/C
Reduces overall width by .75"
 3M Micro Grip [50ACC006] _____ £91
Not available on Plastic Coated or Projection Handrims

Axles

No selection needed with One Arm Drive or Drum Brake

Quick Release [50AXO003]	N/C
Quad Quick Release [50AXO001]	£127

Not available with Mag or Spinergy wheels

Axle Position

Rear Wheel Camber

2° and 4° Not available with Heavy Duty option

0° [50WC00]	2° [50WC02]	4° [50WC04]	N/C
-------------	-------------	-------------	-----

Center of Gravity

Rear [50COG003]	Center [50COG002]	Forward [50COG001]	N/C
-----------------	-------------------	--------------------	-----

No selection needed with Amputee Axle Plate

Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15")	Desk [50ARM039]	N/C
	Full [50ARM040]	
Height Adjustable T-Arm (8.5"-13.5")	Desk [50ARM003]	£248
	Full [50ARM004]	
Height Adjustable T-Arm Low (6.5"-10.5") ¹	Desk [50ARM022]	£248
	Full [50ARM023]	
Height Adjustable Flip Back T-Arm Tall (11.5"-15") ²	Desk [50ARM041]	£230
	Full [50ARM042]	
Height Adjustable Flip Back T-Arm (9.5"-13.5") ²	Desk [50ARM029]	£230
	Full [50ARM030]	

1 If requested configuration is unavailable, Standard Height T-Arm will be used

2 Available only with Curved Rear Frame

3 Not available with Seat Depths less than 16" or Half-Folding Back

4 Available only with Standard Back Posts

5 Available only with Depth Adjustable and Reclining Back Posts

6 Not available with Arm Pads

7 Not Available with 15" Rear Seat Height and 24" wheels

Armrest Type	Length	Price
Flip back (10.5") ^{3,4,7}	Desk [50ARM009]	N/C
	Full [50ARM008]	
Angle Adjustable Locking Extendable Flip Up ⁵	Desk [50ARM010]	£272
	Full [50ARM011]	
Swing Away (10.5"-12.5") ^{4,6}	[50ARM002]	£248
Swing Away with Side Guard (10.5"-12.5") ^{4,6}	[50ARM018]	£248
Omit Armrests	[50ARM007]	N/C

Arm Pads

Standard Arm Pad [50ARM026]	N/C
Waterfall Arm Pad [50ARM020]	£25
<i>Available only with Height Adjustable T-Arm Tall, Height Adjustable T-Arm, and Height Adjustable T-Arm Low</i>	
Foam Grip [50ARM025]	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	

Accessories

Removable Side Guard

Not available with T-Arms, Flip Back or Swing Away with Side Guard Armrests

Adult Composite Side Guards [50SDG002]	N/C
Pediatric Composite Side Guards [50SDG001]	N/C
Aluminum Straight - Black Matte [50SDG024]	£212
Aluminum Straight - Match Frame Color [50SDG025]	£272
<i>Not available with Hydrographics or Titanium Frame</i>	
Aluminum Adjustable Fender w/Foam Pad - Black Matte [50SDG014]	£303
Aluminum Adjustable Fender w/Foam Pad - Match Frame Color [50SDG015]	£303
<i>Not available with Hydrographics or Titanium Frame</i>	
Carbon Fiber Adjustable Fender [50SDG013]	£393

Anti-tips

Rear Anti-Tips [50ATP001]	£91
Swing-In Anti-Tip Single - Left [50ATP007]	£105
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with left side O2 Holder.</i>	
Swing-In Anti-Tip Single - Right [50ATP008]	£105
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with right side O2 Holder.</i>	
Swing-In Anti-Tip Pair [50ATP009]	£213
<i>Not available with 350 lb Weight Capacity, Reclining Back Post or Depth Adjustable Rear Frame. Not available with O2 Holder.</i>	
Omit [50ATP000]	N/C

Pelvic Positioning Belts

1 1/2" Auto Buckle Pelvic Positioning Belt [50POB001]	£68
1 1/2" Padded Auto Buckle Pelvic Positioning Belt [50POB002]	£151
2" Airline Buckle Pelvic Positioning Belt [50POB003]	£110
Small Flex-Mount Positioning Belt (16-24") [50POB007]	£188
Medium Flex-Mount Positioning Belt (19-29") [50POB008]	£188
Large Flex-Mount Positioning Belt (25-37") [50POB009]	£188
Belt Mounting Band Clamp [50ACC054]	£103

O2 Holder

Not Available with Depth Adjustable Rear Frame

O2 Holder - Left Side [50ACC032]	£303
O2 Holder - Right Side [50ACC033]	£303

Other Accessories

Tool Kit [50ACC005]	£20
Removable Underseat Pouch [50ACC015]	£30
Pediatric Ki Mobility Backpack [50ACC026]	£100
Ki Mobility Backpack [50ACC009]	£100
Luggage Carrier [50ACC025]	£73
Black Neoprene Impact Guards [50IMG001]	£86
Cane and Crutch Holder [50ACC007]	£91
IV Pole [50ACC022]	£260

Color

Options	Price	Colors																							
		Gloss														Matte					Metallic				
		Black	Yellow	Light Pink	Dark Pink	Glacier Blue	Electric Green	Radiant Orange	Opalescent Violet	Purple Plum	Rich Teal	Metallic Orange	Granny Smith	Candy Red	Candy Blue	Candy Purple	Tactical Green	Dry Earth Tan	White	Shadow Black	Bright Red	Royal Blue	Charcoal Gray	Iridescent Navy	Iridescent Burgundy
Frame	N/C																								
Mag Color	N/C																								
Color Accent ¹ Package	N/C																								

¹ Includes axle plates and forks. Universal axle plate only. Not available with amputee axle plate, or one-arm drive.

Hydrographics Base Color

Not available with Frame Color. Hydrographics Base Color required with Hydrographics.

Options	Price	Hydrographics					Hydrographics Base Color							
		AOR1 (Desert Digital Camo)	NEXT Camo	Pink Camo	Carbon Fiber	Muddy Girl	Zebra	Black	Blue	Granny Smith	Pink	Purple	Red	
Frame	£725													
Base Color	N/C													

Colored Anodizing Package

Select desired color. If no selection is made, black will be used.

Color	Option	Standard ² Caster Housing	Caster Fork ³	Axle Plate ⁴
		Retail	£61	£61
Color	Black ¹			
	Blue			
	Bright Green			
	Orange			
	Pewter			
	Pink			
	Purple			
Red				

¹ Black is no charge

² Available only if corresponding option selection was made

³ Not available with 2" wide casters, Frog Legs Suspension Fork or Frog Legs II Suspension System

⁴ Color Anodized not available with amputee axle plate

Catalyst Front Seat Height Matrix																															
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster															
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork												
15	13.5	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	16.5	Std	6"	19	17.5	Std	6"	6"											
15.5	14			16.5	15.5	3/4		17.5	16.5	3/4		18.5	17			19.5	18.5	3/4													
15.5 ¹	14 ¹	17.5		16	1 1/2	18.5		17	1 1/2	18.5		17.5	20.5	n/a		1 1/2															
16	14.5	3/4		5"	16	14.5	Std	5"	17	15.5	Std	6"	19	18		3/4	7"	19	17.5		Std	7"									
16.5 ¹	15 ¹				16.5	15			17.5	16			18	16.5				19.5	18				20	18.5							
17	15.5	1 1/2			16.5	15.5	3/4		17.5	16.5			18	17		18.5		17.5	19.5				18	20.5	19	19.5	18.5	20	18.5		
15.5	14	Std			5"	17	16	3/4	5"	17.5	16.5	3/4	6"	18.5		17		Std	7"		20.5		19	Std	7"						
16	14.5					18	16.5			1 1/2	18.5			17.5		19.5					18		20			18.5					
16	15	3/4				6"	16.5	15		Std	6"			18.5		17		1 1/2			7"		19			18	3/4	7"	20.5	n/a	3/4
16.5	15.5						17	15.5	17.5			16	18.5	17		19.5							18	20.5		n/a			21	n/a	
17	15.5	1 1/2	6"				17.5	16	3/4	6"		18	16.5	Std	7"	20				19			1 1/2	7"		21	n/a		1 1/2	7"	
17.5	16						18	16.5			18.5	17	19.5			18		20		18.5	20.5					n/a					
16.5	15	Std					6"	17.5	16		3/4	6"	18			16.5		Std		7"	20		19			1 1/2	7"				20.5
17	15.5			18				16.5	18.5	17			19.5	18	20	18.5	20.5				19										
17	16	3/4		7"				18.5	17	1 1/2	7"		19	17.5	3/4	7"	20				18.5	1 1/2	7"			21			n/a		1 1/2
17.5	16.5							18.5	17			19.5	18	20			18.5	20.5		19											
18	16.5	1 1/2			7"			18.5	17	Std		7"	18.5	17.5			1 1/2	7"	21	n/a	1 1/2	7"			21	n/a			1 1/2		
18.5	17							19	17.5		19.5		18	20	18.5	20.5			19												
18	16.5	Std				7"		18	16.5	Std	7"		19.5	18.5	1 1/2	7"			20	19	1 1/2				7"	21		n/a			1 1/2
18.5	17							19	17.5			19.5	18	20			18.5	20.5	19												
18.5	17.5	3/4	7"					18.5	17.5	3/4		7"	19.5	18			1 1/2	7"	20	18.5	1 1/2			7"		21		n/a	1 1/2	7"	
19	18							18.5	17.5		19.5		18	20	18.5	20.5			19												
19.5	18.5	1 1/2					7"	19	18	3/4	7"		20.5	19	1 1/2	7"			21	n/a	1 1/2					7"	21	n/a			1 1/2
20	19							19.5	18.5			20.5	19	21			n/a														
				19.5				18	1 1/2	7"		6" x 2" caster not available with 5" fork									7"										
				20				18.5			1 1/2	7"					7"														
				20.5	19			7"										7"													

¹ Not available with 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.

Axiom G

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		£85	£85	£85	£85	£85	£85	£85	£85
Depth	14"	[XG1414]							
	15"	[XG1415]	[XG1515]	[XG1615]					
	16"	[XG1416]	[XG1516]	[XG1616]	[XG1716]	[XG1816]	[XG1916]	[XG2016]	
	17"	[XG1417]	[XG1517]	[XG1617]	[XG1717]	[XG1817]	[XG1917]	[XG2017]	
	18"	[XG1418]	[XG1518]	[XG1618]	[XG1718]	[XG1818]	[XG1918]	[XG2018]	[XG2218]
	19"		[XG1519]	[XG1619]	[XG1719]	[XG1819]	[XG1919]	[XG2019]	[XG2219]
	20"		[XG1520]	[XG1620]	[XG1720]	[XG1820]	[XG1920]	[XG2020]	[XG2220]
	21"			[XG1621]		[XG1821]	[XG1921]	[XG2021]	[XG2221]
22"			[XG1622]		[XG1822]	[XG1922]	[XG2022]	[XG2222]	

Axiom G - General Use Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		£125	£125	£125	£125	£125	£125	£125	£125
Depth	14"	[XG1414C2]							
	15"	[XG1415C2]	[XG1515C2]	[XG1615C2]					
	16"	[XG1416C2]	[XG1516C2]	[XG1616C2]	[XG1716C2]	[XG1816C2]	[XG1916C2]	[XG2016C2]	
	17"	[XG1417C2]	[XG1517C2]	[XG1617C2]	[XG1717C2]	[XG1817C2]	[XG1917C2]	[XG2017C2]	
	18"	[XG1418C2]	[XG1518C2]	[XG1618C2]	[XG1718C2]	[XG1818C2]	[XG1918C2]	[XG2018C2]	[XG2218C2]
	19"		[XG1519C2]	[XG1619C2]	[XG1719C2]	[XG1819C2]	[XG1919C2]	[XG2019C2]	[XG2219C2]
	20"		[XG1520C2]	[XG1620C2]	[XG1720C2]	[XG1820C2]	[XG1920C2]	[XG2020C2]	[XG2220C2]
	21"			[XG1621C2]		[XG1821C2]	[XG1921C2]	[XG2021C2]	[XG2221C2]
22"			[XG1622C2]		[XG1822C2]	[XG1922C2]	[XG2022C2]	[XG2222C2]	

Axiom S - Skin Protection Cushion

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
	Retail	£120	£120	£120	£120	£120	£120	£120	£120
Depth	14"	[XS1414]							
	15"	[XS1415]	[XS1515]	[XS1615]					
	16"	[XS1416]	[XS1516]	[XS1616]	[XS1716]	[XS1816]	[XS1916]	[XS2016]	
	17"	[XS1417]	[XS1517]	[XS1617]	[XS1717]	[XS1817]	[XS1917]	[XS2017]	
	18"	[XS1418]	[XS1518]	[XS1618]	[XS1718]	[XS1818]	[XS1918]	[XS2018]	[XS2218]
	19"		[XS1519]	[XS1619]	[XS1719]	[XS1819]	[XS1919]	[XS2019]	[XS2219]
	20"		[XS1520]	[XS1620]	[XS1720]	[XS1820]	[XS1920]	[XS2020]	[XS2220]
	21"			[XS1621]		[XS1821]	[XS1921]	[XS2021]	[XS2221]
22"			[XS1622]		[XS1822]	[XS1922]	[XS2022]	[XS2222]	

Axiom S - Skin Protection Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
	Retail	£160	£160	£160	£160	£160	£160	£160	£160
Depth	14"	[XS1414C2]							
	15"	[XS1415C2]	[XS1515C2]	[XS1615C2]					
	16"	[XS1416C2]	[XS1516C2]	[XS1616C2]	[XS1716C2]	[XS1816C2]	[XS1916C2]	[XS2016C2]	
	17"	[XS1417C2]	[XS1517C2]	[XS1617C2]	[XS1717C2]	[XS1817C2]	[XS1917C2]	[XS2017C2]	
	18"	[XS1418C2]	[XS1518C2]	[XS1618C2]	[XS1718C2]	[XS1818C2]	[XS1918C2]	[XS2018C2]	[XS2218C2]
	19"		[XS1519C2]	[XS1619C2]	[XS1719C2]	[XS1819C2]	[XS1919C2]	[XS2019C2]	[XS2219C2]
	20"		[XS1520C2]	[XS1620C2]	[XS1720C2]	[XS1820C2]	[XS1920C2]	[XS2020C2]	[XS2220C2]
	21"			[XS1621C2]		[XS1821C2]	[XS1921C2]	[XS2021C2]	[XS2221C2]
22"			[XS1622C2]		[XS1822C2]	[XS1922C2]	[XS2022C2]	[XS2222C2]	

Axiom P - Positioning Cushion

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
	Retail	£170	£170	£170	£170	£170	£170	£170	£170
Depth	14"	[XP1414]							
	15"	[XP1415]	[XP1515]	[XP1615]					
	16"	[XP1416]	[XP1516]	[XP1616]	[XP1716]	[XP1816]	[XP1916]	[XP2016]	
	17"	[XP1417]	[XP1517]	[XP1617]	[XP1717]	[XP1817]	[XP1917]	[XP2017]	
	18"	[XP1418]	[XP1518]	[XP1618]	[XP1718]	[XP1818]	[XP1918]	[XP2018]	[XP2218]
	19"		[XP1519]	[XP1619]	[XP1719]	[XP1819]	[XP1919]	[XP2019]	[XP2219]
	20"		[XP1520]	[XP1620]	[XP1720]	[XP1820]	[XP1920]	[XP2020]	[XP2220]
	21"			[XP1621]		[XP1821]	[XP1921]	[XP2021]	[XP2221]
22"			[XP1622]		[XP1822]	[XP1922]	[XP2022]	[XP2222]	

Axiom P - Positioning Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		£210	£210	£210	£210	£210	£210	£210	£210
Depth	14"	[XP1414C2]							
	15"	[XP1415C2]	[XP1515C2]	[XP1615C2]					
	16"	[XP1416C2]	[XP1516C2]	[XP1616C2]	[XP1716C2]	[XP1816C2]	[XP1916C2]	[XP2016C2]	
	17"	[XP1417C2]	[XP1517C2]	[XP1617C2]	[XP1717C2]	[XP1817C2]	[XP1917C2]	[XP2017C2]	
	18"	[XP1418C2]	[XP1518C2]	[XP1618C2]	[XP1718C2]	[XP1818C2]	[XP1918C2]	[XP2018C2]	[XP2218C2]
	19"		[XP1519C2]	[XP1619C2]	[XP1719C2]	[XP1819C2]	[XP1919C2]	[XP2019C2]	[XP2219C2]
	20"		[XP1520C2]	[XP1620C2]	[XP1720C2]	[XP1820C2]	[XP1920C2]	[XP2020C2]	[XP2220C2]
	21"			[XP1621C2]		[XP1821C2]	[XP1921C2]	[XP2021C2]	[XP2221C2]
22"			[XP1622C2]		[XP1822C2]	[XP1922C2]	[XP2022C2]	[XP2222C2]	

Axiom SP Visco - Skin Protection & Positioning Cushion

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		£350	£350	£350	£350	£350	£350	£350	£350	£350
Depth	14"	[XSP1414V]	[XSP1514V]	[XSP1614V]						
	15"	[XSP1415V]	[XSP1515V]	[XSP1615V]	[XSP1715V]	[XSP1815V]				
	16"	[XSP1416V]	[XSP1516V]	[XSP1616V]	[XSP1716V]	[XSP1816V]	[XSP1916V]	[XSP2016V]		[XSP2216V]
	17"	[XSP1417V]	[XSP1517V]	[XSP1617V]	[XSP1717V]	[XSP1817V]	[XSP1917V]	[XSP2017V]		[XSP2217V]
	18"	[XSP1418V]	[XSP1518V]	[XSP1618V]	[XSP1718V]	[XSP1818V]	[XSP1918V]	[XSP2018V]	[XSP2118V]	[XSP2218V]
	19"	[XSP1419V]	[XSP1519V]	[XSP1619V]	[XSP1719V]	[XSP1819V]	[XSP1919V]	[XSP2019V]	[XSP2119V]	[XSP2219V]
	20"	[XSP1420V]	[XSP1520V]	[XSP1620V]	[XSP1720V]	[XSP1820V]	[XSP1920V]	[XSP2020V]	[XSP2120V]	[XSP2220V]
	21"		[XSP1521V]	[XSP1621V]	[XSP1721V]	[XSP1821V]	[XSP1921V]	[XSP2021V]	[XSP2121V]	[XSP2221V]
22"		[XSP1522V]	[XSP1622V]	[XSP1722V]	[XSP1822V]	[XSP1922V]	[XSP2022V]	[XSP2122V]	[XSP2222V]	

Axiom SP Visco - Skin Protection & Positioning Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		£410	£410	£410	£410	£410	£410	£410	£410	£410
Depth	14"	[XSP1414VC2]	[XSP1514VC2]	[XSP1614VC2]						
	15"	[XSP1415VC2]	[XSP1515VC2]	[XSP1615VC2]	[XSP1715VC2]	[XSP1815VC2]				
	16"	[XSP1416VC2]	[XSP1516VC2]	[XSP1616VC2]	[XSP1716VC2]	[XSP1816VC2]	[XSP1916VC2]	[XSP2016VC2]		[XSP2216VC2]
	17"	[XSP1417VC2]	[XSP1517VC2]	[XSP1617VC2]	[XSP1717VC2]	[XSP1817VC2]	[XSP1917VC2]	[XSP2017VC2]		[XSP2217VC2]
	18"	[XSP1418VC2]	[XSP1518VC2]	[XSP1618VC2]	[XSP1718VC2]	[XSP1818VC2]	[XSP1918VC2]	[XSP2018VC2]	[XSP2118VC2]	[XSP2218VC2]
	19"	[XSP1419VC2]	[XSP1519VC2]	[XSP1619VC2]	[XSP1719VC2]	[XSP1819VC2]	[XSP1919VC2]	[XSP2019VC2]	[XSP2119VC2]	[XSP2219VC2]
	20"	[XSP1420VC2]	[XSP1520VC2]	[XSP1620VC2]	[XSP1720VC2]	[XSP1820VC2]	[XSP1920VC2]	[XSP2020VC2]	[XSP2120VC2]	[XSP2220VC2]
	21"		[XSP1521VC2]	[XSP1621VC2]	[XSP1721VC2]	[XSP1821VC2]	[XSP1921VC2]	[XSP2021VC2]	[XSP2121VC2]	[XSP2221VC2]
22"		[XSP1522VC2]	[XSP1622VC2]	[XSP1722VC2]	[XSP1822VC2]	[XSP1922VC2]	[XSP2022VC2]	[XSP2122VC2]	[XSP2222VC2]	

Axiom SP Fluid - Skin Protection & Positioning Cushion

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		£320	£320	£320	£320	£320	£320	£320	£320	£320
Depth	14"	[XSP1414F]	[XSP1514F]	[XSP1614F]						
	15"	[XSP1415F]	[XSP1515F]	[XSP1615F]	[XSP1715F]	[XSP1815F]				
	16"	[XSP1416F]	[XSP1516F]	[XSP1616F]	[XSP1716F]	[XSP1816F]	[XSP1916F]	[XSP2016F]		[XSP2216F]
	17"	[XSP1417F]	[XSP1517F]	[XSP1617F]	[XSP1717F]	[XSP1817F]	[XSP1917F]	[XSP2017F]		[XSP2217F]
	18"	[XSP1418F]	[XSP1518F]	[XSP1618F]	[XSP1718F]	[XSP1818F]	[XSP1918F]	[XSP2018F]	[XSP2118F]	[XSP2218F]
	19"	[XSP1419F]	[XSP1519F]	[XSP1619F]	[XSP1719F]	[XSP1819F]	[XSP1919F]	[XSP2019F]	[XSP2119F]	[XSP2219F]
	20"	[XSP1420F]	[XSP1520F]	[XSP1620F]	[XSP1720F]	[XSP1820F]	[XSP1920F]	[XSP2020F]	[XSP2120F]	[XSP2220F]
	21"		[XSP1521F]	[XSP1621F]	[XSP1721F]	[XSP1821F]	[XSP1921F]	[XSP2021F]	[XSP2121F]	[XSP2221F]
	22"		[XSP1522F]	[XSP1622F]	[XSP1722F]	[XSP1822F]	[XSP1922F]	[XSP2022F]	[XSP2122F]	[XSP2222F]


Axiom SP Fluid - Skin Protection & Positioning Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		£360	£360	£360	£360	£360	£360	£360	£360	£360
Depth	14"	[XSP1414FC2]	[XSP1514FC2]	[XSP1614FC2]						
	15"	[XSP1415FC2]	[XSP1515FC2]	[XSP1615FC2]	[XSP1715FC2]	[XSP1815FC2]				
	16"	[XSP1416FC2]	[XSP1516FC2]	[XSP1616FC2]	[XSP1716FC2]	[XSP1816FC2]	[XSP1916FC2]	[XSP2016FC2]		[XSP2216FC2]
	17"	[XSP1417FC2]	[XSP1517FC2]	[XSP1617FC2]	[XSP1717FC2]	[XSP1817FC2]	[XSP1917FC2]	[XSP2017FC2]		[XSP2217FC2]
	18"	[XSP1418FC2]	[XSP1518FC2]	[XSP1618FC2]	[XSP1718FC2]	[XSP1818FC2]	[XSP1918FC2]	[XSP2018FC2]	[XSP2118FC2]	[XSP2218FC2]
	19"	[XSP1419FC2]	[XSP1519FC2]	[XSP1619FC2]	[XSP1719FC2]	[XSP1819FC2]	[XSP1919FC2]	[XSP2019FC2]	[XSP2119FC2]	[XSP2219FC2]
	20"	[XSP1420FC2]	[XSP1520FC2]	[XSP1620FC2]	[XSP1720FC2]	[XSP1820FC2]	[XSP1920FC2]	[XSP2020FC2]	[XSP2120FC2]	[XSP2220FC2]
	21"		[XSP1521FC2]	[XSP1621FC2]	[XSP1721FC2]	[XSP1821FC2]	[XSP1921FC2]	[XSP2021FC2]	[XSP2121FC2]	[XSP2221FC2]
	22"		[XSP1522FC2]	[XSP1622FC2]	[XSP1722FC2]	[XSP1822FC2]	[XSP1922FC2]	[XSP2022FC2]	[XSP2122FC2]	[XSP2222FC2]

Axiom Back Width	Fits Back Frame Widths	
	Inner Slots	Outer Slots
12" (30.5cm)	11" (28.0cm)	12" -14" (30.5cm - 35.5cm)
14" (35.5cm)	12" -13" (30.5cm - 33.0cm)	14" -16" (35.5cm - 40.5cm)
16" (40.5cm)	14" -15" (35.5cm - 38.0cm)	16" -18" (40.5cm - 45.5cm)
18" (45.5cm)	16" - 17" (40.5cm - 43.0cm)	18" - 20" (45.5cm - 51.0cm)
20" (51.0cm)	18" -19" (45.5cm - 48.5cm)	20" -22" (51.0cm - 56.0cm)
22" (56.0cm)	20" -21" (51.0cm - 53.5cm)	22" -24" (56.0cm - 61.0cm)

Hardware Description	Code
Improved DEX Hardware, removable adjustable hardware (300 lb max user weight)	DEX
New BTX Hardware, removable adjustable hardware not recommended for, or installed on Ultralight Rigid Frame Wheelchairs (400 lb max user weight)	BTX
New Fixed Hardware, non-removeable hardware up to (300 lb max user weight)	FXT



Axiom Back

Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports.

Axiom Posterior Back (PB)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPB1208]	[AXPB1408]	[AXPB1608]	[AXPB1808]	[AXPB2008]	
		10"		[AXPB1410]	[AXPB1610]	[AXPB1810]	[AXPB2010]	[AXPB2210]
		12"	[AXPB1212]	[AXPB1412]	[AXPB1612]	[AXPB1812]	[AXPB2012]	
	Mid	14"		[AXPB1414]	[AXPB1614]	[AXPB1814]	[AXPB2014]	[AXPB2214]
		16"	[AXPB1216]	[AXPB1416]	[AXPB1616]	[AXPB1816]	[AXPB2016]	
		18"		[AXPB1418]	[AXPB1618]	[AXPB1818]	[AXPB2018]	[AXPB2218]

Back Hardware

FXT Hardware is not available with Lateral Support Options

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Lateral Support

Pair	Description	Item Code	Pair
	3 x 4 Lateral Pads with Fixed Mounting Hardware - Pair	PBLF3x4PP	£295
	4 x 5 Lateral Pads with Fixed Mounting Hardware - Pair	PBLF4x5PP	£295
	3 x 4 Lateral Pads with Swing Away Mounting Hardware - Pair	PBLS3X4PP	£628
	4 x 5 Lateral Pads with Swing Away Mounting Hardware - Pair	PBLS4X5PP	£628
	3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair	PBLAS3X4PP	£628
	4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair	PBLAS4X5PP	£628

Axiom Posterior Lateral Back (PL)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPL1208]	[AXPL1408]	[AXPL1608]	[AXPL1808]	[AXPL2008]	
		10"		[AXPL1410]	[AXPL1610]	[AXPL1810]	[AXPL2010]	[AXPL2210]
		12"	[AXPL1212]	[AXPL1412]	[AXPL1612]	[AXPL1812]	[AXPL2012]	
	Mid	14"		[AXPL1414]	[AXPL1614]	[AXPL1814]	[AXPL2014]	[AXPL2214]
		16"	[AXPL1216]	[AXPL1416]	[AXPL1616]	[AXPL1816]	[AXPL2016]	
		18"		[AXPL1418]	[AXPL1618]	[AXPL1818]	[AXPL2018]	[AXPL2218]
	Tail	20"	[AXPL1220]	[AXPL1420]	[AXPL1620]	[AXPL1820]	[AXPL2020]	
		22"		[AXPL1422]	[AXPL1622]	[AXPL1822]	[AXPL2022]	[AXPL2222]

Back Hardware

FXT Hardware is not available with Lateral Support Options

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Lateral Support

Pair	Description	Item Code	Pair
	3 x 4 Lateral Pads with Fixed Mounting Hardware - Pair	PLLF3x4PP	£295
	4 x 5 Lateral Pads with Fixed Mounting Hardware - Pair	PLLF4x5PP	£295
	3 x 4 Lateral Pads with Swing Away Mounting Hardware - Pair	PLLS3x4PP	£628
	4 x 5 Lateral Pads with Swing Away Mounting Hardware - Pair	PLLS4x5PP	£628
	3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair	PLLAS3x4PP	£628
	4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair	PLLAS4x5PP	£628

Axiom Deep Lateral Back (DL)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Mid	14"	[AXDL1214]	[AXDL1414]	[AXDL1614]	[AXDL1814]	[AXDL2014]	[AXDL2214]
		16"		[AXDL1416]	[AXDL1616]	[AXDL1816]	[AXDL2016]	
		18"	[AXDL1218]	[AXDL1418]	[AXDL1618]	[AXDL1818]	[AXDL2018]	[AXDL2218]
	Tail	20"	[AXDL1220]	[AXDL1420]	[AXDL1620]	[AXDL1820]	[AXDL2020]	
		22"		[AXDL1422]	[AXDL1622]	[AXDL1822]	[AXDL2022]	[AXDL2222]

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Headrest

Description	Std. Straight Removable Adjustable Hardware		Multi-Axis Removable Adjustable Hardware		Posalinc Removable Adjustable Hardware		Posalinc Flip-Back Adjustable Hardware	
	Item Code	Retail	Item Code	Retail	Item Code	Retail	Item Code	Retail
Comfort Soft Headrest 8"	[CS08SS]	£268	[CS08MA]	£268	[CS08PL]	£268	[CS08PLF]	N/C
Comfort Soft Headrest 10"	[CS10SS]	£268	[CS10MA]	£268	[CS10PL]	£268	[CS10PLF]	N/C
Comfort Soft Headrest 14"	[CS14SS]	£268	[CS14MA]	£268	[CS14PL]	£268	[CS14PLF]	N/C
Contour Positioning Pad, 8" X 6.25"	[CP08SS]	£303	[CP08MA]	£303	[CP08PL]	£303	[CP08PLF]	N/C
Contour Positioning Pad, 9" X 6.5"	[CP09SS]	£303	[CP09MA]	£303	[CP09PL]	£303	[CP09PLF]	N/C
Contour Positioning Pad, 10" X 7.5"	[CP10SS]	£303	[CP10MA]	£303	[CP10PL]	£303	[CP10PLF]	N/C
Tri-soft Headrest 10"	[TS10SS]	£386	[TS10MA]	£386	[TS10PL]	£386	[TS10PLF]	N/C
Tri-soft Headrest 14"	[TS14SS]	£386	[TS14MA]	£386	[TS14PL]	£386	[TS14PLF]	N/C

Description	Item Code	Retail
Universal Headrest Mounting Plate	[UMP]	£43

Not required for Axiom Headrests

Anterior Trunk Supports

Selection of Anterior Thoracic Mounting Modification required in order to mount posture supports. Not available with Low Backs.

Anterior Trunk Supports		
Description	Item Code	Retail
Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small	[14941]	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium	[14942]	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large	[14943]	£100
Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 18" Width, Small	[14924]	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium	[14927]	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large	[14928]	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small	[14851]	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Small	[14852]	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Medium	[14853]	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Large	[14854]	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small	[14856]	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small	[14857]	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14858]	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large	[14859]	£100
AirLogic Anterior Trunk Support: Standard, Structure, Small	[14282]	£300
AirLogic Anterior Trunk Support: Standard, Structure, Medium	[14283]	£300
AirLogic Anterior Trunk Support: Standard, Structure, Large	[14284]	£300
AirLogic Anterior Trunk Support: Standard, Structure, X-Large	[14285]	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small	[14257]	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14258]	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Large	[14259]	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, X-Large	[14260]	£300

Accessories

Adjustable Shoulder Guides [AXBSG] <i>Includes Anterior Thoracic Mounting Modification. Not available with Low Backs.</i>	£43
Lumbar Pad [LUMPAD] <i>Available only with Axiom Lateral Back and Axiom Deep Lateral back</i>	N/C
Privacy Flap [PRFLAP]	£44

Axiom Back Modifications

Axiom Back Modifications

Extended lead times will apply for all modifications. Makes back non-refundable.

Anterior Thoracic Mounting Modification [AXSMM] <i>Provides holes for Anterior Thoracic Supports and Shoulder Guides</i>	N/C
---	-----

Gloss Finish Back Shell

If no selection is made, standard Shadow Black will be used with no additional lead time

Opalescent Violet [AXMOVI]	N/C
Radiant Orange [AXMROR]	N/C
Electric Green [AXMEGN]	N/C
Glacier Blue [AXMGBL]	N/C
Dark Pink [AXMDPK]	N/C
Light Pink [AXMLPK]	N/C
Rich Teal [AXMRTE]	N/C
Purple Plum [AXMPLM]	N/C
Black [AXMBLK]	N/C
Candy Blue [AXMCBL]	N/C
Candy Purple [AXMPUR]	N/C
Candy Red [AXMRED]	N/C
Granny Smith [AXMGSM]	N/C
Metallic Orange [AXMMOR]	N/C
Yellow [AXMYLW]	N/C

Matte Finish Back Shell

If no selection is made, standard Shadow Black will be used with no additional lead time

Bright Red [AXMBRR]	N/C
White [AXMWHT]	N/C
Dry Earth Tan [AXMTAN]	N/C
Tactical Green [AXMTGR]	N/C
Charcoal Gray [AXMCGM]	N/C
Iridescent Burgundy [AXMIBM]	N/C
Iridescent Navy [AXMINM]	N/C
Royal Blue [AXMRBM]	N/C

Metallic Finish Back Shell

If no selection is made, standard Shadow Black will be used with no additional lead time

Blue Steel [AXMBLS]	N/C
---------------------	-----