



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

**Submit Form**

Effective: 6/17/2024

Rogue ALX	£1999
250 lb Weight Capacity	

## Frame

### Seat Width

See reference A on page 6

14" [R2SWT014]	N/C
15" [R2SWT015]	N/C
16" [R2SWT016]	N/C
17" [R2SWT017]	N/C
18" [R2SWT018]	N/C
19" [R2SWT019]	N/C
20" [R2SWT020]	N/C
21" [R2SWT021]	N/C
22" [R2SWT022]	N/C

### Seat Depth

See reference B on page 6

14" [R2SDE014]	N/C
15" [R2SDE015]	N/C
16" [R2SDE016]	N/C
17" [R2SDE017]	N/C
18" [R2SDE018]	N/C
19" [R2SDE019]	N/C
20" [R2SDE020]	N/C

### Weight Capacity

250 lbs [R2STD001]	N/C
350 lb (Heavy Duty) [R2HVD001]	£303
Requires 6" casters	

### Transit Option

Transit System [R2TRNOPT]	N/C
125 lb maximum weight capacity	
Omit [R2TRNOP2]	N/C

### Sling Position

In Performance Position (0"), sling is at front of frame. The 1" or 2" sling position extends the frame in front of the sling by that amount. See reference C on page 6.

0" (Performance) [R2SP01]	N/C
1" [R2SP03]	N/C
2" (Classic) [R2SP02]	N/C

### Adjustable Depth

Adds the capability to grow seat in depth

No Rear Adjustment [R2AD05]	N/C
1/2" [R2AD06]	N/C
1" [R2AD01]	N/C
2" [R2AD02]	£164
3" [R2AD03]	£164

### Long Sling

Available only with 1", 2" and 3" Adjustable Depth

Long Sling [R2AD04]	N/C
Sling will be maximum 2" longer than seat depth	

### Seat Option

Standard Seat Sling [R2SSS001]	N/C
Aluminum Solid Seat Pan [R2ACCO04]	£260
Seat Pan is 1" longer when Long Sling option is selected	

### Front Frame Bend

Angle between front seat tube and ground. See reference D on page 6.

90° [R2FA01]	N/C
85° [R2FA03]	N/C
80° [R2FA02]	N/C
75° [R2FA04]	£140
70° [R2FA05]	£140

## Backrest

Back Post / Height Range / Back Height							
Back Post	Price	Height Range	Back Height				
Height Adjustable	N/C	Height Adjustable Short Back Post [R2BOP025]	11" <sup>2</sup>	12"	13"	14"	
		Height Adjustable Medium Back Post [R2BOP026]	13"	14"	15"	16"	17"
		Height Adjustable Tall Back Post [R2BOP027]	17"	18"	19"	20"	
Fixed Height <sup>1</sup>	N/C	Fixed Height Straight Back Post [R2BOP011]	9" <sup>2</sup>	10" <sup>2</sup>	11" <sup>2</sup>	12" <sup>3</sup>	13"
			14"	15"	16"	17"	18"

- <sup>1</sup> Not available with back upholstery
- <sup>2</sup> Not available with Fold Down Push Handle
- <sup>3</sup> Fold Down Push Handle requires Height Adjustable Rigidizer Bar

Push Handles	
Omit Push Handle [R2BCK008]	N/C
Standard Push Handle [R2BCK004]	N/C
Low Profile Push Handle [R2BCK007]	£134
<i>2" shorter than standard push handles</i>	
Bolt-On Push Handle [R2BOP005]	£182
<i>Handles will be 1" - 2" above back height</i>	
Fold Down Push Handle [R2BCK006]	£158
<i>Not available with 350 lb Weight Capacity</i>	

Rigidizer Bar	
Adjustable Height Rigidizer Bar [R2BCK001]	N/C
Non-Adjustable Height Rigidizer Bar [R2BCK010]	N/C
<i>Not available with Fixed Height Back Post</i>	

Back Upholstery Options	
<i>Available only with Height Adjustable Back Posts. If selecting both upholstery and Axiom Back, the upholstery will be installed.</i>	
Padded Velcro Adjustable with Standard Cover [R2BUP001]	N/C
<i>Not available with 350 lb Weight Capacity</i>	
Padded Velcro Adjustable with Air Mesh Breathable Cover [R2BUP006]	£73
<i>Not available with 350 lb Weight Capacity</i>	
Tension Adjustable with Standard Cover [R2BUP002]	N/C
Tension Adjustable with Air Mesh Cover [R2BUP007]	£182
Omit [R2BUP003]	N/C

Locking Backrest	
Locking Fold Down Back [R2BCK011]	£43
<i>Not available with 350 lb Weight Capacity</i>	

Seat Back Angle	
<i>Standard is 90° (perpendicular to floor). Minus (-°) closes back angle. Plus (+°) opens back angle. -12° Requires Locking Fold Down Back.</i>	
Standard [R2SBA001]	+3° [R2SBA002]
+9° [R2SBA004]	-3° [R2SBA005]
-9° [R2SBA007]	-12° [R2SBA008]
	+6° [R2SBA003]
	-6° [R2SBA006]
	N/C

## Front Seat Height and Casters

Caster Size / Footrest Taper / Front Seat Height																
<i>Select Front Seat Height and Footrest Taper in row corresponding to Caster Size. 6" casters not available with 85° and 90° Front Frame Bend. Footrest Taper indicates bend of front frame creating footrest width. Taper is measured from outside of seat frame to outside of front tube. 2" Taper not available with 75° Front Frame Bend. Taper is not available with 70° Front Frame Bend. See reference E on page 6. Front Seat Height is measured from floor to top of seat tube at front of seat sling. Adjustable up or down by 1/2". See reference F on page 6.</i>																
Footrest Taper	Caster Size	Front Seat Height														
No Taper (Straight) [TXFT00]	3"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"
	4"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
	5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"		
	6"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"			
1" [TXFT01]	3"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"						
	4"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"							
	5"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"								
	6"	18"	18.5"	19"	19.5"	20"	20.5"									
2" [TXFT02]	3"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"						
	4"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"							
	5"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"								
	6"	18"	18.5"	19"	19.5"	20"	20.5"									

## Casters

		Caster Size				
		Price	3"	4"	5"	6"
Caster Type	.75" Rollerblade	N/C	[R2CAS002]	[R2CAS004]		
	.75" Lighted Rollerblade	£81	[R2CAS001]	[R2CAS005]		
	1" Poly	N/C		[R2CAS006]	[R2CAS010]	[R2CAS014]
	1" Poly Aluminum	£96		[R2CAS007]	[R2CAS011]	
	1" Pneumatic	£103				[R2CAS017]
	1.5" Poly	N/C		[R2CAS032]	[R2CAS025]	[R2CAS026]
	1.5" Soft Roll Aluminum	£96	[R2CAS003]	[R2CAS009]	[R2CAS013]	[R2CAS016]

## Fork Options

Standard Fork [R2CST001]	N/C
Slotted Fork [R2FRK001]	N/C
Single Sided Fork [R2FRK002]	£206
<i>Not available with 6" Caster Size</i>	
Frog Legs II Suspension System [R2CST014]	£424
<i>Lightweight shock absorbing fork. Not available with 350 lb Weight Capacity.</i>	

## Caster Pin Locks

Caster Pin Locks [R2CST005]	£194
<i>Not Available with Single Sided Fork</i>	

## Footrest

### Footplate

Tubular "Open" Footplate [R2FPL004]	N/C
<i>Not available w/footrest widths greater than 17.5". See Reference G on page 6.</i>	
Tubular Footplate with ABS Cover [R2FPL005]	£86
<i>Not available w/footrest widths greater than 15.5". See Reference G on page 6.</i>	
Aluminum Angle Adjustable [R2FPL006]	£182
Aluminum Angle Adjustable with Riser [R2FPL007]	£182
<i>Not available with 350 lb Weight Capacity</i>	
High Mount Aluminum Angle Adjustable [R2FPL008]	£182
<i>Not available with 350 lb Weight Capacity</i>	
Angle Adjustable Flip Under [R2FPL012]	£224
<i>Not available with footrest widths less than 9.5". Not available with 350 lb Weight Capacity</i>	
High Mount Angle Adjustable Flip Under [R2FPL011]	£272
<i>Not available with footrest widths greater than 15.5" or 350 lb Weight Capacity. See Reference G on page 6.</i>	

### Footrest Protector

Footrest Protector [R2ACC020]	N/C
<i>Not available with Angle Adjustable Flip Under, High Mount Angle Adjustable Flip Under or High Mount Aluminum Angle Adjustable Footplates</i>	

### Seat to Footrest Length

*Measured from front edge of seat sling to top, rear of footrest. Footrest length at least 2.5" shorter than Front Seat Height recommended. See reference H on page 6.*

8"	8.5"	9"	9.5"	10"	10.5"
11"	11.5"	12"	12.5"	13"	13.5"
14"	14.5"	15"	15.5"	16"	16.5"
17"	17.5"	18"	18.5"	19"	19.5"

## Rear Seat Height and Wheels

### Wheel Size / Rear Seat Height

*Select Rear Seat Height in row corresponding to desired Wheel Size. For proper setup, selection is required when omitting rear wheels. Measured from floor to top of seat tube at front of back post. Rear Seat Height is custom to needs. See reference I on page 6*

Wheel Size	Rear Seat Height													
22"	13.5" <sup>1</sup>	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
24"	14.5" <sup>1</sup>	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"		
25"	15" <sup>1</sup>	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"			
26"	15.5" <sup>1</sup>	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"				

<sup>1</sup> Available only with Low Profile tires

### Camber Bar Tower

Standard Tower [R2CM01]	N/C
Mx Tower [R2CM02]	N/C

*Standard Tower allows rear seat height adjustment up or down in 1/4" increments. 1" adjustment available; 1/2" in lowest end of range. Mx Tower required with select low seat heights or when large rear seat height adjustment range is required. Provides 3" adjustment range in 1/4" increments.*

## Rear Wheel

		Rear Wheel Size				
		Price	22"	24"	25"	26"
Rear Wheel Type	Spoke	N/C	[R2WHE022]	[R2WHE024]		
	Mag	N/C	[R2WHE122]	[R2WHE124]		
	Maxx Spoke	N/C	[R2WH1122]	[R2WH1124]	[R2WH1125]	[R2WH1126]
	Spinergy Spox	£768	[R2WHE422]	[R2WHE424]	[R2WHE425]	[R2WHE426]
	Spinergy LX	£816	[R2WHE522]	[R2WHE524]	[R2WHE525]	[R2WHE526]

## Spinergy Options

Select spoke color for Spinergy wheels

Black Spokes [R2SSC001]	N/C
Red Spokes [R2SSC002]	Blue Spokes [R2SSC005]
White Spokes [R2SSC003]	£91

## Handrims

### Standard Handrims

Aluminum Anodized [R2HRM001]	N/C
Plastic Coated [R2HRM002]	£121
Projection (8 Vertical Projections) [R2HRM003]	£333
<i>Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels</i>	

### Ergonomic Handrims

Not available with 20" or 26" wheels

Ergo Lite Regular Grip [R2HRM005]	£272
Ergo Standard Regular Grip [R2HRM008]	£272
Ergo Lite Extra Grip [R2HRM006]	£272
Ergo Standard Extra Grip [R2HRM009]	£272
Ergo Lite Omit Grip [R2HRM004]	£272
Ergo Standard Omit Grip [R2HRM007]	£272

## Handrim Options

### Handrim Mount

If no selection is made, Standard Mount will be used

Close Mount [R2TAB001]	N/C
<i>Reduces overall width by .75"</i>	

### Micro-Grip

3M Micro Grip [R2ACC006]	£91
<i>Not available on Plastic Coated or Projection Handrims</i>	

## Center of Gravity Presets

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post. See reference J on page 6.

Amputee Positions				
-1.5" [R2CGA150]	-1" [R2CGA100]	-.5" [R2CGA050]	0" [R2CGS000]	.25" [R2CGS025]
Standard Positions				
.5" [R2CGS050]	.75" [R2CGS075]	1" [R2CGS100]	1.25" [R2CGS125]	1.5" [R2CGS150]
1.75" [R2CGS175]	2" [R2CGS200]	2.25" [R2CGS225]	2.5" [R2CGS250]	2.75" [R2CGS275]
3" [R2CGS300]	3.25" [R2CGS325]	3.5" [R2CGS350]	3.75" [R2CGS375]	4" [R2CGS400]
4.25" [R2CGS425]				

## Rear Tire Type

Full Poly 1 3/8" [R2TIR002]	N/C
<i>Not available on 25" or 26" wheel sizes</i>	
Low Poly 1" [R2TIR003]	N/C
<i>Not available on 25" or 26" wheel sizes or 350 lb Weight Capacity</i>	
Pneumatic with Airless Insert 1 3/8" [R2TIR004]	£110
<i>Not available on 25" wheel size</i>	
Sentinel Solid tire 1" [R2TIR005]	£110
<i>Not available on Superlight Wheels</i>	
Street Pneumatic 1 3/8" [R2TIR001]	£110
<i>Not available on 25" wheel size</i>	
High Pressure 1" [R2TIR006]	N/C
<i>Not available on Mag wheel</i>	
Schwalbe Marathon Plus Evolution 1" [R2TIR010]	£110
<i>Not available on Mag wheel or 22" wheel size</i>	
Mountain 2" Black [R2TIR009]	£151
<i>Available only on 22", 24" and 25" MAXX Spoke and Spinergy wheels</i>	
Schwalbe HS-127 1 3/8" [R2TIR017]	£110
<i>Available only on 24" wheel size. Not available on Mag wheel.</i>	

## Axle

Quick Release [R2AXO003]	N/C
Quad Quick Release [R2AXO001]	£182
<i>Not available with Mag or Spinergy wheels</i>	

## Rear Wheel Spacing

Measured from outside of back to inside of rear tire. Adjustable out 1/2" from setting. 1.5" minimum Rear Wheel Spacing with Height Adjustable T-Arm Armrest. 1" minimum Rear Wheel Spacing with 0° Rear Wheel Camber. See reference K on page 6.

0.75" [R2WS075]	1" [R2WS100]	1.25" [R2WS125]	
1.5" [R2WS150]	1.75" [R2WS175]		N/C

## Rear Wheel Camber

0° [R2WC00]	2° [R2WC02]	4° [R2WC04]	6° [R2WC06]	
8° [R2WC08]				N/C

## Wheel Locks

### Wheel Lock

Push to Lock [R2WLK002]	N/C
Push to Lock Flush Mount [R2WLK001]	£67
Pull to Lock [R2WLK003]	N/C
Under Seat Scissor Lock [R2WLK007]	N/C
Low Profile Scissor Lock [R2WLK012]	N/C
<i>Not available with Full Poly Tires</i>	
Omit [R2WLK008]	N/C

### Wheel Lock Options

Extension Handles [R2WLK004]	£49
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	
Grade Aids [R2WLK009]	£145
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>	

## Armrests

Armrest Type	Length		Price
	Desk	[R2ARM039]	
Height Adjustable T-Arm Tall (11.5"-15") <sup>1</sup>	Full	[R2ARM040]	£212
	Desk	[R2ARM003]	
Height Adjustable T-Arm (8.5"-13.5") <sup>1</sup>	Full	[R2ARM004]	£212
	Desk	[R2ARM022]	
Height Adjustable T-Arm Low (6.5"-10.5") <sup>1</sup>	Desk	[R2ARM022]	£212
	Full	[R2ARM023]	

<sup>1</sup> 1.5" minimum Rear Wheel Spacing

<sup>2</sup> Includes direct mount receiver

Armrest Type		Price
Swing Away (8.5"-12.5")	[R2ARM002]	£140
Swing Away with Side Guard (8.5"-12.5") <sup>2</sup>	[R2ARM018]	£140
Omit Armrests	[R2ARM007]	N/C

## Side Guards

### Direct Mount Side Guards

*Not available with T-Arms or Swing Away with Side Guard Armrests*

Composite Side Guards - Direct Mount Receiver [R2SDG002]	N/C
Aluminum Straight - Black Matte [R2SDG024]	£212
Aluminum Straight - Match Frame Color [R2SDG025]	£272
Aluminum Adjustable Fender with Foam Pad - Black Matte [R2SDG014]	£296
<i>Not available with Adjustable Depth greater than 1/2"</i>	
Aluminum Adjustable Fender with Foam Pad - Match Frame Color [R2SDG015]	£357
<i>Not available with Adjustable Depth greater than 1/2"</i>	
Carbon Fiber Adjustable Fender - Direct Mount Receiver [R2SDG013]	£393
<i>Not available with Adjustable Depth greater than 1/2"</i>	

### Clamp-On Mount Side Guards

*Not available with T-Arms or Swing Away with Side Guard Armrests*

Composite Side Guards - Clamp-On Receiver [R2SDG004]	N/C
Aluminum Straight - Black Matte [R2SDG020]	£212
Aluminum Straight - Match Frame Color [R2SDG023]	£272
Aluminum Adjustable Fender with Foam Pad - Black Matte [R2SDG017]	£296
Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver [R2SDG018]	£357
Carbon Fiber Adjustable Fender - Clamp-On Receiver [R2SDG016]	£393

## Accessories

### Anti - Tips

Rear Anti-Tips [R2ATP001]	£91
User Activated - Pair [R2ATP012]	£320
User Activated - Left [R2ATP010]	£164
User Activated - Right [R2ATP011]	£164
Omit [R2ATP000]	N/C

### Calf Strap

Velcro Adjustable [R2ACC002]	N/C
Padded Velcro Adjustable [R2ACC003]	£67
Padded AirMesh Calf Strap [R2ACC010]	£103

### Pelvic Positioning Belt

1 1/2" Auto Buckle Pelvic Positioning Belt [R2POB001]	£55
1 1/2" Padded Auto Buckle Pelvic Positioning Belt [R2POB002]	£151
Small Flex-Mount Positioning Belt (16-24") [R2POB007]	£188
Medium Flex-Mount Positioning Belt (19-29") [R2POB008]	£188
Large Flex-Mount Positioning Belt (25-37") [R2POB009]	£188

### Other Accessories

Tool Kit [R2ACC005]	£20
Removable Underseat Pouch [R2ACC015]	£30
Ki Mobility Backpack [R2ACC009]	£100
Luggage Carrier [R2ACC025]	£73
Black Neoprene Impact Guards [R2IMG001]	£86
Under Chair Carrier and Removable Bag [R2ACC001]	£30
<i>Not available with scissor locks</i>	
Cane and Crutch Holder [R2ACC007]	£91

### Belt Mounting Kit

Ki Mobility Belt Mounting Clamp [R2ACC023]	£37
Bodypoint Belt Mounting Bracket Kit [R2ACC013]	£103

## Frame Color

Options	Price	Colors																							
		Gloss														Matte						Metallic			
		Black	Yellow	Light Pink	Dark Pink	Glacier Blue	Electric Green	Radiant Orange	Opalescent Violet	Purple Plum	Rich Teal	Metallic Orange	Granny Smith	Candy Red	Candy Blue	Candy Purple	Tactical Green	Dry Earth Tan	White	Shadow Black	Bright Red	Royal Blue	Charcoal Gray	Iridescent Navy	Iridescent Burgundy
Frame	N/C																								
Color Accent <sup>1,2</sup> Package	N/C																								

<sup>1</sup> Includes backplates, forks, and camber tube. Available only with Locking Fold Down Back and Slotted Forks. Not available with Heavy Duty or Transit.

<sup>2</sup> Limited time promotional pricing

## References

### A. Seat Width *See Diagram 1*

Measured from outside of frame tube on one side to outside of frame tube on other side.

### B. Seat Depth *See Diagram 3*

Measured from front of back posts to front edge of seat sling. Seat sling starts at beginning of bend at front of frame.

### C. Sling Position *See Diagram 2*

In Performance Position (0°), sling is at front of frame to keep chair shorter and more maneuverable. The 1" or 2" seat sling position extends the frame in front of the sling by that amount. Adding more frame can improve stability and provide support to aid in transfers.

### D. Front Frame Bend *See Diagram 3*

Angle between front seat tube and ground.

### E. Footrest Taper *See Diagram 1*

Indicates bend of front frame creating footrest. Taper is measured from outside of seat frame to outside of front tube.

### F. Front Seat Height *See Diagram 1*

Measured from floor to top of seat tube at front of seat sling.

### G. Footrest Width *See Diagram 1*

Measured from inside of front tube to inside of front tube on other side. Footrest width is listed below for each seat width and taper.

Seat Width	Straight	1" Taper	2" Taper
Outside Measurement	2" Narrower than Seat Width	4" Narrower than Seat Width	6" Narrower than Seat Width
14"	11.5"	9.5"	
15"	12.5"	10.5"	
16"	13.5"	11.5"	9.5"
17"	14.5"	12.5"	10.5"
18"	15.5"	13.5"	11.5"
19"	16.5"	14.5"	12.5"
20"	17.5"	15.5"	13.5"
21"	18.5"	16.5"	14.5"
22"	19.5"	17.5"	15.5"

### H. Seat to Footrest Length *See Diagram 3*

Measured from front edge of seat sling to top rear of footrest. Footrest length of at least 2.5" shorter than front seat height recommended.

### I. Rear Seat Height *See Diagram 3*

Measured from floor to top of seat tube at front of back post. Rear seat height is custom to needs.

### J. Center of Gravity Preset *See Diagram 3*

Measured from front of back post to center of rear axle. 0° indicates axle will be directly under back post.

### K. Rear Wheel Spacing *See Diagram 1*

Measured from outside of seat back to inside of rear tire. Adjustable out .5" from setting.

Diagram 1




Diagram 2





Diagram 3



## Axiom G

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		£85	£85	£85	£85	£85	£85	£85	£85
Depth	14"	[XG1414]							
	15"	[XG1415]	[XG1515]	[XG1615]					
	16"	[XG1416]	[XG1516]	[XG1616]	[XG1716]	[XG1816]	[XG1916]	[XG2016]	
	17"	[XG1417]	[XG1517]	[XG1617]	[XG1717]	[XG1817]	[XG1917]	[XG2017]	
	18"	[XG1418]	[XG1518]	[XG1618]	[XG1718]	[XG1818]	[XG1918]	[XG2018]	[XG2218]
	19"		[XG1519]	[XG1619]	[XG1719]	[XG1819]	[XG1919]	[XG2019]	[XG2219]
	20"		[XG1520]	[XG1620]	[XG1720]	[XG1820]	[XG1920]	[XG2020]	[XG2220]
	21"			[XG1621]		[XG1821]	[XG1921]	[XG2021]	[XG2221]
22"			[XG1622]		[XG1822]	[XG1922]	[XG2022]	[XG2222]	

## Axiom G - General Use Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		£125	£125	£125	£125	£125	£125	£125	£125
Depth	14"	[XG1414C2]							
	15"	[XG1415C2]	[XG1515C2]	[XG1615C2]					
	16"	[XG1416C2]	[XG1516C2]	[XG1616C2]	[XG1716C2]	[XG1816C2]	[XG1916C2]	[XG2016C2]	
	17"	[XG1417C2]	[XG1517C2]	[XG1617C2]	[XG1717C2]	[XG1817C2]	[XG1917C2]	[XG2017C2]	
	18"	[XG1418C2]	[XG1518C2]	[XG1618C2]	[XG1718C2]	[XG1818C2]	[XG1918C2]	[XG2018C2]	[XG2218C2]
	19"		[XG1519C2]	[XG1619C2]	[XG1719C2]	[XG1819C2]	[XG1919C2]	[XG2019C2]	[XG2219C2]
	20"		[XG1520C2]	[XG1620C2]	[XG1720C2]	[XG1820C2]	[XG1920C2]	[XG2020C2]	[XG2220C2]
	21"			[XG1621C2]		[XG1821C2]	[XG1921C2]	[XG2021C2]	[XG2221C2]
22"			[XG1622C2]		[XG1822C2]	[XG1922C2]	[XG2022C2]	[XG2222C2]	



## Axiom S - Skin Protection Cushion

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		£120	£120	£120	£120	£120	£120	£120	£120
Depth	14"	[XS1414]							
	15"	[XS1415]	[XS1515]	[XS1615]					
	16"	[XS1416]	[XS1516]	[XS1616]	[XS1716]	[XS1816]	[XS1916]	[XS2016]	
	17"	[XS1417]	[XS1517]	[XS1617]	[XS1717]	[XS1817]	[XS1917]	[XS2017]	
	18"	[XS1418]	[XS1518]	[XS1618]	[XS1718]	[XS1818]	[XS1918]	[XS2018]	[XS2218]
	19"		[XS1519]	[XS1619]	[XS1719]	[XS1819]	[XS1919]	[XS2019]	[XS2219]
	20"		[XS1520]	[XS1620]	[XS1720]	[XS1820]	[XS1920]	[XS2020]	[XS2220]
	21"			[XS1621]		[XS1821]	[XS1921]	[XS2021]	[XS2221]
22"			[XS1622]		[XS1822]	[XS1922]	[XS2022]	[XS2222]	

## Axiom S - Skin Protection Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		£160	£160	£160	£160	£160	£160	£160	£160
Depth	14"	[XS1414C2]							
	15"	[XS1415C2]	[XS1515C2]	[XS1615C2]					
	16"	[XS1416C2]	[XS1516C2]	[XS1616C2]	[XS1716C2]	[XS1816C2]	[XS1916C2]	[XS2016C2]	
	17"	[XS1417C2]	[XS1517C2]	[XS1617C2]	[XS1717C2]	[XS1817C2]	[XS1917C2]	[XS2017C2]	
	18"	[XS1418C2]	[XS1518C2]	[XS1618C2]	[XS1718C2]	[XS1818C2]	[XS1918C2]	[XS2018C2]	[XS2218C2]
	19"		[XS1519C2]	[XS1619C2]	[XS1719C2]	[XS1819C2]	[XS1919C2]	[XS2019C2]	[XS2219C2]
	20"		[XS1520C2]	[XS1620C2]	[XS1720C2]	[XS1820C2]	[XS1920C2]	[XS2020C2]	[XS2220C2]
	21"			[XS1621C2]		[XS1821C2]	[XS1921C2]	[XS2021C2]	[XS2221C2]
22"			[XS1622C2]		[XS1822C2]	[XS1922C2]	[XS2022C2]	[XS2222C2]	

## Axiom P - Positioning Cushion

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		£170	£170	£170	£170	£170	£170	£170	£170
Depth	14"	[XP1414]							
	15"	[XP1415]	[XP1515]	[XP1615]					
	16"	[XP1416]	[XP1516]	[XP1616]	[XP1716]	[XP1816]	[XP1916]	[XP2016]	
	17"	[XP1417]	[XP1517]	[XP1617]	[XP1717]	[XP1817]	[XP1917]	[XP2017]	
	18"	[XP1418]	[XP1518]	[XP1618]	[XP1718]	[XP1818]	[XP1918]	[XP2018]	[XP2218]
	19"		[XP1519]	[XP1619]	[XP1719]	[XP1819]	[XP1919]	[XP2019]	[XP2219]
	20"		[XP1520]	[XP1620]	[XP1720]	[XP1820]	[XP1920]	[XP2020]	[XP2220]
	21"			[XP1621]		[XP1821]	[XP1921]	[XP2021]	[XP2221]
22"			[XP1622]		[XP1822]	[XP1922]	[XP2022]	[XP2222]	

## Axiom P - Positioning Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		£210	£210	£210	£210	£210	£210	£210	£210
Depth	14"	[XP1414C2]							
	15"	[XP1415C2]	[XP1515C2]	[XP1615C2]					
	16"	[XP1416C2]	[XP1516C2]	[XP1616C2]	[XP1716C2]	[XP1816C2]	[XP1916C2]	[XP2016C2]	
	17"	[XP1417C2]	[XP1517C2]	[XP1617C2]	[XP1717C2]	[XP1817C2]	[XP1917C2]	[XP2017C2]	
	18"	[XP1418C2]	[XP1518C2]	[XP1618C2]	[XP1718C2]	[XP1818C2]	[XP1918C2]	[XP2018C2]	[XP2218C2]
	19"		[XP1519C2]	[XP1619C2]	[XP1719C2]	[XP1819C2]	[XP1919C2]	[XP2019C2]	[XP2219C2]
	20"		[XP1520C2]	[XP1620C2]	[XP1720C2]	[XP1820C2]	[XP1920C2]	[XP2020C2]	[XP2220C2]
	21"			[XP1621C2]		[XP1821C2]	[XP1921C2]	[XP2021C2]	[XP2221C2]
22"			[XP1622C2]		[XP1822C2]	[XP1922C2]	[XP2022C2]	[XP2222C2]	

## Axiom SP Visco - Skin Protection & Positioning Cushion

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		£350	£350	£350	£350	£350	£350	£350	£350	£350
Depth	14"	[XSP1414V]	[XSP1514V]	[XSP1614V]						
	15"	[XSP1415V]	[XSP1515V]	[XSP1615V]	[XSP1715V]	[XSP1815V]				
	16"	[XSP1416V]	[XSP1516V]	[XSP1616V]	[XSP1716V]	[XSP1816V]	[XSP1916V]	[XSP2016V]		[XSP2216V]
	17"	[XSP1417V]	[XSP1517V]	[XSP1617V]	[XSP1717V]	[XSP1817V]	[XSP1917V]	[XSP2017V]		[XSP2217V]
	18"	[XSP1418V]	[XSP1518V]	[XSP1618V]	[XSP1718V]	[XSP1818V]	[XSP1918V]	[XSP2018V]	[XSP2118V]	[XSP2218V]
	19"	[XSP1419V]	[XSP1519V]	[XSP1619V]	[XSP1719V]	[XSP1819V]	[XSP1919V]	[XSP2019V]	[XSP2119V]	[XSP2219V]
	20"	[XSP1420V]	[XSP1520V]	[XSP1620V]	[XSP1720V]	[XSP1820V]	[XSP1920V]	[XSP2020V]	[XSP2120V]	[XSP2220V]
	21"		[XSP1521V]	[XSP1621V]	[XSP1721V]	[XSP1821V]	[XSP1921V]	[XSP2021V]	[XSP2121V]	[XSP2221V]
22"		[XSP1522V]	[XSP1622V]	[XSP1722V]	[XSP1822V]	[XSP1922V]	[XSP2022V]	[XSP2122V]	[XSP2222V]	

## Axiom SP Visco - Skin Protection & Positioning Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		£410	£410	£410	£410	£410	£410	£410	£410	£410
Depth	14"	[XSP1414VC2]	[XSP1514VC2]	[XSP1614VC2]						
	15"	[XSP1415VC2]	[XSP1515VC2]	[XSP1615VC2]	[XSP1715VC2]	[XSP1815VC2]				
	16"	[XSP1416VC2]	[XSP1516VC2]	[XSP1616VC2]	[XSP1716VC2]	[XSP1816VC2]	[XSP1916VC2]	[XSP2016VC2]		[XSP2216VC2]
	17"	[XSP1417VC2]	[XSP1517VC2]	[XSP1617VC2]	[XSP1717VC2]	[XSP1817VC2]	[XSP1917VC2]	[XSP2017VC2]		[XSP2217VC2]
	18"	[XSP1418VC2]	[XSP1518VC2]	[XSP1618VC2]	[XSP1718VC2]	[XSP1818VC2]	[XSP1918VC2]	[XSP2018VC2]	[XSP2118VC2]	[XSP2218VC2]
	19"	[XSP1419VC2]	[XSP1519VC2]	[XSP1619VC2]	[XSP1719VC2]	[XSP1819VC2]	[XSP1919VC2]	[XSP2019VC2]	[XSP2119VC2]	[XSP2219VC2]
	20"	[XSP1420VC2]	[XSP1520VC2]	[XSP1620VC2]	[XSP1720VC2]	[XSP1820VC2]	[XSP1920VC2]	[XSP2020VC2]	[XSP2120VC2]	[XSP2220VC2]
	21"		[XSP1521VC2]	[XSP1621VC2]	[XSP1721VC2]	[XSP1821VC2]	[XSP1921VC2]	[XSP2021VC2]	[XSP2121VC2]	[XSP2221VC2]
22"		[XSP1522VC2]	[XSP1622VC2]	[XSP1722VC2]	[XSP1822VC2]	[XSP1922VC2]	[XSP2022VC2]	[XSP2122VC2]	[XSP2222VC2]	

## Axiom SP Fluid - Skin Protection & Positioning Cushion

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		£320	£320	£320	£320	£320	£320	£320	£320	£320
Depth	14"	[XSP1414F]	[XSP1514F]	[XSP1614F]						
	15"	[XSP1415F]	[XSP1515F]	[XSP1615F]	[XSP1715F]	[XSP1815F]				
	16"	[XSP1416F]	[XSP1516F]	[XSP1616F]	[XSP1716F]	[XSP1816F]	[XSP1916F]	[XSP2016F]		[XSP2216F]
	17"	[XSP1417F]	[XSP1517F]	[XSP1617F]	[XSP1717F]	[XSP1817F]	[XSP1917F]	[XSP2017F]		[XSP2217F]
	18"	[XSP1418F]	[XSP1518F]	[XSP1618F]	[XSP1718F]	[XSP1818F]	[XSP1918F]	[XSP2018F]	[XSP2118F]	[XSP2218F]
	19"	[XSP1419F]	[XSP1519F]	[XSP1619F]	[XSP1719F]	[XSP1819F]	[XSP1919F]	[XSP2019F]	[XSP2119F]	[XSP2219F]
	20"	[XSP1420F]	[XSP1520F]	[XSP1620F]	[XSP1720F]	[XSP1820F]	[XSP1920F]	[XSP2020F]	[XSP2120F]	[XSP2220F]
	21"		[XSP1521F]	[XSP1621F]	[XSP1721F]	[XSP1821F]	[XSP1921F]	[XSP2021F]	[XSP2121F]	[XSP2221F]
	22"		[XSP1522F]	[XSP1622F]	[XSP1722F]	[XSP1822F]	[XSP1922F]	[XSP2022F]	[XSP2122F]	[XSP2222F]

## Axiom SP Fluid - Skin Protection & Positioning Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		£360	£360	£360	£360	£360	£360	£360	£360	£360
Depth	14"	[XSP1414FC2]	[XSP1514FC2]	[XSP1614FC2]						
	15"	[XSP1415FC2]	[XSP1515FC2]	[XSP1615FC2]	[XSP1715FC2]	[XSP1815FC2]				
	16"	[XSP1416FC2]	[XSP1516FC2]	[XSP1616FC2]	[XSP1716FC2]	[XSP1816FC2]	[XSP1916FC2]	[XSP2016FC2]		[XSP2216FC2]
	17"	[XSP1417FC2]	[XSP1517FC2]	[XSP1617FC2]	[XSP1717FC2]	[XSP1817FC2]	[XSP1917FC2]	[XSP2017FC2]		[XSP2217FC2]
	18"	[XSP1418FC2]	[XSP1518FC2]	[XSP1618FC2]	[XSP1718FC2]	[XSP1818FC2]	[XSP1918FC2]	[XSP2018FC2]	[XSP2118FC2]	[XSP2218FC2]
	19"	[XSP1419FC2]	[XSP1519FC2]	[XSP1619FC2]	[XSP1719FC2]	[XSP1819FC2]	[XSP1919FC2]	[XSP2019FC2]	[XSP2119FC2]	[XSP2219FC2]
	20"	[XSP1420FC2]	[XSP1520FC2]	[XSP1620FC2]	[XSP1720FC2]	[XSP1820FC2]	[XSP1920FC2]	[XSP2020FC2]	[XSP2120FC2]	[XSP2220FC2]
	21"		[XSP1521FC2]	[XSP1621FC2]	[XSP1721FC2]	[XSP1821FC2]	[XSP1921FC2]	[XSP2021FC2]	[XSP2121FC2]	[XSP2221FC2]
	22"		[XSP1522FC2]	[XSP1622FC2]	[XSP1722FC2]	[XSP1822FC2]	[XSP1922FC2]	[XSP2022FC2]	[XSP2122FC2]	[XSP2222FC2]

Axiom Back Width	Fits Back Frame Widths	
	Inner Slots	Outer Slots
12" (30.5cm)	11" (28.0cm)	12" -14" (30.5cm - 35.5cm)
14" (35.5cm)	12" -13" (30.5cm - 33.0cm)	14" -16" (35.5cm - 40.5cm)
16" (40.5cm)	14" -15" (35.5cm - 38.0cm)	16" -18" (40.5cm - 45.5cm)
18" (45.5cm)	16" - 17" (40.5cm - 43.0cm)	18" - 20" (45.5cm - 51.0cm)
20" (51.0cm)	18" -19" (45.5cm - 48.5cm)	20" -22" (51.0cm - 56.0cm)
22" (56.0cm)	20" -21" (51.0cm - 53.5cm)	22" -24" (56.0cm - 61.0cm)

Hardware Description	Code
Improved DEX Hardware, removable adjustable hardware (300 lb max user weight)	<b>DEX</b>
New BTX Hardware, removable adjustable hardware not recommended for, or installed on Ultralight Rigid Frame Wheelchairs (400 lb max user weight)	<b>BTX</b>
New Fixed Hardware, non-removeable hardware up to (300 lb max user weight)	<b>FXT</b>



## Axiom Back

*Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports.*

### Axiom Posterior Back (PB)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPB1208]	[AXPB1408]	[AXPB1608]	[AXPB1808]	[AXPB2008]	
		10"		[AXPB1410]	[AXPB1610]	[AXPB1810]	[AXPB2010]	[AXPB2210]
		12"	[AXPB1212]	[AXPB1412]	[AXPB1612]	[AXPB1812]	[AXPB2012]	
	Mid	14"		[AXPB1414]	[AXPB1614]	[AXPB1814]	[AXPB2014]	[AXPB2214]
		16"	[AXPB1216]	[AXPB1416]	[AXPB1616]	[AXPB1816]	[AXPB2016]	
		18"		[AXPB1418]	[AXPB1618]	[AXPB1818]	[AXPB2018]	[AXPB2218]

### Back Hardware

*FXT Hardware is not available with Lateral Support Options*

DEX Hardware [DEX]    BTX Hardware [BTX]    FXT Hardware [FXT] \_\_\_\_\_ N/C

### Lateral Support

Pair	Description	Item Code	Pair
	3 x 4 Lateral Pads with Fixed Mounting Hardware - Pair	PBLF3x4PP	£295
	4 x 5 Lateral Pads with Fixed Mounting Hardware - Pair	PBLF4x5PP	£295
	3 x 4 Lateral Pads with Swing Away Mounting Hardware - Pair	PBLS3X4PP	£628
	4 x 5 Lateral Pads with Swing Away Mounting Hardware - Pair	PBLS4X5PP	£628
	3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair	PBLAS3X4PP	£628
	4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair	PBLAS4X5PP	£628

### Axiom Posterior Lateral Back (PL)

		Retail	Width					
			12"	14"	16"	18"	20"	22"
		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPL1208]	[AXPL1408]	[AXPL1608]	[AXPL1808]	[AXPL2008]	
		10"		[AXPL1410]	[AXPL1610]	[AXPL1810]	[AXPL2010]	[AXPL2210]
		12"	[AXPL1212]	[AXPL1412]	[AXPL1612]	[AXPL1812]	[AXPL2012]	
	Mid	14"		[AXPL1414]	[AXPL1614]	[AXPL1814]	[AXPL2014]	[AXPL2214]
		16"	[AXPL1216]	[AXPL1416]	[AXPL1616]	[AXPL1816]	[AXPL2016]	
		18"		[AXPL1418]	[AXPL1618]	[AXPL1818]	[AXPL2018]	[AXPL2218]
	Tail	20"	[AXPL1220]	[AXPL1420]	[AXPL1620]	[AXPL1820]	[AXPL2020]	
		22"		[AXPL1422]	[AXPL1622]	[AXPL1822]	[AXPL2022]	[AXPL2222]

### Back Hardware

FXT Hardware is not available with Lateral Support Options

DEX Hardware [DEX]    BTX Hardware [BTX]    FXT Hardware [FXT] \_\_\_\_\_ N/C

### Lateral Support

Pair	Description	Item Code	Pair
	3 x 4 Lateral Pads with Fixed Mounting Hardware - Pair	PLLF3x4PP	£295
	4 x 5 Lateral Pads with Fixed Mounting Hardware - Pair	PLLF4x5PP	£295
	3 x 4 Lateral Pads with Swing Away Mounting Hardware - Pair	PLLS3x4PP	£628
	4 x 5 Lateral Pads with Swing Away Mounting Hardware - Pair	PLLS4x5PP	£628
	3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair	PLLAS3x4PP	£628
	4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair	PLLAS4x5PP	£628

### Axiom Deep Lateral Back (DL)

		Retail	Width					
			12"	14"	16"	18"	20"	22"
		£550	£550	£550	£550	£550	£550	
Height	Mid	14"	[AXDL1214]	[AXDL1414]	[AXDL1614]	[AXDL1814]	[AXDL2014]	[AXDL2214]
		16"		[AXDL1416]	[AXDL1616]	[AXDL1816]	[AXDL2016]	
		18"	[AXDL1218]	[AXDL1418]	[AXDL1618]	[AXDL1818]	[AXDL2018]	[AXDL2218]
	Tail	20"	[AXDL1220]	[AXDL1420]	[AXDL1620]	[AXDL1820]	[AXDL2020]	
		22"		[AXDL1422]	[AXDL1622]	[AXDL1822]	[AXDL2022]	[AXDL2222]

### Back Hardware

DEX Hardware [DEX]    BTX Hardware [BTX]    FXT Hardware [FXT] \_\_\_\_\_ N/C

## Headrest

Description	Std. Straight Removable Adjustable Hardware		Multi-Axis Removable Adjustable Hardware		Posalinc Removable Adjustable Hardware		Posalinc Flip-Back Adjustable Hardware	
	Item Code	Retail	Item Code	Retail	Item Code	Retail	Item Code	Retail
Comfort Soft Headrest 8"	[CS08SS]	£268	[CS08MA]	£268	[CS08PL]	£268	[CS08PLF]	N/C
Comfort Soft Headrest 10"	[CS10SS]	£268	[CS10MA]	£268	[CS10PL]	£268	[CS10PLF]	N/C
Comfort Soft Headrest 14"	[CS14SS]	£268	[CS14MA]	£268	[CS14PL]	£268	[CS14PLF]	N/C
Contour Positioning Pad, 8" X 6.25"	[CP08SS]	£303	[CP08MA]	£303	[CP08PL]	£303	[CP08PLF]	N/C
Contour Positioning Pad, 9" X 6.5"	[CP09SS]	£303	[CP09MA]	£303	[CP09PL]	£303	[CP09PLF]	N/C
Contour Positioning Pad, 10" X 7.5"	[CP10SS]	£303	[CP10MA]	£303	[CP10PL]	£303	[CP10PLF]	N/C
Tri-soft Headrest 10"	[TS10SS]	£386	[TS10MA]	£386	[TS10PL]	£386	[TS10PLF]	N/C
Tri-soft Headrest 14"	[TS14SS]	£386	[TS14MA]	£386	[TS14PL]	£386	[TS14PLF]	N/C

Description	Item Code	Retail
Universal Headrest Mounting Plate	[UMP]	£43

*Not required for Axiom Headrests*

## Anterior Trunk Supports

*Selection of Anterior Thoracic Mounting Modification required in order to mount posture supports. Not available with Low Backs.*

Anterior Trunk Supports		
Description	Item Code	Retail
Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small	[14941]	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium	[14942]	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large	[14943]	£100
Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 18" Width, Small	[14924]	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium	[14927]	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large	[14928]	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small	[14851]	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Small	[14852]	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Medium	[14853]	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Large	[14854]	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small	[14856]	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small	[14857]	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14858]	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large	[14859]	£100
AirLogic Anterior Trunk Support: Standard, Structure, Small	[14282]	£300
AirLogic Anterior Trunk Support: Standard, Structure, Medium	[14283]	£300
AirLogic Anterior Trunk Support: Standard, Structure, Large	[14284]	£300
AirLogic Anterior Trunk Support: Standard, Structure, X-Large	[14285]	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small	[14257]	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14258]	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Large	[14259]	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, X-Large	[14260]	£300

## Accessories

Adjustable Shoulder Guides [AXBSG] <i>Includes Anterior Thoracic Mounting Modification. Not available with Low Backs.</i>	£43
Lumbar Pad [LUMPAD] <i>Available only with Axiom Lateral Back and Axiom Deep Lateral back</i>	N/C
Privacy Flap [PRFLAP]	£44

## Axiom Back Modifications

### Axiom Back Modifications

*Extended lead times will apply for all modifications. Makes back non-refundable.*

Anterior Thoracic Mounting Modification [AXSMM] <i>Provides holes for Anterior Thoracic Supports and Shoulder Guides</i>	N/C
---	-----

### Gloss Finish Back Shell

*If no selection is made, standard Shadow Black will be used with no additional lead time*

Opalescent Violet [AXMOVI]	N/C
Radiant Orange [AXMROR]	N/C
Electric Green [AXMEGN]	N/C
Glacier Blue [AXMGBL]	N/C
Dark Pink [AXMDPK]	N/C
Light Pink [AXMLPK]	N/C
Rich Teal [AXMRTE]	N/C
Purple Plum [AXMPLM]	N/C
Black [AXMBLK]	N/C
Candy Blue [AXMCBL]	N/C
Candy Purple [AXMPUR]	N/C
Candy Red [AXMRED]	N/C
Granny Smith [AXMGSM]	N/C
Metallic Orange [AXMMOR]	N/C
Yellow [AXMYLW]	N/C

### Matte Finish Back Shell

*If no selection is made, standard Shadow Black will be used with no additional lead time*

Bright Red [AXMBRR]	N/C
White [AXMWHT]	N/C
Dry Earth Tan [AXMTAN]	N/C
Tactical Green [AXMTGR]	N/C
Charcoal Gray [AXMCGM]	N/C
Iridescent Burgundy [AXMIBM]	N/C
Iridescent Navy [AXMINM]	N/C
Royal Blue [AXMRBM]	N/C

### Metallic Finish Back Shell

*If no selection is made, standard Shadow Black will be used with no additional lead time*

Blue Steel [AXMBLS]	N/C
---------------------	-----