



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Submit Form

Effective: 6/17/2024

Liberty Folding - Transit Not Included [LFOPTTRN] _____ £3000
 250 lb weight capacity Tilt range: 0° - 20°

Liberty Rigid, Transit Not Included [LFRGDOMT] _____ £3000
 250 lb weight capacity Tilt range: 0° - 20°

Frame

Seat Width

Select Seat Width within desired Seat Width Range

Frame Width Range	Frame Width				
Small 14" - 18" [LFFRW001]	14 [LFSWT014] N/C	15 [LFSWT015] N/C	16 [LFSWT016] N/C	17 [LFSWT017] N/C	18 [LFSWT018] N/C
Medium 16" - 20" [LFFRW002]	16 [LFSWT016] N/C	17 [LFSWT017] N/C	18 [LFSWT018] N/C	19 [LFSWT019] N/C	20 [LFSWT020] N/C
Large 18" - 22" [LFFRW003]	18 [LFSWT018] N/C	19 [LFSWT019] N/C	20 [LFSWT020] N/C	21 [LFSWT021] £455	22 [LFSWT022] £455

Seat Depth

Select Seat Depth within desired Seat Depth Range

Seat Depth Range	Seat Depth		
Short 14" - 16" [LFFDR001]	14 [LFSDE014] N/C	15 [LFSDE015] N/C	16 [LFSDE016] N/C
Medium 16" - 18" [LFFDR002]	16 [LFSDE016] N/C	17 [LFSDE017] N/C	18 [LFSDE018] N/C
Long 18" - 20" [LFFDR003]	18 [LFSDE018] N/C	19 [LFSDE019] N/C	20 [LFSDE020] N/C

Assistive Gas Strut Set-Up

Client Weight > 200 lbs (Setting 1, 450N Spring) [LFGAS001] _____ £121
 Client Weight 176-200 lbs (Setting 2, 450N Spring) [LFGAS002] _____ £121
 Client Weight 151-175 lbs (Setting 3, 450N Spring) [LFGAS003] _____ £121
Standard on chairs less than 16" deep
 Client Weight 151-175 lbs (Setting 1, 300N Spring) [LFGAS004] _____ N/C
Standard on chairs 16" deep or larger
 Client Weight 125-150 lbs (Setting 2, 300N Spring) [LFGAS005] _____ N/C
 Client Weight < 125 lbs (Setting 3, 300N Spring) [LFGAS006] _____ N/C

Seat Option

Standard Seat Sling [LFGSS001] _____ N/C
 Growing Adjustable Seat Pan [LFGACC004] _____ £173
 Omit Seat [LFCUS003] _____ N/C

Transit Option

Transit System [LFGTRNOPT] _____ N/C
No charge when selected with Liberty - Transit Included
 Omit [LFGTRNOP2] _____ N/C

Backrest

Back Post / Height Range / Back Height

Back Post	Price	Height Range	Back Height				
Stroller Handle Back Post ¹ [LFBOP016]	N/C	Short [LFBHE002]	20 [LFBRH200]				
		Tall [LFBHE004]	24 [LFBRH240]				
Height Adjustable Straight with Push Handle [LFBOP001]	£25	Short [LFBHE002]	14 [LFBRH140]	15 [LFBRH150]	16 [LFBRH160]	17 [LFBRH170]	18 [LFBRH180]
		Medium [LFBHE003]	16 [LFBRH160]	17 [LFBRH170]	18 [LFBRH180]	19 [LFBRH190]	20 [LFBRH200]
		Tall [LFBHE004]	18 [LFBRH180]	19 [LFBRH190]	20 [LFBRH200]	21 [LFBRH210]	22 [LFBRH220]
Height Adjustable 8° Bend ² with Push Handle [LFBOP003]	£25	Short [LFBHE002]	14 [LFBRH140]	15 [LFBRH150]	16 [LFBRH160]	17 [LFBRH170]	18 [LFBRH180]
		Medium [LFBHE003]	16 [LFBRH160]	17 [LFBRH170]	18 [LFBRH180]	19 [LFBRH190]	20 [LFBRH200]
		Tall [LFBHE004]	18 [LFBRH180]	19 [LFBRH190]	20 [LFBRH200]	21 [LFBRH210]	22 [LFBRH220]

¹ Not available with Back Upholstery

² Not available with IV Pole

Stroller Handle Back Post Option

Available only with Stroller Handle Back Post

Removable Stroller Handle Extension [LFBOP012] £73

Backrest Options

Fixed (Non-Folding) [LFFDB001] N/C

Angle Adjustable (Folding) [LFFDB002] £92

Back Upholstery

If selecting both upholstery and Axiom Back, the upholstery will be installed.

Standard [LFBUP001] N/C

Tension Adjustable [LFBUP002] £182

Not available on Seat Widths less than 16" combined with 21" and 22"

Backrest Heights

Omit [LFBUP003] N/C

Footrest

Hanger

Hanger	60°	70°	80°	90°		Retail
Extension Tube Mount	[LFHAN001]	[LFHAN002]	[LFHAN003] ¹			N/C
Front Mount		[LFHAN016]	[LFHAN017] ¹	[LFHAN004] ²		N/C
Pro Elevating Leg Rests					[LFHAN044]	£339
Omit Hangers					[LFHAN046]	N/C

¹ Not available with 8" casters

² Not available with 8" casters in combination with 3" extension

Footplates

Composite [LFFPL001] N/C

Not available with Front Mount Hanger

Composite Angle Adjustable [LFFPL014] £170

Not available with Front Mount Hanger or 14" Seat Width

Aluminum Angle Adjustable [LFFPL002] £184

Aluminum Locking Multi-Angle Adjustable [LFFPL020] £303

Available only with Front Mount Hanger

One Piece Flip Up Angle Adjustable [LFFPL009] £188

Not available with Front Mount Hanger or PRO ELR

PRO ELR Cover

Available only on PRO ELR

Pro ELR Gel Pad Cover (PAIR) [LFELRPAD] £304

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length												
		10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	15"	18"			
Extension Tube Hangers (60°, 70°, 80°)	Short (10" - 13.5")	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"					
	Medium (11.5" - 15.5")	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"				
	Long (14" - 18")	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"				
Front Mount Hangers ¹	Short Front Mount Hanger (70°, 80°, 90°)	Standard Selections									Requires 3" Extension			
		4" ²	4.5" ²	5"	5.5"	6"					8.5"	9"	9.5"	
	Medium Front Mount Hanger (70°, 80°, 90°)	6.5"	7"	7.5"	8"						10"	10.5"	11"	
		4" ²	4.5" ²	5"	5.5"	6"	6.5"	7"			10.5"	11"	11.5"	
	Long Front Mount Hanger (70°, 80°)	7.5"	8"	8.5"	9"	9.5"	10"				12"	12.5"	13"	
		5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"		12.5"	13"	13.5"	
		9"	9.5"	10"	10.5"	11"	11.5"	12"		14"	14.5"	15"		
PRO ELR Low Pivot Setting	PRO ELR Low Pivot Setting Short (10.5" - 15.5")	10.5" - 15.5"												
	PRO ELR Low Pivot Setting Long (13" - 19")	13" - 19"												
PRO ELR High Pivot Setting	PRO ELR High Pivot Setting Short (8.5" - 17")	8.5" - 17"												
	PRO ELR High Pivot Setting Long (11.5" - 17.5")	11.5" - 17.5"												

¹ See Minimum Seat Height chart below

² Available with 90° front mount hangers

Minimum Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
	70°	80/90°	70°	80/90°	70°	80°
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Seat Height with Standard Selection	13"	13"	13"	13.5"	14.5"	15.5"
Minimum Seat Height with 3" Extension	13.5"	14"	15.5"	16"	17"	18"

Heel Loops

Not available with One Piece Flip Up Angle Adjustable Footplate

Heel Loops [LFACC008]	£29
Omit Heel Loops [LFACC011]	N/C

Calf Straps

Velcro Adjustable [LFACC002]	£34
Padded Velcro Adjustable [LFACC003]	£66
Padded Airmesh Calf Strap [LFACC010]	£88

Hanger Extension Tubes

Available only with Front Mount Hangers

3" Extension [LFEXT008]	N/C
3" Extension (Provided as Parts) [LFEXT018]	N/C

Residual Limb Support

	Limb Support	Price	Omit	Cover		Pad Size					
				Smooth [LFRLS006]	Air Mesh [LFRLS007]	5 x 8 [LFRLS008]	7 x 10 [LFRLS009]	9 x 10 [LFRLS010]	9 x 14 [LFRLS011]	11 x 10 [LFRLS012]	11 x 14 [LFRLS013]
Left	Left Limb Support	£182	Omit Left Hanger [LFRLS004]	Smooth [LFRLS006]	Air Mesh [LFRLS007]	5 x 8 [LFRLS008]	7 x 10 [LFRLS009]	9 x 10 [LFRLS010]	9 x 14 [LFRLS011]	11 x 10 [LFRLS012]	11 x 14 [LFRLS013]
Right	Right Limb Support	£182	Omit Right Hanger [LFRLS005]	Smooth [LFRLS020]	Air Mesh [LFRLS021]	5 x 8 [LFRLS014]	7 x 10 [LFRLS015]	9 x 10 [LFRLS016]	9 x 14 [LFRLS017]	11 x 10 [LFRLS018]	11 x 14 [LFRLS019]

Hanger Release

4-Way - In-Line [LFREL002]	N/C
<i>Not available with PRO ELR</i>	
4-Way - Out [LFREL003]	N/C

Seat Height, Casters and Rear Wheels

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size			
	4"	5"	6"	8"
	Seat Height			
12"	13"	14"	15"	
	13.5"	14.5"	15.5"	17"
	14"	15"	16"	17.5"
	14.5"	15.5"	16.5"	18"
	15"	16"	17"	18.5"
	15.5"	16.5"	17.5"	19" ²
	16"	17"	18"	
		17.5"	18.5"	
		18"	19" ²	
		18.5"		
	19" ²			
16"	13.5" ¹	14" ¹	15"	
	14" ¹	14.5"	15.5"	17"
	14.5"	15"	16"	17.5"
	15"	15.5"	16.5"	18"
	15.5"	16"	17"	18.5"
	16"	16.5"	17.5"	19"
		17"	18"	19.5"
		17.5"	18.5"	20"
		18"	19"	
		18.5"		
	19"			
20"	13"	14"	15"	
	13.5"	14.5"	15.5"	17"
	14"	15"	16"	17.5"
	14.5"	15.5"	16.5"	18"
	15"	16"	17"	18.5"
	15.5"	16.5"	17.5"	19"
	16"	17"	18"	19.5"
		17.5"	18.5"	20"
		18"	19"	
		18.5"		
	19"			

Rear Wheel Size	Caster Size			
	4"	5"	6"	8"
	Seat Height			
22"	13.5" ¹	14"	15"	
	14"	14.5"	15.5"	17"
	14.5"	15"	16"	17.5"
	15"	15.5"	16.5"	18"
	15.5"	16"	17"	18.5"
	16"	16.5"	17.5"	19"
		17"	18"	19.5"
		17.5"	18.5"	20"
		18"	19" ²	
		18.5"		
	19" ²			
24"	14.5" ¹	14.5" ¹	15"	
	15"	15"	15.5"	17"
	15.5"	15.5"	16"	17.5"
	16"	16"	16.5"	18"
		16.5"	17"	18.5"
		17"	17.5"	19"
		17.5"	18"	19.5"
		18"	18.5"	20"
		18.5"	19"	
		19"		

¹ Available only with Low Profile tires
² Not available with Full Poly tire

Castors

Caster Type		Caster Size				
		Price	4"	5"	6"	8"
Caster Type	1" Poly	N/C	[LFCAS006]	[LFCAS010]	[LFCAS014]	[LFCAS020]
	1" Poly Aluminum	£106	[LFCAS007]	[LFCAS011]		
	1" Pneumatic	£91			[LFCAS017]	[LFCAS021]
	1.5" Poly	N/C	[LFCAS032]	[LFCAS025]	[LFCAS026]	
	1.5" Soft Roll Aluminum	£182	[LFCAS009]	[LFCAS013]	[LFCAS016]	
	2" Poly	£91			[LFCAS024]	
	2" Pneumatic w/Foam Insert	£218				[LFCAS023]



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



Wheels

		Rear Wheel Size					
		Price	12"	16"	20"	22"	24"
Rear Wheel Type	Mag	N/C	[LFWHE112]	[LFWHE116]	[LFWHE120]	[LFWHE122]	[LFWHE124]
	Spoke	N/C				[LFWH1122]	[LFWH1124]

Rear Tire Type

Full Profile

Street Pneumatic 1 3/8" [LFTIR001]	£110
Full Poly 1 3/8" [LFTIR002]	N/C
Pneumatic with Airless Insert 1 3/8" [LFTIR004]	£110

Low Profile

Low Poly 1" [LFTIR003]	N/C
<i>Available only with 16" wheel size</i>	
Sentinel Solid Tire 1" [LFTIR005]	£110
Schwalbe Marathon Plus Evolution 1" [LFTIR010]	£145
<i>Available only on 24" MAXX Spoke</i>	

Handrims

Standard Handrims

Handrims not available on 12" or 16" wheel sizes

Aluminum Anodized [LFHRM001]	N/C
Plastic Coated [LFHRM002]	£121
Projection (8 Vertical Projections) [LFHRM003]	£268
<i>Available only with full profile tires in combination with 22" and 24" MAXX Spoke wheels, or 22" and 24" Mag wheels</i>	
Omit Handrim [LFHRM012]	N/C

Wheel Locks

Wheel locks

Push to Lock [LFWLK002]	N/C
Pull to Lock [LFWLK003]	N/C
<i>Not available with 12" or 16" Full Poly tires</i>	
Hemi Wheel Lock - Left [LFWLK023]	£191
<i>Extension Handle included</i>	
Hemi Wheel Lock - Right [LFWLK024]	£191
<i>Extension Handle included</i>	

Standard Handrim Mount

If no selection is made, Standard will be used

Close Mount [LFHMT002]	N/C
<i>Reduces overall width by .75"</i>	

Handrim Options

Not available on Plastic Coated or Projection Handrims

3M Micro Grip [LFACC006]	£92
--------------------------	-----

Axles

No selection needed with Drum Brake

Fixed [LFAXO002]	N/C
Quick Release [LFAXO003]	N/C

Drum Brake [LFWLK015]	£532
<i>Available only with MAXX Spoke wheels. Axle is Quick Release.</i>	
Dual - Drum Brake with Push to Lock [LFWLK016]	£653
<i>Available only with MAXX Spoke wheels. Axle is Quick Release.</i>	
Dual - Drum Brake with Pull to Lock [LFWLK017]	£653
<i>Available only with MAXX Spoke wheels. Axle is Quick Release.</i>	
Omit Wheel Locks [LFWLK008]	N/C

Wheel Lock Options

Extension Handles [LFWLK004]	£65
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	
Grade Aids [LFWLK009]	£145
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>	

Armrests

Armrest Type	Length		Price
Height Adjustable T-Arm Tall (11.5"-15")	Desk	[LFARM039]	£212
	Full	[LFARM040]	
Height Adjustable T-Arm (8.5"-13.5")	Desk	[LFARM003]	£212
	Full	[LFARM004]	
Angle Adjustable Locking Extendable Flip Up ¹	Desk	[LFARM010]	£574
	Full	[LFARM011]	
Omit Armrests		[LFARM007]	N/C

¹ Available only with 12" and 16" wheel sizes and Stroller Handle Back Post

Accessories

Removable Side Guard

Composite Side Guards [LFSDG002]	£199
<i>Not available with T-Arms</i>	

Pelvic Positioning Belt

1 1/2" Auto Buckle Pelvic Positioning Belt [LFPOB001]	£97
1 1/2" Padded Auto Buckle Pelvic Positioning Belt [LFPOB002]	£151
2" Airline Buckle Pelvic Positioning Belt [LFPOB003]	£106
Small Flex-Mount Positioning Belt (16-24") [LFPOB007]	£188
Medium Flex-Mount Positioning Belt (19-29") [LFPOB008]	£188
Large Flex-Mount Positioning Belt (25-37") [LFPOB009]	£188

Belt Mounting Kit

Seat Belt Mounting Hardware [LFACC024]	£17
Belt Mounting Band Clamp [LFACC054]	£103

O2 Holder

O2 Holder [LFACC018]	£194
<i>Mounted to the left side. Not available with Fold Down Back.</i>	

Frame Color

Black [LFFCS001]	N/C
Candy Blue [LFFCS012]	N/C
Charcoal Gray [LFFCS031]	N/C
Candy Purple [LFFCS010]	N/C
Iridescent Burgundy [LFFCS021]	N/C
Candy Red [LFFCS003]	N/C
Iridescent Navy [LFFCS016]	N/C
Bright Red [LFFCS062]	N/C
Rich Teal [LFFCS061]	N/C
Dry Earth Tan [LFFCS059]	N/C
Tactical Green [LFFCS058]	N/C
Radiant Orange [LFFCS056]	N/C
Dark Pink [LFFCS053]	N/C

Arm Pads

Standard Arm Pad [LFARM026]	N/C
Waterfall Arm Pad [LFARM020]	£25
<i>Not available with Angle Adjustable Locking Extendable Flip Up Armrests</i>	
Foam Grip [LFARM025]	N/C
<i>Available only with Angle Adjustable Locking Extendable Flip Up Armrest</i>	

Anti-Tip

Rear Anti-Tips [LFATP001]	£110
Omit [LFATP000]	N/C

Transit Wheels

Transit wheels [LFACC035]	£303
<i>Not available with Long Depth Seat Frame, Seat Heights less than 14" or greater than 18", 12" or 16" rear wheels, Cane and Crutch Holder, and O2 Holder</i>	

Other Accessories

Tool Kit [LFACC005]	£20
Removable Underseat Pouch [LFACC015]	£30
Ki Mobility Backpack [LFACC009]	£100
Black Neoprene Impact Guards [LFIMG001]	£82
Cane and Crutch Holder [LFACC007]	£188
IV Pole - Mounted Right Side [LFACC030]	£170
<i>Not available with Transit Option or 8° Bend with Push Handle Back Post</i>	
IV Pole - Mounted Left Side [LFACC031]	£170
<i>Not available with Transit Option or 8° Bend with Push Handle Back Post. Left side mount not available with O2 Holder.</i>	

Axiom Back Width	Fits Back Frame Widths	
	Inner Slots	Outer Slots
12" (30.5cm)	11" (28.0cm)	12" -14" (30.5cm - 35.5cm)
14" (35.5cm)	12" -13" (30.5cm - 33.0cm)	14" -16" (35.5cm - 40.5cm)
16" (40.5cm)	14" -15" (35.5cm - 38.0cm)	16" -18" (40.5cm - 45.5cm)
18" (45.5cm)	16" - 17" (40.5cm - 43.0cm)	18" - 20" (45.5cm - 51.0cm)
20" (51.0cm)	18" -19" (45.5cm - 48.5cm)	20" -22" (51.0cm - 56.0cm)
22" (56.0cm)	20" -21" (51.0cm - 53.5cm)	22" -24" (56.0cm - 61.0cm)

Hardware Description	Code
Improved DEX Hardware, removable adjustable hardware (300 lb max user weight)	DEX
New BTX Hardware, removable adjustable hardware not recommended for, or installed on Ultralight Rigid Frame Wheelchairs (400 lb max user weight)	BTX
New Fixed Hardware, non-removeable hardware up to (300 lb max user weight)	FXT



Axiom Back

Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports.

Axiom Posterior Back (PB)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPB1208]	[AXPB1408]	[AXPB1608]	[AXPB1808]	[AXPB2008]	
		10"		[AXPB1410]	[AXPB1610]	[AXPB1810]	[AXPB2010]	[AXPB2210]
		12"	[AXPB1212]	[AXPB1412]	[AXPB1612]	[AXPB1812]	[AXPB2012]	
	Mid	14"		[AXPB1414]	[AXPB1614]	[AXPB1814]	[AXPB2014]	[AXPB2214]
		16"	[AXPB1216]	[AXPB1416]	[AXPB1616]	[AXPB1816]	[AXPB2016]	
		18"		[AXPB1418]	[AXPB1618]	[AXPB1818]	[AXPB2018]	[AXPB2218]

Back Hardware

FXT Hardware is not available with Lateral Support Options

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



Lateral Support

Pair	Description	Item Code	Pair
	3 x 4 Lateral Pads with Fixed Mounting Hardware - Pair	PBLF3x4PP	£295
	4 x 5 Lateral Pads with Fixed Mounting Hardware - Pair	PBLF4x5PP	£295
	3 x 4 Lateral Pads with Swing Away Mounting Hardware - Pair	PBLS3X4PP	£628
	4 x 5 Lateral Pads with Swing Away Mounting Hardware - Pair	PBLS4X5PP	£628
	3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair	PBLAS3X4PP	£628
	4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair	PBLAS4X5PP	£628

Axiom Posterior Lateral Back (PL)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Low	8"	[AXPL1208]	[AXPL1408]	[AXPL1608]	[AXPL1808]	[AXPL2008]	
		10"		[AXPL1410]	[AXPL1610]	[AXPL1810]	[AXPL2010]	[AXPL2210]
		12"	[AXPL1212]	[AXPL1412]	[AXPL1612]	[AXPL1812]	[AXPL2012]	
	Mid	14"		[AXPL1414]	[AXPL1614]	[AXPL1814]	[AXPL2014]	[AXPL2214]
		16"	[AXPL1216]	[AXPL1416]	[AXPL1616]	[AXPL1816]	[AXPL2016]	
		18"		[AXPL1418]	[AXPL1618]	[AXPL1818]	[AXPL2018]	[AXPL2218]
	Tail	20"	[AXPL1220]	[AXPL1420]	[AXPL1620]	[AXPL1820]	[AXPL2020]	
		22"		[AXPL1422]	[AXPL1622]	[AXPL1822]	[AXPL2022]	[AXPL2222]

Back Hardware

FXT Hardware is not available with Lateral Support Options

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Lateral Support

Pair	Description	Item Code	Pair
	3 x 4 Lateral Pads with Fixed Mounting Hardware - Pair	PLLF3x4PP	£295
	4 x 5 Lateral Pads with Fixed Mounting Hardware - Pair	PLLF4x5PP	£295
	3 x 4 Lateral Pads with Swing Away Mounting Hardware - Pair	PLLS3x4PP	£628
	4 x 5 Lateral Pads with Swing Away Mounting Hardware - Pair	PLLS4x5PP	£628
	3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair	PLLAS3x4PP	£628
	4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair	PLLAS4x5PP	£628

Axiom Deep Lateral Back (DL)

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		£550	£550	£550	£550	£550	£550	
Height	Mid	14"	[AXDL1214]	[AXDL1414]	[AXDL1614]	[AXDL1814]	[AXDL2014]	[AXDL2214]
		16"		[AXDL1416]	[AXDL1616]	[AXDL1816]	[AXDL2016]	
		18"	[AXDL1218]	[AXDL1418]	[AXDL1618]	[AXDL1818]	[AXDL2018]	[AXDL2218]
	Tail	20"	[AXDL1220]	[AXDL1420]	[AXDL1620]	[AXDL1820]	[AXDL2020]	
		22"		[AXDL1422]	[AXDL1622]	[AXDL1822]	[AXDL2022]	[AXDL2222]

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com



Headrest

Description	Std. Straight Removable Adjustable Hardware		Multi-Axis Removable Adjustable Hardware		Posalinc Removable Adjustable Hardware		Posalinc Flip-Back Adjustable Hardware	
	Item Code	Retail	Item Code	Retail	Item Code	Retail	Item Code	Retail
Comfort Soft Headrest 8"	[CS08SS]	£268	[CS08MA]	£268	[CS08PL]	£268	[CS08PLF]	N/C
Comfort Soft Headrest 10"	[CS10SS]	£268	[CS10MA]	£268	[CS10PL]	£268	[CS10PLF]	N/C
Comfort Soft Headrest 14"	[CS14SS]	£268	[CS14MA]	£268	[CS14PL]	£268	[CS14PLF]	N/C
Contour Positioning Pad, 8" X 6.25"	[CP08SS]	£303	[CP08MA]	£303	[CP08PL]	£303	[CP08PLF]	N/C
Contour Positioning Pad, 9" X 6.5"	[CP09SS]	£303	[CP09MA]	£303	[CP09PL]	£303	[CP09PLF]	N/C
Contour Positioning Pad, 10" X 7.5"	[CP10SS]	£303	[CP10MA]	£303	[CP10PL]	£303	[CP10PLF]	N/C
Tri-soft Headrest 10"	[TS10SS]	£386	[TS10MA]	£386	[TS10PL]	£386	[TS10PLF]	N/C
Tri-soft Headrest 14"	[TS14SS]	£386	[TS14MA]	£386	[TS14PL]	£386	[TS14PLF]	N/C

Description	Item Code	Retail
Universal Headrest Mounting Plate	[UMP]	£43

Not required for Axiom Headrests

Anterior Trunk Supports

Selection of Anterior Thoracic Mounting Modification required in order to mount posture supports. Not available with Low Backs.

Anterior Trunk Supports		
Description	Item Code	Retail
Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small	[14941]	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium	[14942]	£100
Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large	[14943]	£100
Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 18" Width, Small	[14924]	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium	[14927]	£100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large	[14928]	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small	[14851]	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Small	[14852]	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Medium	[14853]	£100
Dynaform Anterior Trunk Support: Standard, Dynamic, Large	[14854]	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small	[14856]	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small	[14857]	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14858]	£100
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large	[14859]	£100
AirLogic Anterior Trunk Support: Standard, Structure, Small	[14282]	£300
AirLogic Anterior Trunk Support: Standard, Structure, Medium	[14283]	£300
AirLogic Anterior Trunk Support: Standard, Structure, Large	[14284]	£300
AirLogic Anterior Trunk Support: Standard, Structure, X-Large	[14285]	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small	[14257]	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium	[14258]	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Large	[14259]	£300
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, X-Large	[14260]	£300

Accessories

Adjustable Shoulder Guides [AXBSG] <i>Includes Anterior Thoracic Mounting Modification. Not available with Low Backs.</i>	£43
Lumbar Pad [LUMPAD] <i>Available only with Axiom Lateral Back and Axiom Deep Lateral back</i>	N/C
Privacy Flap [PRFLAP]	£44

Axiom Back Modifications

Axiom Back Modifications

Extended lead times will apply for all modifications. Makes back non-refundable.

Anterior Thoracic Mounting Modification [AXSMM] <i>Provides holes for Anterior Thoracic Supports and Shoulder Guides</i>	N/C
---	-----

Gloss Finish Back Shell

If no selection is made, standard Shadow Black will be used with no additional lead time

Opalescent Violet [AXMOVI]	N/C
Radiant Orange [AXMROR]	N/C
Electric Green [AXMEGN]	N/C
Glacier Blue [AXMGBL]	N/C
Dark Pink [AXMDPK]	N/C
Light Pink [AXMLPK]	N/C
Rich Teal [AXMRTE]	N/C
Purple Plum [AXMPLM]	N/C
Black [AXMBLK]	N/C
Candy Blue [AXMCBL]	N/C
Candy Purple [AXMPUR]	N/C
Candy Red [AXMRED]	N/C
Granny Smith [AXMGSM]	N/C
Metallic Orange [AXMMOR]	N/C
Yellow [AXMYLW]	N/C

Matte Finish Back Shell

If no selection is made, standard Shadow Black will be used with no additional lead time

Bright Red [AXMBRR]	N/C
White [AXMWHT]	N/C
Dry Earth Tan [AXMTAN]	N/C
Tactical Green [AXMTGR]	N/C
Charcoal Gray [AXMCGM]	N/C
Iridescent Burgundy [AXMIBM]	N/C
Iridescent Navy [AXMINM]	N/C
Royal Blue [AXMRBM]	N/C

Metallic Finish Back Shell

If no selection is made, standard Shadow Black will be used with no additional lead time

Blue Steel [AXMBLS]	N/C
---------------------	-----

Axiom G

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		£85	£85	£85	£85	£85	£85	£85	£85
Depth	14"	[XG1414]							
	15"	[XG1415]	[XG1515]	[XG1615]					
	16"	[XG1416]	[XG1516]	[XG1616]	[XG1716]	[XG1816]	[XG1916]	[XG2016]	
	17"	[XG1417]	[XG1517]	[XG1617]	[XG1717]	[XG1817]	[XG1917]	[XG2017]	
	18"	[XG1418]	[XG1518]	[XG1618]	[XG1718]	[XG1818]	[XG1918]	[XG2018]	[XG2218]
	19"		[XG1519]	[XG1619]	[XG1719]	[XG1819]	[XG1919]	[XG2019]	[XG2219]
	20"		[XG1520]	[XG1620]	[XG1720]	[XG1820]	[XG1920]	[XG2020]	[XG2220]
	21"			[XG1621]		[XG1821]	[XG1921]	[XG2021]	[XG2221]
22"			[XG1622]		[XG1822]	[XG1922]	[XG2022]	[XG2222]	

Axiom G - General Use Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		£125	£125	£125	£125	£125	£125	£125	£125
Depth	14"	[XG1414C2]							
	15"	[XG1415C2]	[XG1515C2]	[XG1615C2]					
	16"	[XG1416C2]	[XG1516C2]	[XG1616C2]	[XG1716C2]	[XG1816C2]	[XG1916C2]	[XG2016C2]	
	17"	[XG1417C2]	[XG1517C2]	[XG1617C2]	[XG1717C2]	[XG1817C2]	[XG1917C2]	[XG2017C2]	
	18"	[XG1418C2]	[XG1518C2]	[XG1618C2]	[XG1718C2]	[XG1818C2]	[XG1918C2]	[XG2018C2]	[XG2218C2]
	19"		[XG1519C2]	[XG1619C2]	[XG1719C2]	[XG1819C2]	[XG1919C2]	[XG2019C2]	[XG2219C2]
	20"		[XG1520C2]	[XG1620C2]	[XG1720C2]	[XG1820C2]	[XG1920C2]	[XG2020C2]	[XG2220C2]
	21"			[XG1621C2]		[XG1821C2]	[XG1921C2]	[XG2021C2]	[XG2221C2]
22"			[XG1622C2]		[XG1822C2]	[XG1922C2]	[XG2022C2]	[XG2222C2]	

Axiom S - Skin Protection Cushion

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
	Retail	£120	£120	£120	£120	£120	£120	£120	£120
Depth	14"	[XS1414]							
	15"	[XS1415]	[XS1515]	[XS1615]					
	16"	[XS1416]	[XS1516]	[XS1616]	[XS1716]	[XS1816]	[XS1916]	[XS2016]	
	17"	[XS1417]	[XS1517]	[XS1617]	[XS1717]	[XS1817]	[XS1917]	[XS2017]	
	18"	[XS1418]	[XS1518]	[XS1618]	[XS1718]	[XS1818]	[XS1918]	[XS2018]	[XS2218]
	19"		[XS1519]	[XS1619]	[XS1719]	[XS1819]	[XS1919]	[XS2019]	[XS2219]
	20"		[XS1520]	[XS1620]	[XS1720]	[XS1820]	[XS1920]	[XS2020]	[XS2220]
	21"			[XS1621]		[XS1821]	[XS1921]	[XS2021]	[XS2221]
22"			[XS1622]		[XS1822]	[XS1922]	[XS2022]	[XS2222]	

Axiom S - Skin Protection Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
	Retail	£160	£160	£160	£160	£160	£160	£160	£160
Depth	14"	[XS1414C2]							
	15"	[XS1415C2]	[XS1515C2]	[XS1615C2]					
	16"	[XS1416C2]	[XS1516C2]	[XS1616C2]	[XS1716C2]	[XS1816C2]	[XS1916C2]	[XS2016C2]	
	17"	[XS1417C2]	[XS1517C2]	[XS1617C2]	[XS1717C2]	[XS1817C2]	[XS1917C2]	[XS2017C2]	
	18"	[XS1418C2]	[XS1518C2]	[XS1618C2]	[XS1718C2]	[XS1818C2]	[XS1918C2]	[XS2018C2]	[XS2218C2]
	19"		[XS1519C2]	[XS1619C2]	[XS1719C2]	[XS1819C2]	[XS1919C2]	[XS2019C2]	[XS2219C2]
	20"		[XS1520C2]	[XS1620C2]	[XS1720C2]	[XS1820C2]	[XS1920C2]	[XS2020C2]	[XS2220C2]
	21"			[XS1621C2]		[XS1821C2]	[XS1921C2]	[XS2021C2]	[XS2221C2]
22"			[XS1622C2]		[XS1822C2]	[XS1922C2]	[XS2022C2]	[XS2222C2]	

Axiom P - Positioning Cushion

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
	Retail	£170	£170	£170	£170	£170	£170	£170	£170
Depth	14"	[XP1414]							
	15"	[XP1415]	[XP1515]	[XP1615]					
	16"	[XP1416]	[XP1516]	[XP1616]	[XP1716]	[XP1816]	[XP1916]	[XP2016]	
	17"	[XP1417]	[XP1517]	[XP1617]	[XP1717]	[XP1817]	[XP1917]	[XP2017]	
	18"	[XP1418]	[XP1518]	[XP1618]	[XP1718]	[XP1818]	[XP1918]	[XP2018]	[XP2218]
	19"		[XP1519]	[XP1619]	[XP1719]	[XP1819]	[XP1919]	[XP2019]	[XP2219]
	20"		[XP1520]	[XP1620]	[XP1720]	[XP1820]	[XP1920]	[XP2020]	[XP2220]
	21"			[XP1621]		[XP1821]	[XP1921]	[XP2021]	[XP2221]
22"			[XP1622]		[XP1822]	[XP1922]	[XP2022]	[XP2222]	

Axiom P - Positioning Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		£210	£210	£210	£210	£210	£210	£210	£210
Depth	14"	[XP1414C2]							
	15"	[XP1415C2]	[XP1515C2]	[XP1615C2]					
	16"	[XP1416C2]	[XP1516C2]	[XP1616C2]	[XP1716C2]	[XP1816C2]	[XP1916C2]	[XP2016C2]	
	17"	[XP1417C2]	[XP1517C2]	[XP1617C2]	[XP1717C2]	[XP1817C2]	[XP1917C2]	[XP2017C2]	
	18"	[XP1418C2]	[XP1518C2]	[XP1618C2]	[XP1718C2]	[XP1818C2]	[XP1918C2]	[XP2018C2]	[XP2218C2]
	19"		[XP1519C2]	[XP1619C2]	[XP1719C2]	[XP1819C2]	[XP1919C2]	[XP2019C2]	[XP2219C2]
	20"		[XP1520C2]	[XP1620C2]	[XP1720C2]	[XP1820C2]	[XP1920C2]	[XP2020C2]	[XP2220C2]
	21"			[XP1621C2]		[XP1821C2]	[XP1921C2]	[XP2021C2]	[XP2221C2]
22"			[XP1622C2]		[XP1822C2]	[XP1922C2]	[XP2022C2]	[XP2222C2]	

Axiom SP Visco - Skin Protection & Positioning Cushion

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		£350	£350	£350	£350	£350	£350	£350	£350	£350
Depth	14"	[XSP1414V]	[XSP1514V]	[XSP1614V]						
	15"	[XSP1415V]	[XSP1515V]	[XSP1615V]	[XSP1715V]	[XSP1815V]				
	16"	[XSP1416V]	[XSP1516V]	[XSP1616V]	[XSP1716V]	[XSP1816V]	[XSP1916V]	[XSP2016V]		[XSP2216V]
	17"	[XSP1417V]	[XSP1517V]	[XSP1617V]	[XSP1717V]	[XSP1817V]	[XSP1917V]	[XSP2017V]		[XSP2217V]
	18"	[XSP1418V]	[XSP1518V]	[XSP1618V]	[XSP1718V]	[XSP1818V]	[XSP1918V]	[XSP2018V]	[XSP2118V]	[XSP2218V]
	19"	[XSP1419V]	[XSP1519V]	[XSP1619V]	[XSP1719V]	[XSP1819V]	[XSP1919V]	[XSP2019V]	[XSP2119V]	[XSP2219V]
	20"	[XSP1420V]	[XSP1520V]	[XSP1620V]	[XSP1720V]	[XSP1820V]	[XSP1920V]	[XSP2020V]	[XSP2120V]	[XSP2220V]
	21"		[XSP1521V]	[XSP1621V]	[XSP1721V]	[XSP1821V]	[XSP1921V]	[XSP2021V]	[XSP2121V]	[XSP2221V]
22"		[XSP1522V]	[XSP1622V]	[XSP1722V]	[XSP1822V]	[XSP1922V]	[XSP2022V]	[XSP2122V]	[XSP2222V]	

Axiom SP Visco - Skin Protection & Positioning Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		£410	£410	£410	£410	£410	£410	£410	£410	£410
Depth	14"	[XSP1414VC2]	[XSP1514VC2]	[XSP1614VC2]						
	15"	[XSP1415VC2]	[XSP1515VC2]	[XSP1615VC2]	[XSP1715VC2]	[XSP1815VC2]				
	16"	[XSP1416VC2]	[XSP1516VC2]	[XSP1616VC2]	[XSP1716VC2]	[XSP1816VC2]	[XSP1916VC2]	[XSP2016VC2]		[XSP2216VC2]
	17"	[XSP1417VC2]	[XSP1517VC2]	[XSP1617VC2]	[XSP1717VC2]	[XSP1817VC2]	[XSP1917VC2]	[XSP2017VC2]		[XSP2217VC2]
	18"	[XSP1418VC2]	[XSP1518VC2]	[XSP1618VC2]	[XSP1718VC2]	[XSP1818VC2]	[XSP1918VC2]	[XSP2018VC2]	[XSP2118VC2]	[XSP2218VC2]
	19"	[XSP1419VC2]	[XSP1519VC2]	[XSP1619VC2]	[XSP1719VC2]	[XSP1819VC2]	[XSP1919VC2]	[XSP2019VC2]	[XSP2119VC2]	[XSP2219VC2]
	20"	[XSP1420VC2]	[XSP1520VC2]	[XSP1620VC2]	[XSP1720VC2]	[XSP1820VC2]	[XSP1920VC2]	[XSP2020VC2]	[XSP2120VC2]	[XSP2220VC2]
	21"		[XSP1521VC2]	[XSP1621VC2]	[XSP1721VC2]	[XSP1821VC2]	[XSP1921VC2]	[XSP2021VC2]	[XSP2121VC2]	[XSP2221VC2]
22"		[XSP1522VC2]	[XSP1622VC2]	[XSP1722VC2]	[XSP1822VC2]	[XSP1922VC2]	[XSP2022VC2]	[XSP2122VC2]	[XSP2222VC2]	

Axiom SP Fluid - Skin Protection & Positioning Cushion

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		£320	£320	£320	£320	£320	£320	£320	£320	£320
Depth	14"	[XSP1414F]	[XSP1514F]	[XSP1614F]						
	15"	[XSP1415F]	[XSP1515F]	[XSP1615F]	[XSP1715F]	[XSP1815F]				
	16"	[XSP1416F]	[XSP1516F]	[XSP1616F]	[XSP1716F]	[XSP1816F]	[XSP1916F]	[XSP2016F]		[XSP2216F]
	17"	[XSP1417F]	[XSP1517F]	[XSP1617F]	[XSP1717F]	[XSP1817F]	[XSP1917F]	[XSP2017F]		[XSP2217F]
	18"	[XSP1418F]	[XSP1518F]	[XSP1618F]	[XSP1718F]	[XSP1818F]	[XSP1918F]	[XSP2018F]	[XSP2118F]	[XSP2218F]
	19"	[XSP1419F]	[XSP1519F]	[XSP1619F]	[XSP1719F]	[XSP1819F]	[XSP1919F]	[XSP2019F]	[XSP2119F]	[XSP2219F]
	20"	[XSP1420F]	[XSP1520F]	[XSP1620F]	[XSP1720F]	[XSP1820F]	[XSP1920F]	[XSP2020F]	[XSP2120F]	[XSP2220F]
	21"		[XSP1521F]	[XSP1621F]	[XSP1721F]	[XSP1821F]	[XSP1921F]	[XSP2021F]	[XSP2121F]	[XSP2221F]
	22"		[XSP1522F]	[XSP1622F]	[XSP1722F]	[XSP1822F]	[XSP1922F]	[XSP2022F]	[XSP2122F]	[XSP2222F]

Axiom SP Fluid - Skin Protection & Positioning Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		£360	£360	£360	£360	£360	£360	£360	£360	£360
Depth	14"	[XSP1414FC2]	[XSP1514FC2]	[XSP1614FC2]						
	15"	[XSP1415FC2]	[XSP1515FC2]	[XSP1615FC2]	[XSP1715FC2]	[XSP1815FC2]				
	16"	[XSP1416FC2]	[XSP1516FC2]	[XSP1616FC2]	[XSP1716FC2]	[XSP1816FC2]	[XSP1916FC2]	[XSP2016FC2]		[XSP2216FC2]
	17"	[XSP1417FC2]	[XSP1517FC2]	[XSP1617FC2]	[XSP1717FC2]	[XSP1817FC2]	[XSP1917FC2]	[XSP2017FC2]		[XSP2217FC2]
	18"	[XSP1418FC2]	[XSP1518FC2]	[XSP1618FC2]	[XSP1718FC2]	[XSP1818FC2]	[XSP1918FC2]	[XSP2018FC2]	[XSP2118FC2]	[XSP2218FC2]
	19"	[XSP1419FC2]	[XSP1519FC2]	[XSP1619FC2]	[XSP1719FC2]	[XSP1819FC2]	[XSP1919FC2]	[XSP2019FC2]	[XSP2119FC2]	[XSP2219FC2]
	20"	[XSP1420FC2]	[XSP1520FC2]	[XSP1620FC2]	[XSP1720FC2]	[XSP1820FC2]	[XSP1920FC2]	[XSP2020FC2]	[XSP2120FC2]	[XSP2220FC2]
	21"		[XSP1521FC2]	[XSP1621FC2]	[XSP1721FC2]	[XSP1821FC2]	[XSP1921FC2]	[XSP2021FC2]	[XSP2121FC2]	[XSP2221FC2]
	22"		[XSP1522FC2]	[XSP1622FC2]	[XSP1722FC2]	[XSP1822FC2]	[XSP1922FC2]	[XSP2022FC2]	[XSP2122FC2]	[XSP2222FC2]