

BENOIT SYSTEMES

Autonomy is power




LIGHT RANGE 2.1

WWW.BENOITSYSTEMES.COM

POWER

Using the inherent strengths of manual wheelchairs, our power add-ons give people the freedom to tackle many terrains, remove the logistical problems and turn the journeys into more about the destination than the challenge of getting there.

THE FEATURES OF OUR POWER ADD-ONS



*See compatibility guide.


“ Robust and proven technology with strong backup. ”

BENOIT SYSTEMES, since 1994




HOW DOES IT WORK ?

You will be surprised by the lightness and simplicity of use of our power add-ons. They will not change either the overall dimensions of the chair, the folding ability nor the center of balance.



Set-up in 15 seconds



— LIGHT DRIVE 2.1 —
LIGHT DRIVE TURNS YOUR MANUAL WHEELCHAIR INTO AN ULTRA-LIGHT POWERCHAIR IN SECONDS.



— LIGHT ASSIST 2.1 —
ATTENDANT CONTROLLED LIGHT ASSIST TAKES ALL THE EFFORT OUT OF PUSHING THE WHEELCHAIR.





THE LIGHT DRIVE 2.1



USER
operated

Turn your manual wheelchair into an ultra-light power chair in seconds

Transports easily in any vehicle

Free wheel when you need

Independence to go where you want



THE LIGHT ASSIST 2.1



COMPANION
operated

Effortless pushing of the wheelchair for the attendant

Efficient and safe maneuvers

Speed control and braking aid

Journeys together become easier

WHEELCHAIR POWER ASSISTANCE MADE EASY



NEW IN 2021

New anti-tip



NEW IN 2021

More ergonomic clutch system



NEW IN 2021

Extended colour range



NEW IN 2021

New joystick extension

A POWERFUL & COMPACT ATTENDANT CONTROLLED POWER ADD-ON



NEW IN 2021

New storage for the handles



NEW IN 2021

More ergonomic handle locking system



NEW IN 2021

Extended colour range



NEW IN 2021

New anti-tip



NEW IN 2021

More ergonomic clutch system

LIGHT DRIVE 2.1 Mini

A SMALLER VERSION FOR CHILDREN



LIGHT DRIVE 2.1 Mini is instinctive, simple and easy to drive.

It encourages children's independence helping with the social development.

LIGHT DRIVE 2.1 Mini can be refitted as the child grows allowing to change chair without changing power add-on.

TECHNICAL DETAILS

Our power add-ons can be fitted on almost any wheelchair starting with 32 cm seat width and 40 cm between the wheels.

	Maximum user weight	Power add-on weight
LIGHT DRIVE 2.1	120 kg (265 lbs)	8,8 kg (19,4 lbs)
LIGHT ASSIST 2.1	120 kg (265 lbs)	8,1 kg (17,8 lbs)
LIGHT DRIVE 2.1 Mini	120 kg (265 lbs)	8,8 kg (19,4 lbs)
LIGHT DRIVE 2.1 Plus	240 kg (529 lbs)	16 kg (35,2 lbs)

- Max speed** ————— 10 km/h programmable
- Motor** ————— 2 x 130 W motors with brakes
- Motor LIGHT DRIVE 2.1 Plus** ————— 4 x 130 W motors with brakes
- Climbing ability** ————— 20 % gradient based on 70 kg person
- Manufacturer's warranty** ————— 24 months for the power add-on and the battery

HEAVY DUTY VERSION FOR USERS ABOVE 120 KG

LIGHT DRIVE 2.1 Plus

The extra-solid **LIGHT DRIVE 2.1 Plus** supports user's weight up to 240 kg.

Rugged and robust it enables users to enjoy the advantages of four powerful motors while keeping the lightness and mobility of a manual wheelchair.

Autonomy of at least 20 km is provided by a 18 Ah battery.



12AH LITHIUM BATTERY IS COMPATIBLE WITH APPLICABLE IATA REGULATIONS FOR AIR TRAVEL

Model	Weight	Range*
— Lithium Battery	— 3,85 kg (8,5 lbs)	— 15 to 18 km (9,3 to 11,2 m)

*Range depends on terrain and weight of the user



NEW IN 2021

Extended choice of colours (central tube)

- | CLASSICALS | METALLICS | TRENDS |
|------------|------------------|------------|
| ● BLACK | ● NIGHT BLUE | ● SKY BLUE |
| ● WHITE | ● GREY | ● PINK |
| ● BLUE | ● SILVER | ● GREEN |
| ● RED | ● COPPER | ● YELLOW |
| ● KHAKI | ● CHAMPAGNE GOLD | ● ORANGE |

For more information visit our website www.benoitsystemes.com

DISTRIBUTORS

ARGENTINA

Outlandish
www.outlandish.com.ar

AUSTRALIA

Mobility Care
www.mobilitycare.net.au

BELGIUM

Mobility By Olivier
www.mobilitybyolivier.be

CANADA

49 Bespoke
www.49bespoke.com

CYPRUS

Medicom Services
www.medicomstore.com

DENMARK

Sunrise Medical
www.sunrisemedical.dk

FINLAND

Movetta
www.movetta.fi

IRELAND

MMS Medical
www.mmsmedical.ie

ISRAËL

Sagol Direct
Rehabilitation Ltd
www.sagol-center.co.il

ITALY

Progettiamo Autonomia
www.my-pa.it

LUXEMBOURG

Megamed
www.megamed.lu

MALTE

Unicare Ltd
www.unicare.com.mt

NETHERLANDS

Mobility & You
www.mobility-you.nl

PORTUGAL

Rehapoint
www.rehapoint.pt

SINGAPORE

Teng Chiang
Holdings
www.tchold.com

SOUTH AFRICA

Sitwell
www.sitwell.co.za

SPAIN

Eo Prim
<https://www.primclinicassortopedicas.es/>

SWITZERLAND

Degonda Rehab
www.turbo-twist.ch

SWITZERLAND

Gloor Rehabilitation
www.gloorrehab.com

UNITED ARAB EMIRATES

Gulf Drug LLC
www.gulfdrug.com

UNITED KINGDOM

Recare
www.recare.co.uk



WWW.BENOITSYSTEMES.COM



Follow us!

📍 7 Rue du Pont,
21450 Billy Les Chanceaux FRANCE

☎️ (+33) 3 80 96 51 25

@ export@benoitssystemes.com

CE LIGHT RANGE 2.1 is compliant with EU
Regulation 2017/745 on medical devices.
These products are Class I medical devices.

YOUR LOCAL DEALER