

<p>Base Price - £1556.00</p> <p>Valid from 01.01.24 to 31.12.24 Please ensure you are using the latest order form</p>	<p>Order no.:</p> <p>-----</p>
--	--------------------------------

Lisa Rehab folding buggy

Size 1 for children from approx. 3 to 8 years
Size 2 for children from approx. 6 to 16 years



<input type="checkbox"/> Provide estimate only		<input type="checkbox"/> Order	
Billing address/customer number			
Delivery address			
Sales person		Name	
		
		Date	Signature
<p>By placing the order, you agree to comply with the instructions for use/service manuals for accessories and products of other manufacturers when combining them with an Ottobock product, in addition to the instructions for use/service manuals of the respective Ottobock products. Also you confirm to have read, understood and accepted, our general terms and conditions. In addition, you agree, upon placing an order, to conduct your own conformity assessment or approval in accordance with national requirements if and when you implement modifications or attachments that do not meet the specifications of Ottobock instructions for use/service manuals.</p>			

Standard equipment

Compact folding buggy in a very stable, light-weight aluminium design, back support adjustment for one-hand operation

- Seat upholstery and reinforced back support upholstery
- Adjustable and divisible leg support
- Front and rear wheels with ball bearings
- Front wheels with integrated swivel lock
- Wheel lock
- Screw connections guided by sleeves
- Standard lap belt
- Pre-mounted tip-assist, right
- Upholstery colour blue
- Frame colour: silver

Standard tyres:

- EVA wheels, front, 180 mm
- EVA wheels, rear, 250 mm

Standard model			
HR32130001-027	<input type="checkbox"/>	Lisa rehab folding buggy, size 1 with seat upholstery seat depth presetting 30 cm/12"	£1556.00
HR32130002-027	<input type="checkbox"/>	Lisa rehab folding buggy, size 1 with seat upholstery seat depth presetting 35 cm/14"	£1556.00
HR32140001-027	<input type="checkbox"/>	Lisa rehab folding buggy, size 2 with seat upholstery seat depth presetting 35 cm/14"	£1556.00
HR32140002-027	<input type="checkbox"/>	Lisa rehab folding buggy, size 2 with seat upholstery seat depth presetting 40 cm/16"	£1556.00

Accessories			
HR32137100-027	<input type="checkbox"/>	Arm support, piece (please order two pieces)	£76.00
HR16001102	<input type="checkbox"/>	Buggy board	£238.00
HR32135300-343	<input type="checkbox"/>	Simplified padding, size 1 (coro-blue)	£140.00
HR32145300-343	<input type="checkbox"/>	Simplified padding, size 2 (coro-blue)	£119.00
HR32135200-343	<input type="checkbox"/>	Seat pad with multi-functional back, size 1 (coro-blue/dark blue)	£355.00
HR32145200-343	<input type="checkbox"/>	Seat pad with multi-functional back, size 2 (coro-blue/dark blue)	£355.00

HR32137700-031	<input type="checkbox"/>	Head support, size 1, matching multifunctional back support, size 1	£227.00
HR32147700-031	<input type="checkbox"/>	Head support, size 2, matching multifunctional back support, size 2	£227.00
HR32070060-031	<input type="checkbox"/>	Head support size 1, padded (not available in combin. with seat pad with multifunctional back support)	£227.00
HR32086500-031	<input type="checkbox"/>	Head support size 2, padded (not available in combin. with seat pad with multifunctional back support)	£227.00
HR32136100	<input type="checkbox"/>	Five-point belt with neck pad	£214.00
HR32135500-031	<input type="checkbox"/>	Lateral guides with variable seat width reduction (4-6-8 cm) size 1 (only available in combination with pad for multifunctional back support)	£188.00
HR32145500-031	<input type="checkbox"/>	Lateral guides with variable seat width reduction (4-6-8 cm) size 2 (only available in combination with pad for multifunctional back support)	£147.00
HR32137200-031	<input type="checkbox"/>	Flip-down pommel size 1	£303.00
HR32147200-031	<input type="checkbox"/>	Flip-down pommel size 2	£311.00
HR32075800-031	<input type="checkbox"/>	Groin strap for multifunctional back support size 1	£190.00
HR32085800-031	<input type="checkbox"/>	Groin strap for multifunctional back support size 2	£179.00
HR32137800	<input type="checkbox"/>	Grab rail size 1	£122.00
HR32147800	<input type="checkbox"/>	Grab rail size 2	£134.00
HR32135700-031	<input type="checkbox"/>	Padding for grab rail size 1	£62.00
HR32145700-031	<input type="checkbox"/>	Padding for grab rail size 2	£62.00
HR32010080	<input type="checkbox"/>	Hook-and-loop foot straps with instep pad size 1	£36.00
HR32020080	<input type="checkbox"/>	Hook-and-loop foot straps with instep pad size 2	£35.00
HR32200120-343	<input type="checkbox"/>	Winter slip sack coro blue with detachable front part, with slits for straps and accessories, size 2 (seat depth up to approx. 34 cm)	£650.00
HR32400120-343	<input type="checkbox"/>	Winter slip sack coro blue with detachable front part, with slits for straps and accessories, size 3 (seat depth up to approx. 39 cm)	£852.00
HR32100100-074	<input type="checkbox"/>	Rain cape, dark blue size 1	£84.00
HR32010100-074	<input type="checkbox"/>	Rain cape, dark blue size 2	£84.00
481S64=ST005	<input type="checkbox"/>	Bckpack with Otto Bock logo	£69.00
HR32135800	<input type="checkbox"/>	Folding storage bag size 1	£102.00
HR32145800	<input type="checkbox"/>	Folding storage bag size 2	£102.00
HR32136200-094	<input type="checkbox"/>	Canopy in upholstery colour with transparent rain cover	£244.00

Information

Technical data		
	Size 1	Size 2
Seat width (with/without padding)	32/36 cm	37/41 cm
Seat depth	30 or 35 cm	35 or 40 cm
Back support height	65 cm	70 cm
Back support adjustment	30° (in steps 4 x 7.5°)	30° (in steps 4 x 7.5°)
Lower leg length	27 – 40 cm	27 – 40 cm
Max. load capacity	50 kg	60 kg
Weight	13.6 kg	14.6 kg