

## ***Lisa Buggy***

In town and when travelling, Lisa is the ideal choice for children aged three to sixteen, who have mild postural weaknesses. Lisa is compact, lightweight, collapsible, and very reliable. It is easy to manoeuvre even in small spaces thanks to the front swivel wheels.

### **What is special about me?**

- Higher load capacity compared to conventional buggies, and with a footrest
- Larger footrest for improved perception and stability for children with minor postural weakness



## Lisa Buggy

### Standard Model

- 2 sizes
- Colour: blue
- Seat upholstery and reinforced back upholstery
- Height-adjustable and divisible footrest
- Wheel lock
- Sturdy nylon cover
- Standard lap belt
- Folding storage bag
- Back angle adjustment
- Front wheels with integrated caster swivel lock
- EVA wheel set

### Options

- Headrest
- Armrests
- 5-point harness with neck pad
- Lateral guides with seat width reduction
- Pommel
- Groin strap
- Grip bar with padding
- Summer/Winter slip sack
- Rain cover and canopy
- Backpack and storage bag
- Buggy board
- Accessories bag
- For transportation in a wheelchair accessible vehicle



Reclining back



Transportation brackets



Umbrella folding



Underseat storage

### Technical Data

	Size 1	Size 2
Age	3 – 8 years	6 – 16 years
Max user weight	50kg	60kg
Seat depth	30 – 35cm	35 – 40cm
Seat width	32 – 36cm	37 – 41cm
Lower leg length	27 – 40cm	27 – 40cm
Back rest height	65cm	70cm
Folding size (LxWxH)	103 x 45 x 43cm	103 x 45 x 43cm
Wheel Ø front / rear	18 – 25cm	18 – 5cm
Total weight	13.5kg	14kg
Back angle	90° to 110°	90° to 120°