



## Ultra lightweight activity wheelchair in carbon design, suitable for (very) small and/or weak children (already from 12 months)

Included in the base price:

- Complete Carbon package: Curved back plate, seat plate, side guards and foot plate
- Frame versions Retro or Taurus
- Super light weight wheels
- 2 different V-shaped abductions for maximized leg space
- growable without additional parts
- adjustable back angle (80° -100°)
- Activity side guards with integrated bowden cable brake.
- Also included: stabilizer bar; Caster wheels with aluminum fork; wheel cover

*Not suitable as a seat for the transport in a passenger vehicle!*

**Quotation**  
**Order**

Weight capacity: 50 kg

Base price **£5,245**

Prices excl. VAT, valid from 01.03.2021

Delivery address / Company stamp		Billing address: <input type="checkbox"/> same as delivery address	
Name:		Name:	
Address:		Address:	
Ordering information:			
Reference/ internal description:		Date:	
Contact person:		Signature:	

Frame measures				cf. example on p.2
Seat width 10..	Seat depth 11..	Back height 12..	Seat height front in cm	Lower leg length
+ 2 cm extension each <input type="checkbox"/> 18 <input type="checkbox"/> 20 <input type="checkbox"/> 22 <input type="checkbox"/> 24 <input type="checkbox"/> 26 <input type="checkbox"/> 28 <input type="checkbox"/> 30	+ 4 cm extension each <input type="checkbox"/> 18 <input type="checkbox"/> 20 <input type="checkbox"/> 22 <input type="checkbox"/> 24 <input type="checkbox"/> 26 <input type="checkbox"/> 28 <input type="checkbox"/> 30	+ 2.5/5 cm extension each <input type="checkbox"/> 17.5 <input type="checkbox"/> 20 <input type="checkbox"/> 22.5 <input type="checkbox"/> 25 <input type="checkbox"/> 27.5 <input type="checkbox"/> 30 <input type="checkbox"/> 32.5 <input type="checkbox"/> 35	<b>At wheel size 18"</b> <input type="checkbox"/> 32 <input type="checkbox"/> 33,5 <input type="checkbox"/> 35 <b>At wheel size 20"</b> <input type="checkbox"/> 34 <input type="checkbox"/> 35,5 <input type="checkbox"/> 37 <input type="checkbox"/> 38,5 <input type="checkbox"/> 40 <b>At wheel size 22"</b> <input type="checkbox"/> 36 <input type="checkbox"/> 37,5 <input type="checkbox"/> 39 <input type="checkbox"/> 40,5 <input type="checkbox"/> 42 <b>At wheel size 24"</b> <input type="checkbox"/> 38,5 <input type="checkbox"/> 40 <input type="checkbox"/> 41,5 <input type="checkbox"/> 43 <input type="checkbox"/> 44,5 cf. table on p.7	<b>1. Real lower leg length of the user (LLL):</b> _____cm <b>2. Cushion thickness (CT)</b> _____cm

Frame versions				Retro-Frame
Retro-Frame N° 597301 <i>recommended for LLL 15 cm-30 cm*</i>		Taurus-Frame N° 597304 <i>recommended for LLL 20 cm-37 cm*</i>		Taurus-Frame
<input type="checkbox"/> <b>Strong abduction</b>  casters on one line with rear wheels	<input type="checkbox"/> <b>Slight abduction</b>  casters narrower than rear wheels	<input type="checkbox"/> <b>Strong abduction</b>  casters on one line with rear wheels	<input type="checkbox"/> <b>Slight abduction</b>  casters narrower than rear wheels	
N° 597305 <b>£ 0</b>	N° 597306 <b>£ 0</b>	N° 597305 <b>£ 0</b>	N° 597306 <b>£ 0</b>	

\*For other configurations please contact us!

Colours		
Finish alternatives	Moderate metallic effect (except for Black & King's Blue)	Special colours (without additional charge, except for Nightblue)
<input type="checkbox"/> Shiny <b>£ 0</b>	<input type="checkbox"/> Zinc Yellow (02) <input type="checkbox"/> Dark Green (15) <input type="checkbox"/> Pure Orange (21) <input type="checkbox"/> Crimson (04) <input type="checkbox"/> Pastel Blue (20) <input type="checkbox"/> Dark Blue (16)	<b>Strong metallic effect:</b> <input type="checkbox"/> Silver Metallic (36) <input type="checkbox"/> Dormant colours: <input type="checkbox"/> Anthracite Metallic (22) <input type="checkbox"/> Apple Green Dormant (18) <input type="checkbox"/> Fire Red Dormant (45) <input type="checkbox"/> Violet Dormant (50)
<input type="checkbox"/> Matt <b>Not possible with special colours!</b>	not metallic: <input type="checkbox"/> Black (31) <input type="checkbox"/> King's Blue (32)	<b>Coarse structured surface:</b> <input type="checkbox"/> Coarse Black (27) <input type="checkbox"/> Dormant col. with additional charge: <input type="checkbox"/> Nightblue Dorm. (47) <b>£65</b>
N° 1561 <b>£75</b>		



Carbon Elements	Design package color anodization	
Complete Carbon package	Colour examples for anodized colours can be found on our website in the area equipment	
<input checked="" type="checkbox"/> Curved back plate, seat plate, side guards and foot plate made of lightweight and highly stable carbon  <i>cf. configuration on the following pages</i>	<input type="checkbox"/> anodized Blue (01) <input type="checkbox"/> anodized Red (02) <input type="checkbox"/> anodized Green (03) <input type="checkbox"/> anodized Gold Orange (04) <input type="checkbox"/> anodized Violet (05)	The following parts are anodized: Assembly group caster (forks, caps), hand rims and push handle bars  <b>The design package ist only possible in the listed setting.</b>
N° 167302 <b>£ 0</b>	N° 1569.. <b>£195</b>	<b>Restrictions:</b> The hand rims are only available in the design package anodized colours by choosing aluminum hand rims. (not possible with stainless steel, Maxgrepp, PVC, double or nob hand rims. Also not possible with profile wheels).

Abduction		Curved back plate	
Strong abduction	Slight abduction	Model 1	Model 3
Strong abduction for maximized leg space and large LLL	Slight abduction with straight frame tubes for short LLL	Straight upper edge, 2 cm higher than back height	Curved upper edge, 5 cm higher than back height


Measuring the lower leg length (LLL)			
Manufacturer's data	User data	Retailer data	Please indicate both measures
Our data for the possible lower leg length are: The maximum or rather minimal possible distance between the <b>top edge of the seat to the top edge of the foot plate MEASURED VERTICALLY (90°)</b> .	You are measuring at the user: The distance between the <b>hollow of the knee to the floor (=LLL)</b> , measured at a <b>vertical leg position of 90°</b>	When ordering, please indicate <b>TWO</b> measures:  <ul style="list-style-type: none"> <li>• <b>1. the real LLL of the user AND</b></li> <li>• <b>2. the thickness of the cushion used by you</b></li> </ul>	<p><b>1. Real lower leg length of the user (LLL):</b></p> <p>_____ cm</p> <p><b>2. Cushion thickness (CT)</b></p> <p>_____ cm</p>


Example		
Our data	Your data	We are assembling
In the order form, we are indicating e.g.: <b>„possible distance from top edge seat to top edge foot plate at internal mounting of leg support = min. 3.5 up to max. 30 cm“.</b>	The real lower leg length of the user amounts e.g.: <b>„LLL = 31 cm“</b> ; The thickness of the cushion used by you amounts e.g.: <b>„CT = 5 cm“.</b>	Using your data, our technicians are calculating a distance of 26 cm between the top edge of the seat to the top edge of the foot plate. Therefore, the foot plate will be mounted to the position of 26 cm, <b>measured vertically</b> . <i>At a cushion thickness of only 3 cm, our technicians would mount the foot plate to a distance of 28 cm from the top edge of the seat to the top edge of the footplate instead.</i>

Seat system		Curved back plate		
Seat plate	Model 1	Model 3	Back angle	
<input checked="" type="checkbox"/> Seat plate made of carbon, seat height, -depth and -inclination adjustable by seat support frame  Required configuration: <input type="checkbox"/> no seat inclination <input type="checkbox"/> slight seat inclination <input type="checkbox"/> high seat inclination	Curved back plate made of carbon, BH + 2 cm;  <input type="checkbox"/> 20-30 mm curve depth, <b>growable!</b>  <i>curve depth depends on the selected seat width</i>	Curved back plate made of carbon, BH + 5 cm;  <input type="checkbox"/> 20-30 mm curve depth, <b>growable!</b>	<input type="checkbox"/> Adjustable from 80° - 100° (in 5°-steps) incl. stabilizer bar <b>Standard configuration: 90°</b>  Desired back angle: <input type="checkbox"/> ____°	
N° 193079	N° 697901 <b>£ 0</b>	N° 697902 <b>£166</b>	N° 697800 <b>£ 0</b>	




Seat cushions (all covers with zippers and changeable)				
	Standard cushion	Anatomical cushion	Cover design	Incontinence cover
<input type="checkbox"/> <b>No seat cushion</b> (please tick if no cushion is required)	Synthetic foam for firm and comfortable seating <input type="checkbox"/> 3 cm (..03) <input type="checkbox"/> 5 cm (..05)  	<input type="checkbox"/> Anatomically formed, with slight abduction, changeable cover (Thickness of the front upholstery, depending on cushion size, is 5-6 cm)  	in selected cover colour (colour choice and colour code see below) -> Choice of cover colour	<input type="checkbox"/> For <b>standard cushion</b> as chosen on the left
	N° 80.. <b>£90</b>	N° 891002 <b>£419</b>	N° 80..	N° 86.. <b>£89</b>
<b>Additional covers</b> for cushions as chosen on the right <input type="checkbox"/> 1 cover (..01) <input type="checkbox"/> 2 covers (..02)	Pressure relief cushion for medium-soft seating <input type="checkbox"/> 3 cm (..03) <input type="checkbox"/> 5 cm (..05)	<input type="checkbox"/> <b>Additional cover</b> for anatomical cushion as chosen	<input type="checkbox"/> Seat cushion as chosen, <b>without cover</b>	<input type="checkbox"/> For <b>anatomical cushion</b> as chosen on the left
N° 82.. <b>£75/pcs</b>	N° 81.. <b>£119</b>	N° 891003 <b>£166</b>	N° 998001 <b>£30</b>	N° 891004 <b>£166</b>




Back cushions (all covers with zippers and changeable)			
No cushion	Standard cushion	Anat. back cushion	Cover design
<input type="checkbox"/> <b>No back cushion</b> (please tick if no cushion is required)	<input type="checkbox"/> Medium-soft cushion, foam plastic, 3 cm, with zipper, changeable cover	<input type="checkbox"/> Anatomically formed, incl. changeable cover, cavity 3 cm  <b>Only with curved back plate</b> <b>Customized product</b>	in selected cover colour (colour choice and colour code see below) -> Choice of cover colour
		N° 891104 <b>£245</b>	N° 85..
<b>Additional covers</b> for cushions as chosen on the right <input type="checkbox"/> 1 cover (..01) <input type="checkbox"/> 2 covers (..02)		<input type="checkbox"/> <b>Additional cover</b> for anatomical back cushion as chosen	<input type="checkbox"/> Back cushion as chosen, <b>without cover</b>
N° 85.. <b>£75/pcs</b>	N° 8517 <b>£139</b>	N° 891110 <b>£135</b>	N° 998501 <b>£30</b>




Choice of Trevira cover colour	
Note: Color views may differ from reality depending on the screen view / printout and are only for initial orientation.	
Please choose your colour of the selected covers (back / seat): <input type="checkbox"/> black (30) <input type="checkbox"/> grey (31) <input type="checkbox"/> nightblue (32) <input type="checkbox"/> ultramarine blue (33) <input type="checkbox"/> petrol (34) <input type="checkbox"/> applegreen (35) <input type="checkbox"/> burgundy (36) <input type="checkbox"/> red (37) <input type="checkbox"/> pink (38)	
	
N° 80.. (seat cushion/anat. cushion)	N° 85.. (back cushion/anat. back cushion)

Leg support	
<b>Model 1</b>	
<input checked="" type="checkbox"/> Foot plate made of carbon, undivided, adjustable in height, depth and angle	
N° 402679	<b>£ 0</b>


To see the latest versions of our order forms, please visit [www.sorgrollstuhltechnik.de/download-portal](http://www.sorgrollstuhltechnik.de/download-portal)

Push handles		
Extendable	Single-hand push handle	Push bow
Short handle (5 cm) 	Foldable 	Adjustable in angle and foldable 

User brakes		Companion brake
Bowden cable brake	Lockable user brake	Locking brake
		

Tyres		Handrim
Schwalbe Right Run	Schwalbe Marathon Plus	Ergonomic handrim
		

To see the latest versions of our order forms, please visit [www.sorgrollstuhltechnik.de/download-portal](http://www.sorgrollstuhltechnik.de/download-portal)

Push handles				Side guard
Standard	Single-hand push handle	Push bow	Standard	
Push handles extendable by <b>straight</b> : <input type="checkbox"/> 25 cm (.02) <input type="checkbox"/> 35 cm (.03)  <b>cranked</b> : <input type="checkbox"/> 25 cm (.06) <input type="checkbox"/> 35 cm (.04)	<input type="checkbox"/> Extendable push handle with <b>LONG</b> handle ( <b>configure extension on the left</b> ) N° 65.. <b>£230</b> <input type="checkbox"/> Extendable push handle with <b>SHORT</b> handle ( <b>configure extension on the left</b> ) <b>Not in combination with drum brake!</b> N° 65..+ N°623702 <b>£265/P</b>	<input type="checkbox"/> Single-hand push handle, cranked, extendable by 22 cm, foldable and adjustable in angle by steel joint <b>Permitted weight capacity max. 40kg</b> <b>Item only available until December 31, 2021.</b> N° 660202 <b>£277</b>	<input type="checkbox"/> extendable by 20 cm, adjustable in angle <b>made of plastic up to SW 30cm</b> at SW20/22 and BH25/30 only limited height adjustment possible N° 653101 <b>£369</b>	<input checked="" type="checkbox"/> Carbon side guard with integrated wheel cover  N° 502579 <b>£ 0/P</b>

User brake	Companion brake	
Bowden cable brake	Lockable user brake	
<input checked="" type="checkbox"/> Bowden cable brake, integrated in activity side guards, with ball knob N° 353601 <b>£ 0</b>	<input type="checkbox"/> Lockable user brake, <b>lockable in activated or deactivated position</b> , integrated in activity side guards, with ball knob N° 353603 <b>£349</b>	
	<th>Locking brake</th>	Locking brake
	<input type="checkbox"/> Locking brake at the back of the rear wheels <b>To lock in the resting state (No driving brake)</b> N° 353602 <b>£289</b>	

Wheel size and wheel camber	Wheels	
Wheel size	Wheel camber	
Size in inches <input type="checkbox"/> 18" (only for lightweight-wheels with standard hub, only up to seat depth 22) <input type="checkbox"/> 20" <input type="checkbox"/> 22" <input type="checkbox"/> 24" <b>£ 0</b>	<input type="checkbox"/> Wheel camber, firmly welded <input type="checkbox"/> 9° <b>£ 0</b> <input checked="" type="checkbox"/> Assembly set camber 7°/11°, preinstalled for <input type="checkbox"/> 7° <input type="checkbox"/> 11° <b>£69</b>	
	<th>Super lightweight wheel</th>	Super lightweight wheel
	<input checked="" type="checkbox"/> 200g/pair lighter than a lightweight wheel, incl. black aluminium handrims N° 21..78 <b>£ 0</b>	

Tyres	Handrim	
Standard	Performance	
Pneumatic: <input type="checkbox"/> Schwalbe <i>Right Run 1"</i> (.13) N° 20.. <b>£ 0</b>	Pneumatic: <input type="checkbox"/> Schwalbe <i>Marathon Plus</i> , pneumatic flat free tyres, 1" <b>Only for 22" and 24" wheels</b> N° 20..12 <b>£65</b>	
	<th>Ergonomic</th>	Ergonomic
	<input type="checkbox"/> tight <input type="checkbox"/> wide handrim with ergonomic handrim profile, black anodized with built-in rubber grip <b>only for 24" wheels</b> N° 342424 <b>£369</b>	



Spoke cover	Design spoke cover decors / Design curved back		
Transparent	Design spoke covers	Design curved back	
<input type="checkbox"/> transparent (.00) N° 25.. <b>£149</b>	<input type="checkbox"/> <a href="http://www.speichen-schutz.de">www.speichen-schutz.de</a> N° 25..26 <b>TBD</b>	<input type="checkbox"/> <a href="http://www.wheelchairgadgets.com">www.wheelchairgadgets.com</a> N° 25..27 <b>TBD</b>	<input type="checkbox"/> Optional: child's name on spoke cover. Name: _____ <b>Only for Design Decor</b> <b>£ 0</b>
			<input type="checkbox"/> Foil for curved back plate with the same design as the specified spoke cover <b>Only with curved back plate</b> N° 2527 <b>TBD</b>

Caster wheels			
Standard PU	Standard solid rubber	Design PU	Design solid rubber

Belt systems			
2-point pelvic belt	4-point pelvic belt	Chest-shoulder belt	Protection cap
with plastic fastener 	with metal fastener 		only for metal fastener 

Standard head rests <i>Measures in mm (all measures ± 10mm)</i>	Lightweight head rest <i>Measures in mm (all measures ± 10mm)</i>												
<b>Model 1</b>  	<div style="border: 1px solid black; padding: 5px; margin-bottom: 10px;">  Not approved as a headrest for transport in a motor vehicle.                     </div> Pad sizes: <table border="1" style="margin: 5px 0;"> <thead> <tr> <th></th> <th>size 1</th> <th>size 2</th> </tr> </thead> <tbody> <tr> <td>height</td> <td>110</td> <td>130</td> </tr> <tr> <td>width</td> <td>200</td> <td>260</td> </tr> <tr> <td>depth</td> <td>30</td> <td>40</td> </tr> </tbody> </table>  		size 1	size 2	height	110	130	width	200	260	depth	30	40
		size 1	size 2										
height		110	130										
width		200	260										
depth	30	40											
<b>Model 3</b>  													
<b>Model 4</b>  													

Frame options			
Anti-tipper	Outdoor Front End	Water bottle	Sun shade
			 Removable sun shade. Easily adjustable in height and position by flexible holder.

Caster wheels				Caster forks	
Standard PU	Standard solid rubber	Design PU	Design solid rubber	Standard Aluminium	Design Aluminium
Polyurethane: <input type="checkbox"/> 4" (2812) <input type="checkbox"/> 5" (2805)  <b>With 18" rear wheels only 4" caster wheels possible</b>	Solid rubber: <input type="checkbox"/> 4" (2804) <input type="checkbox"/> 5" (2809)  <b>With 18" rear wheels only 4" caster wheels possible</b>	Polyurethane, transparent, with LEDs: <input type="checkbox"/> 4" (2820) <input type="checkbox"/> 5" (2834) <b>With 18" rear wheels only 4" caster wheels possible</b>	Solid rubber, black, with alu rim <input type="checkbox"/> 4" (2835) <input type="checkbox"/> 5" (2836) <b>With 18" rear wheels only 4" caster wheels possible</b>	<input type="checkbox"/> Alu fork, black anodized 	<input type="checkbox"/> Alu fork, coated in chosen frame colour 
£ 0	£ 0	£ 0	£ 0	N° 2826 £ 0	N° 2829 £69

### Positioning aids

Here you will find a small selection to mark. You will find more positioning aids, corresponding mounting possibilities and further equipment in the separate order form „Positioning aids“

**Order form „positioning aids“, available in download area.**




**e.g.:**

- belts,
- fixations and corresponding connections,
- abduction wedges,
- side truss pads,
- head rests,
- and more...

**On the order form refer to this new wheelchair order/enquiry.**


2-point pelvic belt	4-point pelvic belt	Standard head rests	Lightweight head rest
<input type="checkbox"/> with plastic fastener <b>£119</b> <input type="checkbox"/> with metal fastener <b>£139</b> <b>The size is determined by us according to the chosen wheelchair measurements.</b> N° 780502	<input type="checkbox"/> with plastic fastener <b>£149</b> <input type="checkbox"/> with metal fastener <b>£179</b> <b>The size is determined by us according to the chosen wheelchair measurements.</b> N° 780504	<b>All head rests including connection</b> <input type="checkbox"/> Model 1, (...21) shell-shaped small <input type="checkbox"/> Model 3, (...23) curved small <input type="checkbox"/> Model 4, (...24) curved large N° 6860..	Lightweight head rest, adjustable with two conical joints, connection and removable cover included. (cover with same Betex-color as back cushion) <input type="checkbox"/> Size 1 (...20) Pad size 110x200mm <input type="checkbox"/> Size 2 (...21) Pad size 130x260mm  <b>Weight reduction approximately 700-1000g to standard head rest systems</b> <b>Not approved as a headrest for transport in a motor vehicle.</b> N° 686009..
<b>Chest-shoulder belt</b> <input type="checkbox"/> Chest-shoulder belt without zipper N° 7810	<b>Protection cap</b> <input type="checkbox"/> protection cap for metal fastener <b>Makes independent opening for the patient difficult.</b> N° 78050210	<b>£177</b>	<b>£15</b>
<b>£266</b>	<b>£269</b>		

### Frame options

Anti-tipper	Outdoor Front End	Splash guard	Water bottle	Sun shade
Swing-away installation  Mounting <input type="checkbox"/> left <input type="checkbox"/> right <input type="checkbox"/> pair (from SW 22)	<input type="checkbox"/> Outdoor Front End, facilitates significantly active and passive driving on soft or rough surfaces, incl. mounting aid	<input type="checkbox"/> Splash guard for Outdoor Front End, made of plastic  	<input type="checkbox"/> Water bottle 750ml light PET-bottle, incl. holder N° 790101 <b>£35</b>  <b>Bumper</b> <input type="checkbox"/> Bumper, made of rubber to protect furniture, doors, etc. frame form „Taurus“ not possible with Restraint system 	<input type="checkbox"/> variable adjustable sun shade (size 1), in black  <b>Usable height (from the top edge of the seat plate/seat cover to the apex) = selected back height + 40 cm</b> Example: BH 30 + 40 cm = 70 cm 
N° 7213 <b>£149 /pcs</b>	N° 599517 <b>£399</b>	N° 59951212 <b>£65</b>	N° 9102 <b>£51</b>	N° 910301 <b>£220</b>

## Seat heights

**All specifications measured WITHOUT cushions!**  
**All specifications ±5 %**



Combination of wheels	Seat height front measured without cushion from seat edge to floor at horizontally seat				
	Height 1	Height 2	Height 3	Height 4	Height 5
18" x 4" (in 18" fork)	32 cm	33,5 cm	35 cm	-	-
20" x 4" (=B2)	34 cm	35,5 cm	37 cm	38,5 cm	40 cm
20" x 5" (=B1)	34 cm	35,5 cm	37 cm	38,5 cm	40 cm
22" x 4" (=B2)	36 cm	37,5 cm	39 cm	40,5 cm	42 cm
22" x 5" (=B1)	36 cm	37,5 cm	39 cm	40,5 cm	42 cm
24" x 4" (=B2)	38,5 cm	40 cm	41,5 cm	43 cm	44,5 cm
24" x 5" (=B1)	38,5 cm	40 cm	41,5 cm	43 cm	44,5 cm

# Terms and Conditions

## 1. Preface and Scope

(1) Our terms and conditions apply exclusively. Any deviating, opposing, or complementing terms and conditions on the part of the business partner do not form integral part of contract unless explicitly approved of by us. There are no exceptions to the consent requirement, which also applies in cases where we conduct deliveries without restrictions while aware of the business partner's terms and conditions.  
(2) Agreements made in the individual case with the purchaser (including side-agreements, supplements, and alterations) always take priority over these terms and conditions. A contract or our confirmation - both in written form - are authoritative for the content of such arrangements.  
(3) Any legally relevant declarations and notifications the purchaser submits to us after conclusion of the contract (e.g. setting of deadlines, notice of defects, declarations of rescission or abatement) require the written form to be valid.  
(4) Any information on the application of legal regulations serves the mere purpose of clarification. All regulations are legally binding without any such clarifications unless explicitly altered or excluded in these terms and conditions.

## 2. Prices

All prices are ex works, excluding insurance and packaging. The prices are net prices; the value-added tax is stated and charged separately according to the legally binding rate at the time of invoicing. The prices valid on the day of delivery are binding unless explicitly arranged otherwise.

## 3. Conclusion of Contract

SORG's offers are invitatio ad offerendum (invitations to treat). Ordering goods counts as an obligatory offer of contract on the part of the purchaser. Unless the order suggests otherwise, we have the right to accept this contract offer within two weeks time after reception. The acceptance of an offer of contract is effected by the confirmation of the order.

## 4. Terms of Payment

(1) Provided that no other agreements exist, payment is to be made by paying in advance - 50% at placing of order, 50% pre-delivery of goods. The purchaser bears any costs arising from transactions or discounting.  
(2) If a 30-day payment deadline is exceeded, a delay rate of 9% on top of the respective base rate of the European Central Bank falls due without the need for prior notices. Beyond that, the purchaser is charged 3,- € for any reminder sent. We reserve the right to make use of the commercial due rate according to § 353 HGB. A higher damage for delay may be enforced via a suitable verification.

## 5. Offers and Validity

All current price lists are generally valid until December 31 of the respective year, or until new price lists are published. The latest version of the price list will always be provided on our homepage. Offers older than 3 months at that point expire and must be replaced with new offers and prices.

## 6. Time of Delivery and Delay

We state the terms of delivery when accepting the order. The delivery date is non-binding and has mere guideline value - especially with regard to customised products. The terms for delay in delivery are defined by legal regulations. In any case, a writing notice issued by the purchaser is necessary for them to take effect.

## 7. Force Majeure (Act of God)

(1) SORG is not to be held responsible for the following delays of delivery. This also applies if these causes occur with respect to legal representatives, auxiliary persons, or deliverers of SORG: Circumstances of force majeure as well as other extraordinary events that are unpredictable and unavoidable for, and not caused by, SORG, and that only occur after ordering or reserving, and that SORG could not have known about. Moreover, this applies in case of strikes, fire, floods, labour disputes, business disruptions, altered magisterial approbation and legal situation, and magisterial orders which are not due to operational risks.  
(2) In the cases of paragraph 1, SORG is eligible to delay delivery/performance for the time of hindrance as well as an appropriate start-up time. If SORG thus causes unacceptable complications for the purchaser, both SORG and the purchaser are eligible to withdraw from the contract. Circumstances that only lead to a temporal and hence acceptable delay of delivery and performance are excluded from this right. The legal rights of the client remain unaffected.

## 8. Dispatch

We generally dispatch uncovered by insurance and at the recipient's risk. We transfer the risk of coincidental loss or deterioration of the goods as well as delay on to the client at the moment of passing the goods on to the forwarder, bearer, or any other person/institute responsible for delivery. SORG will make any effort to respect the client's wishes in terms of path and kind of delivery. Any costs caused by this (including agreed carriage free delivery) are to be covered by the client. In the case of carriage free delivery conducted by SORG, the recipient is required to put forward any carriage costs for a refund. If the client wishes the dispatch to be delayed, any risk of loss or deterioration of the goods is carried by the recipient from the first day of readiness for dispatch. The goods dispatched are only insured upon special demand and at the cost of the recipient.

9. Damage through Delivery  
Angelieferte Gegenstände sind, auch wenn sie unwesentliche Mängel  
Upon delivery, the goods need to be checked for obvious damage in the presence of the bearer. Any damage is to be documented in a written form and co-signed by the bearer. This protocol is to be forwarded to us immediately. Delayed reclamations, or reclamations not signed by the bearer, cannot be taken into account.

## 10. Warranted Features, Receipt, and Technical Alterations

Delivered objects, even with minor damage, are to be accepted by the client, even if §433 paragraph 1 BGB is not fulfilled. Partial deliveries are permissible. Features only count as warranted if explicitly designated as such. A warranty of condition can only be asserted explicitly and in a written form. SORG reserves the right for technical alterations, especially for solutions of equal or higher value.

## 11. Warranty for Defects

(1) SORG warrants for any defects of delivered new goods for 12 months as long as the purchaser is not the user (§§ 474, 13 BGB). Should the purchaser be the user then the warranty period is two years for new goods, for used goods 12 months. Should the purchaser also be an entrepreneur, then the liability for used items is void, as long as there is no claim for compensation according to N° 12.  
(2) The purchaser's claim for defects presuppose that he has attended his obligation to give notice of defects on the legal grounds for his reasons of inquiry (§§ 377, 381 HGB). If defects are found during the inspection or at a later point, we are to be notified in writing without delay.

The notice counts as not delayed if issued within two weeks. Posting the notice on time suffices for the keeping of the term. Independently from this obligation to inspect and give notice, the purchaser is to report obvious defects (including wrong and short delivery) within two weeks upon delivery (again, posting the notice on time suffices). If the purchaser misses the proper inspection and/or notice of defects, we do not warrant for the reported defects. If the delivered good is detected, we may choose to conduct supplementary performance either by repairing the defect (subsequent improvement) or by delivering flawless goods (compensation delivery). Our right to refuse the chosen method of supplementary performance on the respective legal grounds remains untouched.  
(3) Our liability for defects expires when alteration or maintenance works are done by the client or a third-person-party without our prior written permission. Consequences of excessive strain, negligent or incorrect handling, and violent damage are excluded from the warranty for defects.

## 12. Further Liability

(1) With respect to compensation - regardless of the legal cause - SORG is only liable in case of malice and act of gross negligence. In case of slight negligence, SORG is only liable for damage of life, body, health, and for damage caused by the negligence of an essential contractual obligation (an obligation whose fulfilment make the proper realisation of the contract possible in the first place and in whose fulfilment the contract partner regularly trusts and can trust). In this case, however, our liability does not exceed the covering of predictable, typically occurring damage.  
(2) A free right to termination on the part of the purchaser (acc. to § 651, 649 BGB in particular) is precluded. Otherwise, the usual prerequisites and legal consequences apply.

## 13. Costs of Supplementary Performance

Costs of supplementary performance cannot be charged on SORG's account.

## 14. Voluntary Withdrawal of Goods

SORG is not obliged to take back any delivered end products or parts thereof. If SORG voluntarily decides otherwise in individual cases (especially for therapeutically altered requirements or deaths), SORG refunds a maximum of 75% of the delivery costs for unused parts whose delivery does not date back more than three months. There is no legal obligation concerning withdrawal in this respect. Deliveries which date back more than three months as well as sewed parts and customised products are not returnable. A copy of the delivery note or the bill serves as proof of age of delivery.

## 15. Reservation of Ownership

Goods delivered by SORG remain the property of SORG until any obligations from the parties' business contract on the part of the purchaser have been redeemed. This also applies if SORG has adapted individual or all claims into a running account which has been balanced and acknowledged. The purchaser has the right to dispose of the goods in the proper course of business. The purchaser hereby yields all claims to SORG which arise from the resale against the purchaser or third parties. The purchaser is authorised to withdraw these claims even upon conveyance. SORG's right to withdraw the claims remains untouched. However, SORG is obliged not to withdraw the claims as long as the purchaser complies with his duty of payment and any other duties. SORG has the right to demand that the purchaser makes known the yielded claims and their debtors, that he provides all information necessary for conveyance, and that he hands in all associated documents. As long and insofar the reservation of ownership is intact, the purchaser must not transfer any goods for securing or bailment without SORG's permission. The purchaser may process or alter the goods for the benefit of SORG exclusively. If processed with other goods which do not belong to SORG, SORG obtains co-ownership of these new goods in relation to the invoice value to the cost price of the other processed goods. If a third party interferes with SORG's ownership, the purchaser is to notify SORG immediately. The purchaser must not make any agreements with other purchasers which may exclude or limit SORG's rights. In particular, the purchaser must not make any agreements which destroy or limit the advance assignments to SORG. If the estimated value of the existing guarantees exceeds the claims to be ensured by more than 20%, SORG is insofar obliged to release by request of the purchaser.

## 16. Copyright Protection

If SORG conveys any images, drawings, indications of weight or measurements, technical documentations, or other product descriptions and documents (including electronic ones), the respective indications only apply if explicitly stated that they do. SORG reserves ownership and copyright of all these documents. They may only be made available to third parties upon SORG's explicit permission.

## 17. Place of Fulfilment, Place of Jurisdiction, and Severability Clause

(1) The law of the Federal Republic of Germany applies for these Terms and Conditions and all legal relations between us and our business partners. All international and supranational (contractual) legal systems, the UN Convention on Contracts for the International Sale of Goods in particular, are hereby excluded. Prerequisites and effects of the reservation of ownership acc. to no. 15 are subject to the right of the respective storage location of the good if the choice of law made in favour of the German law is illegitimate or ineffective. The only place of jurisdiction (also on an international basis) for any disputes arising directly or indirectly from the contractual relationship is, according to German Commercial Code, our office in Oberhausen-Rheinhausen, Germany. However, we also have the right to file a suit concerning the general jurisdiction of the purchaser.  
(2) If any clause in these Terms of Conditions should turn out to be ineffective, this does not affect the effectiveness of all other clauses.

Technical alterations and misprints reserved. As of: January 2017  
SORG Rollstuhltechnik GmbH + Co. KG  
Benzstr. 3-5, 68794 Oberhausen-Rheinhausen / Germany  
Phone +49 7254 / 9279-0  
Fax +49 7254 / 9279-10  
info@sorgrollstuhltechnik.de, www.sorgrollstuhltechnik.de