

Activity wheelchair with mechanical tilting, suitable for children and adolescents

Included in the base price:

- Frame versions Retro or Taurus
- 2 different V-shaped abductions for maximized leg space
- Growable back height + 5 cm (optional: Seat width +2 cm, Seat depth + 3 cm)
- mechanical seat tilting from -5° to +35° in the physiological turning point
- back adjustable in angle of 78° -126° and folds down to the front
- Also included: seat plate - tilting in physiological turning point; stabilizer bar; Standard-, Profile- or Lightweight wheels; Caster wheels with aluminum fork; anti-tipper left and right; Cable brake built in to the side guard



Quotation Order

Weight capacity: 50 kg
German aid index N°: 18.50.03.2005

Base price £3,349

Prices excl. VAT, valid from 01.03.2021

Delivery address / Company stamp:

Name:

Address:

Billing address: same as delivery address

Name:

Address:

Ordering information:

Reference/
internal description:

Date:

Contact person:

Signature:

| Frame versions | | | | |
|--|--|--|--|-------------------------|
| Retro frame N° 597301 | | Taurus frame N° 597304 | | Retro frame |
| <input type="checkbox"/> Strong abduction LLL 15-34 cm | <input type="checkbox"/> Slight abduction LLL 15-31 cm | <input type="checkbox"/> Strong abduction LLL 15-34 cm | <input type="checkbox"/> Slight abduction LLL 15-31 cm | Taurus frame |
| casters on one line with rear wheels | casters narrower than rear wheels | casters on one line with rear wheels | casters narrower than rear wheels | |
| N° 597305 £ 0 | N° 597306 £ 0 | N° 597305 £ 0 | N° 597306 £ 0 | |

| Frame- and user measures | | | cf. example on p.6 |
|--|---|---|--|
| Seat width in cm (10..) | Seat depth in cm (11..) | Back height in cm (12..) | |
| <input type="checkbox"/> 20 <input type="checkbox"/> 22 <input type="checkbox"/> 24 <input type="checkbox"/> 26 <input type="checkbox"/> 28 <input type="checkbox"/> 30 <input type="checkbox"/> 32 <input type="checkbox"/> 34 | <input type="checkbox"/> 18 <input type="checkbox"/> 20 <input type="checkbox"/> 22 <input type="checkbox"/> 24 <input type="checkbox"/> 26 <input type="checkbox"/> 28 <input type="checkbox"/> 30 <input type="checkbox"/> 32 <input type="checkbox"/> 34 <input type="checkbox"/> 36 | 04 53 03 52 54 £160 (Back height growable by +5 cm each) we measure: top edge seat board to top edge back rod | 1. Real lower leg length of the user (LLL): _____ cm |
| we measure: inner distance side guard to side guard | we measure: front edge seat to front edge back rod | | 2. Cushion thickness (CT) _____ cm |

| Colours | | |
|---|---|---|
| Finish alternatives | Moderate metallic effect (except for Black and King's Blue) | Special colours (without additional charge, except for Nightblue) |
| <input type="checkbox"/> Shiny £ 0 | <input type="checkbox"/> Zinc Yellow (02) <input type="checkbox"/> Dark Green (15) <input type="checkbox"/> Pure Orange (21) <input type="checkbox"/> Crimson (04) <input type="checkbox"/> Pastel Blue (20) <input type="checkbox"/> Dark Blue (16) <input type="checkbox"/> Brilliant Blue (08) <input type="checkbox"/> Cream White (17) not metallic: <input type="checkbox"/> Girly Pink (19) <input type="checkbox"/> Black (31) <input type="checkbox"/> Chrome Yellow (06) <input type="checkbox"/> King's Blue (32) | Strong metallic effect: Silver Metallic (36) A Anthracite Metallic (22) Coarse structured surface: C Coarse Black (27) Dormant colours: Apple Green Dormant (18) Fire Red Dormant (45) Violet Dormant (50) Dormant col. with additional charge: Nightblue Dorm. (47) £65 |
| N° 1561 £75 | | |
| Design package color anodization (Colour examples for anodized colours can be found on our website in the area equipment) | | |
| <input type="checkbox"/> anodized Blue (01) <input type="checkbox"/> anodized Red (02) <input type="checkbox"/> anodized Green (03) <input type="checkbox"/> anodized Gold Orange (04) <input type="checkbox"/> anodized Violet (05) N° 1570.. £195 | The following parts are anodized: Caster forks, supporting plate for side guards, hand rims The design package ist only possible in the listed setting. | Restrictions: The hand rims are only available in the design package anodized colours by choosing aluminum hand rims. (not possible with stainless steel, Maxgrepp, PVC, double or nob hand rims. Also not possible with profile wheels) |

| Frame options | | | |
|--|---|---|--|
| Adaption set | Wheel base extension | Anti-tipper | Step tube |
| <input type="checkbox"/> Modification set for growable seat width + 2 cm and seat depth up to 3 cm | <input type="checkbox"/> by 4 cm to the back, and adjustment of axis, side guards and brake Attention: Return to standard configuration only possible by new side guards! Not in combination with double handrim. | <input checked="" type="checkbox"/> left and right, swing away installation, adjustable in height | <input type="checkbox"/> Step tube to make tilting more easy |
| N° 597799 195,- € | N° 760177 250,- €/P | N° 7214 £ 0 | N° 7408 70,- € |

| Seat system | | Back system | | |
|--|--|--|---|--|
| Alu seat plate | Back angle | Model 1 | Model 2 | Model 3 |
| <input checked="" type="checkbox"/> screwed on, in the physiological turning point tiltable from -5° to +35° N° 192800 £ 0 | <input checked="" type="checkbox"/> Adjustable in angle by rafter from 81°, 90°, 102°, 114° oder 126°, folds down forward, incl. stabilising rod N° 612105 £ 0 | BH+2 cm; curve depth: <input type="checkbox"/> 15 mm (..015) <input type="checkbox"/> 30 mm (..030) <input type="checkbox"/> 50 mm (..050) Attention: curve depth 50mm only from seat width 30cm Straight upper edge, 2 cm higher than back height N° 693.. £275 | BH+10 cm; curve depth: <input type="checkbox"/> 15 mm (..115) <input type="checkbox"/> 30 mm (..130) <input type="checkbox"/> 50 mm (..150) BH+20 cm; curve depth: <input type="checkbox"/> 15 mm (..215) <input type="checkbox"/> 30 mm (..230) <input type="checkbox"/> 50 mm (..250) Attention: curve depth 50mm only from seat width 30cm Straight upper edge, 10 cm higher than back height N° 693... £355 | BH+5 cm; curve depth: <input type="checkbox"/> 15 mm (..315) <input type="checkbox"/> 30 mm (..330) <input type="checkbox"/> 50 mm (..350) BH+10cm; curve depth: <input type="checkbox"/> 15 mm (..415) <input type="checkbox"/> 30 mm (..430) <input type="checkbox"/> 50 mm (..450) Attention: curve depth 50mm only from seat width 30cm Curved upper edge, 5, 10 or 15 cm higher than back height N° 693... £355 |
| Blocking set | | | | |
| Mounting set to limit tilting: towards the front to max.: <input type="checkbox"/> 0° tilting towards the back to max.: <input type="checkbox"/> +30°, <input type="checkbox"/> +25°, <input type="checkbox"/> +20°, <input type="checkbox"/> +15° tilting N° 192810 £55 | | | | |

| Back extension | Restrictions | Back cover | Fitting set |
|---|---|--|---|
| for curved back | | Adjustable | for solid back |
| <input type="checkbox"/> extension up to 17.5 cm, retractable, with detachable head cushion from SW 22 cm, weight reduction of approx. 1000g compared to a head rest, not possible with head rest (N° 6860..) (Back view) N° 686077 £260 | -Model 1: from BH 35 cm -Model 2: Not possible -Model 3: from BH 30 cm only with version BH+5 cm -Fitting set N° 692598: from SW 30 cm | <input type="checkbox"/> incl. back cushion (thickness 2 cm) Breathable cover, colour selection on page 3 Possible from SW24 N° 6701.. £ 0 | <input type="checkbox"/> one set of clamps, required for mounting a solid back that is manufactured by the specialist trade Curve depth 50mm only from seat width 30cm N° 692598 £355 |

| Seat cushions (All covers with zippers and changeable) | | | | |
|---|---|--|---|--|
| | Standard cushion | Anatomical cushion | Cover design | Incontinence cover |
| <input type="checkbox"/> No seat cushion (please tick if no cushion is required) | Synthetic foam for firm and comfortable seating <input type="checkbox"/> 3 cm (..03) <input type="checkbox"/> 5 cm (..05) N° 80.. £90 | <input type="checkbox"/> Anatomically formed, with slight abduction, changeable cover (Thickness of the front upholstery, depending on cushion size, is 5-6 cm) N° 891002 £420 | The colour of the chosen covers can be selected on page 3. -> choice of Trevira cover colour | <input type="checkbox"/> For standard cushion as chosen on the left |
| Additional covers for cushions as chosen on the right <input type="checkbox"/> 1 cover (..01) <input type="checkbox"/> 2 covers (..02) | Pressure relief cushion for medium-soft seating <input type="checkbox"/> 3 cm (..03) <input type="checkbox"/> 5 cm (..05) | <input type="checkbox"/> Additional cover for anatomical cushion as configured Customized product N° 891003 £160 | <input type="checkbox"/> Seat cushion as chosen, without cover N° 998001 £30 | <input type="checkbox"/> For anatomical cushion as chosen on the left N° 86.. £89 |
| N° 82.. £75/pcs | N° 81.. £120 | N° 891003 £160 | N° 998001 £30 | N° 891004 £160 |

| Back cushions (All covers with zippers and changeable) | | |
|--|---|--|
| Without cushion | Standard cushion | Anatomical back cushion |
| <input type="checkbox"/> No back cushion (please tick if no cushion is required) | <input type="checkbox"/> Standard cushion for curved back plate , 3 cm thick, with zipper, breathable cover (choice of Trevira cover colour below) | Anatomically formed, including cover, changeable, depr. depth: <input type="checkbox"/> 3 cm (..04) <input type="checkbox"/> 5 cm (..06) <i>Only in association with curved back.</i> |
| <input type="checkbox"/> Back cushion as chosen, without cover | N° 8517 £140 Additional covers for cushions as chosen <input type="checkbox"/> 1 cover (..01) <input type="checkbox"/> 2 covers (..02) | N° 8911.. £250 <input type="checkbox"/> Additional cover for anatomical back cushion as configured Customized product |
| N° 998501 £30 | N° 85.. £75/pcs | N° 891110 £135 |

| Side guard |
|--|
| Standard |
| <input checked="" type="checkbox"/> Component of synthetic and aluminum parts with wheel cover |
| |
| N° 5025 £ 0 |

| Choice of Trevira cover colour | |
|---|---|
| Please choose your colour of the chosen covers: <input type="checkbox"/> black (30) <input type="checkbox"/> grey (31) <input type="checkbox"/> night blue (32) <input type="checkbox"/> ultramarine blue (33) <input type="checkbox"/> petrol (34) <input type="checkbox"/> applegreen (35) <input type="checkbox"/> burgundy (36) <input type="checkbox"/> red (37) <input type="checkbox"/> pink (38) | |
| | |
| Note: Color views may differ from reality depending on the screen view / printout and are only for initial orientation. | |
| N° 80.. (seat cushion/anat. cushion) / | N° 85.. (back cushion/anat. back cushion) / N° 6701.. (adjustable back cover) |

| Leg support with foot rest | | | |
|--|--|---|--|
| Model 1 | Model 2 | Calf rest | Heel cups / Heel edge |
| <input type="checkbox"/> Leg support with undivided foot rest. Footplate can be folded up to the rear, angle and height adjustable | <input type="checkbox"/> Leg support with divided foot rest, steplessly adjustable up to 180°, Foot plates can be folded up to the rear, adjustable in height Possible from SW 22 cm and LLL 23 cm - 35 cm | for leg support 401905, Mounting calf rest <input type="checkbox"/> left <input type="checkbox"/> right <input type="checkbox"/> double-sided Possible from LLL 25 cm, foot plates are no longer hingeable backwards. From SW 26 cm | Heel cups , self installation, black <input type="checkbox"/> Size 1 (..01) 7 x 11.5 cm <input type="checkbox"/> Size 2 (..02) 9 x 12.5 cm <input type="checkbox"/> Size 3 (..03) 10 x 15.5 cm |
| N° 401901 £ 0 | N° 401905 £95 | N° 403106 £125 | N° 4020.. £50/P |
| <input type="checkbox"/> Locking mechanism at 90° for chosen foot plate | | | Heel edge <input type="checkbox"/> for foot plate/s |
| N° 4009 £55 | | | N° 401904 £55/pcs |

| Push handles | | | |
|--|--|--|--|
| Standard | Single hand push handle | Push bow made of plastic or metal | |
| extendable, cranked: <input type="checkbox"/> 25 cm (..06) <input type="checkbox"/> 35 cm (..04) | <input type="checkbox"/> Foldable with steel joint, extendable by 20 cm only up to SW30 cm and nominal load max. 40kg Item only available until December 31, 2021. | <input type="checkbox"/> extendable by 20 cm, adjustable in angle made of plastic up to SW 30cm at SW20/22 and BH25/30 only limited height adjustment possible | <input type="checkbox"/> extendable by 20 cm, adjustable in angle made of metal, from SW 32cm at BH25/30 only limited height adjustment possible |
| N° 65.. £230 | N° 660202 £275 | N° 653101 £370 | N° 653102 £370 |

| User brakes | | Companion brakes | |
|---|--|---|--|
| Bowden cable brake | Lockable user brake | Drum brake ¹ | Locking brake ² |
| <input type="checkbox"/> Bowden cable brake, integrated in activity side guards | <input type="checkbox"/> Bowden cable brake, lockable in activated or deactivated position , integrated in activity side guards | <input type="checkbox"/> Drum brake for standard-, profile- or lightweight wheels (N° 21..30/31/82), incl. brake lever for the accompanying person If Bowden cable brake for the user is additionally required please tick also Article N° 3532. Not with 7° wheel camber / 9° wheel camber recommended | <input type="checkbox"/> Locking brake at the back of the rear wheels To lock in the resting state (No driving brake) Not in combination with drum brake (N° 220036) |
| N° 3532 £ 0 | N° 353603 £335 | N° 220036 £285 | N° 353602 £285 |

¹ (Drum brake) When mounting the drum brake in combination with an 11° camber, wheel sizes 22" and 24" may be restricted when the seat width grows.

² (Locking brake) When mounting the locking brake in combination with lateral supports, the back angle adjustment or tilting can be restricted.

Specification of the lower leg length (LLL)

| Measurements | Manufacturer information | Patient information |
|--|--|--|
| When ordering, please state these TWO dimensions: 1. Real lower leg length of the user(LLL): _____ cm 2. Cushion thickness (CT): _____ cm | <p>Our information on the possible lower leg length are: The maximum or minimum possible distance between the upper edge of the seat and the upper edge of the footplate MEASURED vertically (90 °).</p> <p><i>In the order form we give e.g. on:</i> „Possible distance upper edge seat to upper edge footplate = min. 23.5 to max. 32cm“.</p> | <p>You measure on the patient: The distance between the hollow of the knee and the floor (= LLL) with the legs standing upright.</p> <p>Your details: The real lower leg length of the patient is for example: „LLL = 31 cm“; The thickness of the seat cushion you are using is for example: „CT = 5 cm“.</p> |

How we assemble: Our fitters use your information to calculate a distance between the upper seat and the upper base of 26cm. The footplate is therefore mounted **VERTICALLY MEASURED** at a position of 26cm. With a seat cushion thickness of only 3cm, you would instead mount the footplate at a distance of 29cm from the upper edge seat to the upper edge footplate.

Selection of Seat height / Wheel size

| Wheel size | Seat height (in cm) | | | | | | Seat height measured WITHOUT CUSHION from seat edge to floor with horizontal seat. |
|---------------------------------|--------------------------|---|--------------------------|---|--------------------------|---|---|
| | 37,5 | 39 | 39,5 | 41 | 42 | 43,5 | |
| 20" (recommended for SD18 - 28) | <input type="checkbox"/> | <input checked="" type="checkbox"/> <small>N° 19280700 AC 69,- €</small> | | | | | <input type="checkbox"/> |
| 22" (recommended for SD22 - 32) | | | <input type="checkbox"/> | <input checked="" type="checkbox"/> <small>N° 19280700 AC 69,- €</small> | | | <input checked="" type="checkbox"/> |
| 24" (recommended for SD26 - 36) | | | | | <input type="checkbox"/> | <input checked="" type="checkbox"/> <small>N° 19280700 AC 69,- €</small> | <input type="checkbox"/> |
| Possible Caster wheel sizes | 4"+ 5" | 4"+ 5" | 4"+ 5" | 4"+ 5" | 4"+ 5" | 4"+ 5" | |

= Standard seat height for chosen wheel size.

= Seat height for chosen wheel size with increase set. (69,- € ADDITIONAL CHARGE for increase set N° 19280700)

= Not possible




Wheels





| Standard wheel | Profile wheel | Lightweight wheel | Lightweight wheel | Wheel camber |
|---|---|---|--|---|
| <input type="checkbox"/> All-round wheel Additional weight of approx. 900g/wheel compared to lightweight wheels | <input type="checkbox"/> Unity of handrim and rim, reduces the risk of finger injuries Additional weight of approx. 800g/wheel compared to lightweight wheels | <input type="checkbox"/> with standard hub , very lightweight, incl. aluminium handrim Weight reduction of approx. 900g/wheel compared to standard wheel | <input type="checkbox"/> with drum brake hub , very lightweight, incl. aluminium handrim Weight reduction of up to 1,1 kg/wheel compared to standard drum brake wheels. | <input type="checkbox"/> 9° (Standard) £ 0 Mounting set wheel camber 7°/11°, preassembled on <input type="checkbox"/> 7° <input type="checkbox"/> 11° |
| N° 21..30 £ 0 | N° 21..31 £ 0 | N° 21..77 £ 0 | N° 21..82 £ 0 | N° 201518 £ 70 |

¹ (Wheel camber) When mounting the drum brake in combination with an 11° camber, wheel sizes 22" and 24" may be restricted when the seat width grows.

Wheels with double handrim / One-hand drive

| Hub | Opposite side | Distance | Silicone cover for double handrim | PVC coating for double handrim |
|---|--|--|---|---|
| with standard hub and double handrim, <i>only up to 9° camber</i> Mounting <input type="checkbox"/> left <input type="checkbox"/> right | <input type="checkbox"/> Handrim for the opposite side | Distance of double handrim <input type="checkbox"/> standard 18mm (..11) <input type="checkbox"/> tight 9mm (..12) <i>Distance tight: Not in combination with hand rim cover or coated hand rim.</i> | Detachable silicone cover, black <input type="checkbox"/> for inside rim <input type="checkbox"/> for outside rim <input type="checkbox"/> for opposing side | <input type="checkbox"/> Handrim coating for double handrim, both handrims (inside and outside), coated with black PVC N° 262402 £315 <input type="checkbox"/> Handrim coating for opposing side of the double handrim, coated with black PVC |
| N° 2614 £949 with drum brake hub and double handrim Mounting <input type="checkbox"/> left <input type="checkbox"/> right | | | | |
| N° 2615 £1050 | N° 261435 £140 | N° 26.. £ 0 | N° 3205 £70/pcs | N° 262401 £175 |



| Tyres 20"-24" | | | | | |
|--|-----|---|-----|--|---------|
| Standard pneumatic | | Standard solid | | Performance | |
| Pneumatic: <input type="checkbox"/> Schwalbe <i>Right Run 1"</i> , black | | Solid: <input type="checkbox"/> Standard tyres 1", black | | Pneumatic flat free tyres 1", black <input type="checkbox"/> Schwalbe <i>Marathon Plus 1"</i> | |
|  <p>Similar to illustration</p> | |  <p>Similar to illustration</p> | |  <p>Similar to illustration</p> | |
| N° 20..13 | £ 0 | N° 23..11 | £ 0 | N° 20..12 | £70 / P |




| Standard handrims and Mounting <small>Please always indicate the desired type of installation</small> | | | | | |
|---|--|---|--|---|---|
| at standard wheel | at lightweight wheel | Standard / tight | Standard / wide | High / tight | High / wide |
| Mounting: <input type="checkbox"/> Standard / tight <input type="checkbox"/> Standard / wide <input type="checkbox"/> High / tight <input type="checkbox"/> High / wide Material: <input type="checkbox"/> Stainless steel (..01) <input type="checkbox"/> Aluminum (..75) 20" with stainless steel only possible with mounting „high“. | <input type="checkbox"/> Standard / tight <input type="checkbox"/> Standard / wide <input type="checkbox"/> High / tight <input type="checkbox"/> High / wide Material: Aluminum High mounting is not possible at lightweight wheels with drum brakes. |  <p>13mm</p> |  <p>23mm</p> |  <p>13mm</p> |  <p>23mm</p> |
| N° 30.. | £ 0 | N° 30..75 | £ 0 | | |

| More handrims | | | | Handrim options | |
|---|--|--|---|--|--|
| Ergonomic handrims | Design handrims | PVC | Maxgrepp handrims | Silicone | Ergo rubber inlay |
| <input type="checkbox"/> handrim with ergonomic handrim profile, black anodized with built-in rubber grip Only for 24" wheels | <input type="checkbox"/> Aluminium handrims coated in frame colour | <input type="checkbox"/> Handrim coated with black PVC, (<i>Tetra Grip</i>) Not for profile wheels | <input type="checkbox"/> Allows pushing with the palm of hand (soft surface) Not for profile wheels | Detachable silicone cover, black (soft surface) <input type="checkbox"/> for standard wheel (20-24) <input type="checkbox"/> for profile wheel (..02) <input type="checkbox"/> for lightweight wheel (..05) | <input type="checkbox"/> Rubber inlay, 13mm wide Only for profile wheels |
| N° 342424 | £370 | N° 303000 | £48 | N° 33.. | £325 |
| | | N° 34.. | £400 | N° 32.. | £140 |
| | | N° 32..01 | £180 | | |



| Spoke cover | Design Decors (Net prices are not discountable) | | | | |
|---|---|---|---|-----------|-----|
| Transparent | Please specify the desired model of the corresponding manufacturer: | | | | |
| <input type="checkbox"/> transparent (..00) | <input type="checkbox"/> www.speichen-schutz.de _____ child's name on spoke cover | <input type="checkbox"/> www.wheelchairgadgets.com _____ child's name on spoke cover | <input type="checkbox"/> Foil for curved back plate with the corresponding design as chosen Only with curved back plate | | |
| N° 25.. | £146 | N° 25..26 | TBD | N° 25..27 | TBD |
| | | N° 2527 | £110 | | |

| Caster wheels | | | |
|--|---|---|---|
| Standard PU | Standard solid rubber | Design PU | Design solid rubber |
| Polyurethane, black: <input type="checkbox"/> 4" (2812) <input type="checkbox"/> 5" (2805) | Solid rubber, grey: <input type="checkbox"/> 4" (2804) <input type="checkbox"/> 5" (2809) | Polyurethane, transparent, with LEDs: <input type="checkbox"/> 4" (2820) <input type="checkbox"/> 5" (2834) | Solid rubber, black, with alu rim <input type="checkbox"/> 4" (2835) <input type="checkbox"/> 5" (2836) |

| Caster forks | | Caster wheel options | |
|--|---|--|--|
| Standard Aluminum | Design Aluminum | Track fixation | Locking |
| <input type="checkbox"/> Alu fork, black anodized  | <input type="checkbox"/> Alu fork, coated in chosen frame colour  | <input type="checkbox"/> Caster wheel track fixation (straight) (Not in combination with shock protection N°9102) | <input type="checkbox"/> Locking mechanism for track fixation (Not in combination with shock protection N°9102) |
| N° 2826 | £ 0 | N° 2829 | £42 |
| | | N° 2851 | £145 |
| | | N° 2850 | £180 |

| Frame options | | | |
|--|---|---|---|
| Outdoor Front End | Splash guard | Bumper | Sun shade |
| <input type="checkbox"/> facilitates significantly active and passive driving on soft or rough surfaces, incl. jack up aid.  | <input type="checkbox"/> Splash guard for Outdoor Front End, made of black plastic  | <input type="checkbox"/> Bumper, made of rubber to protect furniture, doors, etc.  Frame form „Taurus“ not possible with Restraint system. | <input type="checkbox"/> variable adjustable sun shade, in black Usable height (from the top edge of the seat plate/seat cover to the apex) = selected back height + 40 cm Example: BH 30 + 40 cm = 70 cm N° 910301 £200 <input type="checkbox"/> Extension set Increasing the usable amount of the Sun shade by 14 cm. |
| N° 599516 | £395 | N° 59951212 | £79 |
| | | N° 9102 | £59 |
| | | N° 910305 | £54/ Set |


| Frame options |
|---|
| Water bottle |
| <input type="checkbox"/> Water bottle 750ml light PET-bottle, incl. holder  |
| N° 790101 £35 |

| Therapy tables | |
|--|--|
| Transparent | Design |
| Central installation, plugged under seat plate, adjustable in height and distance <input type="checkbox"/> Size 0 (..10) (39 cm width, recommended for SW 18-24) <input type="checkbox"/> Size 1 (..11) (46 cm width, recommended for SW 26-32) <input type="checkbox"/> Size 2 (..12) (59 cm width, recommended for SW 34-36) Restrictions in combination of therapy table and abduction wedge: Height adjustment of therapy table is limited and therapy table can only be used in combination with abduction wedge holders  | <input type="checkbox"/> Number of design for plexiglas therapy table according to www.wheelchairgadgets.com  |
| N° 9209.. | £300 |
| N° 920099 £85 | |

Lightweight head rests Measures in mm (all measures ± 10mm)

Lightweight head rest, adjustable with two conical joints, connection and removable cover included. (cover with same color as back cushion)


Weight reduction approximately 700-1000g to standard head rest systems.



Pad sizes:

| | size 1 | size 2 |
|--------|--------|--------|
| height | 110 | 130 |
| width | 200 | 260 |
| depth | 30 | 40 |


Not approved as a headrest for transport in a motor vehicle.



| Article name | Article-N° | Price |
|--|-------------|-------|
| <input type="checkbox"/> Lightweight head rest Size 1 (Pad size 110x200mm) | N° 68600920 | £275 |
| <input type="checkbox"/> Lightweight head rest Size 2 (Pad size 130x260mm) | N° 68600921 | £275 |

Standard head rests Measures in mm (all measures ± 10mm)

| Article name | Article-N° | Price |
|--|------------|-------|
| <input type="checkbox"/> Standard head rest model 1 (shell-shaped small) | N° 686021 | £265 |
| <input type="checkbox"/> Standard head rest model 3 (curved small) | N° 686023 | £265 |
| <input type="checkbox"/> Standard head rest model 4 (curved large) | N° 686024 | £265 |





Positioning aids Here you will find a small selection to mark. You will find more positioning aids, corresponding mounting possibilities and further equipment in the separate order form „Positioning aids“

| Separate order form | 2-point pelvic belt | 4-point pelvic belt | Chest-shoulder belt | Protection cap |
|--|--|--|---|--|
| <p>Order form „positioning aids“, available in download area.</p> <p>More products for positioning, e.g.:</p> <ul style="list-style-type: none"> • belts, • fixations and corresponding connections, • abduction wedges, • side truss pads, • head rests and more... <p>On the order form refer to this new wheelchair order/enquiry</p> | <input type="checkbox"/> with plastic fastener £120 <input type="checkbox"/> with metal fastener £145 <p>The size is determined by us according to the chosen wheelchair measurements.</p> <p>Example with plastic fastener</p> | <input type="checkbox"/> with plastic fastener £145 <input type="checkbox"/> with metal fastener £175 <p>The size is determined by us according to the chosen wheelchair measurements.</p> <p>Example with metal fastener</p> | <input type="checkbox"/> Chest-shoulder belt without zipper | <input type="checkbox"/> protection cap for metal fastener <p>Makes independent opening for the patient difficult.</p> <p>Only for metal fastener.</p> |
| | N° 780502 | N° 780504 | N° 7810 189,- € | N° 78050210 £15 |

| Lateral support | | | |
|--|--|---|--|
| From BH 30: attachment possible without restrictions. With BH 25: loss of 5 ° tilting path in combination with tilting blocker (N° 192810). | | | |
| Holder | Pad | | Extensions |
| Holder adjustable, size 1, incl. linkage, without pad <input type="checkbox"/> left side (..01) <input type="checkbox"/> right side (..02) N° 685901.. £175/pcs | for truss pad size 1 Attention: Cushions for size 1 can only be mounted on truss pads size 1. 1 pcs. 2 pcs. <input type="checkbox"/> <input type="checkbox"/> 10x7cm (..04) <input type="checkbox"/> <input type="checkbox"/> 10x9 cm (..05) | for truss pad size 2 Attention: Cushions for size 2 can only be mounted on truss pads size 2. 1 pcs. 2 pcs. <input type="checkbox"/> <input type="checkbox"/> 15 x 10 cm (..04) <input type="checkbox"/> <input type="checkbox"/> 18 x 13 cm (..05) | Additional extension of joint for size 1 (3,8 cm) (..0103) 1 St. 2 St. <input type="checkbox"/> <input type="checkbox"/> for size 2 (5,5 cm) (..0203) 1 pos. 2 pos. <input type="checkbox"/> <input type="checkbox"/> |
| Holder adjustable, size 2, incl. linkage, without pad <input type="checkbox"/> left side (..01) <input type="checkbox"/> right side (..02) N° 685902.. £175/pcs | N° 685901.. £55/pcs | N° 685902.. £55/pcs | N° 6859.. £55/pcs |

1. Connection to back
2. Bottom to swing away
3. Adjustable joint
4. Support for pads
5. Pads



| Transport brackets | |
|---|-------------|
| Transport brackets <input type="checkbox"/> Transport brackets (4 pieces) for transportation in a handicapped accessible passenger car  | |
| N° 58507710 | £329 |

| Electric systems |
|--|
| For all electric aids we recommend casters of 5" or higher. Please regard the possible mounting specifications of the manufacturer. AAT Systems www.aat-online.de Alber Systems www.alber.de |

Terms and conditions

1. Preface and Scope

(1) Our terms and conditions apply exclusively. Any deviating, opposing, or complementing terms and conditions on the part of the business partner do not form integral part of contract unless explicitly approved of by us. There are no exceptions to the consent requirement, which also applies in cases where we conduct deliveries without restrictions while aware of the business partner's terms and conditions.

(2) Agreements made in the individual case with the purchaser (including side-agreements, supplements, and alterations) always take priority over these terms and conditions. A contract or our confirmation - both in written form - are authoritative for the content of such arrangements.

(3) Any legally relevant declarations and notifications the purchaser submits to us after conclusion of the contract (e.g. setting of deadlines, notice of defects, declarations of rescission or abatement) require the written form to be valid.

(4) Any information on the application of legal regulations serves the mere purpose of clarification. All regulations are legally binding without any such clarifications unless explicitly altered or excluded in these terms and conditions.

2. Prices

All prices are ex works, excluding insurance and packaging. The prices are net prices; the value-added tax is stated and charged separately according to the legally binding rate at the time of invoicing. The prices valid on the day of delivery are binding unless explicitly arranged otherwise.

3. Conclusion of Contract

SORG's offers are invitatio ad offerendum (invitations to treat). Ordering goods counts as an obligatory offer of contract on the part of the purchaser. Unless the order suggests otherwise, we have the right to accept this contract offer within two weeks time after reception. The acceptance of an offer of contract is effected by the confirmation of the order.

4. Terms of Payment

(1) Provided that no other agreements exist, payment is to be made by paying in advance - 50% at placing of order, 50% pre-delivery of goods. The purchaser bears any costs arising from transactions or discounting.

(2) If a 30-day payment deadline is exceeded, a delay rate of 9% on top of the respective base rate of the European Central Bank falls due without the need for prior notices. Beyond that, the purchaser is charged 3,- € for any reminder sent. We reserve the right to make use of the commercial due rate according to § 353 HGB. A higher damage for delay may be enforced via a suitable verification.

5. Offers and Validity

All current price lists are generally valid until December 31 of the respective year, or until new price lists are published. The latest version of the price list will always be provided on our homepage. Offers older than 3 months at that point expire and must be replaced with new offers and prices.

6. Time of Delivery and Delay

We state the terms of delivery when accepting the order. The delivery date is non-binding and has mere guideline value - especially with regard to customised products. The terms for delay in delivery are defined by legal regulations. In any case, a writing notice issued by the purchaser is necessary for them to take effect.

7. Force Majeure (Act of God)

(1) SORG is not to be held responsible for the following delays of delivery. This also applies if these causes occur with respect to legal representatives, auxiliary persons, or deliverers of SORG: Circumstances of force majeure as well as other extraordinary events that are unpredictable and unavoidable for, and not caused by, SORG, and that only occur after ordering or reserving, and that SORG could not have known about. Moreover, this applies in case of strikes, fire, floods, labour disputes, business disruptions, altered magisterial approbation and legal situation, and magisterial orders which are not due to operational risks.

(2) In the cases of paragraph 1, SORG is eligible to delay delivery/performance for the time of hindrance as well as an appropriate start-up time. If SORG thus causes unacceptable complications for the purchaser, both SORG and the purchaser are eligible to withdraw from the contract. Circumstances that only lead to a temporal and hence acceptable delay of delivery and performance are excluded from this right. The legal rights of the client remain unaffected.

8. Dispatch

We generally dispatch uncovered by insurance and at the recipient's risk. We transfer the risk of accidental loss or deterioration of the goods as well as delay on to the client at the moment of passing the goods on to the forwarder, bearer, or any other person/institute responsible for delivery. SORG will make any effort to respect the client's wishes in terms of path and kind of delivery. Any costs caused by this (including agreed carriage free delivery) are to be covered by the client. In the case of carriage free delivery conducted by SORG, the recipient is required to put forward any carriage costs for a refund. If the client wishes the dispatch to be delayed, any risk of loss or deterioration of the goods is carried by the recipient from the first day of readiness for dispatch. The goods dispatched are only insured upon special demand and at the cost of the recipient.

9. Damage through Delivery/Angelieferte Gegenstände sind, auch wenn sie unwesentliche Mängel Upon delivery, the goods need to be checked for obvious damage in the presence of the bearer. Any damage is to be documented in a written form and co-signed by the bearer. This protocol is to be forwarded to us immediately. Delayed reclamations, or reclamations not signed by the bearer, cannot be taken into account.

10. Warranted Features, Receipt, and Technical Alterations

Delivered objects, even with minor damage, are to be accepted by the client, even if §433 paragraph 1 BGB is not fulfilled. Partial deliveries are permissible. Features only count as warranted if explicitly designated as such. A warranty of condition can only be asserted explicitly and in a written form. SORG reserves the right for technical alterations, especially for solutions of equal or higher value.

11. Warranty for Defects

(1) SORG warrants for any defects of delivered new goods for 12 months as long as the purchaser is not the user (§§ 474, 13 BGB). Should the purchaser be the user then the warranty period is two years for new goods, for used goods 12 months. Should the purchaser also be an entrepreneur, then the liability for used items is void, as long as there is no claim for compensation according to N°12.

(2) The purchaser's claim for defects presuppose that he has attended his obligation to give notice of defects on the legal grounds for his reasons of inquiry (§§ 377, 381 HGB). If defects are found during the inspection or at a later point, we are to be notified in writing without delay. The notice counts as not delayed if issued within two weeks. Posting the notice on time suffices for the keeping of the term. Independently from this obligation to inspect and give notice, the purchaser

is to report obvious defects (including wrong and short delivery) within two weeks upon delivery (again, posting the notice on time suffices). If the purchaser misses the proper inspection and/or notice of defects, we do not warrant for the reported defects. If the delivered good is detected, we may choose to conduct supplementary performance either by repairing the defect (subsequent improvement) or by delivering flawless goods (compensation delivery). Our right to refuse the chosen method of supplementary performance on the respective legal grounds remains untouched.

(3) Our liability for defects expires when alteration or maintenance works are done by the client or a third-person-party without our prior written permission. Consequences of excessive strain, negligent or incorrect handling, and violent damage are excluded from the warranty for defects.

12. Further Liability

(1) With respect to compensation - regardless of the legal cause - SORG is only liable in case of malice and act of gross negligence. In case of slight negligence, SORG is only liable for damage of life, body, health, and for damage caused by the negligence of an essential contractual obligation (an obligation whose fulfilment make the proper realisation of the contract possible in the first place and in whose fulfilment the contract partner regularly trusts and can trust). In this case, however, our liability does not exceed the covering of predictable, typically occurring damage.

(2) A free right to termination on the part of the purchaser (acc. to § 651, 649 BGB in particular) is precluded. Otherwise, the usual prerequisites and legal consequences apply.

13. Costs of Supplementary Performance

Costs of supplementary performance cannot be charged on SORG's account.

14. Voluntary Withdrawal of Goods

SORG is not obliged to take back any delivered end products or parts thereof. If SORG voluntarily decides otherwise in individual cases (especially for therapeutically altered requirements or deaths), SORG refunds a maximum of 75% of the delivery costs for unused parts whose delivery does not date back more than three months. There is no legal obligation concerning withdrawal in this respect. Deliveries which date back more than three months as well as sewed parts and customised products are not returnable. A copy of the delivery note or the bill serves as proof of age of delivery.

15. Reservation of Ownership

Goods delivered by SORG remain the property of SORG until any obligations from the parties' business contract on the part of the purchaser have been redeemed. This also applies if SORG has adapted individual or all claims into a running account which has been balanced and acknowledged. The purchaser has the right to dispose of the goods in the proper course of business. The purchaser hereby yields all claims to SORG which arise from the resale against the purchaser or third parties. The purchaser is authorised to withdraw these claims even upon conveyance. SORG's right to withdraw the claims remains untouched. However, SORG is obliged not to withdraw the claims as long as the purchaser complies with his duty of payment and any other duties. SORG has the right to demand that the purchaser makes known the yielded claims and their debtors, that he provides all information necessary for conveyance, and that he hands in all associated documents. As long and insofar the reservation of ownership is intact, the purchaser must not transfer any goods for securing or bailment without SORG's permission. The purchaser may process or alter the goods for the benefit of SORG exclusively. If processed with other goods which do not belong to SORG, SORG obtains co-ownership of these new goods in relation to the invoice value to the cost price of the other processed goods. If a third party interferes with SORG's ownership, the purchaser is to notify SORG immediately. The purchaser must not make any agreements with other purchasers which may exclude or limit SORG's rights. In particular, the purchaser must not make any agreements which destroy or limit the advance assignments to SORG. If the estimated value of the existing guarantees exceeds the claims to be ensured by more than 20%, SORG is insofar obliged to release by request of the purchaser.

16. Copyright Protection

If SORG conveys any images, drawings, indications of weight or measurements, technical documentations, or other product descriptions and documents (including electronic ones), the respective indications only apply if explicitly stated that they do. SORG reserves ownership and copyright of all these documents. They may only be made available to third parties upon SORG's explicit permission.

17. Place of Fulfilment, Place of Jurisdiction, and Severability Clause

(1) The law of the Federal Republic of Germany applies for these Terms and Conditions and all legal relations between us and our business partners. All international and supranational (contractual) legal systems, the UN Convention on Contracts for the International Sale of Goods in particular, are hereby excluded. Prerequisites and effects of the reservation of ownership acc. to no. 15 are subject to the right of the respective storage location of the good if the choice of law made in favour of the German law is illegitimate or ineffective. The only place of jurisdiction (also on an international basis) for any disputes arising directly or indirectly from the contractual relationship is, according to German Commercial Code, our office in Oberhausen-Rheinhausen, Germany. However, we also have the right to file a suit concerning the general jurisdiction of the purchaser.

(2) If any clause in these Terms of Conditions should turn out to be ineffective, this does not affect the effectiveness of all other clauses.

Technical alterations and misprints reserved. As of: January 2017

SORG Rollstuhltechnik GmbH + Co. KG
Benzstr. 3-5, 68794 Oberhausen-Rheinhausen / Germany
Phone +49 7254 / 9279-0
Fax +49 7254 / 9279-10
info@sorgrollstuhltechnik.de, www.sorgrollstuhltechnik.de