

Mio Family



Active Wheelchairs
for children and adolescents
Light, extra light or tiltable



Lightweight and agile
active wheelchair



The extra light
active wheelchair



The tiltable
active wheelchair

Mio (Design 2018)

Child active wheelchair

German aid index N° 18.99.02.1023



seat width
18-34 cm
(+2 cm)



back height
17,5-40 cm
(+5 cm)



seat depth
18-36 cm
(+4 cm)



weight capacity
50 kg



wheel sizes
18"- 24"



caster sizes
4" und 5"



Overview

- active wheelchair for children
- for early support (from 12 months old)
- innovative lightweight and compact design
- easy handling
- maximum stability
- suitable as a seat for transport in a vehicle according to ISO 7176-19

Mio Carbon (Design 2018)

Small, agile light weight wheelchair



seat width
18-30 cm
(+2 cm)



back height
17,5-35 cm
(+5 cm)



seat depth
18-30 cm
(+4 cm)



weight capacity
50 kg



wheel sizes
20" - 24"



caster sizes
4" und 5"



Overview

- small, agile light weight wheelchair for children
- for early support (from 12 months old)
- seat plate, firm curved back plate, active side guards and foot plate made of ultra-light carbon
- screws and other parts made mostly from sturdy aluminum rather than steel

Mio Move

The tiltable active wheelchair

German aid index N° 18.99.02.1020



seat width
20-34 cm
(+2 cm)



back height
25-45 cm
(+5 cm)



seat depth
18-36 cm
(+4 cm)



weight capacity
50 kg



wheel sizes
20"- 24"



tilting angle
4" und 5"y



Overview

- tiltable wheelchair for children and adolescents
- for early support (from 12 months old)
- active wheelchair and rehab stroller function in one
- allows easy changes between active and passive phases
- suitable as a seat for transport in a vehicle according to ISO 7176-19



Unique people, unique solutions.

One name - three models

Family overview:

	<i>Mio</i> (Design 2018)	<i>Mio Carbon</i> (Design 2018)	<i>Mio Move</i>
seat width	18 - 34 cm	18 - 30 cm	20 - 34 cm
seat depth	18 - 36 cm	18 - 30 cm	18 - 36 cm
back height	18 - 34 cm	17,5 - 35 cm	24 - 45 cm
weight capacity	50 kg	50 kg	50 kg
wheel sizes	18", 20", 22", 24"	20", 22", 24"	20", 22", 24"
lower leg length	15 - 37 cm	15 - 30 cm	till max. 37 cm
grow in three dimensions (seat width, seat depth, seat height)	✓	✓	✓
abduction (slight or strong)	✓	✓	✓
two frame shapes	✓	✓	✓
suitable as a seat for transport in a vehicle according to ISO 7176-19	✓		✓
tilting			✓

Equipment examples

All equipment options are on the order form of the corresponding wheelchair



outdoor front end

brake integrated into the clothing guard side parts

double hand rim

back unit without back tubes

two frame versions taurus and retro

leg support folds back

one hand push handles

double anti-tipper

therapy table

And much more ...