

A large format extremely durable proportional joystick



ca. 650 grams (1.43 lbf)



40 mm (1.57 in) throw



Hand



Durable joystick for excessive force

Developed for users with excessive force (spasticity, dystonia...). This joystick is fully proportional and characterized by an increased throw and force. It can be used as a hand operated joystick.

Adjustable to individual possibilities

Excessive muscle strength asks for a case-by-case joystick configuration. Our Configurator Software allows tweaking parameters to individual possibilities and needs.

Widely compatible


Often, alternative controls are made for a single model or brand. mo-vis products are compatible with the most common wheelchair operating systems (R-net, DX, LiNX, Q-Logic).



Accessibility



Innovation



Practical solutions



High quality