

Motus CV and CS

The Motus CV is a configurable, folding adaptive wheelchair with a modular design with swing-away removable leg supports. This latter feature gives the user even more flexibility in everyday situations. For example, allowing the drive to move up closer to objects.

The Motus CS is a configurable, folding adaptive wheelchair with a closed front frame and a single-panel foot plate. The single-panel foot plate adds stability to the wheelchair and weighs less than a divided foot plate.

Both the CV and CS versions also have the feature of a half folding back, allowing for easy storage and transportation especially in smaller vehicles. The Motus has also had changes to lower the weight of the chair by 2kg, making this chair attractive to a wide range of individuals.

Motus CV and CS Standard Model

- Aluminium tube frame with cross brace
- Folding
- Leg supports – CV swing away / CS integrated
- Footplate – CV individual angle adj / CS single panel angle adj
- Back support with different height setting options
- Clothing guard side panels
- Drive wheels with hollow rims and quick release axles
- Knee lever wheel lock
- Choice of 8 frame colours
- Suitable for use in transport



Motus CV and CS

Motus CV and CS Standard Model

- Aluminium tube frame with cross brace
- Folding
- Leg supports – CV swing away / CS integrated
- Footplate – CV individual angle adj / CS single panel angle adj
- Back support with different height setting options
- Clothing guard side panels
- Drive wheels with hollow rims and quick release axles
- Knee lever wheel lock
- Choice of 8 frame colours
- Suitable for use in transport

Options

- Angle adjustable back support 30°
- Tray
- Scissor wheel lock
- Crutch holder
- Anti-tipper
- Wide range of cushions



Individual swing away footplates on CV



Optional fold down push handles



Side panel options



Single panel footplate on CS

Technical Data

	Motus CV	Motus CS
Overall Width	505 – 680mm	505 – 680mm
Overall Length	805 – 1155mm	805 – 1155mm
Overall Height	670 – 1100mm	670 – 1100mm
Seat Depth	360 – 540mm	360 – 540mm
Seat Width	355 – 480mm	355 – 480mm
Front Seat Height	380 – 550mm	380 – 550mm
Rear Seat Height	360 – 520mm	360 – 520mm
Seat Angle	10 – 100mm	10 – 100mm
Drive Wheel Camber	0° – 4°	0° – 4°
Back Angle Support	0° – 30°	0° – 30°
Back Support Height	300 – 500mm	300 – 500mm
Lower Leg Length	340 – 550mm	340 – 550mm
Folding Size	270 – 330mm	270 – 330mm
Maximum Load Capacity	125kg	125kg
Total Weight of Chair	11.9 – 20.2kg	11.9 – 19.9kg