

## **Motus XXL**

The Motus XXL has outstanding stability with its reinforced Aluminium frame with double cross brace, which gives this chair strength and reliability in any situation. The high quality materials used to produce this chair keeps the weight down as low as 12.8kg in its lightest configuration.

The Motus XXL is designed for people who require the use of a heavy duty chair, the quality and functionality of the Motus XXL makes this chair a stable folding wheelchair that can take a maximum weight of 140kg and seat width of 55.5cm.

The Motus XXL can also be fitted with a half folding back allowing easy storage and ease of transportation in smaller vehicles.

### **Motus XXL Standard Model**

- Reinforced aluminium tube frame with cross brace
- Folding
- Leg supports swing away or integrated
- Footplate individual angle adj or single panel angle adj
- Back support with different height setting options
- Clothing guard side panels
- Drive wheels with hollow rims and quick release axles
- Knee lever wheel lock
- Choice of 8 frame colours



## Motus XXL

### Options

- Angle adjustable back support 30°
- Positioning belt
- Various push handles
- Tray
- Brake for attendant (drum brake)
- Crutch holder with tip assist
- Anti-tipper
- Wide range of cushions

### Technical Data

	Motus XXL
Overall Width	655 – 755mm
Overall Length	865 – 1155mm
Overall Height	670 – 1100mm
Seat Depth	420 – 540mm
Seat Width	505 – 555mm
Front Seat Height	380 – 550mm
Rear Seat Height	360 – 520mm
Seat Angle	10 – 100mm
Drive Wheel Camber	0°
Back Angle Support	0° – 30°
Back Support Height	300 – 500mm
Lower Leg Length	340 – 550mm
Folding Size	270 – 330mm
Maximum Load Capacity	140kg
Total Weight of Chair	12.8 – 20.7kg



Attendant brake option



Swing away single panel footplate



Swing away anti tipper