



P.O. number: \_\_\_\_\_ Date: \_\_\_\_\_ Ship to: \_\_\_\_\_

Company Name: \_\_\_\_\_ Phone: \_\_\_\_\_ Name: \_\_\_\_\_

Purchasing contact:\* \_\_\_\_\_ Address: \_\_\_\_\_

Clinician/Sales rep/ATP: \_\_\_\_\_ City: \_\_\_\_\_ Country: \_\_\_\_\_

Email: \_\_\_\_\_ Postcode: \_\_\_\_\_

Marked for: \_\_\_\_\_ \* Required Phone: \_\_\_\_\_ Taxes not included

MA083 English Manual (STD)

MODEL												
1	WC01	Ultralight MOVE wheelchair aluminum. (120 kg/265 lb weight capacity)								£1 295	<input type="checkbox"/>	
2	WC09	Heavy duty kit (159 kg/350 lb weight capacity)								£350	<input type="checkbox"/>	
SEAT WIDTH <sup>3</sup>												
3	SW	Seat width:	35,6 (14") <input type="checkbox"/>	38,1 (15") <input type="checkbox"/>	40,6 (16") <input type="checkbox"/>	43,2 (17") <sup>1</sup> <input type="checkbox"/>	45,7 (18") <sup>1</sup> <input type="checkbox"/>	48,3 (19") <sup>1</sup> <input type="checkbox"/>		NC	<input type="checkbox"/>	
4	SW	Seat width:	50,8 (20") <sup>1</sup> <input type="checkbox"/>		55,9 (22") <sup>2</sup> <input type="checkbox"/>					£200	<input type="checkbox"/>	
SEAT DEPTH												
5	SD	Seat depth:	40,6 (16") <input type="checkbox"/>		43,2 (17") <input type="checkbox"/>		45,7 (18") <input type="checkbox"/>			NC	<input type="checkbox"/>	
6	SD	Seat depth:	50,8 (20") <input type="checkbox"/>								£175	<input type="checkbox"/>

1 - Available with HD option

2 - Only available with HD option

[3. Width calculator](#)

### Move FRONT Seat-to-floor Height - Please select your choice

		10,16cm (4") caster			12,7cm (5") caster			15,24cm (6") caster			17,78cm (7") caster			20,32cm (8") caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
Fork/Position		STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"
4 in 10,2cm	P1	31,8cm 12,5"	34,3cm 13,5"	36,8cm 14,5"												
	P2	33cm 13"	35,6cm 14"	38,1cm 15"	34,3cm 13,5"	36,8cm 14,5"	39,4cm 15,5"									
5 in 12,7cm	P1															
	P2	33cm 13"	35,6cm 14"	38,1cm 15"	34,3cm 13,5"	36,8cm 14,5"	39,4cm 15,5"									
	P3	34,3cm 13,5"	36,8cm 14,5"	39,4cm 15,5"	35,6cm 14"	38,1cm 15"	40,6cm 16"	36,8cm 14,5"	39,4cm 15,5"	41,9cm 16,5"						
	P4	35,6cm 14"	38,1cm 15"	40,6cm 16"	36,8cm 14,5"	39,4cm 15,5"	41,9cm 16,5"	38,1cm 15"	40,6cm 16"	43,2cm 17"	39,4cm 15,5"	41,9cm 16,5"	44,5cm 17,5"			
7 in 17,8cm	P1							39,4cm 15,5"	41,9cm 16,5"	44,5cm 17,5"	40,6cm 16"	43,2cm 17"	45,7cm 18"			
	P2				39,4cm 15,5"	41,9cm 16,5"	44,5cm 17,5"	40,6cm 16"	43,2cm 17"	45,7cm 18"	41,9cm 16,5"	44,5cm 17,5"	47cm 18,5"			
	P3				40,6cm 16"	43,2cm 17"	45,7cm 18"	41,9cm 16,5"	44,5cm 17,5"	47cm 18,5"	43,2cm 17"	45,7cm 18"	48,3cm 19"	44,5cm 17,5"	47cm 18,5"	49,5cm 19,5"
	P4	40,6cm 16"	43,2cm 17"	45,7cm 18"	41,9cm 16,5"	44,5cm 17,5"	47cm 18,5"	43,2cm 17"	45,7cm 18"	48,3cm 19"	44,5cm 17,5"	47cm 18,5"	49,5cm 19,5"	45,7cm 18"	48,3cm 19"	50,8cm 20"
	P5	41,9cm 16,5"	44,5cm 17,5"	47cm 18,5"	43,2cm 17"	45,7cm 18"	48,3cm 19"	44,5cm 17,5"	47cm 18,5"	49,5cm 19,5"	45,7cm 18"	48,3cm 19"	50,8cm 20"	47cm 18,5"	49,5cm 19,5"	52,1cm 20,5"

Height available +/- 0,65cm (1/4") (see stem bolt and fork combination at the end of the order form)

CASTERS											
7	CA01	Composite wheel w/ PU tire (1" width)	4" <input type="checkbox"/>	5" <input type="checkbox"/>	6" <input type="checkbox"/>	7" <input type="checkbox"/>	8" <input type="checkbox"/>			NC <input type="checkbox"/>	
8	CA10	Pneumatic casters 6" x 1 1/4"								£140 <input type="checkbox"/>	
9	CA11	Pneumatic casters 8" x 1 1/4"								£140 <input type="checkbox"/>	
10	CA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 1/2"								£70 <input type="checkbox"/>
11	CA15		5" x 1 1/2"								£70 <input type="checkbox"/>
12	CA16		6" x 1 1/2"								£70 <input type="checkbox"/>
13	CA17		4" x 1 1/2"								£120 <input type="checkbox"/>
14	CA18	Newton UltraCaster <i>Soft-Roll</i> Aluminum	5" x 1 1/2"								£125 <input type="checkbox"/>
15	CA19		6" x 1 1/2"								£130 <input type="checkbox"/>

### Move - REAR Seat-to-floor Height - Please select your choice

20" (451)	22" (501)	24" (540)	25" (559)	26" (590)
31,1cm / 12in	33,7cm / 13in	35,6cm / 14in	36,8cm / 14,5in	38,1cm / 15in
33,7cm / 13in	36,2cm / 14in	38,1cm / 15in	39,4cm / 15,5in	40,6cm / 16in
36,2cm / 14in	38,7cm / 15in	40,6cm / 16in	41,9cm / 16,5in	43,2cm / 17in
38,7cm / 15in	41,3cm / 16in	43,2cm / 17in	44,5cm / 17,5in	45,7cm / 18in
41,3cm / 16in	43,8cm / 17in	45,7cm / 18in	47cm / 18,5in	48,3cm / 19in
43,8cm / 17in	46,4cm / 18in	48,3cm / 19in	49,5cm / 19,5in	50,8cm / 20in

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1,25 cm (1/2") of selected seat to floor height"

Footrest									
SWING-AWAY FOOTREST									
16	FR01	60° swing away footrest	34,3 - 43,2cm/13,5-17" <input type="checkbox"/>	38,1 - 48,26cm/15,5-19,5" <input type="checkbox"/>	43,18 - 53,34cm/17-21" <input type="checkbox"/>			Set at:*	NC <input type="checkbox"/>
17	FR02	70° swing away footrest	34,3 - 43,2cm/13,5-17" <input type="checkbox"/>	38,1 - 48,26cm/15,5-19,5" <input type="checkbox"/>	43,18 - 53,34cm/17-21" <input type="checkbox"/>			Set at:*	STD <input type="checkbox"/>
FOOTREST OPTIONS									
18	FR09	High mount (footplate attached to the hanger with clamps)	12,7 - 34,3cm/5-13,5" <input type="checkbox"/>					Set at:*	£66 <input type="checkbox"/>
19	FR10	Elevating legrest with calf pad, 35,5 to 55,88 cm (14-22") length adjustable <sup>1</sup>	Left <input type="checkbox"/>	Right <input type="checkbox"/>			Set at:	£200 <input type="checkbox"/>	

1- The elevating legrest replace the footrest side selected

\* +/- 1/2"

FOOTPLATES									
20	FP01	Standard - Composite							STD <input type="checkbox"/>
21	FP02	Newton adjustable angle - Composite							£183 <input type="checkbox"/>
22	FP03	Newton adjustable angle - Aluminum							£200 <input type="checkbox"/>

FOOTPLATE ACCESSORIES									
23	FP08	Length adjustable heel loop							£35 <input type="checkbox"/>
24	FP09	Calf strap							£30 <input type="checkbox"/>
25	FP30	Bodypoint Aeromesh Padded calf strap							£95 <input type="checkbox"/>

SWING AWAY RESIDUAL LIMB SUPPORT									
26	FR27	Residual limb support - flat 15,24x20,32cm (6"x8")	Left <input type="checkbox"/>	Right <input type="checkbox"/>					£355 <input type="checkbox"/>
27	FR28	Residual limb support - flat 20,32x24,4cm (8"x10")	Left <input type="checkbox"/>	Right <input type="checkbox"/>					£355 <input type="checkbox"/>
28	FR29	Residual limb support - flat 25,4x25,4cm (10"x10")	Left <input type="checkbox"/>	Right <input type="checkbox"/>					£355 <input type="checkbox"/>
29	FR30	Residual limb support - flat 25,4x35,56cm (10"x14")	Left <input type="checkbox"/>	Right <input type="checkbox"/>					£355 <input type="checkbox"/>
30	FR32	Residual limb support - contour 15,24x20,32cm (6"x8")	Left <input type="checkbox"/>	Right <input type="checkbox"/>					£355 <input type="checkbox"/>
31	FR33	Residual limb support - contour 20,32x24,4cm (8"x10")	Left <input type="checkbox"/>	Right <input type="checkbox"/>					£355 <input type="checkbox"/>

Backrest									
FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLE									
32	BC01	8° bend aluminum back cane	40,6 cm (16") <input type="checkbox"/>	45,7 cm (18") <input type="checkbox"/>	50,8 cm (20") <input type="checkbox"/>	STD <input type="checkbox"/>			NC <input type="checkbox"/>
ADJUSTABLE ANGLE BACK									
33	BC02	8° bend aluminum back cane	40,6 cm (16") <input type="checkbox"/>		50,8 cm (20") <input type="checkbox"/>	STD <input type="checkbox"/>			£260 <input type="checkbox"/>
34		Select back cane angle:	90° <input type="checkbox"/>	95° <b>STD</b> <input type="checkbox"/>	100° <input type="checkbox"/>	105° <input type="checkbox"/>	110° <input type="checkbox"/>		

STROLLER HANDLE AND STABILIZER BAR <sup>1</sup>										
35	PH04	Stroller handle (not available with the HD option)							£140	<input type="checkbox"/>
36	PH03	Folding stabilizer bar							£220	<input type="checkbox"/>

1 - Stroller handle and stabilizer bar can't be combined

BACK UPHOLSTERY										
37	BH05	Slip-on Nylon back upholstery							NC	<input type="checkbox"/>
38	BH06	Adjustable tension Nylon back upholstery							£145	<input type="checkbox"/>
39	BH99	Omit back upholstery							STD	<input type="checkbox"/>

SEAT CUSHIONS										
40	SC01	5,08cm (2") Cushion							£145	<input type="checkbox"/>
41	SC02	7,62cm (3") Cushion							£162	<input type="checkbox"/>
42	SC99	Omit seat cushion							STD	<input type="checkbox"/>

### Rear wheels

ONE ARM DRIVE <sup>1</sup> Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims										
43	HR22	One arm drive <sup>2,3</sup>		24"	left <input type="checkbox"/>	right <input type="checkbox"/>			£1 200	<input type="checkbox"/>
44	HR98	Omit handrim on opposite side of the One Arm Drive (OAD)							NC	<input type="checkbox"/>

REAR WHEELS										
45	RW01	Mag <sup>4</sup>	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>				£45	<input type="checkbox"/>
46	RW02	Newton One - spoke wheel <sup>5</sup>	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>				NC	<input type="checkbox"/>
47	RW03	Newton Gravity - ultralight	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>		£275	<input type="checkbox"/>

TIRES Available										
48	TI01	Soft Urethane 3,17cm (1 1/4")				20", 22" 24"			STD	<input type="checkbox"/>
49	TI02	Soft Urethane knobby 3,49cm (1 3/8")				22", 24", 25", 26"			£85	<input type="checkbox"/>
50	TI04	Soft Urethane - Shox 2,54cm (1")				22", 24", 25", 26"			£99	<input type="checkbox"/>
51	TI05	Pneumatic 3,49cm (1 3/8")				20", 22", 24", 25", 26"			NC	<input type="checkbox"/>
52	TI08	Pneumatic with Airless insert 3,49cm (1 3/8") <sup>7</sup>				22", 24", 26"			£70	<input type="checkbox"/>
53	TI06	Pneumatic HP - 2,54cm (1")				20", 22", 24", 25", 26"			£85	<input type="checkbox"/>
54	TI07	Pneumatic HP - Schwalbe Marathon Plus 2,54cm (1")				24", 25", 26"			£95	<input type="checkbox"/>

1- Offered with handrims on each side of the wheelchair

2- Only standard handrim position available

3 - Not available with HD option

4 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position

5- Not compatible with Shox tires

7- Only available with Newton One and Newton Gravity rear wheels

### Handrims

HANDRIM										
55	HR01	Aluminum anodized							STD	<input type="checkbox"/>
56	HR61	Aluminum black hard anodized <sup>1,3</sup>					<b>New</b>		£55	<input type="checkbox"/>
57	HR02	Plastic coated <sup>2</sup>							£85	<input type="checkbox"/>
58	HR03	Q-Grip High friction Neoprene coated <sup>2,3,4</sup>							£295	<input type="checkbox"/>
59	HR07	Newton Air Grip <sup>3</sup>		Black <input type="checkbox"/>					£85	<input type="checkbox"/>
60	HR	Newton Air Grip <sup>3</sup>		(08) Sunkissed Orange <input type="checkbox"/>	(09) Burgundy <input type="checkbox"/>	(10) Sapphire Blue <input type="checkbox"/>			£175	<input type="checkbox"/>
				(11) Monster Green <input type="checkbox"/>	(12) Charcoal <input type="checkbox"/>	(13) Fuchsia <input type="checkbox"/>	(14) Steel Blue <input type="checkbox"/>			
				(15) Ferrari Red <input type="checkbox"/>	(16) Acid Green <input type="checkbox"/>	(17) White <input type="checkbox"/>	(18) Sky Blue <input type="checkbox"/>			
61	HR28	Surge <sup>2,3,4</sup>		(1-1/2" Oval handrim w/ Gripton high friction strip)					£329	<input type="checkbox"/>
62	HR05	Surge LT <sup>2,3,4</sup>		(1-1/8" Oval handrim w/ Gripton high friction strip)					£355	<input type="checkbox"/>
63	HR04	Natural Fit <sup>2,3,4</sup>							£270	<input type="checkbox"/>
64	HR99	Omit handrim							NC	<input type="checkbox"/>

1 - Not available in 26"

2- Only standard handrim position available

3 - Not available with Mag wheels or OAD

4 - Not available in 20"

### Handrims option

Assembly position for aluminum and Newton AirGrip handrims										
65	HR29	Standard Handrim position		Full tab					STD	<input type="checkbox"/>
66	HR30	Narrow Handrim position		Short tab 2 positions					£20	<input type="checkbox"/>
67	HR31	Super Narrow Handrim		Short tab 2 positions STD <input type="checkbox"/>		Short tab 1 position <input type="checkbox"/>			£20	<input type="checkbox"/>
Thumb piece <sup>4</sup> selection for Natural Fit handrims										
68	HR23	Standard grip							£75	<input type="checkbox"/>
69	HR27	Super grip							£120	<input type="checkbox"/>
70	HR099	Omit thumb piece							NC	<input type="checkbox"/>


4- Compatible with Surge, Surge LT handrims

AXLES (0° camber) <sup>1</sup>						
71	AX01	Quick release axle			STD	<input type="checkbox"/>
72	AX02	Permanent axle			NC	<input type="checkbox"/>
73	AX03	Quad release axle			£108	<input type="checkbox"/>
AXLE PLATE						
74	AX04	Amputee adjustable axle plate <sup>2</sup>			£170	<input type="checkbox"/>

1- Not required with one arm drive

2 - May limit seat to floor height available

**CENTER OF GRAVITY**



**Choose position of rear wheel axle plate:**

Measure distance from front of back post to center of axle receiver.  
 Default setting is at 1,9 cm (3/4")

Most stable	Most responsive
AMPUTEE	5,1 cm - 2"
<input type="checkbox"/>	<input type="checkbox"/>

WHEEL LOCKS						
75	WL08	Newton Wheel lock - push to lock			STD	<input type="checkbox"/>
76	WL02	Pull to lock			NC	<input type="checkbox"/>
77	WL04	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)			£135	<input type="checkbox"/>
78	WL13 WL14	15,24cm (6") removable extension handle			£66	<input type="checkbox"/>
79	WL06	Attendant Lock (additional set of wheel lock mounted to the back)			£175	<input type="checkbox"/>
80	WL99	Omit wheel locks			NC	<input type="checkbox"/>

### Armrests

FLIP BACK ARMREST (Convertible to "T" armrest) <sup>1,2</sup>		Height Adjustment		Removable function		
81	AR01	25,4cm (10") Desk length arm pad	20 - 30cm (8 - 12") STD <input type="checkbox"/>	25 - 36cm (10 - 14") <input type="checkbox"/>	Locked <input type="checkbox"/>	Unlocked <input type="checkbox"/> STD £310 <input type="checkbox"/>
82	AR02	35,56cm (14") Full length arm pad	20 - 30cm (8 - 12") STD <input type="checkbox"/>	25 - 36cm (10 - 14") <input type="checkbox"/>	Locked <input type="checkbox"/>	Unlocked <input type="checkbox"/> STD £310 <input type="checkbox"/>
"T" ARMREST <sup>1</sup>		Height Adjustment				
83	AR03	25,4cm (10") Desk length arm pad	20 - 30cm (8 - 12") STD <input type="checkbox"/>	25 - 36cm (10 - 14") <input type="checkbox"/>		£285 <input type="checkbox"/>
84	AR04	35,56cm (14") full length arm pad	20 - 30cm (8 - 12") STD <input type="checkbox"/>	25 - 36cm (10 - 14") <input type="checkbox"/>		£285 <input type="checkbox"/>
GEL ARM PAD						
85	AR12	35,56cm (14") Full Length Arm Pad - Gel Ovation				£133 <input type="checkbox"/>
86	AR99	Omit armrest				NC <input type="checkbox"/>

1 - May have a 1/2" adjustment restriction on both limits of the range

2 - Not available with attendant lock

POSITIONING BELTS						
87	SB01	Velcro adjustable	Standard 88,9 cm/35" <input type="checkbox"/>	Long 152,40 cm/60" <input type="checkbox"/>	£40	<input type="checkbox"/>
88	SB02	Auto buckle	Standard 88,9 cm/35" <input type="checkbox"/>	Long 152,40 cm/60" <input type="checkbox"/>	£75	<input type="checkbox"/>
89	SB07	Hardware attachment for 2 point seatbelt (1 pair)			NC	<input type="checkbox"/>
90	SB03	Additional hardware for 2 extra attachment points (needed for the 4 point seat belt) (1 pair)			£15	<input type="checkbox"/>
91	SB04	BodyPoint Hip Belt - 2 point not padded		M (142,24 cm/56") <input type="checkbox"/>	L (154,94 cm/61") <input type="checkbox"/>	£60 <input type="checkbox"/>
92	SB05	BodyPoint Hip Belt - 2 point padded	S (129,54 cm/51") <input type="checkbox"/>	M (142,24 cm/56") <input type="checkbox"/>	L (154,94 cm/61") <input type="checkbox"/>	£150 <input type="checkbox"/>
93	SB06	BodyPoint Hip Belt - 4 point padded	S (129,54 cm/51") <input type="checkbox"/>	M (142,24 cm/56") <input type="checkbox"/>	L (154,94 cm/61") <input type="checkbox"/>	£165 <input type="checkbox"/>
94	SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps		M (73,66 cm/18") <input type="checkbox"/>	L (93,98 cm/24") <input type="checkbox"/>	£250 <input type="checkbox"/>
95	SB99	Omit positioning belt				NC <input type="checkbox"/>

SIDE GUARDS (not available with Flip Back and "T" Armrest)						
96	SG06	Plastic fender side guards (removable) <sup>1</sup>			£235	<input type="checkbox"/>
SIDE GUARDS OPTION						
97	SG07	Back cane mounted sideguard clip			NC	<input type="checkbox"/>

1- Not compatible with 20" and 22" rear wheels

HEADREST										
98	AC01	Profiled Headrest	S (8,8x15,2 cm/3,5"x6") <input type="checkbox"/>	M (11,4x17,8 cm/4,5"x7") <input type="checkbox"/>	L (14x2,3 cm/5,5"x8") <input type="checkbox"/>				£158	<input type="checkbox"/>
99	AC28	Headrest mounting hardware with stabiliser bar <sup>1</sup>							£258	<input type="checkbox"/>

1 - Not available with a stroller handle or a folding stabiliser bar

ACCESSORIES										
100	AC03	Anti tippers							£120	<input type="checkbox"/>
101	AC29	Swing-away anti-tippers <sup>1</sup>	Left <input type="checkbox"/>	Right <input type="checkbox"/>			New		£116	<input type="checkbox"/>
102	AC27	Tip assist <sup>2</sup>							£54	<input type="checkbox"/>
103	AC04	Transit Tie Down (Unoccupied)							NC	<input type="checkbox"/>
104	AC05	WC-19 Transit Tie Down w/Q'straint belt <sup>3</sup>	Short belt 20-41" <input type="checkbox"/>	Long belt 20-46" <input type="checkbox"/>					£280	<input type="checkbox"/>
105	AC06	Cane carrier	Left <input type="checkbox"/>	Right <input type="checkbox"/>					£104	<input type="checkbox"/>
106	AC07	Oxygen tank holder - E Type	Left						£270	<input type="checkbox"/>
107	AC08	Spoke guards <sup>4</sup>	Available: 20", 22", 24", 26"						£100	<input type="checkbox"/>
108	AC09	Footrest neoprene impact guard (pair) - Black							£79	<input type="checkbox"/>
109	AC10	Seat pouch							£65	<input type="checkbox"/>
110	AC26	BodyPoint Mobility Bag							£108	<input type="checkbox"/>
111	AC23	Newton Backpack							£100	<input type="checkbox"/>
112	AC12	Velcro Closing Strap							£45	<input type="checkbox"/>
113	AC13	Tool kit							£39	<input type="checkbox"/>
114	A7AC16	Motion Composites Polo shirt	W/S <input type="checkbox"/>	W/M <input type="checkbox"/>	W/L <input type="checkbox"/>				£29	<input type="checkbox"/>
			M/S <input type="checkbox"/>	M/M <input type="checkbox"/>	M/L <input type="checkbox"/>	M/XL <input type="checkbox"/>	M/XXL <input type="checkbox"/>	M/XXXL <input type="checkbox"/>		

1 - Must select the pair (left and right) with HD option


2 - Cannot be combined with Anti Tippers

3 - May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

4 - Not available with Mag wheels and Spinergy

FRAME COLOR (with black aluminum crossbrace)											
115	FC01	Black <input type="checkbox"/>							NC	<input type="checkbox"/>	
116	FC	(03) Sunkissed Orange <input type="checkbox"/>	(04) Burgundy <input type="checkbox"/>	(05) Sapphire Blue <input type="checkbox"/>	(06) Monster Green <input type="checkbox"/>				£95	<input type="checkbox"/>	
		(07) Charcoal <input type="checkbox"/>	(08) Fuchsia <input type="checkbox"/>	(09) Steel Blue <input type="checkbox"/>	(17) Ferrari Red <input type="checkbox"/>	(18) Acid Green <input type="checkbox"/>					
117	FC99	Custom Color	<a href="#">Specify PPG Paintit color:</a>						Additional 10 days lead-time	£300	<input type="checkbox"/>

<b>TOTAL</b>										
--------------	--	--	--	--	--	--	--	--	--	--



**Notes :**