




TREKKING
& TREKKING POWER


TREKKING
& TREKKING POWER


by Gianni Conte

WWW.TRIRIDEITALIA.COM    


by Gianni Conte

Multi Traction Wheelchair

Discover a new mobility


TREKKING
& TREKKING POWER

MTW

REAR MOTORIZATION KIT FOR MANUAL WHEELCHAIRS

Born from Triride technology, it allows you to transform your manual wheelchair into an advanced technology electronic wheelchair. The MTW Trekking / Trekking Power kit attaches and removes in a few seconds, returning the wheelchair to its original state without any structural changes.

Indoor, Trekking or Trekking Power

Combined with any Triride, Tribike Hybrid or Triwheel model, it allows you to have a high-performance electric wheelchair with one, two or three-wheel drive (3X3). Thanks to the special control you can choose the type of drive you want. By removing the front power traction it is possible to use the wheelchair in indoor mode using the practical joystick.



- indoor -

- trekking -

- trekking power -

GENERAL FEATURES*

- Maximum power: from 600 W to 2400 W (Rh+Lh)
- Single or double battery: 36 V or 48 V
- Traction selector: Front (I), Rear (II), Multi Traction 3x3 (III)
- Range*: up to 50 Km
- Programmable electronics via Triride software for a complete customisation
- Rear LED light with stop function, flashing when braking
- Anti-tip arm
- Wheel and tyre diameter of your choice
- Wide range of accessories

*For all features see the technical data sheet



UNIVERSAL

Compatible and installable on most wheelchairs on the market, with both rigid and foldable frame.

EASY TO INSTALL

The wheelchair is equipped with patented Triride attachments to which it is possible to hook and unhook the MTW rear axle without making any structural changes.

GRAVEL, SAND, GRASS, DIRT, SNOW

The motorized wheels supplied provide rear-wheel drive that allows you to tackle any type of off-road route or very demanding climb.

