

NOIR 2.0



Seat width	33 - 36 - 39 - 42 - 45 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 cm
Seat height	Made-to-measure
Footplate distance	Made-to-measure (standard) Adjustable (on request)
Backrest height	From 24,5 to 47 cm (adjustable)
Front frame angle	85° - 90° - 95°
Total length	85 cm approx. (std conf.)
Total width	Seat width + 17 cm (with 0° camber)
Backrest angle	From 74° to 94° to the seat (adjustable)
Camber	0° - 2° - 4°
Weight without rear wheels	Starting from 4,0 kg approx.
Wheelchair weight	Starting from 6,6 kg approx.
Maximum load	125 kg