

# permobil Order Form - Progeo noir 2.0

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts, AL7 1AN United Kingdom - Phone: +44 1484 722888

Date .....	UK (English)	
User .....	Company stamp	Rigid frame wheelchair
End user weight (kg) .....		<b>£ 8313,00</b>
Language user manual .....	Responsible .....	<b>Retail price - VAT not included</b>

**Carbon fibre frame made-to-measure.**

**NOIR is a wheelchair without possibilities of adjustments after its manufacturing.**

See reference tables

Take care to write all the requested measures (in centimetres)

**Frame Angle**

85°  90°  95°

**Seat depth (PS)**

35     37,5  
 40     42,5  
 45     47,5

**See note (1)**

Backrest rake adjustable from -4° to +16°

**Seat width (LS)**

33     36  
 39     42  
 45

**camber**

2°     0°     4°

**Setting C**

Extremely prudential 7,5  
 Prudential 9  
 Standard 10,5  
 Active 12  
 Extreme 13,5

**Front height (min 44-max 54)**

**Rear height (min 37-max 48)**  
*calculated with rear wheels 24"*

Std. Value 33,5 (min 26-max 47) adjustment steps 1.5cm

**Footplate option**

Carbon fibre frame with integrated footplate, not adjustable  
*Lightest configuration*

Aluminium footplate adjustable distance +/- 1,5cm (rake and depth adjustable)  
**Minimum footplate distance 35cm.**

Internal   
  2/3 Internal   
  2/3 External   
  External

11531

**Seat**

Seat cloth with pocket (under seat)

Seat with tension straps 08883    **231,00**

Carbon rigid seat 08336    **489,00**

Pocket (under seat) [25x11x4,5] 08281    **76,00**

Pag. 1 / 6

UK\_202404

Rev. 2024.1\_Apr


customer\_support\_uk@permobil.com - https://www.permobil.com/en-gb

Date	Order	User	Address
------	-------	------	---------

Since that the paint is sensitive to the surface material and texture, the colour selected for the wheelchair frame can also differ considerably from what is displayed on a RAL colour chart or also compared to what is displayed on a monitor.

### CARBON

☆ Polished transparent painting



Opaque transparent painting 08925 643,00

### ESSENCE

#### Polished transparent painting

<input type="radio"/> ANTHRACITE Essence	10351	902,00
<input type="radio"/> BLUE Essence	08928	902,00
<input type="radio"/> GREEN Essence	11163	902,00
<input type="radio"/> LIME Essence	08926	902,00
<input type="radio"/> ORANGE Essence	10349	902,00
<input type="radio"/> RED Essence	08927	902,00
<input type="radio"/> SAND Essence	11165	902,00
<input type="radio"/> VIOLET Essence	11167	902,00
<input type="radio"/> WHITE Essence	10350	902,00



#### Opaque transparent painting

<input type="radio"/> ANTHRACITE Opaque Essence	08929	1.155,00
<input type="radio"/> BLUE Opaque Essence	10354	1.155,00
<input type="radio"/> GREEN Opaque Essence	11162	1.155,00
<input type="radio"/> LIME Opaque Essence	10352	1.155,00
<input type="radio"/> ORANGE Opaque Essence	10355	1.155,00
<input type="radio"/> RED Opaque Essence	10353	1.155,00
<input type="radio"/> SAND Opaque Essence	11164	1.155,00
<input type="radio"/> VIOLET Opaque Essence	11166	1.155,00
<input type="radio"/> WHITE Opaque Essence	10356	1.155,00

my ESSENCE colour RAL \_\_\_\_\_

Glossy

Note


Opaque

11186 1.332,00

### RACE

#### Polished transparent painting

<input type="radio"/> ANTHRACITE Race	11168	902,00
<input type="radio"/> BLUE Race	11170	902,00
<input type="radio"/> GREEN Race	11172	902,00
<input type="radio"/> LIME Race	11174	902,00
<input type="radio"/> ORANGE Race	11176	902,00
<input type="radio"/> RED Race	11178	902,00
<input type="radio"/> SAND Race	11180	902,00
<input type="radio"/> VIOLET Race	11182	902,00
<input type="radio"/> WHITE Race	11184	902,00



#### Opaque transparent painting

<input type="radio"/> ANTHRACITE Opaque Race	11169	1.155,00
<input type="radio"/> BLUE Opaque Race	11171	1.155,00
<input type="radio"/> GREEN Opaque Race	11173	1.155,00
<input type="radio"/> LIME Opaque Race	11175	1.155,00
<input type="radio"/> ORANGE Opaque Race	11177	1.155,00
<input type="radio"/> RED Opaque Race	11179	1.155,00
<input type="radio"/> SAND Opaque Race	11181	1.155,00
<input type="radio"/> VIOLET Opaque Race	11183	1.155,00
<input type="radio"/> WHITE Opaque Race	11185	1.155,00

my RACE colour RAL \_\_\_\_\_

Glossy

Note

Opaque

11187 1.332,00

### FLASH

#### Polished transparent painting

<input type="radio"/> BLUE Flash	11190	902,00
<input type="radio"/> WHITE Flash	11192	902,00
<input type="radio"/> RED Flash	11188	902,00



#### Opaque transparent painting

<input type="radio"/> BLUE Opaque Flash	11191	1.155,00
<input type="radio"/> WHITE Opaque Flash	11193	1.155,00
<input type="radio"/> RED Opaque Flash	11189	1.155,00

#### Colour set: includes rear wheel hub and rim and front forks

Only standard Progeo 22"/24" wheels have coloured hubs and rims


☆ Black	22" 24" 25" 26"		
<input type="radio"/> Bright silver	22" 24"	00833	130,00
<input type="radio"/> Bright lime	22" 24"	00833	130,00
<input type="radio"/> Bright blue	22" 24"	00833	130,00
<input type="radio"/> Bright red	22" 24"	00833	130,00
<input type="radio"/> Bright fuchsia	22" 24"	00833	130,00
<input type="radio"/> Pearl white	22" 24"	00833	130,00
<input type="radio"/> Tekno anthracite	22" 24"	00833	130,00
<input type="radio"/> Matt Midnight blue	22" 24"	00833	130,00
<input type="radio"/> Matt grey	22" 24"	00833	130,00
<input type="radio"/> Matt red	22" 24"	00833	130,00

Date	Order	User	Address
------	-------	------	---------

### Backrest adjustment and fixing system

- Backrest with double locking (upright and folded) with rake adjustable from -4° to +16°
- Backrest with upright locking system ONLY (no locking when folded) with rake adjustable from -4° to +16°
- Fixed back(non-folding back) with rake adjustable from -4° to +16°


- Integrated backrest (Banded, not folding back)
  - Exclude installation of the titanium backrest
  - Exclude installation of the Tubular armrests



10841    **1.285,00**

### Backrest

- Adjustable velcro straps and padded transpiring upholstery
- Carbon fibre backrest
- Titanium backrest  
Not available with: Integrated backrest (Banded, not folding back)



- AIR** Integrated posture backrest (NOT REMOVABLE)
  - Carbon shell - Lateral retaining 7,5cm
  - Depth adjustment ±1cm. - Angle adjustment ±10°
  - not available for seat width 48

**1.078,00**

<b>H 28</b>		<b>H 34</b>		<b>H 40</b>	
Backrest Height (TS) min 34		Backrest Height (TS) min 38		Backrest Height (TS) min 44	
(.....)		(.....)		(.....)	
Lateral retaining wing A=high - R=regular					
<b>A</b>	<b>R</b>	<b>A</b>	<b>R</b>	<b>A</b>	<b>R</b>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

#### Installation on:

- Aluminium backrest tubes  
WITHOUT backrest rigidizer  
(Carbon tubes with integrated backrest)



- Carbon fibre backrest



- Titanium backrest  
Not available with: Integrated backrest (Banded, not folding back)



### Handle tube

- Tubes with short handle
- SPORT tubes Without push handles *Lightest configuration*
- Folding handles **Minimum backrest height 30,5cm**



02958    **260,00**

### Casters

- 3" Sport *Lightest configuration*
- 4" Soft
- 3" Active Soft (wide) (Aluminium alloy rim)
- 4" Active Sport (narrow) (Aluminium alloy rim)
- 4" Active Soft (wide) (Aluminium alloy rim)
- 5" Active Soft (wide) (Aluminium alloy rim)
- 5" Pneumatic

### Rear wheels - See chart N°3

- 24" ULTRA-LIGHT PROGEO KEY SPORT *Lightest configuration*
- 24" Ultra-light "Spinergy LX"
- 24" Ultra-light "Spinergy LXL"  Black hub and white spokes      Hub and spokes black
- 24" Ultra-light "Spinergy Flex rim" Integrated aluminium/rubber handrim
- 25" Ultra-light "Spinergy LX"
- 25" Ultra-light "Spinergy LXL"  Black hub and white spokes
- 26" Ultra-light "Spinergy LX"

04353    **543,00**


Date	Order	User	Address
------	-------	------	---------

<b>Rear wheel tyre</b>			
<input checked="" type="checkbox"/> High pressure tyre	(24"-25"-26")	<i>Lightest configuration</i>	
<input type="checkbox"/> High pressure profile tyre	(24")		
<input type="checkbox"/> Marathon Plus	(24"-25")	<b>with anti-puncture protection</b>	<b>01669 65,00</b>
<input type="checkbox"/> Anti-puncture - Colour Black	(24")		<b>11083 130,00</b>
<b>Not available with Rear wheels : 24" ULTRA-LIGHT PROGEO KEY SPORT, 24" Ultra-light "Spinergy LXL"</b>			



<b>Handrim</b>		<b>Wheel distance</b> <input checked="" type="checkbox"/> <b>3 cm</b>	<input type="checkbox"/> <b>4 cm</b>		
<input checked="" type="checkbox"/> Aluminium	(24"-25"-26")	<i>Lightest configuration</i>			
<input type="checkbox"/> Aluminium black SOFT-TOUCH	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)		<b>09170 130,00</b>	
<input type="checkbox"/> Chromed Steel	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)		<b>01684 152,00</b>	
<input type="checkbox"/> Titanium	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)		<b>00016 283,00</b>	
<input type="checkbox"/> Ultra Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)		<b>00017 324,00</b>	
<input type="checkbox"/> Max Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)		<b>00834 359,00</b>	
<input type="checkbox"/> TOP FRICTION	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)		<b>11023 411,00</b>	
<i>Oval countour with a rubber insert</i>					
<input type="checkbox"/> NOVA H	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)		<b>11025 271,00</b>	
<i>Ergonomic black profile</i>					
<input type="checkbox"/> NOVA D	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)		<b>11026 295,00</b>	
<i>Ergonomic profile with rubber coating (ideal for quadriplegics)</i>					
<input type="checkbox"/> SURGE LT	(25")	Oval countour with a rubber insert		<b>06427 371,00</b>	
<input type="checkbox"/> Quad	(24")			<b>00693 206,00</b>	
<input type="checkbox"/> Silicon cover for standard aluminium handrim		Available for measures (22"- 24" - 26")		<b>04050 165,00</b>	

<b>Wheel locks</b>			
<input checked="" type="checkbox"/> Sport brake			
<input type="checkbox"/> Bent lever push to lock			
<input type="checkbox"/> Straight lever push to lock			
<input type="checkbox"/> Straight lever pull to lock			
<input type="checkbox"/> Brake with extension lever push to lock			
<input type="checkbox"/> Brake with extension lever pull to lock			

<b>Mudguard</b>			
<input checked="" type="checkbox"/> Carbon fibre mudguard			
<input type="checkbox"/> Carbon fibre Straight sideguard	<input checked="" type="checkbox"/> Aligned wheel	<input type="checkbox"/> Under wheel ..... cm	<input type="checkbox"/> Above wheel ..... cm
<input type="checkbox"/> Removable carbon fibre mudguard			<b>04109 260,00</b>
<input type="checkbox"/> Removable Carbon fibre straight sideguard	<input checked="" type="checkbox"/> Aligned wheel	<input type="checkbox"/> Under wheel ..... cm	<input type="checkbox"/> Above wheel ..... cm
			<b>04109 260,00</b>

<b>Armrests</b>			
<input type="checkbox"/> Tubular armrests	(Swing away, removable and height adjustable armrest)		<b>00694 260,00</b>
<i>Not available with Integrated backrest (Banded, not folding back)</i>			

<b>Accessories</b>			
<input type="checkbox"/> Pad for front frame protectors			<b>07058 54,00</b>
<input type="checkbox"/> Backrest-seat connection cloth			<b>10241 18,00</b>
<input type="checkbox"/> Removable storage bag (with quick coupling)			<b>09811 130,00</b>

<b>Optional 1</b>			
<input type="checkbox"/> Transit wheels			<b>03290 152,00</b>
<input type="checkbox"/> Anti-tipper with rapid hookup	(inside rotation to reduce the encumbrance)		
	<input type="checkbox"/> Right		<b>06821 130,00</b>
	<input type="checkbox"/> Left		<b>06822 130,00</b>

<b>Optional 2</b>			
<input type="checkbox"/> Safety pelvic belt			<b>03972 65,00</b>
<input type="checkbox"/> Harness			<b>04103 130,00</b>
<input type="checkbox"/> Lumbar support	(to insert in the backrest upholstery)		<b>00013 41,00</b>
<input type="checkbox"/> Tetraclip	(Quick release axle for quadriplegic)		<b>03947 182,00</b>


<b>Additional Kit</b>			
<input type="checkbox"/> 24" Off-road rear wheels	Off-Road tyre - aluminium handrim (black colour)		<b>03378 513,00</b>
<input type="checkbox"/> 24" FAT rear wheels	FAT tyre 3" - aluminium handrim (black colour)		<b>09357 902,00</b>

Date .....	Order .....	User .....	Address .....
------------	-------------	------------	---------------


### PHYSIO - Posture backrest (Std. Backrest upholstery not provided)

Removable with quick release  
Suitable for tubes with diameter 18÷26mm  
Depth adjustment (4cm)  
Angle adjustment (±18°)

- AIR** Carbon shell  
not available for seat width 48

<b>BACKREST HEIGHT</b>	<b>28</b>	<b>34</b>	<b>40</b>	<b>1.078,00</b>	
Shell's shape  	Lateral retaining wing A=high - R=regular				
	<b>A</b>	<b>R</b>	<b>A</b>	<b>R</b>	<b>A</b>
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

- POSTURE** Aluminum shell

<b>BACKREST HEIGHT</b>	<b>34</b>	<b>40</b>	<b>52</b>	<b>1.043,00</b>	
Shell's shape  	Lateral retaining wing R=regular				
	<b>R</b>	<b>R</b>	<b>R</b>	<b>R</b>	<b>R</b>
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Configuration DEEP with side support PAD (height adjustable)


- Aluminum lateral with padded, breathable upholstery.



<p><b>Small (h 13cm)</b></p> <div style="display: flex; justify-content: space-around;"> <div style="text-align: center;"> <input type="checkbox"/> Right <b>10272</b> <b>30,00</b> </div> <div style="text-align: center;"> <input type="checkbox"/> Left <b>10272</b> <b>30,00</b> </div> </div>	<p><b>Large (h 20cm)</b></p> <div style="display: flex; justify-content: space-around;"> <div style="text-align: center;"> <input type="checkbox"/> Right <b>10273</b> <b>41,00</b> </div> <div style="text-align: center;"> <input type="checkbox"/> Left <b>10273</b> <b>41,00</b> </div> </div>
--	--

### SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve  
Aluminum lateral with padded, breathable upholstery.



<input type="radio"/> Large (18x12)		<b>05647 188,00</b>
<input type="radio"/> Small (14x9)		<b>05643 188,00</b>
<input type="radio"/> Lateral support extension	Quantity (number pcs).....	<b>10271 71,00</b>
<input type="radio"/> Lateral support extension	Quantity (number pcs).....	
<input type="radio"/> Small (14x9)		<b>05642 188,00</b>
<input type="radio"/> Large (18x12)		<b>05646 188,00</b>

- Headrest adjustable - Double ball mounting system  
Available only with backrest Physio Posture



<input type="checkbox"/> Large (34x14)		<b>05648 295,00</b>
<input type="checkbox"/> Small (28x14)		

Date	Order	User	Address
------	-------	------	---------

### REFERENCE TABLES

**Progeo (Español)** NOIR is a wheelchair without possibilities of adjustments after its manufacturing.

(1) The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.

1	A Footplate width			
	DP	35÷39	40÷45	46÷50
LS				
33		21,5	21	20,5
36		24,5	24	23,5
39		27,5	27	26,5
42		30,5	30	29,5
45		33,5	33	32,5
48		36,5	36	35,5

2	Setting	C :
	Extremely prudential	7,5
	Prudential	9
	Standard	10,5
	Active	12
	Extreme	13,5

3	Maximum width
Camber	
0°	LS + 18
2°	LS + 21
4°	LS + 24

Maximum width calculated including 3 cm. distance between tyre and handrim

Send this page only if any optional have been chosen	
<input checked="" type="checkbox"/>	standard (supplied if no requests are indicated)
<input type="checkbox"/>	on request (without extra charge)
<input type="checkbox"/>	optional (cost to be added at the basic price)

Wheelchair : Progeo noir 2.0 - Rev. 2024.1

Prices - VAT excluded  
to be added to the price list