

Juvo

The Juvo B5 and B6 power wheelchair have a compact design and have a customisable modular design and is available in both front and rear wheel drive options.

The Juvo features a variety of seating options with the Elite versions coming with the new VAS (Vario Adjust Seat) can be tailored to meet your individual requirements. It comes with a padded back support which is tension adjustable and a low cut seat cushion. The biomechanical Ergo back angle adjustment is free from unwanted shear forces.

The Juvo uses the R-net VR2 control system which can be programmed to suit your requirements. However, Ottobock have worked closely with Curtiss-Wright and have devised a controller which is unique to Ottobock. This is the Ten° controller, which has been designed to have everything at the optimum 10° reach angle ensuring ease and comfort of use, with a built in hand pad it is not only comfortable for the user but also gives stability when using the joystick. An added feature is that you can have it colour coded to match your chair.

What is special about me?

- Good directional stability
- Large range of possibilities to adapt your sitting position
- Relatively small turning radius
- Good cornering ability
- Easy to approach furniture etc.
- Unique features with the VAS Seat and Ten° Controller

Standard Model

- Standard seat
- Seat cushion
- Back upholstery, adaptable
- Seat Inclination 0°
- Mechanical back adjustment
- Lap belt
- Leg support adjustment
- Individual plastic footplate
- 53Ah / 63Ah batteries with fan
- 9" puncture proof casters

- 14" puncture proof drive wheels
- VR2 controller
- Swing away controller

Options

- Front wheel drive
- VAS seat system
- USB charging socket
- Ten° controller
- Mechanical or Gyro track stabilisation
- 15 frame colour options

Quality for life



Technical Data

| | |
|---|--|
| Speed | 4mph / 5mph / 6mph / 8mph |
| Battery size | 53Ah / 63Ah / 73Ah |
| Range under test conditions approx. | 35km |
| Control unit | VR2 or Ottobock Ten° |
| Max user weight | Size 1: 140kg Size 2: 160kg |
| Max. obstacle negotiation height | FWD – 100mm RWD – 60mm optionally 100mm |
| Climbing ability | 10° or 17% |
| Turning radius | Size 1: 600mm Size 2: 640mm |
| Seat height (depending on seating solution) | 400mm – 550mm |
| Seat width | 340mm – 560mm |
| Seat depth | 340mm – 500mm |
| Back angle | 0° to +30° Electric back angle adjustment -5 to 30° on VAS |
| Seat tilt | 45° electric |
| Back height | 450mm – 550mm |
| Armrest height | 225mm – 350mm |
| Seat inclination | -3° / 0° / 3° / 6° |
| Max. charging time | 12hr |
| Weight from | 107kg |
| Total width | 600mm – B5 640mm – B6 |
| Overall Length | 1050mm – 1250mm |