

## Kimba Buggy

The Kimba combines maximum functionality while meeting the highest aesthetic criteria and providing numerous individual options. When in the buggy the child can actively experience the environment and train their motor skills. Parents can easily adapt the Kimba to the needs of their child with a few simple adjustments with the tilt and recline functions.

The Kimba also features suspension which allows the child to have a comfortable ride even over uneven ground. The seat unit can be either forward facing so the child can take in the surroundings while out, or it can be rear facing so the parent or guardian can see the child should that be required.

### Standard Model

#### Outdoor Mobility Base

- Available in 2 sizes
- Seat tilt adjustability up to 30°
- Height adjustable push bar with ratchet joint
- Rear axle suspension
- Swivelling front wheels with swivel lock
- EVA wheels with wheel lock
- Integrated anchor points for mounting in a wheelchair accessible vehicle
- Fixed storage tray with carrying capacity of 15kg



## Kimba Buggy

### Kimba Multifunctional Seating Unit

- Available in 2 sizes, size 1: 1-6years and size 2: 4-10 years approx.
- Seating unit made of aluminium
- 2x Fabric choices Spacer Fabric and Microfiber
- Removable, washable covers
- Adjustable seat depth, width and back height
- Lateral support pads, height adjustable
- Knee angle and height of footrest assembly can be adjusted

### Options


- Thoracic supports
- Fixation vest, chest and shoulder harnesses
- Lumbar Support padding
- Pommel
- Various positioning belts
- Footrest block 5cm for shorter lower leg length
- Seat depth reduction
- Grab rail, tray and armrest
- Sun canopy in upholstery colour, including rain cover



Back recline adjustment



Buggy board



Integrated basket



Seat tilt

### Technical Data

| Multifunctional Seating Unit         | Size 1          | Size 2          |
|--------------------------------------|-----------------|-----------------|
| Age                                  | 1 – 6 years     | 4 -10 years     |
| Max. user weight                     | 40 kg           | 40 kg           |
| Seat width                           | 19 – 31 cm      | 24 – 40 cm      |
| Seat depth                           | 20 – 30 cm      | 26 – 35 cm      |
| Lower leg length                     | 19 – 31 cm      | 20 – 37 cm      |
| Backrest height                      | 41 – 61 cm      | 56 – 70 cm      |
| Back angle                           | 80° to 180°     | 80° to 180°     |
| Folding size (LxWxH)                 | 58 x 37 x 42 cm | 59 x 43 x 58 cm |
| Seat unit weight                     | 7 kg            | 8.5 kg          |
| Seat depth with seat depth reduction | From 14 cm      | From 14 cm      |

### Technical Data

| Outdoor Mobility Base             | Size 1          | Size 2          |
|-----------------------------------|-----------------|-----------------|
| Max. user weight                  | 55 kg           | 55 kg           |
| Total width                       | 60 cm           | 70 cm           |
| Wheel diameter front/rear         | 17 / 28 cm      | 17 / 28 cm      |
| Push handle                       | 71 - 116 cm     | 71 - 116 cm     |
| Folding size without seat (LxWxH) | 90 x 60 x 53 cm | 97 x 67 x 59 cm |
| Tilt in Space                     | -10° to +35°    | -10° to +35°    |
| Base weight                       | 10.5 kg         | 11.5 kg         |