

## **Kimba Inline**

The Kimba Inline buggy is a mobility base for twins or siblings, who require positioning during transport. It can be combined with various seating solutions, seating shells, modular seating, Kimba seat unit, Leckey Squiggles and Leckey Mygo systems.

The Kimba Inline twin buggy is easy to push and steer and as it is the same width as other buggies it therefore does not have problems with doorways. The positioning of the seats can be arranged to suit you with them either facing the direction travel, facing the parent or guardian or facing each other allowing the siblings to interact.

### **What is special about me?**

- Seats position one behind the other
- Consistently narrow mobility base even with child growth
- Easy to operate in everyday life
- The frame length is adjustable to meet the needs of the children
- Children can sit in the direction of travel, facing the attendant for facing each other, by turning the seating units



## Kimba Inline

### Standard Model

#### Outdoor Mobility Base

- Pneumatic tyres
- Rear wheels with drum brake
- Seat frame tilt-able in both directions

#### Kimba Multifunctional Seating Unit

- Size 1: 1-6years a prox.
- Seating unit made of aluminium
- Powder Coated
- Removable, washable covers
- Adjustable seat depth, width and back height
- Lateral support pads, height adjustable
- Knee angle and height of footrest assembly can be adjusted



### Options

- Footrest assembly for seating shell interface
- Parallel seating shell adapter

### Technical Data

	Outdoor Mobility Base
Push bar height	93 cm -103 cm
Width	64 cm
Seat tilt	-30° to +40°
Wheel diameter front/back	20cm / 25cm
Folding size (LxWxH) cm	90 x 64 x 48
Weight	10.6 kg
Load capacity	50 kg
Maximum weight of child	50 kg