

Kimba Cross

The Kimba Cross with its three wheel mobility base is versatile, practical and also looks good with its aluminium oval tubes. On rugged terrain, the seat suspension absorbs the impact allowing the child to remain comfortable at all times.

The large pneumatic tires, drum brakes and swivelling front wheel allows the buggy to be easy to manoeuvre. The front wheel can be easily removed by the quick release mechanism, which then allows the bicycle bar to be attached and then Kimba Cross can then be towed by a push bike for more adventures.

What is special about me?

- Swivelling wheel with quick lock attachment for confined spaces
- Cycle trailer package (draw bar with quick lock, cycle coupling, pennant, light)



Kimba Cross

Standard Model

Outdoor Mobility Base

- Telescoping frame fits various sized seating units
- Spring system for the seat
- 16" front wheel, pneumatic tyre, rigid quick lock attachment
- 18" rear wheels, quick release axle
- Drum brake with lock option
- Height adjustable pushbar
- Can be folded
- Attractive design

Kimba Multifunctional Seating Unit

- Available in 2 sizes, size 1: 1-6 years and size 2: 4-10 years approx.
- Seating unit made of aluminium
- Powder Coated
- Removable, washable covers
- Adjustable seat depth, width and back height
- Lateral support pads, height adjustable
- Knee angle and height of footrest assembly can be adjusted

Options

- Anti-tipper, clothing protectors
- Spoke protectors
- Storage bag
- Choice of material colours

Technical Data

	Outdoor Mobility Base	Size 1	Size 2
Push bar height	90 cm -115 cm		
Width	66 cm		
Length	140 cm – 170 cm		
Height	75 cm – 120 cm		
Seat depth		19 cm -31 cm	22 cm – 40 cm
Seat width		20 cm – 30 cm	26 cm – 35 cm
Lower leg length		19 cm – 31 cm	20 cm – 37 cm
Backrest height/angle		46 cm – 66 cm / 0° to + 90°	56 cm – 70 cm / 0° to +90°
Wheel diameter front/back	40cm / 45cm (swivelling wheel 30cm / 45cm)		
Seat tilt	-5° to +30°		
Load capacity	50 cm	35 cm	35 cm
Weight	15.8kg (incl. splash guard)	7 cm	8.8 cm
Folding size (LxWxH) cm	94 x 60 x 32	72 x 40 x 22	82 x 44 x 27