

Zenit product range

Lightweight and dynamic




Enhanced dynamics and efficiency

Weight-optimised components and improved rigid-frame properties enhance the driving quality considerably. The new front wheel adapter allows the front wheels to be adjusted for precise directional stability.

Numerous other adjustment options also permit optimal adaptation for even greater comfort. When you get a wheelchair for active use from the Zenit line, you get a product that makes zero compromises.

i You may request the instructions for use with comprehensive information on the product as a PDF file at the email address oa@ottobock.com. The PDF file can also be displayed in a larger size.



Adjustable, patented front wheel connection and folding foot plate



Permanently welded caster wheel connection



Folding, angle-adjustable back support



Permanently welded drive wheel connection

Additional options

- Back support tubes curved to measure
- Option of using other back support systems, e.g. Baxx
- Back support upholstery, e.g. breathable or ultra light
- Push handles, e.g. fold-down
- Side panels made of aluminium or carbon
- Foot plates and foot plate bar, e.g. carbon
- Various wheel lock systems, e.g. scissors or knee lever wheel lock 4.0
- Drive wheel rims, e.g. Infinity Ultralight
- Handrims, e.g. stainless steel or rubber-coated
- Caster wheels, e.g. soft casters
- Various colour packages for drive wheel rims and caster fork

Only with Zenit / Zenit R

- Aluminium or carbon frame tube

Only with Zenit / Zenit CLT

- Lock for leg support (patented fixation)

Only with Zenit R CLT

- Lower leg length adjustable or permanently welded

Technical data	Zenit	Zenit CLT	Zenit R	Zenit R CLT
Overall length	720–1,040 mm	640–1,040 mm	640–1,040 mm	640–1,040 mm
Overall width	490–725 mm	495–710 mm	495–710 mm	495–710 mm
Overall height	720–1,040 mm	620–1,050 mm	620–1,050 mm	620–1,050 mm
Seat width	320–460 mm	320–460 mm	320–460 mm	320–460 mm
Seat depth	360–500 mm	360–500 mm	360–500 mm	360–500 mm
Front seat height	380–540 mm	380–540 mm	380–540 mm	380–540 mm
Rear seat height	370–500 mm	370–500 mm	370–500 mm	370–500 mm
Back support height	250–500 mm	250–500 mm	250–500 mm	250–500 mm
Lower leg length	280–510 mm	280–510 mm	280–510 mm	280–510 mm
Anterior frame angle	75°/85°	75°/85°	75°/85°	75°/85°
Max. load	120 kg	100 kg	120 kg	100 kg
Weight from approx.	8.9 kg	8.9 kg	8.8 kg/8.5 kg	7.4 kg
Transport weight (without drive wheels) from approx.	–	6.3 kg	6.3 kg/6.0 kg	4.9 kg

Material



Aluminium



Carbon

The style package includes coloured drive wheel rims and front wheel adapters in the colours Sparkle light orange, Candy red, Marine blue and Sparkle Granny Smith.



Zenit

Article no. 480F160=1_AA02_C

Product features

- Folding wheelchair with the look and driving characteristics of a rigid-frame wheelchair
- Adjustable, patented front wheel connection
- Adjustable seat height and lower leg length
- Centre of gravity adjusted to the millimetre
- Folding, angle-adjustable back support
- Easy-Go folding system
- Folding foot plate
- Open frame design
- Choice of frame offsets




Zenit CLT

Article no. 480F160=1_AA03_C

Product features

- Folding wheelchair with the look and driving characteristics of a rigid-frame wheelchair
- Permanently welded components:
 - Back support connection
 - Drive wheel connection
 - Wheel lock connection
 - Front wheel connection
 - Side panel mounting
 - Easy-Go folding system
 - Folding foot plate
- Open frame design
- Choice of frame offsets
- Adjustable lower leg length (optional)
- Folding back support (optional)

RAL colours (powder coating)



White RAL 9016



Jet black RAL 9005



Jet black RAL 9005 matte (texture)

Effect colours (powder coating)



Sparkle light orange



Candy red



Vinho sparkle



Ice blue



Sky blue



Marine blue



Sparkle granny-smith



Brilliant silver



Anthracite metallic



Shimano matte

Anodised colours (Not for Zenit CLT and Zenit R CLT available.)



Gold, matte anodised



Orange, matte anodised



Red, matte anodised



Violet, matte anodised



Blue, matte anodised



Green, matte anodised



Anthracite, matte anodised



Black, matte anodised



Zenit R

Article no. 480F160=1_AA04_C

Product features

- Rigid-frame wheelchair
- Ergonomic oval axle
- Adjustable, patented front wheel connection
- Adjustable seat height and lower leg length
- Centre of gravity adjusted to the millimetre
- Folding, angle-adjustable back support
- Open frame design
- Choice of frame offsets



Zenit R CLT

Article no. 480F160=1_AA05_C

Product features

- Rigid-frame wheelchair
- Ergonomic oval axle
- Permanently welded components:
 - Back support connection
 - Drive wheel connection
 - Wheel lock connection
 - Front wheel connection
 - Side panel mounting
- Folding back support (optional)
- Open frame design
- Choice of frame offsets

