

Permobil

# Panthera S3 Series



permobil

# Designed to fit you and your active lifestyle

The Panthera S3 series is designed to accommodate users of different builds. Although there are variations in dimensions, every S3 offers the comfort of a full-width frame, a low transport weight, and a complete range of accessories. Available in seven models.

## Key benefits:

- Made to maneuver
- Low weight and robust
- Adjusts for comfort
- Integrated anti-tip device
- Easy to manage



## Specifications

	S3	S3 Short	S3 Short Low
<b>Seat width</b>	33 - 36 - 39 - 42 - 45 - 50 cm	30 - 33 - 36 - 39 - 42 - 45 cm	33 - 36 - 39 cm
<b>Seat depth</b>	40 cm	35 cm	27 - 33 cm, 35 cm
<b>Seat height</b>	Front 47 cm (24" rear wheels) Back 43 cm (24" rear wheels)	Front 47 cm (24" rear wheels) Back 43 cm (24" rear wheels)	Front 47 cm (24" rear wheels) Back 43 cm (24" rear wheels)
<b>Seat angle</b>	7°	7°	7°
<b>Total length</b>	78 - 90 cm	73 - 85 cm	73 - 85 cm
<b>Total width</b>	Seat width +21 cm	Seat width +21 cm	Seat width +21 cm
<b>Front frame angle</b>	105°	105°	105°
<b>Backrest height</b>	20 - 45 cm	20 - 45 cm	20 - 45 cm
<b>Backrest angle</b>	13° - (-8°)	13° - (-8°)	13° - (-8°)
<b>Footplate distance</b>	35 - 43 cm	35 - 43 cm	35 - 43 cm
<b>Camber</b>	2.2°, 4° (seat width 33 - 45 cm) 6° (seat width 33 cm) 8° (seat width 33 - 39 cm)	2.2°, 4° (seat width 30 - 45 cm) 6° (seat width 30 - 33 cm) 8° (seat width 30 - 39 cm)	2.2°, 4°, 6° (seat width 33 cm), 8°
<b>Weight without rear wheels</b>	4.2 - 5.3 kg	4.0 - 4.7 kg	4.5 - 4.8 kg
<b>Wheelchair weight</b>	7.8 - 9.0 kg	7.6 - 8.3 kg	8.0 - 8.3 kg
<b>Maximum load</b>	100 kg, 150 kg (seat width 45- 50 cm)	100 kg, 150 kg (seat width 45 cm)	100 kg
<b>Material (chassis)</b>	Chrome molybdenum tubing	Chrome molybdenum tubing	Chrome molybdenum tubing
<b>Material (upholstery)</b>	Polyurethane covered polyester	Polyurethane covered polyester	Polyurethane covered polyester



Personalise your ride\*

ANODIC BLACK

PEARL WHITE

GREY METALLIC

BLUE CLASSIC

PINE GREEN

Actual product colors may vary.

S3 Large	S3 Long	S3 0°	S3 Short Abd
39 - 42 - 45 - 50 cm	33 - 36 - 39 - 42 - 45 - 50 cm	33 - 36 - 39 - 42 - 45 - 50 cm	27 - 30 - 33 cm
45 cm	45 cm	35 - 42.5 cm	27 - 33 cm
Front 49.5 cm (24" rear wheels) Back 45.5 cm (24" rear wheels)	Front 47 cm (24" rear wheels) Back 43 cm (24" rear wheels)	Front 45 cm (24" rear wheels) Back 45 cm (24" rear wheels)	Front 47 cm (24" rear wheels) Back 43 cm (24" rear wheels)
7°	7°	0°	4°
83 - 95 cm	83 - 95 cm	84 cm	73 - 85 cm
Seat width +21 cm	Seat width +21 cm	Seat width +21 cm	Seat width +28 cm (6° camber)
105°	105°		105°
20 - 45 cm	20 - 45 cm	20 - 45 cm	20 - 45 cm (25 - 40 cm seat width 27 cm)
13° - (-8°)	13° - (-8°)	13° - (-8°)	13° - (-8°)
41 - 51 cm	35 - 43 cm	35 - 43 cm	35 - 43 cm
2.2°, 4° (seat width 39 - 45 cm)	2.2°, 4° (seat width 33 - 45 cm) 6° (seat width 33 cm) 8° (seat width 33 - 39 cm)	2.2°, 4° (seat width 33 - 45 cm) 6° (seat width 33 cm) 8° (seat width 33 - 39 cm)	2° - 4° - 6° - 8°
4.8 - 5.2 kg	4.7 - 5.3 kg	4.5 - 5.4 kg	4410 - 4570 grams
8.4 - 8.8 kg	8.3 - 8.8 kg	8.0 - 8.9 kg	7970 - 8130 grams
100 kg, 150 kg (seat width 45 - 50 cm)	100 kg, 150 kg (seat width 45 - 50 cm)	100 kg, 150 kg (seat width 45 - 50 cm)	100 kg
Chrome molybdenum tubing	Chrome molybdenum tubing	Chrome molybdenum tubing	Chrome molybdenum tubing
Polyurethane covered polyester	Polyurethane covered polyester	Polyurethane covered polyester	Polyurethane covered polyester



Permobil

## Panthera S3 Series

Designed to fit you  
and your **active lifestyle**