

Adjusted by:

Tel:

Save the form to the desktop and  
open it with Adobe Reader.

E-mail:

Date:

The price changes depending on  
the choices made in the form £:

Standard configuration includes; chassis, seat, calf strap, backrest lock, standard colour option, backrest, carbon rear axle, castor and fork, titanium foot rest or foldable footplates, push to lock S3 brake, Panthera crosslaced steel spoke rear wheel, titanium pushrims, and Schwalbe right Run tires.

**Choice required:**Max. user weights : Seat width **27-42 cm 100 kg**, Seat width **45-50 cm 150 kg**, (S3 Large) Seat width **39-50 cm 150 kg**

Optional

Upcharge /Add-price

Marked yellow

Click on "[blue underlined text](#)" for more info or image**Original design****1) Chassis**[S3](#)Seat width **33-50** cm, Seat depth 40 cm[S3 Short](#)Seat width **30-45** cm, Seat depth 27-33 cm (seat width 30-33), 35 cm (seat width 36-45)[S3 Short Low](#)Seat width **33-39** cm, Seat depth 35 cm[S3 Large](#)Seat width **39-50** cm, Seat depth 45 cm[S3 Long](#)Seat width **33-50** cm, Seat depth 45 cm[S3 0°](#)Seat width **33-50** cm, Seat depth 35, 37.5, 40, 42.5 cm[S3 Short Abd](#)Seat width **27-33** cm, Seat depth 27-33 cm

Clear

**Alternative options****Other options****2) Chassis colour**[Black](#) - Anodic Black[Gray](#) - Metallic gray

Clear

Optional colour - A delivery time of 20 days  
for optional colour. Fill in art.no. in the square**3) Seat width**

Check allowed seat width at "1) Chassis"

Clear	27 cm	30 cm	33 cm	36 cm	39 cm	42 cm	45 cm	50 cm
-------	-------	-------	-------	-------	-------	-------	-------	-------

**4) Seat depth**

Clear	35 cm	37.5 cm	40 cm	42.5 cm	45 cm
-------	-------	---------	-------	---------	-------

**5) Back rest**Height (HH) [INFO](#)[Standard](#)

20 cm

35 cm

[Tapered](#) 525BBHH VER2[INFO](#)[Tapered 6 cm](#)

523BBHH VER2

25 cm

40 cm

[Curved back](#) 524BBHH VER2[Flared 3 cm](#)

Clear

30 cm

45 cm

[For fixed backrest](#)[Curved-tapered](#)

Clear

**6) Rear axle**[Carbon std](#)

33023BB

Clear

[Carbon +15 mm](#)  
33023BC\*

\* seat width BB+1

Use this rear axle when choosing arm support  
at Accessories and/or increased seat height[Other](#) rear axle or seat height[Increase the seat height](#)**7) Castor and fork**[art.no.](#)[90 mm](#)[120 mm](#)[150 mm](#)

Clear

**8) Footrest**[Footrest Titanium](#)

36001BB

[Footrest U-shape, small heel bar](#)

3600134 (width 33)

Clear

[Titanium, extended, reinforced](#) \*

36111BB

\* std for S3 Large

[Foldable footplates;](#) [art.no.](#)[Heel strap for footplates](#)[Plastic plate for footrest](#)

36010BB

[Carbon footplate](#)

65040BB

**9) Brake**[Brake S3](#)

4630000

Clear

[One hand brake, right](#)[art.no.](#)[Assistant brake \\*](#)[One hand brake, left](#)

Art.no.

\* Not for short / low chassis

All prices for rearwheels are for pair

### 10) Rear wheel - without pushrim

[Rear wheel sizes](#)

[24" std](#)

[art.no.](#)

Clear

22" [standard](#)

24" [Spinergy X](#) black spokes

25" [Spinergy X](#) black spokes

26" [Spinergy X](#) black spokes

24" [Spinergy X](#) white spokes

25" [Spinergy X](#) white spokes

26" [Spinergy X](#) white spokes

[One arm drive](#)

[Rear wheel 24x2" MTB](#)  
Titanium pushrims

### Tyre - choose tyre

All prices for tyres are for pair

[Schwalbe Right Run](#) 22", 24", 25", 26"

Clear

[Rugged 24"](#)

[Schwalbe Marathon Plus](#), 24", 25", 26"

[Schwalbe One](#) 24", 25"

[PU 24" \(massive\)](#)

reTyre does not fit PU tyres

[reTyre skins 24"](#)

[reTyre skins 25"](#)

[reTyre skins 26"](#)

Clear

### 11) Pushrim

All prices for pushrims are for pair

[Titanium](#)

190DD11 (x2 / wheelchair)

Clear

[Paragrip](#)

19051DDBT x 2

[Tetragrip](#)

19052DDBT x 2

[Maxgrip 22", 24"](#)

1903060 x 2 (60=24", 80=22"

(contains latex)

[Aluminum X](#)

1904404 x2 (04=24", 05=25", 06=26")

## Accessories

### Side guard [art.no.](#)

[Side guard S3](#)

[video](#)

Clear

### Arm support [art.no.](#)

[Arm support S3](#) \*

[video](#)

\* 26" std for 24" wheels for Short Low  
26" std for Assistant brake

### Side guard Carbon [art.no.](#)

[Without fender](#)

[Small fender](#)

[Large fender](#)

### Push handles [Video](#)

[Height adjustable S3](#)

4252060 VER2

[Height adjustable S3 +10 cm](#)

4252063 VER2

Max. userweight 100 kg

Clear

### Push handle foldable

[Foldable S3](#)

4245000

[Video](#)

### Push bar

[Push bar S3](#)

42400BB VER2

Max. userweight 100 kg

### Anti-Tip [Video](#)

[Anti-tip](#) std seat height

4353000

[Anti-tip](#) S3 Short Low

4355000

[Anti-tip](#) S3 Large

4353080

Clear

[Other seat height](#)

Art.no.

### Other

[Tube protection, front](#)

9900251

[Tube protection, rear](#)

9900251

[Crutch holder S3](#)

9900301 (x1/ crutch)

[Backrest wedge](#)

4710012

[Freewheel](#)

### Cushion [Wedge cushion](#)

[2,5 cm](#)

572BB40

[Wedge cushion](#)

57710BB

[5 cm](#)

575BB40

Clear

### Spoke protectors

[24" logo](#)

4500000

[22" logo](#)

4500100

[24" no logo](#)

4500011

[22" no logo](#)

4500111

Clear

### Spoke protectors with theme

[22" Splash](#)

4500322

[24" Splash](#)

4500324

[22" Car rim](#)

4500122

[24" Car rim](#)

4500124

[22" Wild animals](#)

4500422

[24" Wild animals](#)

4500424

[22" Smiley](#)

4500222

[24" Smiley](#)

4500224

Clear

Chose art.no. for Freewheel according to the footrest height of the wheelchair  
Read more about Freewheel on the [homepage](#).

SPLASH CAR RIM WILD ANIMALS SMILEY

## **Set-up adjustments**

### **Back rest angle**

Note measure (A) for desired back rest angle.  
14-28 mm. Standard = 20 mm

A= \_\_\_\_\_mm

### **Rear axle adjustment**

Note measure (A) for desired balance  
Standard = 60 mm. at raised rear axle = 51 mm

A= \_\_\_\_\_mm

### **Footrest height - Foot plates**

Note measure (A) and (B) for desired footrest

A= \_\_\_\_\_mm    B= \_\_\_\_\_mm

## **Extra**