

Permobil

# Panthera S3 Swing Series



permobil

# Adapts to the active life you lead

The Panthera S3 Swing offers the same comfort and low transport weight as the Panthera S3 series, the S3 Swing provides the added flexibility of a swing-away foot rest for easier transfers. Also available in S3 Swing Short and Long option to suit smaller or taller builds.

## Key benefits:

- Made to maneuver
- Low weight and robust
- Adjusts for comfort
- Integrated anti-tip device
- Easy to maintain



Personalise your ride\*

ANODIC BLACK

PEARL WHITE

GREY METALLIC

BLUE CLASSIC

PINE GREEN

Actual product colors may vary.

## Specifications

	S3 Swing	S3 Swing Short
<b>Seat width</b>	36 - 39 - 42 - 45 - 50 cm	33 - 36 - 39 cm
<b>Seat depth</b>	40 cm	35 cm
<b>Seat height</b>	Front 45 cm (24" rear wheels) Back 43 cm (24" rear wheels)	Front 45 cm (24" rear wheels) Back 43 cm (24" rear wheels)
<b>Seat angle</b>	4°	4°
<b>Total length</b>	78 - 90 cm	88 cm
<b>Total width</b>	Seat width +21 cm	Seat width +21 cm
<b>Front frame angle</b>	105°	105°
<b>Backrest height</b>	20 - 45 cm	20 - 45 cm
<b>Backrest angle</b>	13° - (-8°)	13° - (-8°)
<b>Footplate distance</b>	38 - 46 cm	38 - 46 cm
<b>Camber</b>	2.2°, 4° (seat width 36 - 45 cm) 8° (seat width 36 - 39 cm)	2.2°, 4°, 6 (seat width 33 cm) 8°
<b>Weight without rear wheels</b>	4.1 - 4.9 kg	4.0 - 4.1 kg
<b>Wheelchair weight</b>	9.0 - 10.0 kg	8.8 - 9.0 kg
<b>Maximum load</b>	100 kg, 150 kg (seat width 45 - 50 cm)	100 kg
<b>Material (chassis)</b>	Chrome molybdenum tubing	Chrome molybdenum tubing
<b>Material (upholstery)</b>	Polyurethane covered polyester	Polyurethane covered polyester

	S3 Swing Low 0°	S3 Swing Long
<b>Seat width</b>	36 - 39 - 42 - 45 - 50 cm	36 - 39 - 42 - 45 - 50 cm
<b>Seat depth</b>	35 - 37.5 - 40 cm	37.5 - 40 - 42.5 - 45 cm
<b>Seat height</b>	Front 39 cm (24" rear wheels) Back 39 cm (24" rear wheels)	Front 45 cm (24" rear wheels) Back 43 cm (24" rear wheels)
<b>Seat angle</b>	4°	4°
<b>Total length</b>	78 - 90cm	83 - 95 cm
<b>Total width</b>	Seat width +21 cm	Seat width +21 cm
<b>Front frame angle</b>	105°	105°
<b>Backrest height</b>	20 - 45 cm	20 - 45 cm
<b>Backrest angle</b>	13° - (-8°)	13° - (-8°)
<b>Footplate distance</b>	38 - 46 cm	38 - 46cm
<b>Camber</b>	2.2° 4° (seat width 33-45cm) 8° (seat width 33-39cm)	2.2° 4° (seat width 33-45cm) 8° (seat width 33-39cm)
<b>Weight without rear wheels</b>	4.6 - 5.4 kg	4.9 - 5.7 kg
<b>Wheelchair weight</b>	9.4 - 10.2 kg	9.7 - 10.5 kg
<b>Maximum load</b>	100 kg, 150 kg (seat width 45- 50 cm)	100 kg, 150 kg (seat width 45- 50 cm)
<b>Material (chassis)</b>	Chrome molybdenum tubing	Chrome molybdenum tubing
<b>Material (upholstery)</b>	Polyurethane covered polyester	Polyurethane covered polyester

# Innovating for individuals