

Adjusted by:

Tel:

Save the form to the desktop and
open it with Adobe Reader.

E-mail:

Date:

The price changes depending on **£:**
the choices made in the form

Standard configuration includes chassis, seat, calfstrap, backrest lock, standard colour option, backrest, carbon rear axle, castor and fork, swing away foot rest with foldable footplates, push to lock S3 brake, Panthera crosslaced steel spoke rear wheel, titanium pushrims, and Schwalbe right Run tires.

Choice required:

Max. user weights : Seat width **33-42 cm 100 kg**, Seat width **45-50 cm 150 kg**,

Optional

Upcharge /Add-price

Marked yellow

Click on "[blue underlined text](#)" for more info or image**Original design****Alternative options****Other options****1 Chassis**[S3 Swing](#)Seat width **36-50** cm, Seat depth 40 cm[S3 Swing Short](#)Seat width **33-39** cm, Seat depth 35 cm[S3 Swing Low 0°](#)Seat width **36-50** cm, Seat depth 35,37.5,40 cm (seat width 45 cm) 37.5,40 cm (seat width 50 cm)[S3 Swing Long](#)Seat width **36-50** cm, Seat depth 37.5,40,42.5, 45 cm

Clear

2 Chassis colour

Black- Anodic black

Optional colour - A delivery time of 20 days
for optional colour. Fill in art.no. in the square

Gray- Metallic gray

Clear

3 Seat width

33 cm 36 cm 39 cm 42 cm 45 cm 50 cm

Clear

Check allowed seat width at "1 Chassis"

4 Seat depth

35 cm 37.5 cm 40 cm 42.5cm 45 cm

Clear

5) Back rest Height (HH) [INFO](#)[Standard](#)

20 cm 35 cm

[Tapered](#) 525BBHH VER2[INFO](#)[Tapered 6 cm](#)

523BBHH VER2

25 cm 40 cm

[Curved back](#) 524BBHH VER2[Flared 3 cm](#)

Clear

30 cm 45 cm

[For fixed backrest](#)[Curved-tapered](#)

Clear

6) Rear axle[Carbon std](#)

33023BB

Clear

[Carbon +15 mm](#)
33023BC*

* seat width BB+1

Use this rear axle when choosing arm support at
Accessories and/or increased seat height[Other](#) rear axle or seat height[Increase the seat height](#)**7) Castor and fork**[art.no.](#)[90 mm](#)[120 mm](#)

Clear

8) Footrest[art.no.](#)[Swing S3 cpl pair](#)[Footplate large](#), extended 80 mm forward

Clear

[Extended tubes. 40 mm](#)[Heel strap for footplates](#)**9) Brake**[Brake S3](#)

4630000

Clear

Art.no.

[Assistant brake *](#)

* Not for short / low chassis

All prices for rearwheels are for pair

10) Rear wheel - without pushrim

[Rear wheel sizes](#)

24" std	art.no.	22" standard	One arm drive
Clear		24" Spinergy X black spokes	
		25" Spinergy X black spokes	
		26" Spinergy X black spokes	
		24" Spinergy X white spokes	
		25" Spinergy X white spokes	
		26" Spinergy X white spokes	

11) Tyre - choose tyre

All prices for tyres are for pair

Schwalbe Right Run 22", 24", 25", 26"	Rugged 24"	PU 24" (massive)	
Clear	Schwalbe Marathon Plus , 24", 25", 26"	reTyre does not fit PU tyres	
	Schwalbe One 24", 25"	reTyre skins 24"	
		reTyre skins 25"	
		reTyre skins 26"	Clear

12) Pushrim

All prices for pushrims are for pair

Titanium 190DD11 (x2 / wheelchair)	Paragrip 19051DDBT x 2	Maxgrip 22", 24" 1903060 x 2 (60=24", 80=22" (contains latex)	
Clear	Tetragrip 19052DDBT x 2	Aluminum X 1904404 x2 (04=24", 05=25", 06=26")	

Accessories

Side guard [art.no.](#)

[Side guard S3](#)
video
Clear

Arm support [art.no.](#)

[Arm support S3](#) *
video
* 26" std for 24" wheels for Short Low
26" std for Assistant brake

Side guard Carbon [art.no.](#)

[Without fender](#)
[Small fender](#)
[Large fender](#)

Push handles [Video](#)

[Height adjustable S3](#)
4252060 VER2
[Height adjustable S3 +10 cm](#)
4252063 VER2
Max. userweight 100 kg
Clear

Push handle foldable

[Foldable S3](#) [Video](#)
4245000

Push bar

[Push bar S3](#)
42400BB VER2
Max. userweight 100 kg

Anti-Tip [Video](#)

[Anti-tip](#) std seat height
4353000
[Anti-tip](#) short low
4355000
[Anti-tip](#) Large
4353080
Clear
[Other seat height](#)

Art.no.

Cushion [Wedge cushion](#)

[2,5 cm](#) [Wedge cushion](#)
572BB40 57710BB
[5 cm](#)
575BB40
Clear

Spoke protectors

[24" logo](#) [22" logo](#)
4500000 4500100
[24" no logo](#) [22" no logo](#)
4500011 4500111
Clear

Spoke protectors with theme

22" Splash 4500322	24" Splash 4500324
22" Car rim 4500122	24" Car rim 4500124
22" Wild animals 4500422	24" Wild animals 4500424
22" Smiley 4500222	24" Smiley 4500224

Clear

Other

[Tube protection, front](#)
9900251
[Tube protection, rear](#)
9900251
[Crutch holder S3](#)
9900301 (x1/ crutch)
[Backrest wedge](#)
4710012

[Table S3](#)

6810020
Only together with arm support

SPLASH CAR RIM WILD ANIMALS SMILEY

Set-up adjustments

Back rest angle

Note measure (A) for desired back rest angle.
14-28 mm. Standard = 20 mm

A= _____ mm

Rear axle adjustment

Note measure (A) for desired balance
Standard = 40 mm.

A= _____ mm

Footrest height

Note measure (A) and (B) for desired height

A= _____ mm B= _____ mm

Extra