

Permobil

# Panthera U3 Light



permobil

# Optimal ergonomics to enable your lifestyle

A low-weight, maneuverable manual wheelchair, the Panthera U3 Light offers stability and comfortable posture support. Designed for optimal ergonomics when seated or driving, this wheelchair can also be easily lifted into a car thanks to its balanced, effortless-grip chassis. Available as U3 Light and U3 Light L.

## Key benefits:

- Ultra-light weight
- Extremely maneuverable
- Easy to transport
- Adaptable to individual needs
- Easy to maintain



## Specifications

	U3 Light	U3 Light Y-front
Seat width	33 - 36 - 39 - 42 - 45 cm	33 - 36 - 39 - 42 - 45 cm
Seat depth	35 - 46 cm	35 - 46 cm
Seat height	Front 47 cm (w. 24" rear wheels) Back 43 cm (w. 24" rear wheels)	Front 47 cm (w. 24" rear wheels) Back 43 cm (w. 24" rear wheels)
Seat angle	7°	7°
Total length	82 cm	82 cm
Total width	Seat width +21 cm	Seat width +21 cm
Front frame angle	100°	100°
Backrest height	25 - 35 cm	25 - 35 cm
Backrest angle	11.5° - (-4.5°)	11.5° - (-4.5°)
Footplate distance	35 - 43 cm	35 - 43 cm
Camber	2.2°	2.2°
Weight without rear wheels	3.4 - 4.0 kg	3.7 - 3.9 kg
Wheelchair weight	5.7 - 6.3 kg	6.0 - 6.2 kg
Maximum load	100 kg	100 kg
Material (chassis)	Chrome molybdenum tubing	Chrome molybdenum tubing
Material (upholstery)	Polyurethane covered polyester	Polyurethane covered polyester





Personalise your ride\*

ANODIC  
BLACK

PEARL  
WHITE

GREY  
METALLIC

BLUE  
CLASSIC

PINE  
GREEN

Actual product colors may vary.

U3 Light Y-front 90°	U3 Light Y-front High 4°	U3 Light Y-front Low
33 - 36 - 39 - 42 - 45 cm	33 - 36 - 39 - 42 - 45 cm	33 - 36 - 39 - 42 - 45 cm
35 - 46 cm	35 - 46 cm	35 - 46cm
Front 47 cm (24" rear wheels) Back 43 cm (24" rear wheels)	Front 51 cm (24" rear wheels) Back 45 cm (24" rear wheels)	Front 44.5 cm (w. 24" rear wheels) Back 40.5 cm (w. 24" rear wheels)
7°	4°	7°
73 cm	82 cm	82 cm
Seat width +21 cm	Seat width +21 cm	Seat width +21 cm
90°	100°	100°
25 - 35 cm	25 - 35 cm	25 - 35 cm
11.5° - (-4.5°)	11.5° - (-4.5°)	11.5° - (-4.5°)
35 - 43 cm	36 - 44 cm	36 - 44 cm
2.2°	2.2°	2.2°
4.0 - 4.2 kg	3.7 - 3.9 kg	3.7 - 3.9 kg
6.3 - 6.5 kg	6.0 - 6.2 kg	6.0 - 6.2 kg
100 kg	100 kg	100 kg
Chrome molybdenum tubing	Chrome molybdenum tubing	Chrome molybdenum tubing
Polyurethane covered polyester	Polyurethane covered polyester	Polyurethane covered polyester

A woman with blonde hair, wearing a beige puffer jacket with a fur-lined hood and black pants, is smiling while riding a motorized bicycle. She is seated in a wheelchair that is attached to the bicycle frame. The bicycle has a black frame with a red seat and handlebars. A young girl in a purple jacket and white helmet stands next to the bicycle, looking towards the camera. The background shows a house with a grey roof and a green lawn.

**Enjoy  
the ride**

Panthera  
**U3 Light**

For experienced manual wheelchair users who need a low-weight, maneuverable mobility device, the Panthera U3 Light offers stability and comfortable posture support.

Pair it with Mototronik\* power assist device for added convenience.

\*Mototronik - see picture

**permobil**

permobil.com | info@permobil.co.uk | T. +44 (0)1484 722 888



PANTHERAU3LIGHTDATASHEETUK221228