

Adjusted by:

Tel:

Save the form to the desktop and
open it with Adobe Reader.

E-mail:

Date:

The price changes depending on
the choices made in the form £:

Standard configuration includes; chassis, seat, calf strap, backrest lock, standard colour option, backrest, carbon rear axle, castor and fork, titanium foot rest or foldable footplates, push to lock S3 brake, Panthera crosslaced steel spoke rear wheel, titanium pushrims and Schwalbe right Run tires.

Choice required:Max. user weights : Seat width **33-42 cm 100 kg**, seat width **45 cm 150 kg**Click on "[blue underlined text](#)" for more info or image**Original design****1) Chassis**[U3](#)Seat width **33-45** cm Seat depth 40 cm

Clear

2) Chassis colour**Black** - Anodic Black**Gray** -Metallic Gray

Clear

3) Seat width[33 cm](#)[36 cm](#)[39 cm](#)[42 cm](#)[45 cm](#)

Clear

4) Seat depth[40 cm](#)

Clear

5) Back rest Height (HH)[INFO](#)[Standard](#)**20 cm****35 cm**[For fixed backrest](#)

523BBHH VER2

25 cm**40 cm**

Clear

30 cm**45 cm**

Clear

6) Rear axle[Carbon std](#)
33023BB[Carbon +15 mm](#)
33023BC*

* seat width BB+1

[Other](#) rear axle or seat height

Clear

Use this rear axle when choosing arm support
at Accessories and/or increased seat height[Increase the seat height](#)**7) Castor and fork**[art.no.](#)[90 mm](#)[120 mm](#)[150 mm](#)

Clear

8) Footrest[Footrest Titanium](#)
36001BB[Titanium, extended, reinforced *](#)

36111BB

* std for S3 Large

[Carbon footplate](#)

65050BB

[Footrest U-shape, small heel bar](#)
3600134 (width 33)[Foldable footplates;](#) [art.no.](#)[Heel strap for footplates](#)

Clear

[Plastic plate for footrest](#)

36010BB

9) Brake[Brake S3](#)

4630000

Clear

10) Rear wheel - without pushrim

[Rear wheel sizes](#)

[24" std](#)

[art.no.](#)

22" [standard](#)

Clear

Tyre

[Schwalbe Right Run](#) 22", 24", 25", 26"

Clear

11) Pushrim

[Titanium](#)

190DD11 (x2 / wheelchair)

Clear

[Aluminum X](#)

1904404 x2 (04=24", 05=25", 06=26")

DD=Wheel diameter

Accessories

Side guard [art.no.](#)

[Side guard S3](#)

[video](#)

Clear

Anti-Tip [Video](#)

[Anti-tip](#) std seat height
4353000

Clear

[Other seat height](#)

[Art.no.](#)

Set-up adjustments

Back rest angle

Note measure (A) for desired back rest angle.
14-28 mm. Standard = 20 mm

A= _____mm

Rear axle adjustment

Note measure (A) for desired balance
Standard = 60 mm. at raised rear axle = 51 mm

A= _____mm

Footrest height - Foot plates

Note measure (A) and (B) for desired footrest

A= _____mm B= _____mm

Extra