

Permobil

Panthera X3



permobil

Strong and ultralight for exceptional maneuverability

Designed for experienced manual wheelchair users who need a dynamic wheelchair to fit into their active lifestyles, the Panthera X3 is one of the lightest wheelchairs on the market. It offers multiple configurable options and increased durability with a modern design that enhances both the form and function, all while prioritizing users comfort.

Key benefits:

- Ultra-light weight
- Reinforced frame for users up to 125 kg
- Easy to transport
- Increased rear axel configurability
- Extremely maneuverable



Specifications

Seat width	330 - 360 - 390 - 420 - 450 mm
Seat depth	350 - 470 mm
Seat height	Front 455 - 470 mm (w. 24" rear wheels) Back 415 - 430 mm (w. 24" rear wheels)
Seat angle	7°
Total length	830 - 845 mm
Total width	Seat width +220 mm
Front frame angle	100°
Backrest height	220 - 350 mm
Backrest angle	14° - (-4.5°)
Footplate distance	360 - 440 mm
Camber	2.2°
Weight without rear wheels	2.3 kg*
Wheelchair weight	5 kg**
Maximum load	125 kg
Material (chassis)	Carbon fibre
Material (upholstery)	Polyurethane covered polyester

*Weight configuration: SW: 330 mm. BRH: 220 mm. Excluding Rearwheels, Brake, Sideguards, cushions and other accessories.

**Weight configuration: SW: 330 mm. BRH: 220 mm. Excluding Sideguards, cushions and other accessories.

Permobil reserves the right to change specifications and design without prior notice.

Personalise your ride

CARBON
FIBRE

permobil

customer_support_uk@permobil.com | permobil.com



GR5096-0.2UK240808

Actual product color may vary.