

M3 Corpus[®]

Made for more.



permobil

M3 Corpus[®]

Reach for more

The M3 Corpus, is made for more. Its compact base offers maneuverability in tight spaces without sacrificing performance. Even more functionality is now possible with enhanced stability in various seating positions. With powerful LED lights, increased driving range and battery life, the colourful new M3 Corpus aims to give you more.

Active Reach

Tilts the seat forward by up to 20° and raises it to extend functional reach, aid in transfers and your daily activities.

Active Height

Enables you to drive at up to 5 km/h while fully elevated and combines seamlessly with *Active Reach* to give you the ultimate advantage in accessing your environment.

Corpus: the ultimate seat system

Improved comfort, pelvic positioning, lateral stability and immersion are all key objectives of the new Corpus cushions.

Corpus cushions feature softer dual-density foam on the backrest, a re-designed pelvic well and Stretch-Air™, a breathable, water-resistant cushion fabric that is machine-washable at 60°C.

More control

Innovative *Smooth* suspension is designed to keep you stable and in control as you traverse the obstacles in your world

More drivability

The new Permobil battery and Permobil VoltPro charger work together to provide you with up to 25% greater driving distance, optimized charging times and longer battery life.



More freedom

MORE COMFORT

The softer dual-density backrest cushion and re-designed pelvic well increase immersion into the seat, giving more comfort *and* stability

MORE SUPPORT

BodiLink® positioning accessories specially designed for the Corpus seating system (option)

MORE COLOUR

A fresh new aesthetic allows you to mix & match shroud/hub insert colour options to match your style.

MORE VISIBILITY

LED lights front and back are twice as bright for improved visibility.

MORE FREEDOM

to move in a greater variety of drive positions with enhanced stability.



M3 Corpus[®] Specifications

Max. speed	10 km/h (6 km/h option)
Driving Range (G34: 65Ah / G24: 85Ah batteries)	30/40 km*
Max. user weight	150 kg
Obstacle climbing (with run up)	50 mm (75 mm)
Min. turning diameter	1120 mm
Pivot width	1100 mm
Suspension	Smooth
Electronics R-Net	120 A
Crash test ISO 7176-19 – 4-point tie-down	Yes – max. 150 kg user weight
Crash test ISO 7176-19 – Dahl docking system	Yes – max. 136 kg user weight
Tested according to EN12182/EN12184	Yes
Max. Length incl. seat	1080 mm
Base width	615 mm
Width incl. seat and armrest	650–790 mm
Height	960–1170 mm
Stowage length/ height	855 mm/823 mm
Weight incl. batteries 65 Ah	176 kg
Seat to floor height (excl. cushion)	450 mm
Seat elevator height (excl. cushion)	450–750 mm
Seat depth	370–570 mm (25 mm increments)
Seat width	420–570 mm (50 mm increments)
Backrest width	360–510 mm (50 mm increments)
Backrest height (excl. seat cushion)	470; 545–670 mm (25 mm increments)
Armrest height	185–320 mm
Distance between armrests	380–480/480–580 mm
Backrest angle adjustment (manual)	90°–135°
Seat plane adjustment – Anterior/posterior tilt	-20° to 50°

*Range may vary depending on user weight and if the wheelchair is used frequently on slopes, rough ground or to climb kerbs often.



LEG REST 90°–180°



POSTERIOR TILT 0 – 50°



SEAT ELEVATOR



ELECTRIC BACK REST
85°–180°



ANTERIOR TILT



More colour

SKY

AQUA

SPRING

CANARY

SUNSET

CHILI

FLAMINGO

VIOLET

MIDNIGHT

ARCTIC

GOLD

BRONZE

PLATINUM



MyPermobil



Longer battery life