

## Q1450 WITH SYNERGY SEATING

### DEALERS DETAILS

|              |           |
|--------------|-----------|
| Dealers Name |           |
| Order Date   |           |
| Address      |           |
|              | Postcode  |
| Phone        | Email     |
| PO No        | Order Ref |

### DELIVERY DETAILS

|                    |  |
|--------------------|--|
| Dealers Address    |  |
| Postcode           |  |
| Quantum Specialist |  |
| Ⓞ Clients Weight   |  |

### BASE MODEL

**£10,195**

- 5 Colour Options
- Height Adjustable Removable Armrests
- Flip-up Centre Mount Foot Platform
- Angle Adjustable Back Support
- Spectrum Foam SPP Seat Cushion
- 3" x 14" Flat Free Drive Tyres
- 2" x 9" Flat Free Rear Castors



- 1.5" Anti-tip Wheels
- Swing-arm Suspension System
- VR2 Control System
- 2 x 70 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight:  
42 st 8lbs. (272 kg)

Standard = Standard Factory Feature  
FOC = Select Instead of Standard Factory Feature

### DRIVE WHEELS

|                                |          |
|--------------------------------|----------|
| 3" x 14" Flat Free Grey Tyres. | Standard |
| 3" x 14" Pneumatic Grey Tyres. | FOC      |

### REAR CASTORS

|                               |          |
|-------------------------------|----------|
| 2" x 9" Flat Free Grey Tyres. | Standard |
| 2" x 9" Pneumatic Grey Tyres. | FOC      |

### COLOURS

|                  |          |
|------------------|----------|
| Back in Black    | Standard |
| Electric Blue    | FOC      |
| Red Corvette     | FOC      |
| Evolution Orange | FOC      |
| Pretty in Pink   | FOC      |

Please download our free PQ-365 app to view colours and other product features in more detail.



## SEAT PACKAGES - SEAT WIDTH

|               |
|---------------|
| 20"-24"W.     |
| 24"-26"W      |
| 26"-28"W      |
| 28"-30"W      |
| 30"-32"W      |
| Set width at: |

## SEAT PACKAGES - SEAT DEPTH

|               |
|---------------|
| 20"-24"D.     |
| 24"-28"D      |
| Set depth at: |

## STEALTH SEAT CUSHIONS

|   |     |
|---|-----|
| No Seat Cushion Required  | FOC |
| Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i> | STD |
| Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>  | FOC |
| TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>                  | FOC |

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

## SEAT POSITIONING

|  |                   |
|--|-------------------|
| Static Seat.   | Standard          |
| HD Powered Seat Tilt. <i>0° to 45°. Uses 1 power option.</i> | £1600             |
| Static Seat to Floor Height. <i>From seat pan.</i>           | 16.5" 17.5" 18.5" |
| HD Tilt Seat to Floor Height. <i>From seat pan.</i>          | 19.5" 20.5"       |

## SYNERGY BACK SUPPORTS

|   |          |
|---|----------|
| Synergy Visco Padded Back. <i>Height range 17" to 25" tall from seat pan.</i> | Standard |
| HD Sling Back. <i>Height range 14" to 25" tall from seat pan.</i>             | £150     |
| Enter height preference here.   |          |

## BACK SUPPORT POSITIONING

|  |                               |
|--|-------------------------------|
| Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>                      | Standard                      |
| 1" Endo Cane Mounts. <i>Reduces backrest width by 2" compared to seat width.</i> | £40                           |
| 2" Endo Cane Mounts. <i>Reduces backrest width by 4" compared to seat width.</i> | £40                           |
| Back Cane Height.  | 16"H. 17"H. 18"H. 19"H. 20"H. |
| <i>From seat pan.</i>  | 21"H. 22"H. 23"H. 24"H. 25"H. |

## ARM SUPPORTS

|  |   |          |
|--|---|----------|
| HD Single Post Tool Height Adjustable Removable.   | Standard  |          |
| Bariatric Single Post Tool Height Adjustable Removable.  | £150  |          |
| HD & Bariatric Height Package. <i>From seat pan.</i>   | 8"-12" FOC 10"-14" Standard                         |          |
| Cantilever Flip-Up Armrests. <i>Adj. range 11"-17". Not available with HD sling back support. Increases overall width by 7".</i> | £150  |          |
| Arm Support Pads.  | 10" Standard. <i>Not available with Bariatric.</i>  | FOC      |
|  | 14" Standard. <i>Not available with Bariatric.</i>  | Standard |
|  | 10" Waterfall. <i>Not available with Bariatric.</i> | FOC      |
|  | 14" Waterfall. <i>Not available with Bariatric.</i> | FOC      |
|  | Bariatric. <i>Bariatric arm supports only.</i>      | FOC      |

## CENTRE MOUNTED LEG SUPPORTS

|  |   |          |  |         |     |          |
|--|---|----------|--|---------|-----|----------|
| 17"W x 13"D Base Mounted Flip-up Foot Platform. <i>Mounted to chassis. Will not move with power tilt. Ignore other centre mount options if selected.</i> |   |          |  |         |     | FOC      |
| 17"W x 16"D Base Mounted Flip-up Foot Platform. <i>Mounted to chassis. Will not move with power tilt. Ignore other centre mount options if selected.</i> |   |          |  |         |     | Standard |
| Manually Adjustable Centre Mount Foot Platform.  |   |          |  |         |     | £244     |
| Manual Adjustable Lower Leg Extensions.  |   | 10"-16"  | Standard   | 13"-19" | FOC |          |
| <i>From seat pan.</i>  |   |          |  |         |     |          |
| Footplates.  | 10"W x 11"D Split Footplates.   |          |  |         |     | Standard |
|  | 17"W x 13"D One Piece Footplate.  |          |  |         |     | FOC      |
|  | 17"W x 16"D One Piece Footplate.  |          |  |         |     | FOC      |
| Calf Support.  | Individual Calf Pads. 6"W x 6"H.  | Standard | Calf Panel. 18"W x 8"H.  |         | £42 |          |
|  | Reverse Dartex Calf Pad Covers.   | £42      | Reverse Dartex Calf Panel Cover.   |         | £62 |          |
| Footplate Accessories.   | 4" Width Extension Kit. <i>Not available with 2" length kit or 10"W x 6"D footplates.</i> |          |  |         |     | £16      |
|  | 2" Length Extension Kit. <i>Not available with 4" width kit or 10"W x 6"D footplates.</i> |          |  |         |     | £16      |
|  | Footplate Connector.  | £10      | Outer Edge Supports.   |         | £10 |          |
|  | Wheel Bumpers.  | £10      | Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. |         |     |          |

## INDIVIDUALLY MOUNTED LEG SUPPORTS

|                                     |            |          |             |          |             |      |
|-------------------------------------|------------|----------|-------------|----------|-------------|------|
| 70° HD Drop-in Swing-away Legrests. |            |          |             |          |             | £340 |
| 70° HD Lower Extensions.            | 8.5"-12"   | FOC      | 13"-18"     | Standard | 16"-21"     | FOC  |
| <i>From seat pan.</i>               |            |          |             |          |             |      |
| 70° HD Footplates.                  | 6"W x 8"D. | Standard | 6"W x 11"D. | FOC      | 8"W x 11"D. | FOC  |

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops.

|  |             |      |              |      |              |
|--|-------------|------|--------------|------|--------------|
| Stealth Residual Limb Support (each). <i>Includes a 60-90° Swing-away leg support hanger for mounting.</i> | <b>LEFT</b> | £300 | Limb Support | £300 | <b>RIGHT</b> |
|--|-------------|------|--------------|------|--------------|

## STEALTH HEAD SUPPORTS

|   |      |
|---|------|
| UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>           | £219 |
| UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>          | £219 |
| UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>          | £219 |
| UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i> | £327 |
| Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.          | £88  |

Please visit [stealthproducts.com](http://stealthproducts.com) for alternative head supports.

## PELVIC SOLUTIONS

|   |             |     |             |          |             |          |
|---|-------------|-----|-------------|----------|-------------|----------|
| Auto Buckle Lap Belt, Push Button Release.      |             |     |             |          |             | Standard |
| Auto Buckle Lap Belt Length.                    | 38"         | FOC | 50"         | FOC      | 60"         | Standard |
|   | 70"         | FOC | 80"         | FOC      | 90"         | FOC      |
| Retractable Lap Belt, Push Button Release.      |             |     |             |          |             | £156     |
| Stealth Positioning 2-Point Padded Pelvic Belt. | Medium      | £78 | Single Pull | Standard | Push Button | Standard |
|   | Large       | £78 | Dual Pull   | FOC      |             |          |
|   | Extra Large | £78 | Rear Pull   | FOC      |             |          |

## NON-EXPANDABLE CONTROLS

|   |   |          |
|---|---|----------|
| VR2 Controls. <i>With no power options.</i> | Standard  |          |
| VR2 Controls. <i>2 power option ready.</i>  | £125  |          |
| Control Method.                             | Hand Control. <i>Select position in Hand Control Mounts below.</i>          | Standard |
|   | Attendant Control. <i>Hand control mounted at rear.</i>                     | £125     |
|   | Dual Controls. <i>Hand control and secondary attendant control at rear.</i> | £430     |

## EXPANDABLE CONTROLS

Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth. Q-Logic 3 Enhanced Display is supplied with built-in Bluetooth and Infrared.

|  |       |
|--|-------|
| Q-Logic 3e Controls. <i>No power options.</i>  | £225  |
| Q-Logic 3e Controls. <i>2 power option ready.</i>  | £300  |
| Q-Logic 3 EX Controls. <i>No power options.</i>  | £275  |
| Q-Logic 3 EX Controls. <i>2 power option ready.</i>  | £350  |
| Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches &amp; flexible display mount.</i>   | £1074 |
| IR Blaster. <i>Infrared module for Q-Logic 3 EX hand control.</i>  | £50   |
| SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i> | £316  |

iAccess OLED 4 Button Switch Module.  
*Mounted to opposite armrest unless specified.*

£300

Push Button iAccess  
Toggle Lever iAccess

If iAccess is selected choose either  
Push Button or Toggle lever option.

|                 |  |          |
|-----------------|--|----------|
| Control Method. | Hand Control. <i>Select position in Hand Control Mounts below.</i>                             | Standard |
|                 | Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i> | -£100    |
|                 | Attendant Control. <i>Hand control mounted at rear.</i>  | £125     |
|                 | Dual Controls. <i>Hand control and secondary attendant control at rear.</i>                    | £430     |

## HAND CONTROL MOUNTS

|   |                 |          |             |          |                   |     |
|---|-----------------|----------|-------------|----------|-------------------|-----|
| Fixed.  | Standard        |          |             |          |                   |     |
| Swing-away. <i>Swings in &amp; out with height, tension and pivot adjustment.</i> | £75             |          |             |          |                   |     |
| Retract 4 Swing-away. <i>Full 180° Movement.</i>                                  | £150            |          |             |          |                   |     |
| Mounting Position.  | Left            | FOC      | Right       | Standard |                   |     |
| Hand Control Joystick Handle.   | Standard Carrot | Standard | T-Bar       | £35      | Extended I-Handle | £35 |
|   | Mushroom        | £35      | Medium Ball | £35      | 3" Goal Post      | £50 |
|   | Chin Cup        | £35      | Soft Ball   | £35      | 4" Goal Post      | £50 |

## SPEED PROGRAM & LIGHTS PACKAGE

|   |          |
|---|----------|
| 4 Mph Program. <i>Lights and indicators not included.</i>                             | Standard |
| 4 Mph Program with Lights and indicators.   | £450     |
| Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i> | £260     |

## ACCESSORIES

|  |      |
|--|------|
| Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i> | £379 |
| Wooden Trim Kit for Tray. <i>Only with swing-away tray table.</i>  | £100 |

## NOTES & SPECIAL REQUESTS

## PRICE

**Q1450 WITH  
SYNERGY SEATING**



**BASE PRICE** **£10,195**

**COST OPTIONS**

**TOTAL**

5



**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk  
**RETAIL UK - 24/06/23 V1**



V1: 14/06/2022