



Q4

The Quantum Q4 is a compact Mid-Wheel 6 Drive power chair, designed for both indoor and outdoor use. Featuring a tight turning radius of 508 mm (20"), the Q4 is easily manoeuvrable and can adapt to any environment.

QUANTUM[®]

Specifications

WEIGHT CAPACITY	21 stone 6 lbs. (136 kg)
MAXIMUM SPEED ¹	Up to 6 mph
GROUND CLEARANCE ²	63 mm (2.5") at suspension arm 76 mm (3") at center of frame
TURNING RADIUS ²	508 mm (20")
OVERALL SIZE ^{2, 6}	Length: 901 mm (35.5") Width: 610 mm (24")
DRIVE WHEELS	355 mm (14") pneumatic or solid available
CASTER WHEELS	Front: 152 mm (6") solid Rear: 152 mm (6") solid 203 mm (8") solid
SUSPENSION	ATX Suspension (Active-Trac® with long trail arm)
COMPONENT WEIGHTS ^{5, 7}	Base: 137 lbs (62 kg) Batteries: 38 lbs (17.24 kg) each
DRIVETRAIN	Two motor, Mid-Wheel 6 Drive
BRAKES	"Intelligent Braking" electronic regenerative, disc park brake
BATTERIES ³	Two 12-volt, deep-cycle, 50 Ah batteries
BATTERY CHARGER	8-amp, Off-board
MAXIMUM OBSTACLE CLIMBING ABILITY	76 mm (3")

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability

Features

- 508 mm (20") turning radius
- Maximum speed up to 6 mph
- Crash tested
- Mid-Wheel 6 design
- Improved ATX Suspension
- ION Seat
- TB Flex Seat

Options

- ION Seating Options:
16" x 16", 406 mm x 406 mm
18" x 18", 457 mm x 457 mm
- TB Flex Seating Options:
15" x 15", 381 mm x 381 mm
21" x 20", 533 mm x 508 mm
- Power and Manual Options:
Fixed Tilt (tool adjustable) | Manual Recline (folding).
Power Tilt | Power Recline | Power Elevating Leg Rests.

Colours



Blue



Red



Black