



Q6 edge[®] 2.0

With features like standard 4 mph motors and advanced responsiveness, the Q6 Edge[®] 2.0 is engineered to meet the performance needs of the most active user. The Q6 Edge 2.0 accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value, and is available with our iLevel[®] technology.

QUANTUM[®]

Specifications

WEIGHT CAPACITY ⁴	21 stone 6 lbs. (136 kg)
MAXIMUM SPEED ¹	Up to 6 mph Up to 3.5 mph w/iLevel activated
GROUND CLEARANCE ²	69 mm (2.7")
TURNING RADIUS ²	521 mm (20.5")
OVERALL SIZE ^{2, 6}	Length: 902 mm (35.5") without front rigging Width: 622 mm (24.4") 735 mm (28.9") with armrests
DRIVE WHEELS	356 mm (14") pneumatic or solid available
CASTER WHEELS	Front: 152 mm (6") solid Rear: 152 mm (6") solid (8" solid front & rear optional)
SUSPENSION	ATX Suspension (Active-Trac® with extra stability)
COMPONENT WEIGHTS ^{5, 7}	Base: 138.6 lbs (63 kg) Batteries: 38 lbs (17.24 kg) each
DRIVETRAIN	Two motor, Mid-wheel 6@ drive
BRAKES	"Intelligent Braking" electronic regenerative, disc park brake
BATTERIES ³	Two 12-volt, deep-cycle 62 Ah batteries
BATTERY CHARGER	8-amp, Off-board
MAXIMUM OBSTACLE CLIMBING ABILITY	90mm (3.5")

SEATING SIZES

TRU-Balance® 3 Static Seat	Width: 12-22" Depth: 12-22"
TRU-Balance® 3 Power Tilt	Width: 12-22" Depth: 12-22"
TRU-Balance® 3 Tilt and iLevel®	Width: 12-22" Depth: 12-22"

SEAT-TO-FLOOR HEIGHTS

TRU-Balance® 3 Static Seat ⁵	15.875-17.875"
TRU-Balance® 3 Power Tilt ⁵	17.5-19.5"
TRU-Balance® 3 Tilt & iLevel® ⁵	17.5-19.5"
TRU-Balance® Synergy Seating	15.875-19.875"

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability

Features

- Standard 4 mph motors (6 mph optional) offer advanced performance, precise control and responsiveness via innovative single-stage gearbox
- Bold aesthetics
- iLevel® option available with up to 12" of seat elevation at walking speed (up to 3.5 mph)
- iLevel is retrofittable on the Q6 Edge® 2.0

Options

- iLevel® option allows 12" of seat elevation while at walking speed (3.5 mph)
- TRU-Balance® 3 Seating
- Synergy Seating
- TRU-Balance 3 Power Tilt
- TRU-Balance 3 Power Recline
- LED fender lights (optional)
- 8" solid castors
- Powered Elevating Leg Rests or Footboard
- NE+
- Q-Logic 3 EX
- VR2
- R-Net

Colours*



* Please contact your representative for the full range of options and colours.