

Q6 EDGE 2.0 WITH SYNERGY SEATING

DEALERS DETAILS

Dealers Name _____
 Order Date _____
 Address _____
 _____ Postcode _____
 Phone _____ Email _____
 PO No _____ Order Ref _____

DELIVERY DETAILS

Dealers Address _____
 _____ Postcode _____
 Quantum Specialist _____
 Ⓢ Clients Weight _____

BASE MODEL

£4,825

- 13 Colour Options
- Height Adjustable Removable Armrests
- Centre Mount Foot Platform
- Synergy Seating
- 14" Flat Free Drive Tyres
- 2" x 6" Flat Free Front Castors



- 2" x 6" Flat Free Rear Castors
- ATX Suspension System
- NE+ Control System
- 2 x 62 ah Batteries
- 4 Mph (6 Km/h) Motors
- Max User Weight: 21 st 6lbs. (136 kg)

Standard = Standard Factory Feature
 FOC = Select Instead of Standard Factory Feature

DRIVE WHEELS

3" x 14" Flat Free Black Tyres.	Standard
3" x 14" Pneumatic Black Tyres.	FOC
4" x 14" Pneumatic Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£220
4" x 14" Flat Free Black Tyre, Knobby, Extended Traction. <i>Increases overall base width by 1.5".</i>	£320

FRONT CASTORS

2" x 6" Flat Free Black Tyres.	Standard
2" x 8" Flat Free Black Tyres.	£364

COLOURS

Back in Black	Standard	Lemon Crush*	FOC
Green Machine	FOC	Night Train*	FOC
Midnight Blue	FOC	Evolution Orange*	FOC
Electric Blue	FOC	Military Green*	FOC
Deep Purple	FOC	Grey Street*	FOC
Red Corvette	FOC	White Satin*	FOC
Pretty in Pink	FOC	Custom Hydro Dip*	£600

* Extended lead-time.

Please download our free PQ-365 app to view colours and other product features in more detail.



SEAT PACKAGES - SEAT WIDTH

14"-18"W.
16"-20"W
20"-22"W
Set width at:

SEAT PACKAGES - SEAT DEPTH

14"-18"D.
15"-20"D
20"-22"D
Set depth at:

STEALTH SEAT CUSHIONS

No Seat Cushion Required	FOC
Simplicity G: General Use, Fluid Resistant Stealth-Tek Cover. <i>158.7 Kg weight limit.</i>	Standard
Zen SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>160 Kg weight limit.</i>	£75
TRU-Comfort 2 SPP: Skin Protection, Positioning, 4-way Stretch Coolcore Cover. <i>158.7 Kg weight limit.</i>	£100
Glacial SP: Skin Protection, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£115
Solution SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£115
Essence SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover. <i>158.7 Kg weight limit.</i>	£135
Premiere P: Positioning, 4-way stretch, Coolcore cover. <i>158.7 Kg weight limit.</i>	£135
Spectrum Foam SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190
Spectrum Gel SPP: Skin Protection, Positioning, 4-way Stretch, Coolcore Cover, Urethane Liner. <i>158.7 Kg weight limit.</i>	£190

Please note that seat cushions may not be available in all width and depth combinations. Refer to our Stealth Cushion Order Forms for size availability.

SEAT POSITIONING

Static Seat.				Standard
Static Seat to Floor Height. <i>From seat pan.</i>	15.875"	16.875"	17.875"	

BACKREST SPECIFICATIONS

Back Canes System.							Standard
Back Cane Height.	16"H.	17"H.	18"H.	19"H.	20"H.		
<i>From seat pan.</i>	21"H.	22"H.	23"H.	24"H.	25"H.		

BACK SUPPORT

Active Contour Back Support with TRU-Comfort 4 Back Cushion and Wedge Inserts.				Standard
Tension Adjustable Back Upholstery (TAB). <i>May increase lead-time.</i>				£120
Back Support Height.	16"H.	17"H.	18"H.	
<i>From seat pan. Back support height must be similar to back cane height. Ignore section if using an ADI back support.</i>	19"H.	20"H.	21"H.	
	22"H.	23"H.	24"H.	
	25"H.			
ADI Back Support. <i>Skin protection, 4 way stretch, Coolcore cover. Mounts to back canes with cane clamps.</i>				POA
Enter ADI Q-Back 2 part number here for quotation.				
Privacy Flap. <i>12.75"W. Closes gap between back support and seat pan. Not available with a TAB.</i>				£50

BACK SUPPORT POSITIONING

Tool Angle Adjustable. <i>85° to 120° in 5° increments.</i>	Standard
Flip-forward Back Canes. <i>Not suitable for occupied transport or with Endo mounts.</i>	£200
2" Endo Cane Mounts. <i>Not suitable for occupied transport. 1" each side. Reduces back width by 2". Minimum seat width 12".</i>	£40
4" Endo Cane Mounts. <i>Not suitable for occupied transport 2" each side. Reduces back width by 4". Minimum seat width 14".</i>	£40

ARM SUPPORTS

Single Post Drop-in Removable. <i>Adjustment range 11"-17" from seat pan. Minimum seat depth 14".</i>	Standard
Singe Post HD Tool Adjustable Removable. <i>Adjustment range 10"-14" from seat pan. Minimum seat depth 14".</i>	£220
2-Post Flip-up Removable. <i>Adjustment range 9" to 13" from seat pan. Minimum seat depth 16".</i>	£145
Paediatric Single Post Drop-in Removable. <i>Adjustment range 7.75"-9.75" from seat pan. 10" pads only.</i>	FOC
Cantilever Flip-Up. <i>Not available with a TAB. Increases overall width by 7". Height adj. 11"-17". Not available with Flip Forward Backrest.</i>	£280
Arm Support Pads.	10" Standard Standard 14" Standard FOC
	10" Waterfall. FOC 14" Waterfall. FOC
	10" Falcon Gel. <i>Includes mounting plate.</i> £80 14" Falcon Gel. <i>Includes mounting plate.</i> £80

CENTRE MOUNTED LEG SUPPORTS

Base Mounted Flip-up Foot Platform. <i>Mounted to chassis. Will not move with powered lift or tilt.</i>	FOC		
Manually Adjustable Centre Mount Foot Platform. <i>Minimum seat width 14".</i>	Standard		
Paediatric Adjustable Centre Mount Foot Platform. <i>12"-14" seat widths only.</i>	FOC		
Power Elevating Articulating Centre Mount Foot Platform. <i>Minimum seat width 14". Uses 1 power option.</i>	£700		
Power Elevating Articulating & Lower To Floor Centre Mount Foot Platform. <i>Minimum seat width 14". Lowers to the floor to aid standing transfers. Uses 2 power options.</i>	£1112		
Lower Leg Extensions. <i>From seat pan.</i>	8"-12". <i>Paediatric centre mount only.</i> FOC	10"-16". <i>Manually adj. & powered elevating only</i> FOC	13"-19". <i>Manually adj. & powered elevating only</i> STD
Footplates.	10"W x 6"D Split Footplates. <i>Paediatric centre mount only.</i>		FOC
	10"W x 8"D Split Footplates.	FOC	10"W x 11"D Split Footplates. Standard
Calf Support.	Individual Calf Pads. <i>6"W x 6"H.</i>	Standard	Calf Panel. <i>18"W x 8"H.</i> £42
	Reverse Dartex Calf Pad Covers.	£42	Reverse Dartex Calf Panel Cover. £62
Footplate Accessories.	4" Width Extension Kit. <i>Not available with 2" length kit or 10"W x 6"D footplates.</i>		£16
	2" Length Extension Kit. <i>Not available with 4" width kit or 10"W x 6"D footplates.</i>		£16
	Footplate Connector.	£10	Outer Edge Supports. £10
	Wheel Bumpers.	£10	Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information.

INDIVIDUALLY MOUNTED LEG SUPPORTS

70° Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 12" to 15". Minimum seat width 15".</i>	FOC
90° Paediatric Swing-Away Leg Supports. <i>4"W x 6"D angle & depth adj. footplates. Adjustment range 7" to 11".</i>	£110
Manual Elevating Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 15" to 20". Min. seat width 15".</i>	£130
Power Elevating Swing-Away Leg Supports. <i>Flip-up footplates. Adj. range 14.5" to 19.5". Min. seat width 16". Uses 2 power options.</i>	£650

INDIVIDUALLY MOUNTED CONFIGURABLE LEG SUPPORTS

60° to 90° Tool Adjustable Swing-Away Leg Supports. *Select lowers and footplates below. Minimum seat width 15"* £110

60° TO 90°	Lower Leg Extensions. <i>From seat pan.</i>	11"-13"	FOC	13.5"-16.5"	Standard	15.5"-19.5"	FOC
	Footplate Type.	Flip-up Extruded. <i>8.5"W x 6"D.</i>	STD	Angle & Depth Adjustable. <i>Select size below.</i>	£115	Multi-axis. <i>Select size below.</i>	£263
	Footplate Size. <i>Angle adj. & Multi-axis only.</i>	4"W x 6"D	FOC	5"W x 6"D	FOC	5"W x 8"D	FOC
		6"W x 8"D	STD	6"W x 11"D	FOC	8"W x 11"D	FOC
Accessories.	Wheel Bumpers						£10

Power Elevating Articulating Swing-Away Leg Supports. *Inc. calf pads. Adj. range 15" to 20". Min. seat width 15". Uses 2 power options.* £1250

POWER ELEVATING ARTICULATING	Lower Leg Extensions. <i>From seat pan.</i>	11"-13"	FOC	14.5"-18"	Standard	17.5"-21"	FOC
	Footplate Type.	Angle & Depth Adjustable. <i>Select size below.</i>					Standard
	Footplate Size.	4"W x 6"D	FOC	5"W x 6"D	FOC	5"W x 8"D	FOC
		6"W x 8"D	Standard	6"W x 11"D	FOC	8"W x 11"D	FOC
Accessories.	Small Calf Pads		Standard	Wheel Bumpers		£10	
	Large Calf Pads		FOC	Gel Padded Knee Adductors		£60	

Certain configurations may result in reduced ground clearance below footplates. Contact Quantum Technical for information. Adjustment ranges above are from seat pan to footplates. Individual legrests are supplied with heel loops.

Stealth Residual Limb Support (each). *Includes a 60-90° Swing-away leg support hanger for mounting.*

LEFT

£300

Limb Support

£300

RIGHT

STEALTH HEAD SUPPORTS

UniLink 8" Comfort Plus. <i>HMO475-17M / CP08CCFPM.</i>	£219
UniLink 10" Comfort Plus. <i>HMO475-17M / CP10CCFPM.</i>	£219
UniLink 14" Comfort Plus. <i>HMO475-17M / CP14CCFPM.</i>	£219
UniLink 14" Adjustable Comfort Plus. <i>HMO475-17M / ACP150M.</i>	£327
Upgrade Mounting Hardware from HMO475-17M to UHB3D33200.	£88

Please visit stealthproducts.com for alternative head supports.

STEALTH TRUNK SOLUTIONS

Thoracic Lateral Supports.	LEFT	£150	Lateral Thoracic Swing-away Hardware.		£150	RIGHT		
		£40	3" H x 5" L Curved Lateral Pad.		£40			
		£40	4" H x 6" L Curved Lateral Pad.		£40			
Thoracic Lateral Support Extras.	Additional 1.5" Link (each).					£52		
	Additional 2" Link (each).					£52		
	Additional 3" Link (each).					£52		
4-Point Harness.	Extra Small	£137	Classic	STD*	Dynamic	STD*	Vinyl Face	FOC
	Small	£137		Contoured		FOC		Structured
	Medium	£137						
	Large	£137						
	Extra Large	£137						
Chest Strap.	Extra Small	£87	Small	£87	Side Release		Standard	
	Medium	£87	Large	£87	Velcro Release		£19	
	Extra Large	£87						

*STD = Standard option **Classic Only

STEALTH PELVIC & THIGH SOLUTIONS

UniLink Thigh Supports.	LEFT	£176	UniLink Hip Guide Long Removable.	£176	RIGHT	
		£176	Unilink Hip Guide Short Removable.	£176		
		£48	4x4 Pelvic Thigh Pad.	£48		
		£48	4x8 Pelvic Thigh Pad.	£48		
		£48	4x12 Pelvic Thigh Pad.	£48		
Knee Adductors.	LEFT	£156	Fixed Cane Mounted.	£156	RIGHT	
		£170	Swing-Away Track Mounted.	£170		
		£36	3" Round Knee Pad.	£36		
		£36	4" Round Knee Pad.	£36		
Pelvic Belts.	Auto Buckle Lap Belt, Push Button Release.			Standard		
	Retractable Lap Belt, Push Button Release.			£156		
Positioning 2-Point Padded Pelvic Belt.	Extra Small	£78	Single Pull	Standard	Side Release	Standard
	Small	£78		FOC		FOC
	Medium	£78	Dual Pull	FOC	Push Button	FOC
	Large	£78				
	Extra Large	£78				
Stealth Positioning 4-Point Linear Padded Pelvic Belt, Single Pull.	Extra Small	£120	Small	£120	Push Button	Standard
	Medium	£120	Large	£120		
Medial Supports (Pommels).	Fixed Flip-down Hardware. <i>MTHW-660.</i>		£180	Small Wedge. <i>MTMD-12.</i>		£103
	Medium Wedge. <i>MTMD-13.</i>		£103	Large Wedge. <i>MTMD-14.</i>		£103
	Adjustable Wedge. <i>MTMD-42.</i>		£199	Please visit stealthproducts.com for alternative lateral, pelvic & lower extremity solutions.		

NON-EXPANDABLE CONTROLS

NE + Controls. <i>With no power options.</i>		Standard
NE + Controls. <i>2 power option ready.</i>		£125
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard
	Attendant Control. <i>Hand control mounted at rear.</i>	£100
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430

EXPANDABLE CONTROLS

Q-Logic 3e and 3 EX hand controls are supplied with built-in Bluetooth. Q-Logic 3 Enhanced Display is supplied with built-in Bluetooth and Infrared.

Q-Logic 3e Controls. <i>No power options.</i>	£225
Q-Logic 3e Controls. <i>2 power option ready.</i>	£300
Q-Logic 3 EX Controls. <i>No power options.</i>	£275
Q-Logic 3 EX Controls. <i>2 power option ready.</i>	£350
Q-Logic 3 Enhanced Display. <i>Includes 2 egg switches & flexible display mount.</i>	£1074
IR Blaster. <i>Infrared module for Q-Logic 3 EX hand control.</i>	£50
SCIM Module Kit. <i>Allows Q-Logic 3 EX joystick display to be used for speciality controls. Not available with 3e controls.</i>	£316

iAccess OLED 4 Button Switch Module. <i>Mounted to opposite armrest unless specified.</i>	£300	<table border="0"> <tr> <td>Push Button iAccess</td> <td rowspan="2">If iAccess is selected choose either Push Button or Toggle lever option.</td> </tr> <tr> <td>Toggle Lever iAccess</td> </tr> </table>	Push Button iAccess	If iAccess is selected choose either Push Button or Toggle lever option.	Toggle Lever iAccess
Push Button iAccess	If iAccess is selected choose either Push Button or Toggle lever option.				
Toggle Lever iAccess					
Control Method.	Hand Control. <i>Select position in Hand Control Mounts below.</i>	Standard			
	Omit Hand Control. <i>Often selected when ED is used as master remote with special device.</i>	-£100			
	Attendant Control. <i>Hand control mounted at rear.</i>	£125			
	Dual Controls. <i>Hand control and secondary attendant control at rear.</i>	£430			

HAND CONTROL MOUNTS

Fixed.						Standard
Swing-away. <i>Swings in & out with height, tension and pivot adjustment.</i>						£75
Retract 4 Swing-away. <i>Full 180° Movement.</i>						£150
Mounting Position.	Left	FOC	Right	Standard		
Hand Control Joystick Handle.	Standard Carrot	Standard	T-Bar	£35	Extended I-Handle	£35
	Mushroom	£35	Medium Ball	£35	3" Goal Post	£50
	Chin Cup	£35	Soft Ball	£35	4" Goal Post	£50

SPEED PROGRAM & LIGHTS PACKAGE

4 Mph Program. <i>Lights and indicators not included.</i>	Standard
4 Mph Program with Lights and indicators.	£450
6 Mph Program with Lights and indicators.	£450
Lights & Indicator Surcharge. <i>Required if no power options have been selected.</i>	£260
Accu-Trac Motors. <i>Requires Q-Logic 3 EX. Delivers enhanced traction and tracking, especially with switched drive controls. Only available with 6 mph speed program and lights package.</i>	£450
LED Fender Lights. <i>Includes single switch button mounted to side rail. Not road legal.</i>	£190

ACCESSORIES

USB Mobile Device Charger. <i>Please specify mounting position. i.e arm support, backrest, seat rail. Not available with VR2 or R-Net.</i>	£110
Glove Box.	£30
Hydration System.	£125
Oxygen Holder.	£50
Stealth Tablet Holder.	£342
Magnetic to XLR Charger Adapter. <i>Only with Q-Logic 3e.</i>	£49
Swing-Away Tray Table. <i>Supplied with relevant matching cut-out for joystick, or no cut out if attendant control only.</i>	£379
Wooden Trim Kit for Tray. <i>Only with swing-away tray table.</i>	£100
Prepare Chassis for Dahl Docking System. <i>Dahl system not included.</i>	£100

NOTES & SPECIAL REQUESTS

PRICE

Q6 EDGE 2.0 WITH SYNERGY SEATING



BASE PRICE **£4,825**

COST OPTIONS

TOTAL



Quantum®
 32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
RETAIL UK - 15/06/23 V1



Crash Tested



V1: 14/06/2022