



Q6 edge[®] Z

The Q6 Edge[®] Z utilizes powerful 4-pole motors for maximum torque and Mid-wheel 6[®] design for aggressive outdoor performance and precision indoor maneuverability. The Q6 Edge[®] Z accepts a broad range of exclusive seating and electronic options for unmatched adaptability and rehab capability, and is available with our iLevel[®] technology.

QUANTUM[®]

Specifications

WEIGHT CAPACITY ⁴	21 stone 6 lbs. (136 kg)
MAXIMUM SPEED ¹	Up to 6 mph (standard) Up to 4 mph (optional) Up to 8 mph (optional)
GROUND CLEARANCE ²	76 mm (3")
TURNING RADIUS ²	501 mm (19.75")
OVERALL SIZE ^{2, 6}	Length: 904 mm (35.6") without front riggings 1127 mm (44.4") with front riggings Width: 666 mm (26.25")
DRIVE WHEELS	356 mm (14") pneumatic or solid available
CASTER WHEELS	Front: 152 mm (6") solid Rear: 152 mm (6") solid
SUSPENSION	ATX Suspension (Active-Trac® with extra stability)
COMPONENT WEIGHTS ^{5, 7}	Base: 143.5 lbs (65 kg) Batteries: 52 lbs (23.58 kg) each
DRIVETRAIN	4-pole, Mid-wheel 6@ drive
BRAKES	"Intelligent Braking" electronic regenerative, disc park brake
BATTERIES ³	Two 12-volt, deep-cycle 73 Ah (standard)
BATTERY CHARGER	8-amp, Off-board
SEATING SIZES	
TRU-Balance® 3 Static Seat	Width: 12-24" Depth: 12-24"
TRU-Balance® 3 Power Tilt	Width: 12-24" Depth: 14-24"
TRU-Balance® 3 Tilt and iLevel®	Width: 12-24" Depth: 12-24"
SEAT-TO-FLOOR HEIGHTS	
TRU-Balance® 3 Static Seat ⁵	16.5-18.5"
TRU-Balance® 3 Power Tilt ⁵	18.125-20.125"
TRU-Balance® 3 Tilt & iLevel® ⁵	19.75-21.375"

1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.

2) Without front riggings

3) Battery capacity may change when adding power positioning systems.

4) Weight capacity is dictated by the selection of seating.

5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.

6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.

7) Battery weight may vary +/- 3 lbs. based on manufacturer.

8) All Quantum Rehab® units can be ordered with specialty controls.

9) Call for availability

Features

- Mid-Wheel 6@ allows six wheels on the ground for maximum stability
- Low impact OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- iLevel® option available with up to 12" of seat elevation at walking speed (up to 3.5 mph)
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters and semi-independent rear casters beam for enhanced performance over more varied terrain

Options

- iLevel® 12" power adjustable seat height
- Q-Logic 3 EX Drive Controls
- TRU-Balance® 3 Seating
- 8 mph motor package
- TRU-Balance® 3 Power Tilt
- TRU-Balance® 3 Power Recline

Colours*



Pretty
in Pink

Lime
Light

Yellow
Blaze

Deep
Purple

Orange
Crush



Electric
Blue



Red
Corvette



Midnight
Blue



White
Satin



Grey
Street



Back
in Black

* Please contact your representative for the full range of options and colours.