



RETAIL

SPECTRUM GEL SPP CUSHION

DEALERS DETAILS

Dealers Name _____

Order Date _____

Address _____

Postcode _____

Email _____

Phone _____ Fax _____

DELIVERY DETAILS

Dealers Address _____

Postcode _____

PO No _____

Order Type _____ Quote _____ Order _____

DETAILS

This cushion offers skin protection and positioning with a Coolcore, 4-way stretch, breathable cover and an urthane liner.

SunMate FRG foam top (SunMate FRG XXSoft (blue); Medsoft (green)) with multi-layered foam insert or twin cell ultra-lightweight viscous gel insert in ischial well, sacral relief



Weight capacities
 158.75 kg: 16"-20" width,
 358.33 kg: 22" width+
 Overall weight: 1.3 kg
 Overall height 3"
 14x18" up to 32x28"
 Odd sizes in 17" up to 19"
 Custom sizes available

CUSHION

| | | ADD WIDTH | ADD DEPTH | |
|-----------------------------|---------------------|----------------------|----------------------|------|
| Spectrum Gel SPP. 16"-20"W. | ST-SPCGSPP | <input type="text"/> | <input type="text"/> | £265 |
| Spectrum Gel SPP. 22"-24" W | ST-SPCGSPP | <input type="text"/> | <input type="text"/> | £310 |
| Cushion Cover. | Coolcore. | SPCGSPP-NBLK | <input type="text"/> | STD |
| | Mesh. | SPCGSPP-MBLK | <input type="text"/> | FOC |
| | Stealth Tek. | SPCGSPP-TBLK | <input type="text"/> | FOC |
| Fluid Proof Liner. | Urethane. | URL1709070 | <input type="text"/> | STD |
| | Polyurethane. | PUL1709030 | <input type="text"/> | £48 |
| Extra Covers And Liners. | Coolcore Cover. | SPCGSPP-NBLK | <input type="text"/> | £70 |
| | Mesh Cover. | SPCGSPP-MBLK | <input type="text"/> | £44 |
| | Stealth Tek Cover. | SPCGSPP-TBLK | <input type="text"/> | £50 |
| | Urethane Liner. | URL1709070 | <input type="text"/> | £40 |
| | Polyurethane Liner. | PUL1709030 | <input type="text"/> | £48 |

| | Width | | | | | | | |
|----|-------|----|----|----|----|----|----|----|
| | 14 | 16 | 17 | 18 | 19 | 20 | 22 | 24 |
| 15 | | | | | | | | |
| 16 | | | | | | | | |
| 17 | | | | | | | | |
| 18 | | | | | | | | |
| 19 | | | | | | | | |
| 20 | | | | | | | | |
| 21 | | | | | | | | |
| 22 | | | | | | | | |
| 24 | | | | | | | | |
| 26 | | | | | | | | |

Select your size, then complete the Part Number by adding the Width first then Depth.

ADD SOLID BASE

ABS FLAT RIGIDISER

| | | | | | |
|--------------------------------------|-----|--------------------------------------|-----|--------------------------------------|-----|
| 14"W x 18"D. <i>DWR1010E1910I.</i> | £34 | 19"W x 17"D. <i>MEC1804152-1917.</i> | £34 | 24"W x 21"D. <i>DWR1010E236K.</i> | £34 |
| 16"W x 15"D. <i>MEC1804152-1615.</i> | £34 | 19"W x 18"D. <i>MEC1804152-1917.</i> | £34 | 24"W x 22"D. <i>DWR1010E236M.</i> | £34 |
| 16"W x 16"D. <i>DWR1010E191G.</i> | £34 | 19"W x 19"D. <i>MEC1804152-1919.</i> | £34 | 24"W x 24"D. <i>DWR1010E236O.</i> | £34 |
| 16"W x 17"D. <i>DWR1010E191G.</i> | £34 | 19"W x 20"D. <i>MEC1804152-1919.</i> | £34 | 24"W x 26"D. <i>DWR1010E236Q.</i> | £34 |
| 16"W x 18"D. <i>DWR1010E191I.</i> | £34 | 20"W x 16"D. <i>DWR1010E193G.</i> | £34 | 26"W x 18"D. <i>DWR1010E240I.</i> | £34 |
| 16"W x 20"D. <i>DWR1010E191I.</i> | £34 | 20"W x 18"D. <i>DWR1010E193I.</i> | £34 | 26"W x 20"D. <i>DWR1010E240K.</i> | £34 |
| 17"W x 16"D. <i>MEC1804152-1715.</i> | £34 | 20"W x 19"D. <i>DWR1010E193I.</i> | £34 | 26"W x 22"D. <i>DWR1010E240M.</i> | £34 |
| 17"W x 17"D. <i>DWR1010E234H.</i> | £34 | 20"W x 20"D. <i>DWR1010E193K.</i> | £34 | 26"W x 24"D. <i>DWR1010E240O.</i> | £34 |
| 17"W x 18"D. <i>DWR1010E234H.</i> | £34 | 20"W x 21"D. <i>DWR1010E193K.</i> | £34 | 28"W x 18"D. <i>DWR1010E244I.</i> | £34 |
| 17"W x 19"D. <i>DWR1010E234J.</i> | £34 | 20"W x 22"D. <i>DWR1010E193M.</i> | £34 | 28"W x 20"D. <i>DWR1010E244K.</i> | £34 |
| 17"W x 20"D. <i>DWR1010E234J.</i> | £34 | 20"W x 24"D. <i>DWR1010E19O.</i> | £34 | 28"W x 22"D. <i>DWR1010E244M.</i> | £34 |
| 18"W x 16"D. <i>DWR1010E192G.</i> | £34 | 22"W x 16"D. <i>DWR1010E194G.</i> | £34 | 28"W x 24"D. <i>DWR1010E244O.</i> | £34 |
| 18"W x 17"D. <i>DWR1010E192G.</i> | £34 | 22"W x 18"D. <i>DWR1010E194I.</i> | £34 | 28"W x 28"D. <i>DWR1010E244S.</i> | £34 |
| 18"W x 18"D. <i>DWR1010E192I.</i> | £34 | 22"W x 19"D. <i>DWR1010E194I.</i> | £34 | 30"W x 18"D. <i>DWR1010E341II.</i> | £34 |
| 18"W x 19"D. <i>DWR1010E192I.</i> | £34 | 22"W x 20"D. <i>DWR1010E194K.</i> | £34 | 30"W x 20"D. <i>DWR1010E341IK.</i> | £34 |
| 18"W x 20"D. <i>DWR1010E192K.</i> | £34 | 22"W x 22"D. <i>DWR1010E194M.</i> | £34 | 30"W x 22"D. <i>DWR1010E341IM.</i> | £34 |
| 18"W x 21"D. <i>DWR1010E192K.</i> | £34 | 22"W x 24"D. <i>DWR1010E194O.</i> | £34 | 30"W x 24"D. <i>DWR1010E341IO.</i> | £34 |
| 18"W x 22"D. <i>DWR1010E192M.</i> | £34 | 24"W x 16"D. <i>DWR1010E194G.</i> | £34 | 30"W x 28"D. <i>DWR1010E341IS.</i> | £34 |
| 18"W x 24"D. <i>DWR1010E192O.</i> | £34 | 24"W x 18"D. <i>DWR1010E236I.</i> | £34 | 32"W x 26"D. <i>MEC1804152-3226.</i> | £34 |
| 19"W x 16"D. <i>MEC1804152-1915.</i> | £34 | 24"W x 20"D. <i>DWR1010E236K.</i> | £34 | 32"W x 28"D. <i>MEC1804152-3228.</i> | £34 |

22" width and greater cushions include a flat rigidizer as standard.

ACM CONTOURED RIGIDISER SUITABLE FOR SLING SEATS

| | | | | | |
|--------------------------------|-----|--------------------------------|-----|--------------------------------|-----|
| 16"W x 16"D. <i>CUS142467.</i> | £56 | 18"W x 18"D. <i>CUS142466.</i> | £56 | 20"W x 19"D. <i>CUS142465.</i> | £56 |
| 16"W x 17"D. <i>CUS142467.</i> | £56 | 18"W x 19"D. <i>CUS142466.</i> | £56 | 20"W x 20"D. <i>CUS142465.</i> | £56 |
| 16"W x 18"D. <i>CUS142467.</i> | £56 | 18"W x 20"D. <i>CUS142466.</i> | £56 | | |
| 16"W x 16"D. <i>CUS142466.</i> | £56 | 20"W x 18"D. <i>CUS142465.</i> | £56 | | |

For use on manual wheelchairs only.


NOTES & SPECIAL REQUESTS

Special Order Spectrum Gel SPP Cushion.

£450

Please detail your special order cushion request. All special order cushion requests must be approved with the Quantum Custom Department. Please contact Quantum Inside Sales or Custom Department with questions. Extended lead time will apply. Cushion part numbers involved will be the cushion part number listed in previous sections modified to meet special order request stated in details below. Drawings are available for sketch details.

TOP VIEW



FRONT VIEW



RIGHT SIDE VIEW



PRICE

TOTAL

2



Stealth Products

32 Wedgwood Road, Bicester,

Oxon, OX26 4UL

Phone: 01869 324600

Fax: 01869 323070

Email: sales@quantumrehab.co.uk

Web: www.quantumrehab.co.uk

SPECTRUM GEL SPP - 15/01/24 V1



V1: 14/06/2022