

Light, Portable Easy to Fold



Pride
QUEST



With the Pride Quest Folding Scooter, the world is at your fingertips! With a unique folding design, this scooter folds in just three simple steps and fits easily in a closet, corner or vehicle trunk! Enjoy full lighting and speeds up to 4 mph.

PRIDE[®]
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST[®]

www.pride-mobility.co.uk

Specifications

Weight Capacity	18 st 1 lbs (115 kg)
Maximum Speed ¹	Up to 4 mph
Ground Clearance ³	60 mm (2.3") at motor
Turning Radius ³	820 mm (32.2")
Overall Length ³	930 mm (36.6")
Overall Width ³	485 mm (19.09")
Max. Climb Angle	7 degrees
Drive Tyres	180 mm (7")
Suspension	Indego Suspension
Total Weight Without Batteries ^{3, 4}	22 kg (48.5 lbs)
Heaviest Piece When Disassembled ³	N/A
Drivetrain	4-pole 270W
Dual Breaking System	-
Battery Requirements ³	2 x 11.5 Ah battery
Battery Weight ⁴	3.6 kg (7.5 lbs)
Battery Charger	2.5 Amp off-board
Colour	Red, Black, White

Features

- Height adjustable tiller for comfort
- Front and rear LED lighting
- Easy access levers for quick, secure folding
- Easy to transport
- Available in a variety of colours

Also available in:

- Red
- Black
- White



¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

² Includes standard seat weight.

³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

⁴ Battery weight may vary based on battery manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.