



## TECHNICAL DATA

<b>Max. user weight</b>	140kg	<b>Frame versions</b>	Swing-Away / Fixed Front
<b>Seat width</b>	340 - 540mm	<b>Frame angle</b>	75 / 85 (FF) & 70 / 80 / 90
<b>Seat depth</b>	250 - 500mm	<b>Front wheel</b>	4", 5", 6", 7"
<b>Front seat height</b>	360 - 540mm	<b>Rear wheel</b>	22", 24", 25", 26"
<b>Rear seat height</b>	350 - 520mm	<b>Wheelchair weight</b>	9.9kg
<b>Backrest height</b>	250 - 450mm	<b>Dismanteled weight</b>	7.3kg
<b>Backrest angle</b>	from -15° to +10°	<b>Mass of heaviest part</b>	2.2kg
<b>Total width</b>	SW +170 (0 Camber & 25" Wheels)	<b>Removable parts</b>	Rear wheels (for transportation)
<b>Total length</b>	1113mm	<b>Frame colours</b>	24 colours
<b>Total height</b>	1020mm	<b>Crash tested</b>	Yes (approved for transport in a vehicle according ISO 7176-19)
<b>Max. pivot width</b>	840mm	<b>Transport without user</b>	Fix the wheelchair firmly and securely in the vehicle (e.g. car or plane)
<b>Max. safe slope for use of parking brakes</b>	7°	<b>Intended use and environment</b>	For people unable to walk and with limited mobility; inside and outside. Operation by user or assistant.

Sunrise Medical cares about your safety. Please read your owner's manual for safe operating instructions and maintenance guidelines. In line with our policy of continual product improvement, Sunrise Medical reserves the right to make changes at any time, and without notice. For the latest information, please visit [sunrisemedical.co.uk](http://sunrisemedical.co.uk). If you are visually impaired, this document can be viewed in PDF format at [sunrisemedical.co.uk](http://sunrisemedical.co.uk)



## QS5 X FIXED SPEC GUIDE

To make the QS5 X as easy to work with as possible, we have created a number of specifications that offer a variety of benefits as detailed inside. All you need to do is choose your seat width, seat depth and the colour. All configurations are available in both fixed front & swing-away.



Use our online visualiser to create the QUICKIE QS5 X of your dreams at: [sunrisemedical.co.uk](http://sunrisemedical.co.uk)



Thorns Road, Brierley Hill, West Midlands, DY5 2LD  
t: +44 (0) 845 605 6688 | e: [marketing@sunmed.co.uk](mailto:marketing@sunmed.co.uk)

[sunrisemedical.co.uk](http://sunrisemedical.co.uk)

## QS5 X SPECIFICATIONS

To make the QS5 X as easy to work with as possible, we have created a number of specifications that offer a variety of benefits as detailed below. All you need to do is choose your seat width, seat depth and the colour. All configurations are available in both fixed front & swing-away.



### CORE - £1,850

An entry-level specification that offers all the benefits of the QS5 X in an affordable package.

#### Specification:

- 50mm Seat Cushion
- Fixed Backrest
- Medium Push Handles
- EXO Pro Upholstery
- Composite Removeable Sideguards
- Anti-Tip, Swingaway, Lightweight
- Composite Divided Footplates
- Aluminium Fork & Castors
- Lightweight Wheels
- Puncture Proof Tyres
- Aluminium Anodised Handrims
- Standard Brakes

#### Swing-Away Frame as per the above including:

- Single Post Sideguard & Armpad
- 80 Degree Alu Swing-Away Hanger



### LIGHTWEIGHT - £2,550

The lightest configuration in our range, for users where weight is key.

#### Specification:

- 50mm Seat Cushion
- Fixed Backrest
- Omit Push Handles
- EXO Pro Upholstery
- Carbon Fibre Fender Sideguards
- Composite Divided Footplates
- One Arm Fork & Aluminium Castors
- Proton Wheels
- Scwalbe One Tyres
- Aluminium Anodised Handrims
- Compact Brake EVO

#### Swing-Away Frame as per the above including:

- Tubular Armrest
- 80 Degree Alu Swing-Away Hanger



### PERFORMANCE - £2,550

For users who want the most rigid ride, focused on drive performance and useability.

#### Specification:

- 50mm Seat Cushion
- Angle Adjustable Backrest
- Fold Down Push Handles
- EXO Pro Upholstery
- Aluminium Fender Sideguards
- Auto Folding Fixation
- Aluminium Platform Footplates
- Aluminium Fork & Castors
- Lightweight Wheels
- Marathon Plus Tyres
- Ellipse 3R Handrims
- Compact Brake EVO

#### Swing-Away Frame as per the above including:

- 80 Degree Alu Swing-Away Hanger



### CARBON - £2,945

Our most stylish configuration, carbon sideguards and footplate for a premium look and finish.

#### Specification:

- 50mm Seat Cushion
- Angle Adjustable Backrest
- Fold Down Push Handles
- EXO Pro Upholstery
- Carbon Fibre Fender Sideguards
- Auto Folding Fixation
- Carbon Fibre Platform Footplates
- One Arm Fork & Aluminium Castors
- Lightweight Wheels
- Marathon Plus Tyres
- Ellipse 3R Handrims
- Compact Brake EVO

#### Swing-Away Frame as per the above including:

- 80 Degree Alu Swing-Away Hanger