

# QS5 X PERFORMANCE SPEC (SWING-AWAY)

£2,550

## Order Information

Account Code	
Purchase Order	

## Delivery Address

--

## Standard Specification

The following options are standard on all QS5 X Performance Swing-Away Specifications:

	QS5 X	QSC050001	Aluminium Hanger 80
QSC020014	Seat sling, nylon, 1 utility bag	QSC050050	Aluminium Divided Footplates
QSC020102	Seat cushion, foam (medium-soft), black, 50 mm	QSC080029 / 64	Aluminium Castor Fork
QSC030012	Angle Adjustable Backrest	QSC080306	Castor, solid (soft), rim: aluminium
QSC030203	Fold Down Push Handles	QSC070003	Rear wheel, Lightweight, 24 spokes (silver)
QSC030317	EXO PRO Backrest sling	QSC070104	Marathon Plus Tyres
QSC040103	Aluminium Fender	QSC070212	Ellipse 3R Handrims
QSC010015	Swing-Away Frame	QSC070317	Handrim mount, narrow (17 mm)
QSC010056	Auto folding fixation	QSC070602	2 ° camber
QSC090102	with Crash Test option	QSC060023	Compact Brake Evo

## Standard Dimensions

COG50	Centre of gravity 5 cm	24" Rear Wheel
	5" Castor Wheel, 115mm Castor Fork	40cm Rear Seat Height
	46cm Front Seat Height	40cm Lower Leg Length
	Backrest Height Range 40-45 cm	Backrest Height 40 cm
	90° Seat Angle	

## Configured Options - Please choose based on client

## Frame Dimensions

<b>Seat Width</b>	34	36	38	40	42	44	46	48	50	52	54
<b>Seat Depth</b>	40	42	44	46	48	50					

## Colour Options

<b>Frame Colour</b>	Brilliant Marine	Brilliant Purple	Brilliant Red	Marsred (Orange)	Mamba Green
	Matt Black	Signal White	Brilliant Silver	Blue	Pink
<b>Frame Stickers</b>	Orange	Silver	Red	Blue	Black
<b>Upholstery Colour</b>	Orange	Silver	Red	Blue	Black
<b>Lightweight Wheels</b>	Orange	Black	Red	Blue	
<b>Aluminium Fork</b>	Marsred	Silver	Red Ruby	White	Black
<b>Castor Wheel</b>	Orange	Silver	Red	Blue	Black

## Optional Upgrades

QSC110156	Twin LED Lights	£150	QSC070005	18 Spoke Spinergy	£245
-----------	-----------------	------	-----------	-------------------	------