

QS5 X PERFORMANCE SPEC (FIXED FRONT)

£2,550

Order Information		Delivery Address	
Account Code			
Purchase Order			

Standard Specification

The following options are standard on all QS5 X Performance Fixed Front Specifications:

	QS5 X	QSC050102	Platform Aluminium Footplate
QSC020014	Seat sling, nylon, 1 utility bag	QSC080029 / 64	Aluminium Castor Fork
QSC020102	Seat cushion, foam (medium-soft), black, 50 mm	QSC080306	Castor, solid (soft), rim: aluminium
QSC030012	Angle Adjustable Backrest	QSC070003	Rear wheel, Lightweight, 24 spokes (silver)
QSC030203	Fold Down Push Handles	QSC070104	Marathon Plus Tyres
QSC030317	EXO PRO Backrest sling	QSC070212	Ellipse 3R Handrims
QSC040103	Aluminium Fender	QSC070602	2 ° camber
QSC010016	Frame FF 85° (integrated footrest) & 2cm Inset	QSC060023	Compact Brake Evo
QSC010056	Auto folding fixation	QSC070317	Handrim mount, narrow (17 mm)
QSC090102	with Crash Test option	QSC010057	Touch Points

Standard Dimensions

COG50	Centre of gravity 5 cm	24" Rear Wheel
	5" Castor Wheel, 115mm Castor Fork	40cm Rear Seat Height
	46cm Front Seat Height	40cm Lower Leg Length
	Backrest Height Range 40-45 cm	Backrest Height 40 cm
	90 ° Seat Angle	

Configured Options - Please choose based on client

Frame Dimensions

Seat Width	34	36	38	40	42	44	46	48	50	52	54
Seat Depth	40	42	44	46	48	50					

Colour Options

Frame Colour	Brilliant Marine	Brilliant Purple	Brilliant Red	Marsred (Orange)	Mamba Green
	Matt Black	Signal White	Brilliant Silver	Blue	Pink
Frame Stickers	Orange	Silver	Red	Blue	Black
Upholstery Colour	Orange	Silver	Red	Blue	Black
Lightweight Wheels	Orange	Black	Red	Blue	
Aluminium Fork	Marsred	Silver	Red Ruby	White	Black
Castor Wheel	Orange	Silver	Red	Blue	Black

Optional Upgrades

QSC110156	Twin LED Lights	£150	QSC070005	18 Spoke Spinergy	£245
-----------	-----------------	------	-----------	-------------------	------