

QS5 X CORE SPEC (SWING-AWAY)

£1,950

Order Information

Account Code	
Purchase Order	

Delivery Address

--

Standard Specification

The following options are standard on all QS5 X Core Swing-Away Specifications:

QSC020014	Seat sling, nylon, 1 utility bag	QSC050001	Aluminium Hanger, 80°
QSC020102	Seat cushion, foam (medium-soft), black, 50 mm	QSC050020	Footboard, divided, composite, angle-adjustable
QSC030010	Backrest, fixed (90°)	QSC080029	Castor fork, Aluminium
QSC030201	Push handles, Medium	QSC080306	Castor, solid (soft), rim: aluminium
QSC030317	EXO PRO Backrest sling	QSC070003	Rear wheel, Lightweight, 24 spokes (silver)
QSC040052	Single-post sideguard, armpad	QSC070107	Tyre, pneumatic, slick, Schwalbe One
QSC010015	Swing-Away Frame	QSC070202	Handrims, aluminium, anodised, silver
QSC090036	Anti-tip tubes, aluminium, lightweight	QSC070317	Handrim mount, narrow (17 mm)
QSC010056	Auto folding fixation	QSC070602	2 ° camber
QSC090102	with Crash Test option	QSC060001	Standard brake

Standard Dimensions

COG50	Centre of gravity 5 cm	24" Rear Wheel
	5" Castor Wheel, 115mm Castor Fork	40cm Rear Seat Height
	46cm Front Seat Height	40cm Lower Leg Length
	Backrest Height Range 40-45 cm	Backrest Height 40 cm

Configured Options - Please choose based on client

Frame Dimensions

Seat Width	36	38	40	42	44	46	48	50	52	54
Seat Depth	40	42	44	46	48	50				

Colour Options

Frame Colour	Brilliant Marine	Brilliant Purple	Brilliant Red	Marsred (Orange)	Mamba Green
	Matt Black	Signal White	Brilliant Silver	Blue	Pink
Frame Stickers	Orange	Silver	Red	Blue	Black
Upholstery Colour	Orange	Silver	Red	Blue	Black
Lightweight Whels	Orange	Black	Red	Blue	
Aluminium Fork	Marsred	Silver	Red Ruby	White	Black
Castor Wheel	Orange	Silver	Red	Blue	Black

Optional Upgrades

QSC110156	Twin LED Lights	£150	QSC070005	18 Spoke Spinergy	£405
-----------	-----------------	------	-----------	-------------------	------