



QUICKIE
Live without limits.

SALSA R²
COMPACT AND
POWERFUL

SALSA R²

Compact and powerful



ENHANCES YOUR MOBILITY. INSIDE AND OUTSIDE.



The Salsa R² incorporates a great indoor and outdoor performance. Indoors the compact base and low seat to floor height (starting from 42cm) ensures easy access under tables. Outdoors the adjustable suspension and powerful motors (2-pole or 4-pole) always offer a comfortable and active driving characteristic up to a speed of 10kph. With simple access to the batteries from the rear they can be easily changed while the user is in the chair.

- > **Active swing-away hangers**
For easy transfer. Also available as elevating legrest (manual or powered)

Simple battery and controller access while user is in the chair

Batteries available with 60Ah or 70Ah

12.5" and 14" drive wheels

9" and 10" castor wheels



FULLY ADJUSTABLE SEATING SYSTEM

The Salsa R² is equipped with a fully adjustable seating system which can be quickly adjusted throughout the full range of width and depth, back angle and height without additional parts. The seat is designed for perfect adaptation to the user and enables prompt servicing and refurbishment. Alternative backs and cushions can be easily fitted.



Powered tilt module
(0° to 30°) with centre of
gravity adjustment to suit
correct users balance.



Configure the wheelchair to your individual needs - The Salsa R² offers a wide range of options. For further information on product options please check the Salsa R² order form.



// POWERED TILT MODULE WITH CENTRE OF GRAVITY ADJUSTMENT

Offers a tilt from 0° up to 30° with an additional weight shift mechanism to correct users balance.



// POWERED SINGLE PIVOT TILT OPTION

For basic tilt from 3° to 20° to quickly change between a driving and a relaxed seat position



// POWERED RECLINE BACK

This recline backrest offers a range from -5° to 30°. Easily retrofittable if developing needs requires higher support.



// POWERED ARTICULATING LEGREST

For advanced leg support. Also available as powered elevating legrest



// TRUNK SUPPORT LATERALS

Multi-adjustable laterals offer advanced positioning of the trunk.



// FULL RANGE OF CONTROL DEVICES AVAILABLE

- VR2 standard control (up to two actuators)
- Up to four actuator controls through a separate switch box
- Dual control
- R-net optional controller
- Full range of special controls



> **Easy transfer**

The smart flip-up armrest design enables easy transfer. It can be infinitely adjusted in height, depth and angle.

// FULL RANGE OF POWERED SEATING AND SPECIAL CONTROL OPTIONS

Change your seating position with three different types of tilt modules – the powered or mechanical single pivot tilt option for a basic tilt, the powered tilt module with integrated centre of gravity adjustment and the lift and tilt combi module.

The backrest can be optionally reclined with an actuator or a mechanical gas strut.

Finally the full range of special controls or input devices can be added to the chair to find the best operating performance for the user.



New base integrated LED lights and indicators (optional)

4 point tie-down system for transportation in a vehicle. Optional Dahl car docking system available.

Adjustable drive wheel suspension



> **Stylish and compact base**
Available in 3 colours - blue, red and white

// **SALSA R² IS APPROVED FOR TRANSPORT IN A VEHICLE**

The Salsa R² is transit tested according to ISO 7176-19 and therefore approved for transportation in a vehicle.

The tests are carried out with the 4 point tie-down system and the optional Dahl docking system for an even easier fixation.

The Salsa R² can be ordered in a configuration which is ready to fit the Dahl docking system.

DAHL
ENGINEERING

TECHNICAL DATA

SALSA R² // Powerchair

Maximum user weight:	Std: 140kg	Battery size:	60Ah and 70Ah
Castor wheel options:	9" and 10"	Maximum dynamic stability:	8° (14%)
Drive wheel options:	12.5" and 14"	Overall width:	60cm (12.5" wheels); 62cm (14" wheels)
Seat width:	41cm – 51cm	Turning radius:	Maximum 120cm
Seat depth:	41cm – 51cm	Maximum kerb climbing:	10cm (with kerb climber)
Maximum speed:	6kph and 10kph	Colour choices:	3 colours - red, blue and white
Maximum range:	30km (35km with 70Ah)	Transport (ISO 7176-19):	4 point tie-down and Dahl docking system

Find out more today!

Visit www.SunriseMedical.co.uk for further information on product options, or contact your local dealer or representative.



Get SOCIAL with QUICKIE



Keep up-to-date with the latest news on Quickie products through our Facebook page and YouTube channel:

www.facebook.com/QuickieWheelchairs

www.youtube.com/quickiewheelchairseu

Sunrise Medical Limited

Thorns Road, Brierley Hill, West Midlands, DY5 2LD

Tel: +44 (0) 845 605 66 88 Fax: +44 (0) 845 605 66 89

www.SunriseMedical.co.uk

