

# Rugged, Durable Made for Adventure



Feel the crisp, fresh air, smell the blossoms on the trees and enjoy the peaceful serenity of the outdoors. Your adventure awaits, and nothing stands in your way. Forge a new path with the Ranger®, an aggressive, outdoor scooter designed for extreme and rugged conditions. With powerful suspension, bright lighting and dual hydraulic brakes, this beast of a scooter is a powerful force.



## Specifications

Weight Capacity	25 st (159 kg)
Maximum Speed <sup>1</sup>	Up to 8 mph
Ground Clearance <sup>3</sup>	140 mm (5.5")
Turning Radius <sup>3</sup>	2407 mm (94.75")
Overall Length <sup>3</sup>	1556 mm (61.25")
Overall Width <sup>3</sup>	829 mm (32.5")
Max. Obstacle Climb	100 mm (3.9")
Drive Tyres	368 mm (14.5") - pneumatic
Suspension	Enhanced front and rear suspension
Total Weight Without Batteries <sup>3, 4</sup>	124.25 kg (273 lbs)
Heaviest Piece When Disassembled <sup>3</sup>	99.25 kg (219 lbs)
Drivetrain	Rear-wheel drive, dual inline, 24 volt DC motor
Dual Breaking System	Regenerative and electromechanical
Battery Requirements <sup>3</sup>	2 x 75 Ah batteries
Battery Weight <sup>4</sup>	21.75 kg (48 lbs)
Battery Charger	8 Amp off-board
Colour	Stone Grey, Desert Storm

## Features

- LCD Delta tiller display
- Full LED lights package
- Reversing sensor with on/off switch
- 14.5" chunky drive tyres
- USB charger located on the tiller
- Enhanced front and rear suspension

Available in:

- Desert Storm
- Stone Grey



<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

<sup>2</sup> Includes standard seat weight.

<sup>3</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

<sup>4</sup> Battery weight may vary based on battery manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.