

# Rascal Rhythm SL

Mid-range powerchairs



Unique & Reliable

Your new route to independence

**RASCAL** 

# Rhythm Seat Lift

## Rhythm SL Specification:

Rear Wheel Drive (RWD) and Front Wheel Drive (FWD) variations

Type Class RWD	B / 2
Type Class FWD	A / 1
Overall length	105.2cm / 41.4"
Overall width	61.6cm / 24.2"
Floor to seat height	57- 69.5cm (min-max)
Seat width	48.3cm / 19"
Seat depth	45.7cm / 18"
Maximum carrying capacity	159kgs / 25 stone
Battery voltage & capacity	12V / 50Ah x 2
Maximum speed RWD	6.4km/h / 4mph
Maximum speed FWD	5km/h / 3mph
Ground clearance	4cm / 1.6"
Maximum kerb height RWD	5cm / 2"
Range*	up to 25km / 15.6 miles
Maximum safe slope RWD	6°
Maximum safe slope FWD	3°
Turn-around width	1.18m / 3'9"
Colour	White

\*Range on full charge and flat ground

See full Tech Page on Website & Owner's Manual for details

## Rascal Rhythm SL Profile:

The Rascal Rhythm with Seat Lift (SL) offers graceful simplicity to the user including the flexibility of dual-directional drive with a rotating seat giving you the option of Front-Wheel Drive for superb manoeuvrability indoors or Rear-Wheel Drive for extra traction outdoors.

Experience a refined drive with smooth acceleration and turning as the Rhythm SL moves in harmony with you. Well defined driving positions, a swing-away joystick and a fully adjustable seat and headrest combine for proven all-day comfort.

Local Rascal Mobility Retailer:

## Main Features:

- 5 inch Column Seat Lift
- Front-Wheel or Rear-Wheel Drive
- Flip-up width adjustable armrests
- Convenient joystick charging port
- Swing-away joystick
- High back adjustable seat
- Fold-away footplates front and rear

OPTIONAL  
ACCESSORY

## Optional Extras:

- Option of swing-away legrests
- Option of elevating legrests
- Attendant Control



Specifications subject to change without prior notice  
Images are for illustrative purposes only and may differ  
from the actual product. E. and O.E - v 2 Final-Mar 2023