

# Vippi **LiFe**

Transportable scooters



**LiFetime Guarantee**

on Lithium Ferro-Phosphate Battery\*



The Lithium Ferro-Phosphate Battery is **'built to last'** & the **SAFEST** Lithium Battery in the Mobility sector



**RASCAL**



## Vippi Specification:

Type Class	B / 2
Overall length	100cm / 39.4"
Overall width	49.5cm / 19.5"
Floor to seat height	34-39cm (min-max)
Seat width	45cm / 17.7"
Seat depth	35cm / 13.7"
Maximum carrying capacity	130kgs / 20.4 stone
Battery voltage & capacity	25.6V / 5Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	3.5cm / 1.3"
Maximum kerb height	5cm / 2"
Range*	up to 20.3km / 13.8 miles
The maximum safe slope	6°
Turn-around width	1.35m / 4'4"
Colour	Blue Moon / Green Lightning / Red Sky

Single Heaviest component 15.4kgs / 33.9lbs

\*Range on full charge and flat ground

See full Tech Page on Website & Owner's Manual for details

## Rascal Vippi **LiFe** Profile:

Trendy, transportable & fuss free, the Vippi LiFe is lightweight, versatile & compact. Offering peace of mind with a LiFetime Guarantee on the battery, it's all about the State-of-The-Art Lithium Ferro-Phosphate Battery weighing just 4.4kg. Designed & engineered to make every day easy & spontaneous, whether you are off on an adventure or running an essential errand. White stitching on the soft padded seat provides a premium finish complemented by a stylish black tiller with a fantastic pop of colour. Re-discover your power with this super scooter.

It dismantles in a flash with an easy take-apart mechanism deconstructing into 5 parts for storage & transportation with the single heaviest part weighing only 15.4kgs.

Local Rascal Mobility Retailer:

† LiFetime Battery Guarantee T&Cs apply



## Main Features:

- British Made / Medical /Military grade **Lithium Ferro-Phosphate Battery**
- Car transportable
- Fold down backrest
- Flip-up width adjustable armrests
- Removable front basket
- Tiller adjustment
- Solid tyres all-round
- Swivel seat

## Optional:

- Rascal arm rest bag



## 10 Ah weighs just 4.4kg:

Relax & enjoy the ride in the knowledge that the **Lithium Ferro-Phosphate** (LiFePO4) Battery pack is the safest, **greenest** with an increased efficiency & performance to stand up to everyday life.



Specifications subject to change without prior notice  
Images are for illustrative purposes only  
E. and O.E. v2 March 2024