

# ROBOOTER

R E B O O T Y O U R L I F E

## ROBOOTER X40

TECHNOLOGY MEETS  
FUNCTIONAL ELEGANCE  
A TRUE ORIGINAL



DESIGN  
AWARD  
2022



[www.roobooter.com](http://www.roobooter.com)

# ROBOOTER X40

## BUILT FOR YOU

Ergonomic & Aesthetically Pleasing Design

- Independent shock suspension system
- Memory foam seating
- Auto-folding mechanism
- Smart control with app integration
- Monitoring with fault alerts



Speed Settings



Cruise Control



Remote Drive




Voice Settings



Battery Monitoring



Gyroscope for Stability



## DESIGNED TO GO

Everyday Functionality



Tight turning radius

Omni-wheel built for rough terrain

Auto-folds to fit in small places

Removable lithium battery

Built-in lighting for safety

## PRODUCT PARAMETERS

|   |                             |                                 |                           |
|---|-----------------------------|---------------------------------|---------------------------|
| Folding                                 | Electric Folding            | Battery capacity                | 20Ah(Optional 10ah)       |
| Maximum speed                           | 3.75mph                     | Charging Time                   | 4-6H                      |
| Maximum obstacle clearance              | 1.77"                       | Overall dimension               | (L) 40"x (W) 23"x (H) 36" |
| Maximum gap clearance                   | 4"                          | Folding shape dimension         | (L) 33"x (W) 23"x (H) 15" |
| Maximum climbing angle                  | 10°                         | Seat wide                       | 16.5"                     |
| Weight without battery                  | 92lbs.                      | Maximum range                   | 15miles                   |
| Battery number                          | 1                           | Side Arms                       | Angle Adjustable          |
| Battery weight (20ah)                   | 7.5lbs.                     | Maximum Load                    | 120kg                     |
| Specifications of front and rear wheels | Front: 10-in. / Rear 10-in. | Recommended Ambient Temperature | -10°C ~ 50°C              |
| Maximum capacity                        | ≤265lbs.                    | Waterproof Grade                | IPX4                      |