



Etac UK 0(121)561-2222 | orders.uk@etac.com |
www.kimobility.com | www.etac.com



Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Submit Form

Effective: 1/10/2024

Rogue XP £2796
Rogue Width Growable Frame Price includes one time free growth kit 200 lb weight capacity Maximum Seat Width of 18"

Frame

Seat Width

See reference A on page 8.

| | |
|----------------|-----|
| 10" [RXSWT010] | N/C |
| 11" [RXSWT011] | N/C |
| 12" [RXSWT012] | N/C |
| 13" [RXSWT013] | N/C |
| 14" [RXSWT014] | N/C |
| 15" [RXSWT015] | N/C |
| 16" [RXSWT016] | N/C |
| 17" [RXSWT017] | N/C |
| 18" [RXSWT018] | N/C |

Seat Depth

See reference B on page 8.

| | |
|----------------|------|
| 10" [RXSDE010] | N/C |
| 11" [RXSDE011] | N/C |
| 12" [RXSDE012] | N/C |
| 13" [RXSDE013] | N/C |
| 14" [RXSDE014] | N/C |
| 15" [RXSDE015] | N/C |
| 16" [RXSDE016] | N/C |
| 17" [RXSDE017] | N/C |
| 18" [RXSDE018] | N/C |
| 19" [RXSDE019] | N/C |
| 20" [RXSDE020] | £303 |

Transit Option

| | |
|-------------------------------|-----|
| Transit System [RXTRNOPT] | N/C |
| <i>Cannot be retrofitted.</i> | |
| Omit [RXTRNOP2] | N/C |

Seat Option

| | |
|-------------------------------|------|
| Standard Seat Sling [RXSS001] | N/C |
| Solid Seat Pan [RXACC004] | £260 |
| Omit Seat [RXCUS003] | N/C |

Sling Position

In Performance Position (0"), sling is at front of frame. The 1" or 2" sling position extends the frame in front of the sling by that amount. See reference C on page 8.

| | |
|---------------------------|-----|
| 0" (Performance) [RXSP01] | N/C |
| 1" [RXSP03] | N/C |
| 2" (Classic) [RXSP02] | N/C |

Adjustable Depth

Adds the capability to grow seat in depth. Seat depth plus Adjustable Depth cannot exceed 20.5".

| | |
|-----------------------------|-----|
| No Rear Adjustment [RXAD05] | N/C |
| 1/2" [RXAD06] | N/C |
| 1" [RXAD01] | N/C |
| 2" [RXAD02] | N/C |
| 3" [RXAD03] | N/C |

Front Frame Bend

Angle between front seat tube and ground. See reference D on page 8.

| | |
|--------------|------|
| 90° [RXFA01] | £145 |
| 85° [RXFA03] | £145 |
| 80° [RXFA02] | N/C |
| 75° [RXFA04] | £145 |
| 70° [RXFA05] | £145 |

Backrest

| Backrest | | | | | | |
|---|--|------------------|------------------|------------------|------------------|-----|
| Back Post | Height Range | Back Height | | | | |
| Height Adjustable (Standard or 1" Inward Offset) | Extra Short ^{1,2} [RXBOP024] | 9" | 10" | 11" | 12" | |
| | Short 2 [RXBOP025] | 11" | 12" | 13" | 14" ¹ | |
| | Medium ² [RXBOP026] | 13" | 14" | 15" | 16" | |
| | Tall 2 [RXBOP027] | 15" | 16" | 17" | 18" | |
| Fixed Height (Straight or 1" Inward Offset) | Fixed ³ [RXBOP011] | 9" ¹ | 10" ¹ | 11" ¹ | 12" | 13" |
| | | 14" | 15" | 16" | 17" | 18" |
| Fixed Height 2" Offset | Fixed Offset 2" ³ [RXBOP022] | 11" ¹ | 12" ¹ | 13" ¹ | 14" ¹ | |
| | | 15" | 16" | 17" | | |

¹ Not available with Fold Down Push Handle

² Each Offset Back Post is offset 1/2" creating 1" overall offset. Only Available with Standard or 1" Inward Offset.

³ Uses 1" tubes compatible with after-market backs. Not available with back upholstery. ² Inward or Outward Offset Back Posts are offset 1" creating 2" overall offset. ¹ Offset Back Post is offset 1/2" creating 1" overall offset

Push Handle

| | |
|--|------|
| Omit Push Handle [RXBCK008] | N/C |
| Standard Push Handle [RXBCK004] | N/C |
| Low Profile Push Handle [RXBCK007] | £130 |
| <i>2" shorter than standard push handles.</i> | |
| Fold Down Push Handle [RXBCK006] | £158 |
| <i>Not available on Height Adjustable Back Posts with 9"-11" Back Heights or Offset Back Posts with 11"-14" Back Heights</i> | |
| Bolt-On Push Handle [RXBOP005] | £170 |
| <i>Handles will be 1" - 2" above back height.</i> | |

Stroller Handle

Stroller Handles are removable and height adjustable.

| | |
|---|------|
| Ergo Stroller Handle [RXBCK009] | £212 |
| Single Stroller Handle [RXBOP007] | £164 |
| Double Stroller Handle [RXBOP006] | £212 |
| <i>Not available with 10" Seat Width.</i> | |

Backrest Release Cable Color

Select desired color. If no selection is made, black will be used.

| | |
|----------------------------------|------------------------------------|
| Black Backrest Cable [RXSBC008] | Red Backrest Cable [RXSBC005] |
| White Backrest Cable [RXSBC007] | Blue Backrest Cable [RXSBC006] |
| Yellow Backrest Cable [RXSBC004] | Orange Backrest Cable [RXSBC003] |
| Green Backrest Cable [RXSBC002] | Pink Backrest Cable [RXSBC001] N/C |

Seat Back Angle

Standard is 90° (perpendicular to floor). Minus (-°) closes back angle. Plus (+°) opens back angle.

| | | |
|---------------------|-----------------|----------------|
| Standard [RXSBA001] | -2° [RXSBA014] | +2° [RXSBA009] |
| -4° [RXSBA015] | +4° [RXSBA010] | -6° [RXSBA016] |
| +6° [RXSBA011] | -8° [RXSBA017] | +8° [RXSBA012] |
| -10° [RXSBA018] | +10° [RXSBA013] | N/C |

Back Upholstery Options

Available only with Height Adjustable Back Posts. If selecting both upholstery and Axiom Back, the upholstery will be installed.

| | |
|--|------|
| Padded Velcro Adjustable with Standard Cover [RXBUP001] | N/C |
| Padded Velcro Adjustable with Air Mesh Breathable Cover [RXBUP006] | N/C |
| <i>Not available on chair widths less than 12".</i> | |
| Tension Adjustable with Standard Cover [RXBUP002] | N/C |
| Tension Adjustable with Air Mesh Cover [RXBUP007] | £182 |
| <i>Not available on chair widths less than 12".</i> | |
| Omit [RXBUP003] | N/C |

Back Post Option

| | |
|--|-----|
| Straight - No Offset [RXOBO004] | N/C |
| 1" Inward Offset [RXOBO003] | N/C |
| 2" Inward Offset [RXOBO001] | N/C |
| <i>Not available with Height Adjustable Back Post.</i> | |
| 2" Outward Offset [RXOBO002] | N/C |
| <i>Not available with Height Adjustable Back Post or Swing Away Armrest.</i> | |

Front Seat Height and Casters

Caster Size / Footrest Taper / Front Seat Height

Select Front Seat Height and Footrest Taper in row corresponding to Caster Size. 6" casters not available with 85° and 90° Front Frame Bend.

Footrest Taper indicates bend of front frame creating footrest width. Taper is measured from outside of seat frame to outside of front tube. Taper is not available with 70° Front Frame Bend. See reference E on page 8.

Front Seat Height is measured from floor to top of seat tube at front of seat sling. Adjustable up or down by 1/2". See reference F on page 8.

| Footrest Taper | Caster Size | Front Seat Height | | | | | | | | | | | | | | |
|------------------------------|-------------|-------------------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| | | 13.5" | 14" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" | 20.5" |
| No Taper (Straight) [RXFT00] | 3" | 13.5" | 14" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" | 20.5" |
| | 4" | 14" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" | 20.5" | |
| | 5" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" | 20.5" | | |
| | 6" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" | 20.5" | | | |
| 1" Taper [RXFT01] | 3" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" | 20.5" | | | | | | |
| | 4" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" | 20.5" | | | | | | | |
| | 5" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" | 20.5" | | | | | | | | |
| | 6" | 18" | 18.5" | 19" | 19.5" | 20" | 20.5" | | | | | | | | | |

Casters

| Caster Type | Caster Description | Caster Size | | | | |
|-------------|--------------------------|-------------|------------|------------|------------|------------|
| | | Price | 3" | 4" | 5" | 6" |
| Caster Type | .75" Rollerblade | N/C | [RXCAS002] | [RXCAS004] | | |
| | .75" Lighted Rollerblade | £81 | [RXCAS001] | [RXCAS005] | | |
| | 1" Poly | N/C | | [RXCAS006] | [RXCAS010] | [RXCAS014] |
| | 1" Poly Aluminum | £96 | | [RXCAS007] | [RXCAS011] | |
| | 1" Pneumatic | N/C | | | | [RXCAS017] |
| | 1.5" Poly | N/C | | [RXCAS032] | [RXCAS025] | [RXCAS026] |
| | 1.5" Soft Roll Aluminum | £96 | [RXCAS003] | [RXCAS009] | [RXCAS013] | [RXCAS016] |

Fork Options

| | |
|---|------|
| Standard Fork [RXFRK001] | N/C |
| Single Sided Fork [RXFRK003] | £206 |
| Frog Legs II Suspension System [RXCST014] | £424 |

Lightweight shock absorbing fork. Not included in Colored Anodizing Package.

Footrest

Footplates

| | |
|--|------|
| Tubular "Open" Footplate [RXFPL004] | N/C |
| <i>Not available with footrest widths less than 9.5". See Reference G on page 8.</i> | |
| Tubular Footplate with ABS Cover [RXFPL005] | £91 |
| <i>Not available with footrest widths less than 9.5". See Reference G on page 8.</i> | |
| Aluminum Angle Adjustable [RXFPL006] | £182 |
| Aluminum Angle Adjustable with Riser [RXFPL007] | £182 |
| High Mount Aluminum Angle Adjustable [RXFPL008] | £182 |
| Angle Adjustable Flip Under [RXFPL012] | £224 |
| <i>Not available with footrest widths less than 9.5". See Reference G on page 8.</i> | |
| High Mount Angle Adjustable Flip Under [RXFPL011] | £272 |
| <i>Not available w/footrest widths greater than 15.5". See Reference G on 8.</i> | |

Footrest Protector

| | |
|--|-----|
| Footrest Protector [RXACC020] | N/C |
| <i>Not available with Angle Adjustable Flip Under, High Mount Angle Adjustable Flip Under or High Mount Aluminum Angle Adjustable Footplates</i> | |

Seat to Footrest Length

Measured from front edge of seat sling to top, rear of footrest. Footrest length at least 2.5" shorter than Front Seat Height recommended. See reference H on page 8.

| 8" | 8.5" | 9" | 9.5" | 10" | 10.5" |
|-----|-------|-----|-------|-----|-------|
| 11" | 11.5" | 12" | 12.5" | 13" | 13.5" |
| 14" | 14.5" | 15" | 15.5" | 16" | 16.5" |
| 17" | 17.5" | 18" | 18.5" | 19" | 19.5" |

Rear Seat Height and Wheels

Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size. For proper setup, selection is required when omitting rear wheels. Measured from floor to top of seat tube at front of back post. Rear Seat Height is custom to needs. See reference 1 on page 8.

| Wheel | Rear Seat Height | | | | | | | | | | | | | | | |
|-------|--------------------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-----|-------|-----|-------|-----|
| 18" | 11.5" ¹ | 12" | 12.5" | 13" | 13.5" | 14" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" |
| 20" | 12.5" ¹ | 13" | 13.5" | 14" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" |
| 22" | 13.5" ¹ | 14" | 14.5" | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" | | |
| 24" | 14.5" ¹ | 15" | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" | | | | |
| 25" | 15" ¹ | 15.5" | 16" | 16.5" | 17" | 17.5" | 18" | 18.5" | 19" | 19.5" | 20" | | | | | |

¹ Available only with Low Profile tires

One Arm Drive

One Arm Drive - Quick Release [RXDRV007] £846

Not available with 10" seat width. Only available with 20", 22" and 24" spoke wheels and Quick Release Axles. Spoke wheel style is exclusive to OAD.

Standard handrim will be Aluminum Anodized. Limitations exist based on rear wheel spacing, camber and armrest selections.

Left Side Drive [RXDRV002] N/C

Right Side Drive [RXDRV003] N/C

Aluminum Drive Handrim [RXDRV004] N/C

Poly Coated Drive Handrim [RXDRV005] £212

Wheels

No selection needed with One Arm Drive.

| | | Rear Wheel Size | | | | | |
|-----------------|--------------------------|-----------------|------------|------------|------------|------------|-------------------------|
| | | Price | 18" | 20" | 22" | 24" | 25" |
| Rear Wheel Type | Maxx Spoke | N/C | [RXWH1118] | [RXWH1120] | [RXWH1122] | [RXWH1124] | [RXWH1125] |
| | Mag | £121 | | [RXWH1220] | [RXWH1222] | [RXWH1224] | |
| | Colored Mag ² | £182 | | [RXWH1120] | [RXWH1322] | [RXWH1324] | |
| | Superlight Spoke | £121 | | | | [RXWHE224] | [RXWHE225] ¹ |
| | Spinergy Spox | £768 | | [RXWHE420] | [RXWHE422] | [RXWHE424] | [RXWHE425] |
| | Spinergy LX | £996 | | | [RXWHE522] | [RXWHE524] | [RXWHE525] |

¹ Available only with Superlight and Plastic Coated Handrims

² Mag color selection is on page 7

Spinergy Options

Select spoke color and handrim tab for Spinergy wheels.

Black Spokes [RXSSC001] N/C

Red Spokes [RXSSC002] White Spokes [RXSSC003]

Blue Spokes [RXSSC005] £91

Rear Tire Type

Full Profile

Street Pneumatic 1 3/8" [RXTIRO01] £182
Not available on 25" wheel size.

Full Poly 1 3/8" [RXTIRO02] N/C
Not available on Superlight wheel or 25" wheel size.

Pneumatic with Airless Insert 1 3/8" [RXTIRO04] £110
Not available on Superlight wheel or 25" wheel size.

Schwalbe HS-127 1 3/8" [RXTIRO17] £110
Available only on 24" wheel size. Not available on Mag wheel.

Mountain 2" Black [RXTIRO09] £151
Available only on 22", 24" and 25" MAXX Spoke and Spinergy wheels

Low Profile

Low Poly 1" [RXTIRO03] £110
Not available on Superlight wheel or 25" wheel size.

Sentinel Solid Tire 1" [RXTIRO05] £110
Not available on Superlight wheel, 18" wheel size or 20" and 22" Mag wheels.

High Pressure 1" [RXTIRO06] N/C
Not available on Mag wheel or 18" wheel size.

Schwalbe One 1" [RXTIRO18] £110
Available only on 24" and 25" wheel sizes. Not available on Mag wheel.

Schwalbe Marathon Plus Evolution 1" [RXTIRO10] £110
Available only on 24" and 25" wheel sizes. Not available on Mag wheel.



Etac UK 0(121)561-2222 | orders.uk@etac.com |
www.kimobility.com | www.etac.com



Handrims

No selection needed with One Arm Drive.

Standard Handrims

| | |
|--|------|
| Aluminum Anodized [RXHRM001] | N/C |
| Superlight Aluminum Anodized [RXHRM010] | N/C |
| <i>Available only on 24" and 25" wheel sizes. Available only on MAXX Spoke, Superlight or Spinergy wheels.</i> | |
| Plastic Coated [RXHRM002] | £121 |
| Projection (8 Vertical Projections) [RXHRM003] | £333 |
| <i>Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels</i> | |
| Omit Handrim [RXHRM012] | N/C |

Ergonomic Handrims

Not available with 20" or 26" wheels. Tab mount only available with 24" wheels.

| | |
|---------------------------------------|------|
| Ergo Lite Regular Grip [RXHRM005] | £272 |
| Ergo Standard Regular Grip [RXHRM008] | £272 |
| Ergo Lite Extra Grip [RXHRM006] | £272 |
| Ergo Standard Extra Grip [RXHRM009] | £272 |
| Ergo Lite Omit Grip [RXHRM004] | £272 |
| Ergo Standard Omit Grip [RXHRM007] | £272 |

Handrim Option

Handrim Mount

If no selection is made, Standard Mount will be used

| | |
|------------------------|-----|
| Close Mount [RXTAB001] | N/C |
|------------------------|-----|

Micro-Grip

| | |
|--|-----|
| 3M Micro Grip [RXACC006] | £91 |
| <i>Not available on Plastic Coated or Projection Handrims.</i> | |

Axles

No selection needed with Drum Brake.

| | |
|---|------|
| Quick Release [RXAXO003] | N/C |
| Quad Quick Release [RXAXO001] | £182 |
| <i>Not available with Mag or Spinergy wheels.</i> | |

Rear Wheel Camber

The following apply if One Arm Drive is selected: 6° and 8° not available. 4° not available with T-Arms. 0° required with the combination of a 5th Wheel and Seat Widths less than 16".

| | | | | |
|-------------|-------------|-------------|-------------|-----|
| 0° [RXWC00] | 2° [RXWC02] | 4° [RXWC04] | 6° [RXWC06] | |
| 8° [RXWC08] | | | | N/C |

Center of Gravity presets

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post. 2" minimum Center of Gravity with 5th Wheel and One Arm Drive. See reference J on page 8.

| Amputee Positions | | | | |
|--------------------|------------------|------------------|------------------|------------------|
| -1.5" [RXCGA150] | -1" [RXCGA100] | -5" [RXCGA050] | 0" [RXCGS000] | .25" [RXCGS025] |
| Standard Positions | | | | |
| .5" [RXCGS050] | .75" [RXCGS075] | 1" [RXCGS100] | 1.25" [RXCGS125] | 1.5" [RXCGS150] |
| 1.75" [RXCGS175] | 2" [RXCGS200] | 2.25" [RXCGS225] | 2.5" [RXCGS250] | 2.75" [RXCGS275] |
| 3" [RXCGS300] | 3.25" [RXCGS325] | 3.5" [RXCGS350] | 3.75" [RXCGS375] | 4" [RXCGS400] |
| 4.25" [RXCGS425] | | | | |

Rear Wheel Spacing

1" minimum with 0° camber. 1.5" minimum with T-Arms. The following apply if One Arm Drive is selected: 1.75" minimum with T-Arms. 1.25 minimum with select Rear Seat Heights. No selection needed with Wheelchair Drive Ready. See Reference K on page 8.

| | | | |
|-----------------|--------------|-----------------|----------------|
| 0.75" [RXWS075] | 1" [RXWS100] | 1.25" [RXWS125] | 1.5" [RXWS150] |
| 1.75" [RXWS175] | | | N/C |

Wheel Locks

Wheel Locks

| | |
|--|------|
| Push to Lock [RXWLK002] | N/C |
| Push to Lock Flush Mount [RXWLK001] | N/C |
| Pull to Lock [RXWLK003] | N/C |
| Omit [RXWLK008] | N/C |
| Under Seat Scissor Lock [RXWLK007] | £37 |
| Low Profile Scissor Lock [RXWLK012] | £121 |
| <i>Not available with Low Poly, Full Poly or Sentinel Solid tires.</i> | |
| Short Thro Scissor Lock [RXWLK006] | £175 |
| Drum Brake [RXWLK015] | £532 |
| <i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Low Profile or Fold Down Push Handles, or Ergo and Single Stroller Handles without integrated push handles. Axle is Quick Release.</i> | |

| | |
|--|------|
| Dual - Drum Brake with Push to Lock [RXWLK016] | £653 |
| <i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Low Profile or Fold Down Push Handles, or Ergo and Single Stroller Handles without integrated push handles. Axle is Quick Release.</i> | |

| | |
|--|------|
| Dual - Drum Brake with Pull to Lock [RXWLK017] | £653 |
| <i>Available only with 20", 22" and 24" MAXX Spoke wheels. Not available with Low Profile or Fold Down Push Handles, or Ergo and Single Stroller Handles without integrated push handles. Axle is Quick Release.</i> | |

Wheel Lock Options

| | |
|--|------|
| Extension Handles [RXWLK004] | £49 |
| <i>Available only with Push to Lock or Pull to Lock Wheel Locks.</i> | |
| Grade Aids [RXWLK009] | £145 |
| <i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i> | |

Armrests

| Armrest Type | Length | | Price |
|--|--------|------------|-------|
| Height Adjustable T-Arm Tall (11.5"-15") | Desk | [RXARM039] | £212 |
| | Full | [RXARM040] | |
| Height Adjustable T-Arm (8.5"-13.5") ¹ | Desk | [RXARM003] | £212 |
| | Full | [RXARM004] | |
| Height Adjustable T-Arm Low (6.5"-10.5") ² | Desk | [RXARM022] | £212 |
| | Full | [RXARM023] | |
| Angle Adjustable Locking Extendable Flip Up ³ | Desk | [RXARM010] | £266 |
| | Full | [RXARM011] | |

¹ Not available with Seat Depths less than 13". 1.5" minimum Rear Wheel Spacing

² Not available with Seat Depths less than 12". 1.5" minimum Rear Wheel Spacing. If requested configuration is unavailable, Standard Height T-Arm will be used.

³ Requires Fixed Height Backrest. Not available with backrest heights less than 12"

⁴ Not available with Arm Pads or Outward Offset Back Posts. If selection is not available with rear wheel and rear seat height configuration, then Tall Swing Away armrest will be used.

⁵ Not available with Arm Pads or Outward Offset Back Posts.

⁶ Includes Composite Side Guard with Direct Mount Receiver. Not available with Arm Pads or Outward Offset Back Posts. If selection is not available with rear wheel and rear seat height configuration, then Tall Swing Away with Side Guard armrest will be used.

⁷ Includes Composite Side Guard with Direct Mount Receiver. Not available with Arm Pads or Outward Offset Back Posts.

| Armrest Type | | Price |
|---|------------|-------|
| Short Swing Away (7" - 9") ⁴ | [RXARM001] | £255 |
| Tall Swing Away (9.25" - 11.25") ⁵ | [RXARM021] | £255 |
| Short Swing Away with Side Guards (7"-9") ⁶ | [RXARM027] | £255 |
| Tall Swing Away with Side Gueards (9.25" - 11.25") ⁷ | [RXARM028] | £255 |
| Omit Armrests | [RXARM007] | N/C |

Arm Pads

Not available with Swing Away Armrests.

Standard Arm Pad [RXARM026]

Foam Grip [RXARM025]

Available only with Angle Adjustable Locking Extendable Flip Up Armrest. N/C

Waterfall Arm Pad [RXARM020]

Not available with Angle Adjustable Locking Extendable Flip Up Armrest. £30

Side Guards

Direct Mount Side Guards

Not available with T-Arms or Swing Away with Side Guard Armrests.

Composite Side Guards - Direct Mount Receiver [RXSDG002] N/C
^{12" minimum Seat Depth.}

Aluminum Straight - Direct Mount Receiver - Black Matte [RXSDG024] £212
^{12" minimum Seat Depth.}

Aluminum Straight - Direct Mount Receiver - Match Frame Color [RXSDG025] £272
^{12" minimum Seat Depth. Not available with Hydrographics}

Aluminum Adjustable Fender w/Foam Pad - Direct Mount Receiver [RXSDG014] £296
^{Not available with Adjustable Depth greater than 1/2".}

Aluminum Adjustable Fender w/Foam Pad - Direct Mount Receiver [RXSDG015] £357
^{Not available with Adjustable Depth greater than 1/2". Not available with Hydrographics.}

Carbon Fiber Adjustable Fender - Direct Mount Receiver [RXSDG013] £393
^{Not available with Adjustable Depth greater than 1/2".}

Clamp-On Mount Side Guards

Not available with T-Arms or Swing Away with Side Guard Armrests.

Composite Side Guards - Clamp-On Receiver [RXSDG004] £199
^{12" minimum Seat Depth.}

Composite Side Guards - Pediatric - Clamp-On Receiver [RXSDG007] £199
Aluminum Straight - Clamp On Receiver - Black Matte [RXSDG020] £212
^{12" minimum Seat Depth.}

Aluminum Straight - Clamp-On Receiver - Match Frame Color [RXSDG023] £272
^{12" minimum Seat Depth. Not available with Hydrographics.}

Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver - Black [RXSDG017] £296

Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver - [RXSDG018] £357
^{Not available with Hydrographics.}

Carbon Fiber Adjustable Fender - Clamp-On Receiver [RXSDG016] £393

Accessories

5th Wheel

Rotating removable rear caster. Not available with 25" rear wheel.

| | |
|--|------|
| Standard 5th Wheel [RXATP002] | £206 |
| Dynamic 5th Wheel [RXATP004] | £206 |
| Soft Elastomer [RXELA001] Firm Elastomer [RXELA002] | N/C |

Both elastomers are included; the selected elastomer will be installed.

Anti-Tip

| | |
|---------------------------|-----|
| Rear Anti-Tips [RXATP001] | £89 |
| Omit [RXATP000] | N/C |

Pelvic Positioning Belts

| | |
|--|------|
| 1" Auto Buckle Pelvic Positioning Belt [RXPOB004] | £55 |
| 1" Padded Auto Buckle Pelvic Positioning Belt [RXPOB005] | £151 |
| 1 1/2" Auto Buckle Pelvic Positioning Belt [RXPOB001] | £55 |
| 1 1/2" Padded Auto Buckle Pelvic Positioning Belt [RXPOB002] | £151 |
| Small Bodypoint Evoflex Pelvic Stabilizer (16-24") [RXPOB007] | £188 |
| Medium Bodypoint Evoflex Pelvic Stabilizer (19-29") [RXPOB008] | £188 |
| Large Bodypoint Evoflex Pelvic Stabilizer (25-37") [RXPOB009] | £188 |

Belt Mounting Kit

| | |
|--|------|
| Ki Mobility Belt Mounting Clamp [RXACC023] | £37 |
| Bodypoint Belt Mounting Bracket Kit [RXACC013] | £103 |

Calf Strap

| | |
|--------------------------------------|------|
| Velcro Adjustable [RXACC002] | £37 |
| Padded Velcro Adjustable [RXACC003] | £67 |
| Bodypoint Padded AeroMesh [RXACC010] | £103 |

Spoke Guards

| | |
|--|------|
| Clear Spoke Guards [RXSPG001] | £148 |
| <i>Not available with Mag wheel.</i> | |
| Black Spoke Guards [RXSPG002] | £148 |
| <i>Not available with 18" or 20" wheel sizes or Mag wheel.</i> | |
| Blue Monster Spoke [RXSPG003] ¹ | £148 |
| Robo Fun Spoke Guards [RXSPG004] ¹ | £148 |
| Tutu Cute Spoke Guards [RXSPG005] ¹ | £148 |
| Tuxedo Time Spoke Guards [RXSPG008] ¹ | £148 |

1 Available only with 20" and 22" wheel sizes. Not available with Mag wheel.

Other Accessories

| | |
|---|------|
| Tool Kit [RXACC005] | £20 |
| Removable Underseat Pouch [RXACC015] | £30 |
| Pediatric Ki Mobility Backpack [RXACC026] | £100 |
| Ki Mobility Backpack [RXACC009] | £100 |
| Luggage Carrier [RXACC025] | £73 |
| Black Neoprene Impact Guards [RXIMG001] | £86 |
| Under Chair Carrier and Removable Bag [RXACC001] | £30 |
| <i>Not available with Seat Widths less than 14" or scissor locks.</i> | |
| Cane and Crutch Holder [RXACC007] | £91 |

Color

Frame Color

*Frame : Select desired frame and accent colors. Accent color is used on backrest rigidizer bar; if no accent color is selected, frame color will be used.
Hydrographics: Available on frame only; please choose accent color to left. Not available with Transit Option.*

| Options | Price | Colors | | | | | | | | | | | | | | | | | | | |
|--------------------------------------|-------|--------|--------------|----------------|------------|--------------|-----------|--------|--------------|-----------------|--------|-------|--------|---------------|---------------|---------------------|-----------------|----------------|-----------|------------|--------------|
| | | Gloss | | | | | | | | | | Matte | | | | | | | | | |
| | | Black | Blue Skies | Bubblegum Pink | Candy Blue | Candy Purple | Candy Red | Copper | Granny Smith | Metallic Orange | Silver | Teal | Yellow | Charcoal Gray | Emerald Green | Iridescent Burgundy | Iridescent Navy | Military Green | Retro Red | Royal Blue | Shadow Black |
| Frame | N/C | | | | | | | | | | | | | | | | | | | | |
| Accent | N/C | | | | | | | | | | | | | | | | | | | | |
| Mag Color | | | | | | | | | | | | | | | | | | | | | |
| Color Anodizing Package ¹ | Black | Blue | Bright Green | Orange | Pewter | Pink | Purple | Red | Gold | | | | | | | | | | | | |
| | N/C | £164 | £164 | £164 | £164 | £164 | £164 | £164 | £164 | | | | | | | | | | | | |

1 Select desired color. If no selection is made, black will be used. Included in Colored Anodizing Package: Backrest plates, camber tube (not available with Wheelchair Drive Ready selection), towers, CG clamps (not available with negative Center of Gravity or Wheelchair Drive Ready selections), side guard receivers, RP tuning housing, front caster forks (standard or Single Sided) and Dynamic 5th Wheel split clamp.

References

A. Seat Width See Diagram 1

Measured from outside of frame tube on one side to outside of frame tube on other side.

B. Seat Depth See Diagram 3

Measured from front of back posts to front edge of seat sling. Seat sling starts at beginning of bend at front of frame.

C. Sling Position See Diagram 2

In Performance Position (0°), sling is at front of frame to keep chair shorter and more maneuverable. The 1" or 2" seat sling position extends the frame in front of the sling by that amount. Adding more frame can improve stability and provide support to aid in transfers.

D. Front Frame Bend See Diagram 3

Angle between front seat tube and ground.

E. Footrest Taper See Diagram 1

Indicates bend of front frame creating footrest. Taper is measured from outside of seat frame to outside of front tube.

F. Front Seat Height See Diagram 1

Measured from floor to top of seat tube at front of seat sling.

G. Footrest Width See Diagram 1

Measured from inside of front tube to inside of front tube on other side. Footrest width is listed below for each seat width and taper.

| Seat Width | Straight | 1" Taper | 2" Taper |
|---------------------|-----------------------------|-----------------------------|-----------------------------|
| Outside Measurement | 2" Narrower than Seat Width | 4" Narrower than Seat Width | 6" Narrower than Seat Width |
| 14" | 11.5" | 9.5" | |
| 15" | 12.5" | 10.5" | |
| 16" | 13.5" | 11.5" | 9.5" |
| 17" | 14.5" | 12.5" | 10.5" |
| 18" | 15.5" | 13.5" | 11.5" |
| 19" | 16.5" | 14.5" | 12.5" |
| 20" | 17.5" | 15.5" | 13.5" |
| 21" | 18.5" | 16.5" | 14.5" |
| 22" | 19.5" | 17.5" | 15.5" |

H. Seat to Footrest Length See Diagram 3

Measured from front edge of seat sling to top rear of footrest. Footrest length of at least 2.5" shorter than front seat height recommended.

I. Rear Seat Height See Diagram 3

Measured from floor to top of seat tube at front of back post. Rear seat height is custom to needs.

J. Center of Gravity Preset See Diagram 3

Measured from front of back post to center of rear axle. 0° indicates axle will be directly under back post.

K. Rear Wheel Spacing See Diagram 1

Measured from outside of seat back to inside of rear tire. Adjustable out .5" from setting.

Diagram 1




Diagram 2





Diagram 3



Axiom G

| | | Width | | | | | | | | | | | | | | | |
|--------|-----|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|
| | | 14" | | 15" | | 16" | | 17" | | 18" | | 19" | | 20" | | 22" | |
| | | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover |
| Retail | | £85 | £70 | £85 | £70 | £85 | £70 | £85 | £70 | £85 | £70 | £85 | £70 | £85 | £70 | £85 | £70 |
| Height | 14" | | | | | | | | | | | | | | | | |
| | 15" | | | | | | | | | | | | | | | | |
| | 16" | | | | | | | | | | | | | | | | |
| | 17" | | | | | | | | | | | | | | | | |
| | 18" | | | | | | | | | | | | | | | | |
| | 19" | | | | | | | | | | | | | | | | |
| | 20" | | | | | | | | | | | | | | | | |
| | 21" | | | | | | | | | | | | | | | | |
| | 22" | | | | | | | | | | | | | | | | |

Axiom S

| | | Width | | | | | | | | | | | | | | | |
|--------|-----|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|
| | | 14" | | 15" | | 16" | | 17" | | 18" | | 19" | | 20" | | 22" | |
| | | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover |
| Retail | | £120 | £70 | £120 | £70 | £120 | £70 | £120 | £70 | £120 | £70 | £120 | £70 | £120 | £70 | £120 | £70 |
| Height | 14" | | | | | | | | | | | | | | | | |
| | 15" | | | | | | | | | | | | | | | | |
| | 16" | | | | | | | | | | | | | | | | |
| | 17" | | | | | | | | | | | | | | | | |
| | 18" | | | | | | | | | | | | | | | | |
| | 19" | | | | | | | | | | | | | | | | |
| | 20" | | | | | | | | | | | | | | | | |
| | 21" | | | | | | | | | | | | | | | | |
| | 22" | | | | | | | | | | | | | | | | |

Axiom P

| | | Width | | | | | | | | | | | | | | | |
|--------|-----|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|---------|-------------|
| | | 14" | | 15" | | 16" | | 17" | | 18" | | 19" | | 20" | | 22" | |
| | | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover | Cushion | Outer Cover |
| Retail | | £170 | £70 | £170 | £70 | £170 | £70 | £170 | £70 | £170 | £70 | £170 | £70 | £170 | £70 | £170 | £70 |
| Height | 14" | | | | | | | | | | | | | | | | |
| | 15" | | | | | | | | | | | | | | | | |
| | 16" | | | | | | | | | | | | | | | | |
| | 17" | | | | | | | | | | | | | | | | |
| | 18" | | | | | | | | | | | | | | | | |
| | 19" | | | | | | | | | | | | | | | | |
| | 20" | | | | | | | | | | | | | | | | |
| | 22" | | | | | | | | | | | | | | | | |

Axiom SP - Visco

| | | Width | | | | | | | | | | | | | | | | | | | | | | | | | | |
|--------|-----|---------|-------|------|---------|-------|-----|---------|-------|------|---------|-------|-----|---------|-------|------|---------|-------|-----|---------|-------|------|---------|-------|--|-----|--|--|
| | | 14" | | | 15" | | | 16" | | | 17" | | | 18" | | | 19" | | | 20" | | | 21" | | | 22" | | |
| | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | | | |
| Retail | | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | £350 | £70 | | | | | |
| Height | 14" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 15" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 16" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 17" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 18" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 19" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 20" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 22" | | | | | | | | | | | | | | | | | | | | | | | | | | | |

Axiom SP - Fluid

| | | Width | | | | | | | | | | | | | | | | | | | | | | | | | | |
|--------|-----|---------|-------|------|---------|-------|-----|---------|-------|------|---------|-------|-----|---------|-------|------|---------|-------|-----|---------|-------|------|---------|-------|--|-----|--|--|
| | | 14" | | | 15" | | | 16" | | | 17" | | | 18" | | | 19" | | | 20" | | | 21" | | | 22" | | |
| | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | Cushion | Cover | | | | |
| Retail | | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | £320 | £70 | | | | | |
| Height | 14" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 15" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 16" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 17" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 18" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 19" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 20" | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | 22" | | | | | | | | | | | | | | | | | | | | | | | | | | | |


Solid Seat Insert

Weight capacity: 350 lbs

| | | Width | | | | | | | | |
|--------|-----|---------|---------|---------|---------|---------|---------|---------|---------|---------|
| | | 14" | 15" | 16" | 17" | 18" | 19" | 20" | 21" | 22" |
| | | Cushion | Cushion | Cushion | Cushion | Cushion | Cushion | Cushion | Cushion | Cushion |
| Retail | | £55 | £55 | £55 | £55 | £55 | £55 | £55 | £55 | £55 |
| Height | 14" | | | | | | | | | |
| | 15" | | | | | | | | | |
| | 16" | | | | | | | | | |
| | 17" | | | | | | | | | |
| | 18" | | | | | | | | | |
| | 19" | | | | | | | | | |
| | 20" | | | | | | | | | |
| | 21" | | | | | | | | | |
| | 22" | | | | | | | | | |

| Axiom Back Width | Fits Back Frame Widths | |
|------------------|-----------------------------|-----------------------------|
| | Inner Slots | Outer Slots |
| 12" (30.5cm) | 11" (28.0cm) | 12" -14" (30.5cm -35.5cm) |
| 14" (35.5cm) | 12" -13" (30.5cm -33.0cm) | 14" -16" (35.5cm -40.5cm) |
| 16" (40.5cm) | 14" -15" (35.5cm -38.0cm) | 16" -18" (40.5cm -45.5cm) |
| 18" (45.5cm) | 16" - 17" (40.5cm - 43.0cm) | 18" - 20" (45.5cm - 50.5cm) |
| 20" (51.0cm) | 18" -19" (46.0cm -48.5cm) | 20" -22" (50.5cm -55.5cm) |
| 22" (56.0cm) | 20" -21" (51.0cm -53.5cm) | 22" -24" (56.0cm -61.0cm) |

| Hardware Description | Code |
|--|------------|
| Improved DEX Hardware, removable adjustable hardware (300 lb max user weight) | DEX |
| New BTX Hardware, removable adjustable hardware not recommended for, or installed on Ultralight Rigid Frame Wheelchairs (400 lb max user weight) | BTX |
| New Fixed Hardware, non-removeable hardware up to (300 lb max user weight) | FXT |



Axiom Back

Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports. Weight capacity: 300 lbs.

Axiom Posterior Back (PB)

| | | Width | | | | | | |
|--------|-----|-------|------------|------------|------------|------------|------------|------------|
| | | 12" | 14" | 16" | 18" | 20" | 22" | |
| Retail | | £550 | £550 | £550 | £550 | £550 | £550 | |
| Height | Low | 8" | [AXPB1208] | [AXPB1408] | [AXPB1608] | [AXPB1808] | [AXPB2008] | |
| | | 10" | | [AXPB1410] | [AXPB1610] | [AXPB1810] | [AXPB2010] | [AXPB2210] |
| | | 12" | [AXPB1212] | [AXPB1412] | [AXPB1612] | [AXPB1812] | [AXPB2012] | |
| | Mid | 14" | | [AXPB1414] | [AXPB1614] | [AXPB1814] | [AXPB2014] | [AXPB2214] |
| | | 16" | [AXPB1216] | [AXPB1416] | [AXPB1616] | [AXPB1816] | [AXPB2016] | |
| | | 18" | | [AXPB1418] | [AXPB1618] | [AXPB1818] | [AXPB2018] | [AXPB2218] |

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Lateral Support for AXPB Backs

| Lateral Trunk Supports for AXPB Backs | | | | | | | |
|---------------------------------------|------|-------|---|-----------|------|------|--|
| Pair | Left | Right | Description | Item Code | Pair | Each | |
| | | | 3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back | PBLF3x4P | £295 | £148 | |
| | | | 4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Back | PBLF4x5P | £295 | £148 | |
| | | | 3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back | PBLS3x4P | £628 | £314 | |
| | | | 4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Back | PBLS4x5P | £628 | £314 | |
| | | | 3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware | PBLAS3x4P | £628 | £314 | |
| | | | 4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware | PBLAS4x5P | £628 | £314 | |

Axiom Posterior Lateral Back (PL)

| | | Width | | | | | | |
|--------|-----|------------|------------|------------|------------|------------|------------|------------|
| | | 12" | 14" | 16" | 18" | 20" | 22" | |
| Retail | | £550 | £550 | £550 | £550 | £550 | £550 | |
| Height | Low | 8" | [AXPL1208] | [AXPL1408] | [AXPL1608] | [AXPL1808] | [AXPL2008] | |
| | | 10" | | [AXPL1410] | [AXPL1610] | [AXPL1810] | [AXPL2010] | [AXPL2210] |
| | | 12" | [AXPL1212] | [AXPL1412] | [AXPL1612] | [AXPL1812] | [AXPL2012] | |
| | Mid | 14" | | [AXPL1414] | [AXPL1614] | [AXPL1814] | [AXPL2014] | [AXPL2214] |
| | | 16" | [AXPL1216] | [AXPL1416] | [AXPL1616] | [AXPL1816] | [AXPL2016] | |
| | | 18" | | [AXPL1418] | [AXPL1618] | [AXPL1818] | [AXPL2018] | [AXPL2218] |
| Tall | 20" | [AXPL1220] | [AXPL1420] | [AXPL1620] | [AXPL1820] | [AXPL2020] | | |
| | 22" | | [AXPL1422] | [AXPL1622] | [AXPL1822] | [AXPL2022] | [AXPL2222] | |

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Lateral Support for AXPL

| Lateral Trunk Supports for AXPL & AXLB Backs | | | | | | |
|--|------|-------|---|-----------|------|------|
| Pair | Left | Right | Description | Item Code | Pair | Each |
| | | | 3 x 4 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back | PLLF3x4P | £295 | £295 |
| | | | 4 x 5 Lateral Pads with Fixed Mounting Hardware for Axiom Posterior Lateral Back | PLLF4x5P | £295 | £295 |
| | | | 3 x 4 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back | PLLS3x4P | £628 | £314 |
| | | | 4 x 5 Lateral Pads with Swing Away Mounting Hardware for Axiom Posterior Lateral Back | PLLS4x5P | £628 | £314 |
| | | | 3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware | PLLAS3x4P | £628 | £314 |
| | | | 4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware | PLLAS4x5P | £628 | £314 |

Axiom Deep Lateral Back (DL)

| | | Width | | | | | | |
|--------|------|-------|------------|------------|------------|------------|------------|------------|
| | | 12" | 14" | 16" | 18" | 20" | 22" | |
| Retail | | £550 | £550 | £550 | £550 | £550 | £550 | |
| HCPC | | E2620 | E2620 | E2620 | E2620 | E2620 | E2621 | |
| Height | Mid | 14" | [AXDL1214] | [AXDL1414] | [AXDL1614] | [AXDL1814] | [AXDL2014] | [AXDL2214] |
| | | 16" | | [AXDL1416] | [AXDL1616] | [AXDL1816] | [AXDL2016] | |
| | | 18" | [AXDL1218] | [AXDL1418] | [AXDL1618] | [AXDL1818] | [AXDL2018] | [AXDL2218] |
| | Tall | 20" | [AXDL1220] | [AXDL1420] | [AXDL1620] | [AXDL1820] | [AXDL2020] | |
| | | 22" | | [AXDL1422] | [AXDL1622] | [AXDL1822] | [AXDL2022] | [AXDL2222] |

Back Hardware

DEX Hardware [DEX] BTX Hardware [BTX] FXT Hardware [FXT] _____ N/C

Headrest

| Description | Std. Straight Removable Adjustable Hardware | | Multi-Axis Removable Adjustable Hardware | | Posalinc Removable Adjustable Hardware | | HCPCS |
|-------------------------------------|---|--------|--|--------|--|--------|---------------|
| | Item Code | Retail | Item Code | Retail | Item Code | Retail | |
| Comfort Soft Headrest 8" | [CS08SS] | £268 | [CS08MA] | £268 | [CS08PL] | £268 | E0955 / E1028 |
| Comfort Soft Headrest 10" | [CS10SS] | £268 | [CS10MA] | £268 | [CS10PL] | £268 | E0955 / E1028 |
| Comfort Soft Headrest 14" | [CS14SS] | £268 | [CS14MA] | £268 | [CS14PL] | £268 | E0955 / E1028 |
| Contour Positioning Pad, 8" X 6.25" | [CP08SS] | £303 | [CP08MA] | £303 | [CP08PL] | £303 | E0955 / E1028 |
| Contour Positioning Pad, 9" X 6.5" | [CP09SS] | £303 | [CP09MA] | £303 | [CP09PL] | £303 | E0955 / E1028 |
| Contour Positioning Pad, 10" X 7.5" | [CP10SS] | £303 | [CP10MA] | £303 | [CP10PL] | £303 | E0955 / E1028 |
| Tri-soft Headrest 10" | [TS10SS] | £386 | [TS10MA] | £386 | [TS10PL] | £386 | E0955 / E1028 |
| Tri-soft Headrest 14" | [TS14SS] | £386 | [TS14MA] | £386 | [TS14PL] | £386 | E0955 / E1028 |

| Description | Item Code | Retail | HCPCS |
|-----------------------------------|-----------|--------|-------|
| Universal Headrest Mounting Plate | [UMP] | £43 | K0108 |

Not required for Axiom Headrests

Anterior Thoracic Supports

Selection of a modification from the Axiom Back Modifications section required in order to mount posture supports. Not available with Low Backs.

| Anterior Thoracic Supports | | | |
|--|-----------|-------|--------|
| Description | Item Code | HCPCS | Retail |
| Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small | [14941] | E0960 | £100 |
| Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium | [14942] | E0960 | £100 |
| Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large | [14943] | E0960 | £100 |
| Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 20" Width, Small | [14924] | E0960 | £100 |
| Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium | [14927] | E0960 | £100 |
| Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large | [14928] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small | [14851] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Standard, Dynamic, Small | [14852] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Standard, Dynamic, Medium | [14853] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Standard, Dynamic, Large | [14854] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small | [14856] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small | [14857] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium | [14858] | E0960 | £100 |
| Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large | [14859] | E0960 | £100 |
| AirLogic Anterior Trunk Support: Standard, Structure, Small | [14282] | E0960 | £300 |
| AirLogic Anterior Trunk Support: Standard, Structure, Medium | [14283] | E0960 | £300 |
| AirLogic Anterior Trunk Support: Standard, Structure, Large | [14284] | E0960 | £300 |
| AirLogic Anterior Trunk Support: Standard, Structure, X-Large | [14285] | E0960 | £300 |
| AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small | [14257] | E0960 | £300 |
| AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium | [14258] | E0960 | £300 |
| AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Large | [14259] | E0960 | £300 |
| AirLogic Anterior Trunk Support: Slim Cut, Dynamic, X-Large | [14260] | E0960 | £300 |

Accessories

| | |
|--|-----|
| Adjustable Shoulder Guides (K0108) [AXBSG] <i>Includes Anterior Thoracic Mounting Modification. Not available with Low Backs.</i> | £43 |
| Extra Lumbar Pad (K0108) [LUMPAD] <i>Available only with Axiom Lateral Back and Axiom Deep Lateral back</i> | N/C |
| Privacy Flap (K0108) [PRFLAP] | £44 |

Axiom Back Modifications

Axiom Back Modifications

Extended lead times will apply for all modifications. Makes back non-refundable.

| | |
|---|-----|
| Anterior Thoracic Mounting Modification [AXSMM] <i>Provides holes for Anterior Thoracic Supports and Shoulder Guides</i> | N/C |
|---|-----|

Gloss Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

| | |
|-------------------------------------|-----|
| Black Back Shell [AXMBLK] | N/C |
| Blue Skies Back Shell [AXMSKY] | N/C |
| Bubblegum Pink Back Shell [AXMPNK] | N/C |
| Candy Blue Back Shell [AXMCBL] | N/C |
| Candy Purple Back Shell [AXMPUR] | N/C |
| Candy Red Back Shell [AXMRED] | N/C |
| Copper Back Shell [AXMCOP] | N/C |
| Granny Smith Back Shell [AXMGSM] | N/C |
| Metallic Orange Back Shell [AXMMOR] | N/C |
| Yellow Back Shell [AXMYLW] | N/C |
| Silver Back Shell [AXMSIL] | N/C |
| Teal Back Shell [AXMTEA] | N/C |

Matte Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

| | |
|---|-----|
| Military Green Back Shell [AXMAGM] | N/C |
| Charcoal Gray Back Shell [AXMCGM] | N/C |
| Emerald Green Back Shell [AXMEGM] | N/C |
| Iridescent Burgundy Back Shell [AXMIBM] | N/C |
| Iridescent Navy Back Shell [AXMINM] | N/C |
| Retro Red Back Shell [AXMRRM] | N/C |
| Royal Blue Back Shell [AXMRBM] | N/C |
| Smooth White Back Shell [AXMSWM] | N/C |