

SANGO xxi

Version: 4
Valid from: 01-01-2021
Valid until: New version or latest 31-12-2021
Price list no.: 1000102

SANGO xxi with SEGO xxi £ 9.795,-

Includes the following options:

- Front Wheel Drive base
- 6 km/h
- Castor wheels 9", air
- Drive Wheels 14", air
- Suspension on drive wheels
- 78Ah batteries
- 12A charger
- R-net 120A controller
- R-net LED joystick module
- Fixed joystick bracket, right hand
- Seat height set at 45 cm
- Anatomic formed seat 60 x 55 cm (w x d)
- Contour standard backrest 60 x 50 cm (w x h)
- Swing away footrests with calf strap
- Detachable armrests
- PU armpads
- Anti-tippers



This picture may deviate from the standard configuration

Type: Order Quotation

Name dealer: _____

Reference no.: _____

Delivery address: _____

Remarks: _____

Specifications

Maximum user weight 250 kg
Transport brackets ISO7176-19
CE-marked

Gross prices in GBP excluding VAT.
Subject to change without prior notice.

Standard = STD

DRIVE BASE			
1000102	<input type="checkbox"/>	SANGO xxi with SEGO xxi	9.795
6000001	<input type="checkbox"/>	Front Wheel Drive (FWD)	STD
Maximum speed			
6000294	<input type="checkbox"/>	6 km/h HT motor	STD
Frame colour			
6000011	<input type="checkbox"/>	Matt blue metallic	STD
6000013	<input type="checkbox"/>	Different colour, please specify RAL colour	357
Batteries & charger			
6000016	<input type="checkbox"/>	Batteries 78Ah	STD
6000295	<input type="checkbox"/>	Batteries 100Ah	254
6000018	<input type="checkbox"/>	Charger 24V - 12A	STD

WHEELS			
Wheel size			
6000020	<input type="checkbox"/>	14" Drive wheels with 9" castor wheels	STD
Wheel type			
6000024	<input type="checkbox"/>	Castor wheels, air	STD
6000026	<input type="checkbox"/>	Castor wheels, solid	0
6000028	<input type="checkbox"/>	Drive wheels, air	STD
6000029	<input type="checkbox"/>	Drive wheels, solid	0
Tyres			
6000030	<input type="checkbox"/>	Black	STD
6000031	<input type="checkbox"/>	Grey	0
Mudguards			
6000032	<input type="checkbox"/>	Mudguards castor wheels 9"	68
6000033	<input type="checkbox"/>	Mudguard drive wheels	88

LIGHTING & POWERED OPTIONS			
Lighting			
6000040	<input type="checkbox"/>	LED lights & indicators	303
Seat adjustments			
6000034	<input type="checkbox"/>	Powered backrest recline 30°	(1 actuator) 587
6000263	<input type="checkbox"/>	Powered seat tilt in space 25°	(1 actuator) 899
6000297	<input type="checkbox"/>	Powered seat lift HD 20cm and tilt in space 25°, MAX 250 kg (Special = extended delivery time)	(2 actuators) 2.884
6000288	<input type="checkbox"/>	Powered "Catapult" seat, 25°, please specify seat depth	(1 actuator) 1.711
Legrest adjustments (incl. calf supports or PU calf plates) (for manual legrests and legrest accessories, see Manual Legrests)			
6000038	<input type="checkbox"/>	Powered centre mount legrest, built-in length compensation	(angle adj footplate) (1 actuator) 1.061
6000039	<input type="checkbox"/>	Powered centre mount Transfer legrest, built-in length compensation	(restrictions for RWD) (footplate can sink to the floor) (2 actuators) 1.321
6000422	<input type="checkbox"/>	Powered centre mount Transfer legrest, independent length and angle adjustment	" " 1.321
6000415	<input type="checkbox"/>	Footplate foam layer	(5 mm high density foam/ cellular rubber) (center mount footplate only) 34
6000037	<input type="checkbox"/>	Powered articulating legrests	(6°-70°, angle adjustable footplates and PU calf plates) (2 actuators) 791
6000135	<input type="checkbox"/>	Calf plates upholstered instead of PU	55
6000136	<input type="checkbox"/>	Kneepad, right	46
6000137	<input type="checkbox"/>	Kneepad, left	46

CONTROLS			
R-net system			
Joystick modules (all underneath LED and LCD joystick modules can operate 6 actuators and lights)			
6000052	<input type="checkbox"/>	R-net LED	STD
6000053	<input type="checkbox"/>	R-net LCD-1	(CJSM1) (2 jack sockets) 102
6000055	<input type="checkbox"/>	R-net LCD-2 Infra-red	(CJSM2) (2 jack sockets) 204
6000057	<input type="checkbox"/>	R-net LCD-2 Infra-red and Bluetooth	(CJSM2) (2 jack sockets) 357
6000056	<input type="checkbox"/>	HD kit for R-net LCD-2	807
Mounting of the Joystick module			
6000041	<input type="checkbox"/>	Fixed bracket, right	STD
6000042	<input type="checkbox"/>	Fixed bracket, left	0
6000043	<input type="checkbox"/>	Swing-away bracket, right	59
6000044	<input type="checkbox"/>	Swing-away bracket, left	59
Power modules			
6000061	<input type="checkbox"/>	120A Power module	STD
6000062	<input type="checkbox"/>	DAC box, direct actuator control. Button	(only with R-net LCD joystick modules) 347
6000362	<input type="checkbox"/>	DAC box, direct actuator control. Toggle	(only with R-net LCD joystick modules) 347
6000063	<input type="checkbox"/>	Gyro module	296
6000065	<input type="checkbox"/>	R-net i-Device Bluetooth module	(for Apple IOS) 250
Attendant control			
6000067	<input type="checkbox"/>	R-net + attendant control, right	(dual) 332
6000068	<input type="checkbox"/>	R-net + attendant control, left	(dual) 332
6000230	<input type="checkbox"/>	R-net joystick module as attendant control, right	61
6000231	<input type="checkbox"/>	R-net joystick module as attendant control, left	61
Omni module (Special controls, see end of order form)			
6000070	<input type="checkbox"/>	Mounted on fixed bracket, right	867
6000071	<input type="checkbox"/>	Mounted on fixed bracket, left	867
6000072	<input type="checkbox"/>	Mounted on Swing-away bracket, right	928
6000073	<input type="checkbox"/>	Mounted on Swing-away bracket, left	928
Push button switches			
6000081	<input type="checkbox"/>	Small Ø 17mm , red	86
6000082	<input type="checkbox"/>	Small Ø 17mm, green	86
6000083	<input type="checkbox"/>	Medium Ø 30mm, red	108
6000084	<input type="checkbox"/>	Medium Ø 30mm, green	108
6000388	<input type="checkbox"/>	Big Ø 50mm, red	108
6000389	<input type="checkbox"/>	Big Ø 50mm, green	108
Other R-net options			
Joystick knobs			
6000074	<input type="checkbox"/>	Mushroom knob	48
6000075	<input type="checkbox"/>	Ball mini, Ø32 mm	48
6000416	<input type="checkbox"/>	Golf ball, Ø43 mm	48
6000387	<input type="checkbox"/>	Foam ball, Ø45 mm	48
6000076	<input type="checkbox"/>	Soft foam ball, Ø70 mm	48
6000077	<input type="checkbox"/>	Stick joystick	48
6000078	<input type="checkbox"/>	T-grip	56
6000443	<input type="checkbox"/>	U-Grip vertical, small	(6 cm wide) 71
6000079	<input type="checkbox"/>	U-Grip vertical, large	(10,5 cm wide) 71
6000444	<input type="checkbox"/>	U-Grip horizontal, small	(8 cm wide) 71
6000080	<input type="checkbox"/>	U-Grip horizontal, large	(10 cm wide) 71

6000255	<input type="checkbox"/>	Single USB charger, 3A	100
6000392	<input type="checkbox"/>	External charging socket	153
6000393	<input type="checkbox"/>	System locking plug	47 (anti-theft)

SEATS			
Seat height (floor to top seat frame in cm and 0° tilt)			
6000395	<input type="checkbox"/>	45,5 cm	STD
6000396	<input type="checkbox"/>	47,5 cm	
6000397	<input type="checkbox"/>	49,5 cm	
Seat cushion size			
6000204	<input type="checkbox"/>	Width 60 cm	STD
6000205	<input type="checkbox"/>	Width 64 cm	0
6000191	<input type="checkbox"/>	Depth 48 cm	0
6000192	<input type="checkbox"/>	Depth 55 cm	STD
6000091	<input type="checkbox"/>	Different seat cushion size, specify desired size	204
Seat cushion type			
6000092	<input type="checkbox"/>	Contour	finish: vinyl STD
6000094	<input type="checkbox"/>	Contour, partially strapped	air mesh 92
6000096	<input type="checkbox"/>	Anatomic Formed	air mesh (partially strapped) 163
6000098	<input type="checkbox"/>	Anatomic Formed Visco	air mesh (partially memory foam and webbed) 189
Seat cushion replacements			
6000101	<input type="checkbox"/>	Seat board	finish: vinyl 0
6000103	<input type="checkbox"/>	Seat board, partially strapped	vinyl 91
6000104	<input type="checkbox"/>	Threshold, front	vinyl (upholstered, H=5cm) 33
6000105	<input type="checkbox"/>	Threshold, rear	vinyl (upholstered, H=5cm) 33
6000106	<input type="checkbox"/>	Acrylic threshold, rear	acryl (depth adjustable, H=5cm) 38
Alternative finish			
6000107	<input type="checkbox"/>	Air mesh instead of vinyl	71
6000108	<input type="checkbox"/>	Stretchrubber instead of air mesh	39
Seat accessories			
6000110	<input type="checkbox"/>	Thigh support upholsterd small, set	(not with catapult seat) (22 x 11 cm LxH) 347
6000111	<input type="checkbox"/>	Thigh support upholstered large, set	(not with catapult seat) (36 x 13cm LxH) 367
6000232	<input type="checkbox"/>	Acrylic thigh support instead of upholstered	0

BACKRESTS			
Backrest cushion size			
6000206	<input type="checkbox"/>	Width 60 cm	STD
6000207	<input type="checkbox"/>	Width 64 cm	0
6000200	<input type="checkbox"/>	Height 50 cm	STD
6000201	<input type="checkbox"/>	Height 57 cm	0
6000115	<input type="checkbox"/>	Different backrest cushion size, specify desired size	213
Backrest cushion type			
6000116	<input type="checkbox"/>	Contour	finish: vinyl STD
6000117	<input type="checkbox"/>	Deep contour	air mesh (width adj. built-in laterals) 105
6000118	<input type="checkbox"/>	Deep contour vertical	air mesh (vertical elastic straps and width adj. built-in laterals) 166
6000119	<input type="checkbox"/>	Deep contour horizontal	air mesh (horizontal tension adj. straps and width adj. built-in laterals) 186
6000120	<input type="checkbox"/>	Deep contour Flex	air mesh (2 wooden strips, horiz.tension adj. straps and width adj. built-in laterals) 216
6000121	<input type="checkbox"/>	Extra Deep contour	air mesh (horizontal tension adj. straps and width adj. built-in laterals) 230
6000122	<input type="checkbox"/>	Extra Deep contour Flex,	air mesh (2 wooden strips, horiz.tension adj. straps and width adj. built-in laterals) 260
Alternative finish			
6000123	<input type="checkbox"/>	Air mesh instead of vinyl	81
6000124	<input type="checkbox"/>	Stretchrubber instead of air mesh	44
Backrest accessories			
6000126	<input type="checkbox"/>	Thoracic support, adjustable, right	(Stealth) 281
6000127	<input type="checkbox"/>	Thoracic support, adjustable, left	(Stealth) 281
6000129	<input type="checkbox"/>	Folding backrest	133

MANUAL LEGRESTS			
Centre mount legrest (for powered legrests see Powered options and lighting)			
6000130	<input type="checkbox"/>	Centre mount legrest, angle and height adjustable	0
Centre mount legrest accessories			
6000131	<input type="checkbox"/>	Calf plates	138
Swing-away legrests			
6000289	<input type="checkbox"/>	Standard footrests, calf strap	(angle adjustable foot plates) STD
6000133	<input type="checkbox"/>	Standard footrests, depth adjustable PU calf plates	(angle adjustable foot plates) 138
6000134	<input type="checkbox"/>	Angle adjustable legrests, depth adjustable PU calf plates	(6°-70° gas spring operated and angle adjustable foot plates) 265
Swing-away legrest accessories			
6000135	<input type="checkbox"/>	Calf plates upholstered instead of PU	55
6000136	<input type="checkbox"/>	Kneepad, right	46
6000137	<input type="checkbox"/>	Kneepad, left	46
6000138	<input type="checkbox"/>	Leather heel support, right	41
6000139	<input type="checkbox"/>	Leather heel support, left	41
6000142	<input type="checkbox"/>	Amputation support, right	122
6000143	<input type="checkbox"/>	Amputation support, left	122

ARMRESTS			
Armrests			
6000144	<input type="checkbox"/>	Armrest detachable	STD
Armpads			
6000146	<input type="checkbox"/>	PU 34 cm, right	STD
6000147	<input type="checkbox"/>	PU 34 cm, left	STD
6000148	<input type="checkbox"/>	Upholstered 30 cm, right	56
6000149	<input type="checkbox"/>	Upholstered 30 cm, left	56
6000150	<input type="checkbox"/>	Upholstered 40 cm, right	61
6000151	<input type="checkbox"/>	Upholstered 40 cm, left	61
6000152	<input type="checkbox"/>	P-shape upholstered 40 cm, right	77
6000153	<input type="checkbox"/>	P-shape upholstered 40 cm, left	77
6000353	<input type="checkbox"/>	Air mesh instead of vinyl	54

HEADRESTS			
Headrest metal bracket			
6000154	<input type="checkbox"/>	Adjustable in height and depth	97
6000155	<input type="checkbox"/>	Adjustable in height, depth and width	117
Headrest cushion			
6000157	<input type="checkbox"/>	PU small	87
6000158	<input type="checkbox"/>	Padded butterfly shape	97
6000159	<input type="checkbox"/>	Padded flexible small	112
6000160	<input type="checkbox"/>	Padded flexible large	122
6000354	<input type="checkbox"/>	Air Mesh instead of vinyl	49
		Finish:	
		PU	
		Vinyl	
		Vinyl	
		Vinyl	

ACCESSORIES			
6000170	<input type="checkbox"/>	Lapbelt long, 2.90 m	STD
6000441	<input type="checkbox"/>	Armrest bag, right	(30 x 20 x 5 cm) 22
6000442	<input type="checkbox"/>	Armrest bag, left	" " 22
6000173	<input type="checkbox"/>	Cane/crutch holder, right	(for 2 canes/crutches) 56
6000352	<input type="checkbox"/>	Cane/crutch holder, left	(for 2 canes/crutches) 56
6000176	<input type="checkbox"/>	Bag/backpack bracket (2 hooks)	(mounted on headrest bracket) 97

REMARKS:

Special controls

Different delivery time / different discount
Additional on chair, as spare part see spare part list

Prices in GBP excluding VAT.
Subject to change without prior notice.

Standard = STD

SPECIAL CONTROLS (R-net only)			
OMNI Special input devices		(operating force)	
8000003	<input type="checkbox"/> Micro joystick	(8,5 gram)	2.357
8000004	<input type="checkbox"/> Multi joystick	(50 gram)	2.130
8000005	<input type="checkbox"/> All-round Light joystick	(120 gram)	1.694
8000006	<input type="checkbox"/> All-round Light joystick , chin	(incl. 2 satellite switches and foam ball) (120 gram)	2.058
8000007	<input type="checkbox"/> All-round joystick	(250 gram)	1.065
8000008	<input type="checkbox"/> All-round joystick , chin	(incl. 2 satellite switches and foam ball) (250 gram)	1.427
8000031	<input type="checkbox"/> All-round heavy duty joystick	(650 gram)	2.308
R-net (direct plug-in) Special input devices		(operating force)	
8000032	<input type="checkbox"/> Micro joystick	(8,5 gram)	2.475
8000033	<input type="checkbox"/> Multi joystick	(50 gram)	2.248
8000034	<input type="checkbox"/> All-round Light joystick	(120 gram)	1.810
8000035	<input type="checkbox"/> All-round Light joystick , chin	(incl. 2 satellite switches and foam ball) (120 gram)	2.176
8000036	<input type="checkbox"/> All-round joystick	(250 gram)	1.186
8000037	<input type="checkbox"/> All-round joystick , chin	(incl. 2 satellite switches and foam ball) (250 gram)	1.549
8000038	<input type="checkbox"/> All-round heavy duty joystick	(650 gram)	2.427
Mounting of special input device			
8000009	<input type="checkbox"/> On armrest, right		84
8000010	<input type="checkbox"/> On armrest, left		84
8000011	<input type="checkbox"/> In /on tray , specify dimensions and required mounting position in a detailed drawing		735
Swing arms for chin control:			
8000012	<input type="checkbox"/> Manual horizontal, right		1.337
8000030	<input type="checkbox"/> Manual horizontal, left		1.337
8000013	<input type="checkbox"/> Powered horizontal, right		2.487
8000014	<input type="checkbox"/> Powered horizontal, left		2.487
8000015	<input type="checkbox"/> Powered vertical , right		2.487
8000016	<input type="checkbox"/> Powered vertical, left		2.487
8000039	<input type="checkbox"/> Switch Ø 17mm to operate swingarm, right	(mounted on headrest bracket)	427
8000040	<input type="checkbox"/> Switch Ø 17mm to operate swingarm, left	(mounted on headrest bracket)	427
Switches			
<input type="checkbox"/>	Please specify required type, location and way of mounting of the switches		POR
REMARKS:			