



WHERE YOU WANT TO BE.

Sapphire[®] 2





Sapphire[®] 2

Take it apart, take it away, take it further!





Detachable

To get your Sapphire 2 ready for transporting, simply lift off the seat, lift out the batteries, and lift the rear handle to separate the front and rear of the scooter. It really is that simple! By the way, it goes back together just as smoothly!



The Sterling Lock system is cable and connector free so there are no fiddly bits to deal with.

Transportable

Now the Sapphire 2 can be loaded into your vehicle. Luckily, as one of the lightest scooters in its class, this is no problem with the Sterling Lift dual handles fitted to the batteries and rear drive unit. These make lifting and handling the parts safer and more comfortable.



With its large batteries carrying you in style and comfort, the Sapphire 2 is the scooter you'll want to go further on!

The tiller can be safely secured in position for transporting.

Comfortable

Now you're ready! Stretch out your legs and take advantage of the Sterling Comfort package: three wheel legroom with four wheel stability, infinitely adjustable tiller, wrap-around handle bars and finger-thumb directional control all fitted as standard.



Sapphire[®] 2




Adj.
42 - 64cm [16.5 - 25"]



Adj.
42 - 55cm
(16.5 - 21.7")




43cm
(17")



52cm
(20")

96cm
(37.8")



121cm (47.5")



59cm (23")

Technical Data Sheet

| | | | |
|-------------------------------|------------------|-------------------------------------|---|
| Version | 4 Wheels | Drive unit weight | 19.5kg (43lbs) |
| Speed | 4mph | Battery weight | unit 12kg (26.4lbs) |
| Batteries | 35Ah | Total weight when folded/dismantled | 22.5kg (49.5lbs) |
| Charger off board | 5Ah | Turning diameter | 143cm (57.2") |
| Maximum range | 30km (19 miles) | Maximum safe slope | 14° - 113kg (17 stone) user weight / 10° - 150kg (23.8 stone) user weight |
| Colour | Carbon Metallic | Ground clearance | 9.5cm (3.7") |
| Maximum user weight | 150kg (23 stone) | Tyres | Pneumatic |
| Total weight | 79kg (174lbs) | Front wheel size | 22.5cm (9") |
| Heaviest part (front section) | 22.5kg (49.6lbs) | Rear wheel size | 26.5cm (10.5") |
| Seat weight | 12.5kg (27.5lbs) | Stylish lights | Yes |

Find out more today!

Visit www.sterlingscooters.co.uk for more information, to find your local dealer or to request a demonstration.

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.co.uk or alternatively is available on request in large text.



Sunrise Medical Limited
Thorns Road
Brierley Hill
West Midlands
DY5 2LD
Tel: +44 (0) 845 605 66 88
Fax: +44 (0) 845 605 66 89
www.SunriseMedical.co.uk